

TECHNE

QUANTICA[®]

Real-time nucleic acid
detection system



real multiplex
capability...

real chemistry
versatility...

The *real* journey begins...



Barloworld Scientific

TECHNE

QUANTICA®

Real-time nucleic acid detection system



The Thermal System:

- 96-well low profile microplate format with heat sealed optical film
- Temperature controlled heated lid between 10°C above ambient and 110°C
- Ramp rates up to 2.6°C/sec
- Temperature range 4°C – 99°C
- Robot friendly CD-type block loading mechanism
- Block uniformity of better than $\pm 0.25^\circ\text{C}$


Quantica® Features:

- Solid state white light source and photon counting photomultiplier tube (PMT) detector giving an impressive excitation range of 470nm to 650nm and a detection range of 500nm to 710nm.
- Multiplex - multiple wavelengths are detectable per sample using up to 4 paired excitation and emission filters housed in individual cartridge systems.
- Wide dynamic range of at least 9 orders of magnitude from starting copy number and high sensitivity detecting down to 1nM Fluorescein and single copy templates.

The Heated Lid:

- A fully flexible optical heated lid system designed to prevent loss of sample and condensation. Flexible temperature control allows further method optimization where different heated lid temperatures would be advantageous during different stages of a program run.

Optical Path:





Quantica® – the *real-time* nucleic acid detection system from Techne has been designed with the advantage of having an open architecture and chemistry format that allows the end user full flexibility of the methods and research they wish to pursue.

Quantica® provides at least comparable sensitivity and reproducibility with greater flexibility than many competitive instruments currently on the market whilst using a standard 96 well format.





Accompanying Quantica® is our unique and fully flexible, intuitive wizard based application software Quansoft.

Quansoft Home Page

By employing a series of user friendly windows Quansoft enables any real-time experiment to be created with ease.



Experiment Editor

Using combinations of plate layout, program and analysis parameters this unique Quansoft feature allows the user even easier management of experiment protocols.





Plate Layout Editor

Within a matter of seconds a 96 well microplate can be left blank or assigned with unknowns, controls, standards, user defined samples or test samples. All are colour coded for ease of identification.




Program Editor

Individual cycles and steps or a ramp read are added quickly to build up and display the thermal program and read points.



Analysis Wizards

Select the analysis wizard of your choice to take you quickly through the setup of the analysis type and its associated parameters for your specific experiment, either pre or post run.



Choose from:

- None
- Baseline corrected
- Quantification
- Dissociation curves
- Plus minus scoring
- Allelic discrimination
- Multi-read




QUANTICA[®]

Real-time nucleic acid detection system

Results Editor

View and easily modify the results parameters of your experiment post run using Quansoft's multiple document interface (MDI).

Export results tables, visually customised graphs or even your raw fluorescence data to most commonly used software packages.



Quantica Reports

Comprehensive, completely user customisable report format can be selected for all the available analysis types.

Quantica® Technical Specifications*

Thermal Cycler

Block format	96 x 0.2ml well
Thermal cycler block specification	8 x Peltier block employing quad circuit technology to enhance performance
Thermal cycler block	Uniformity at 50°C: better than $\pm 0.25^\circ\text{C}$
Thermal cycler max ¹ ramp rate	2.6°C/sec
Thermal cycler temperature range	4°C – 99°C
Sample volume	15 – 50 μl
Heated lid	Adjustable between 10°C above ambient to 110°C
Connectability	Connect up to 4 Quantica to one computer ¹
Plate format	Low profile 96 well plate

Optical Detection System

Excitation source	Solid state white light source
Detector	Photon counting photomultiplier tube
Multiplex dye detection	Up to 4 dyes per reaction tube
User selected filters	Maximum of 4 paired excitation/emission filter cartridge systems suitable with currently used dyes
Fluorescence excitation range	470nm – 650nm
Fluorescence detection range	500nm – 710nm
Dynamic range	At least 9 orders of magnitude from starting copy number
Sensitivity ²	1nM fluorescein in a 20 μl sample. Single initial template copy detection.

Dimensions

Weight	25kg
Size and footprint	41cm x 51cm x 53cm (WxHxD)
Input power	750W
Power supply	100/115/230v 50/60Hz

Ordering information

Description

¹Quantica real-time nucleic acid detection system - please contact your local Techne representative for details or visit www.techne.com

* Specification correct at going to print.

¹ Minimum computer specifications apply.

Computer available as an option in the UK only.

² Assay dependent



Rest of World:

Techne
Duxford Cambridge
CB2 4PZ England
Tel: +44 (0) 1223 832401
Fax: +44 (0) 1223 836838
Email: sales@techne.com
www.techne.com

North and South America:

Techne Incorporated
3 Terri Lane Suite 10 Burlington
NJ 08016 USA
Toll free: 800-225-9243
Tel: 609-589-2560
Fax: 609-589-2571
Email: quantica@techneusa.com
www.techneusa.com

Middle East:

Techne Middle East Ltd
Po Box 27887
Engomi 2433
Nicosia Cyprus
Tel: +357 22 660427
Fax: +357 22 660356
Email: salesme@techne.com

Germany:

Techne AG
Obere Hauptstrasse 67b
D-09235 Burkhardtsdorf
Germany
Tel: +49 (0)3721-271888
Fax: +49 (0)3721-270608
Email: info@techne.de