


CAPP[®]

D E N M A R K

CRR-0316



USER MANUAL

TABLE OF CONTENT

1. Product Description	1
2. Safety Precautions	4
3. Standard Parts Listing	6
4. Installation	7
5. Removing & Installing Disc Plate	8
6. Operation	9
7. Warranty Statement	9
8. Product Disposal	11

1. PRODUCT DESCRIPTION

1.1 INTRODUCTION

This manual provides important safety information for the Disc Rotator. It should be kept near the equipment for quick & easy reference. The rotating motion offers thorough mixing. The disc plate can be easily removed and replaced for cleaning.

1.2 INTENDED USE

This Disc Plate rotator is ideal for biological samples, particularly blood samples. Also suitable for separating and mixing of powdered and liquid samples.

1.3 STANDARD ACCESSORIES

1. Power supply (adaptor)
2. Product user manual & warranty card

NOTE: *Before using the instrument, please read this user manual carefully. This user manual is intended to assist with the operation and care of the unit only and not its repair. For repair please contact the supplier.*

For your reference, make a note of serial number, date of purchase and supplier here.

Serial No.:

Purchase Date:

Supplier:

•1.4 FEATURES

- AC synchronous motor for noiseless operation.
- Accommodate upto 4 ESR Tubes & 12 Blood collection tubes
- Circular rotating action gives gentle but effective mixing
- Fix speed of 30 RPM with an accuracy of +/-2 RPM
- Continuous / infinite mode of operation
- Fix disc plate angle of 38°
- Robust construction with strong metal design
- Maximum operational weight capacity of 1.5kgs
- Compact saves bench space design
- Easy to remove and install disc plates for replacement and cleaning.

1.5 TECHNICAL SPECIFICATIONS

Motor Type	AC Synchronous Motor
Motion Type	Rotating - fix degree at 38°
Number of Disc Plate	1
Speed Range	30 RPM
Speed Accuracy	± 2 RPM
Disc Plate Angle	38°
Operation Type	Continuous Operation
Max. Loading Capacity	1.5 Kg
Timer Range	Continuous mode
Removable Disc Plate	Yes
Weight	1.7 Kg
Permissible Ambient Temp.	5° C to 40 °C
Permissible Relative Moisture	< 80%
Air Pressure	80 to 106 kPa
Dimensions (W x D x H)	278 x 252 x 183 mm
Input Voltage	110 - 220 VAC, 50 / 60 Hz
Power Consumption	115V : 3W, 220V : 6W

2. SAFETY PRECAUTIONS


- Never use the instrument in any manner not specified in this manual.
- Equipment used in any manner not specified in this manual or by the manufacturer might result in the lapse of warranty.
- Do not try to stop the disc plate by blocking / holding it. Do not attempt to remove disc plate during operation.
- Never move the instrument while the equipment is under operation.
- Repairs must only be performed by authorized service technician.
- Only use recommended original spare parts for best result & product safety.
- If liquids gets spilled on the disc plate, the instrument must be cleaned carefully and properly before being used again.
- Prior to mixing, the tubes should be visually inspected for material damage.
- Damaged tubes must not be used as this can result in sample loss and the contamination of the product.
- The maximum capacity of 2 kg must not be exceeded.
- This product must be used for specified applications only. It must not be operated in a hazardous or flammable environment and must not be used to mix explosive or highly reactive substances.
- Do not place the potential hazardous material within the clearance area/envelope.
- Wear your personal protective equipment in accordance with the hazards category of the medium to be processed. There is risk of liquid splashing.
- Place the instrument on a flat, leveled and solid surface. Make sure




- the feet of the instrument are clean and undamaged.
- Load the tubes safely in the clamps provided on the disc plate. Always distribute the samples evenly on the device.
- Do not move the instrument while it is operating or connected to the main power supply.
- Ensure that the disc plate knob is properly tightened before performing any operation.
- Do not apply any pressure to the disc plate knob. It should be used only for locking and unlocking disc plate
- Ensure that only closed tubes are used for mixing.
- Do not lean on the equipment. It may damage the equipment or harm the operator.
- Do not fill tubes near the device. Liquid spillage may harm the device.
- In the event of contamination caused by aggressive agents, the device must be cleaned immediately using a natural cleaning liquid. If any damage is seen, contact the service technician.
- The power adaptor given with the unit is designed to be used for this particular unit. Do not use any other power adaptor.
- Only process media which will not react dangerously to the extra energy produce by mixing.
- Do not operate the appliance in explosive atmosphere with hazardous substance or under water.

3. STANDARD PARTS LISTING

Font Side



Back Side



4. INSTALLATION

Gently remove the upper packaging and take out the Disc Rotator by holding it from the bottom. When this equipment is used for the first time, ensure that all the packaging accessories are removed from the product. Please keep all the packaging in safe storage for at least 2 years for warranty purposes.


4.1 LOCATION & MOUNTING

Place the device on a flat and leveled surface and ensure that all four feet stand on the surface firmly. Avoid installing on a slippery surface or surface prone to vibration.

1. Ideal ambient temperature is $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$; avoid placing the device in direct sunlight.
2. Keep clearance of at least at least 30 cm from all sides to guarantee cooling efficiency.
3. Keep this away from heat or water to avoid sample temperature issues.
4. Do not place the device in such a place that it becomes difficult to operate.

4.2 CONNECTING POWER ADAPTOR

1. Connect one side of the power adaptor to the rear side of device and other to the power supply as shown in the figure .
2. Ensure the power switch is OFF while connecting the power adaptor.




5. REMOVING & INSTALLING DISC PLATE

5.1 REMOVING DISC PLATE

Disc plate can be removed by a simple turning of a knob in counter-clockwise direction. Turn the knob as shown in the below image and then remove disc plate.

Below image shows the process to remove disc plate.


Turn the knob in Counter Clockwise direction



5.2 INSTALLING DISC PLATE

Follow the exact reverse process to insert the disc plate back to position


Turn the knob in Clockwise direction



6. OPERATION

6.1 SWITCH ON THE DEVICE

After connecting the power adaptor to the side of the device as shown in below picture, switch ON the main power supply & then switch ON the device from the rear side. Make sure that tubes are distributed evenly before use.



6.2 START / STOP OPERATION

Once all the pre-requisites procedure are done, switch on the main power supply to START the operation. Speed gradually increases from 0 RPM to 30 RPM. This equipment runs on continuous mode (for infinite time).

Note: *There is no provision to change Speed in this equipment*

7. WARRANTY STATEMENT

This product is warranted to be free from defects in material and workmanship for a period of two (2) years from date of purchase. Your product will be duly repaired upon prompt notification in compliance with the following conditions :

- This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in this instruction manual.

This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship. This warranty does not cover any incidental or consequential damages, commercial loss or any other damages from the use of this product.

- The warranty is invalidated by any non-factory modification, which will immediately terminate all liabilities on us for the products or damages caused by its use. The customer shall be responsible for the product or use of products as well as any supervision required for safety. If requested the products must be returned to the distributor in well packed and insured manner and all shipping charges must be paid.

***Note:** Some states do not allow limitation on the length of implied warranties or the exclusion or limitation of incidental or consequential damages. This warranty gives you specific legal rights. This warranty is given expressly in lieu of all other warranties, expressed or implied.*

The purchaser agrees that there is no warranty of merchantability or of fitness for any intended purpose and there are no other remedies or warranties, implied, which extend beyond the description on the face of the agreement. This warranty is only applicable to the original purchaser.

- Product received without proper authorization will not be processed for warranty or service. All items returned for service should be sent with postage prepaid in the original packaging or other suitable packaging, padded to avoid damage. We will not be responsible for damage incurred by improper packaging.

***Note:** This warranty is valid only if the warranty is registered with the supplier within 30 days from the date of purchase.*

8. PRODUCT DISPOSAL

In case the product is to be disposed of, the relevant legal regulations are to be observed.

Information on the disposal of electrical and electronic devices in the European Community

The disposal of electrical devices is regulated within the European Community by national regulations based on EU Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). According to these regulations, any devices supplied after 13.06.05 in the business to business sphere, to which this product is assigned, may no longer be disposed off in municipal or domestic waste. They are marked with the following symbol to indicate this.



As disposal regulations within the EU may vary from country to country, please contact your supplier if necessary.



Landbrugsvej 10 | DK 5260 | Odense S | Denmark

Tlf +45 6613 6140 | Fax +45 6613 2770

www.capp.dk | info@capp.dk