# Thermo Shakers and BlockThermostats with Easy Control MHR 11 – TH 21

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# Content

<b>1. Basic Data</b>	<b>4</b>
Safety	4
Warranty	5
Disclaimer of Warranty	5
Important Information	5
Environmental conditions	6
Temperature Accuracy	6
<b>2. First Steps</b>	<b>6</b>
Setup	6
Temperatur Setting	7
Operation Shaking	10
<b>3. Assembly of Accessories</b>	<b>10</b>
Anti-condensation plate BA 24 / 96	10
<b>4. Troubleshooting</b>	<b>11</b>
<b>5. Maintenance</b>	<b>11</b>
Cleaning	11
Service	13
<b>6. Technical Data</b>	<b>14</b>
<b>7. Article description blocks &amp; accessories</b>	<b>15</b>
<b>8. Explanations</b>	<b>17</b>

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## 1. Basic Data

The HLC by DITABIS Thermo Shakers and Block Thermostats are used for tempering and mixing of solutions in closed reaction tubes and plates. Please pay attention that the target temperature of the samples should not exceed the boiling point of the solution, which should be heated. Please consider this when setting the target temperature via the software of the Thermo Shaker or Block Thermostat. If there are any uncertainties regarding the liquids, which should be heated or cooled, please contact DITABIS. The HLC by DITABIS devices are meant for indoor usage only. Please use only HLC by DITABIS accessories or accessories recommended by HLC by DITABIS.

### Safety

The systems comply with the standards and directives mentioned in the applicable CE declaration.

**Please take note of the following safety measures:**

- » Do only connect the system to an earthed mains power socket of 230 V, 50 Hz.
- » Do not place hot blocks on inflammable or not heat-resistant surfaces.
- » Make sure to only use containers in the block which are suitable for the desired temperature range.
- » If liquid gets into the system, immediately pull out the plug and contact our service department to ensure complete safety.
- » Please do not touch or transport hot systems.
- » Contact with highly flammable fluids must be avoided.
- » The ventilation slots of the device must not be blocked, to ensure the proper ventilation at all times.
- » Damaged mains cable should only be replaced with equal ones.

## Warranty

All functions of the systems have been tested thoroughly. After that, the system and the accessories were packed carefully in perfect condition.

If, however, any damages or defects are detected during installation or setup, please inform DITABIS immediately (see our Terms and Conditions at [www.ditabis.com](http://www.ditabis.com)).

If used under normal laboratory conditions according to the Operating Manual, DITABIS grants a warranty of one year for the material and workmanship, starting with the day of shipment.

## Disclaimer of warranty

In case of unspecified use of the device, the manufacturer assumes no responsibility.

## Important Information

This manual complies with the standards and directives mentioned in the EN DIN 61010-1 Norm. IQ, OQ and PQ documents are available on request.

The systems are provided with a high-precision self-optimising temperature controller. Due to the self-optimisation it might occur that the set temperature is exceeded in individual cases.

Please note, that by heating your samples dangerous gases may be emitted. In this case, the device must be used with an extractor hood.

If the Thermo Shakers are operated with a high shaking frequency, the vibrations of the system might transmit to the underground. Please do consider this when placing the system.

Contact with highly flammable liquids must be avoided. The exchangeable block may be hot.

The ventilation slots of the device must not be blocked, because the proper ventilation of the device must be assured at all times.

## Environmental conditions

The ideal area of operation of the devices is achieved at 80% relative air humidity at most and lays between +3°C - +50°C ambient air temperature. This range should not be underrun or exceeded.

## Temperature Accuracy

Every HLC by DITABIS system is calibrated with a gauged high-precision temperature measuring device to temperature accuracy and reproducibility. This process is performed at a room temperature of +20°C and a block temperature of +37°C. The absolute calibration of the heating unit can also be performed for a certain temperature range if the customer desires. HLC by DITABIS reaches a very high accuracy by measuring the block temperature directly inside the block. It can, however, not be avoided that the cold room temperature affects the sample temperature in case of high temperature difference between room and block temperature. The anti-condensation plate reduces this effect. We recommend using the anti-condensate plate BA 24 / BA 96 or the tempering tub BT 01 / BT 02, if high temperature accuracy in the sample is required. When setting a process time, please note that the sample temperature is always reached a little later than the displayed block temperature.

## 2. First Steps

### Setup

Plug the mains cable in the IEC socket on the back of the system and connect it with the 230 V, 50 Hz via an earthed mains power socket.

## Changeable blocks

At delivery, the changeable blocks are packed separately. Place the blocks on the tempering plates and take care of the proper alignment of the central threaded bolt and the two small corner bolts ensuring that the block lays accurately on the tempering plate. Then attach the block with the short socket screw tightly at the bottom of the tempering plate to get the block fixed for shaking and to obtain even better temperature accuracy.

Check the tight attachment of the block by pulling it up. If it loosens, please use the longer screw.

Never start the process without the block being tightly attached. Do always use two blocks for two-block systems, even if only one of them is used. This ensures safety as well as a smooth shaking operation.

Switch on the system with the mains switch on the back of the system, the actual temperature value is indicated on the display.

## Temperature Setting


---

Actual  
value:

After power-on, the current actual  
temperature value is displayed.

---


Display set  
value:

Push **Set**. The last entered set value  
is displayed and a hyphen appears  
on the top left as a marking.

---

Back to  
actual  
value:


Push **OK** and the current actual val-  
ue is displayed.

---

Change set  
value:

Push **Set** to get to the set value.  
Push **Set** again and then release the  
button. The cursor stands now at the  
figure in the middle (1. Digit posi-


---

tion).

With **Set** and **Start / Stop** you can navigate to the digit of the set value which should to be changed.


Change values:

The digit of the set value is increased or decreased by pushing **+** or **-**. By exceeding the value 9, the next decade is increased automatically. To activate the hundreds digit, exceed the value 9 on the 2. digit.


By falling below the value 0, the next-lower decade is selected automatically.

Confirm set value:

Push **OK** and the current actual value is displayed again.


Start heating process:

If the system heats, a “:” emerges. The green LED lights. If a process time is preset, the horizontal hyphen flashes as an  to indicate the preset time.

Push **Start / Stop**.to start the heating process. Depending on the heat output it either flashes or is permanently visible.

Stop heating process:

The “:” turns off and the current temperature actual value is displayed. If a process time is preset and if the process is stopped before the expiration of this time, the clock


---

is reset to 0:00. For that push **Start / Stop**.

---


Set time:

To set the time, please return to actual-value (see page 6) and push **OK** for **three seconds**, then release the button. The display indicates :00 (= infinite), to activate the hour digit press **Set**. Up to 19 hours and 59 min can be reset. The flashing cursor stands at the ten-minute-digit, the other digits can be selected with **Set** and **Start / Stop**.

---


Change values:

The selected time digit is increased or decreased by pushing **+** or **-**. By exceeding the value 9, the next decade is increased automatically. By falling below the value 0, the next-lower decade is selected automatically.

---


Confirm process time:

Push **OK** and the current actual value is displayed.

---

Start process:


To start a process push **Start / Stop**. The preset process time only starts if the set value is reached. If a process time is preset, the horizontal hyphen flashes. After the expiration of the time, the heating switches off.

---

---

Displayed  
remaining  
time:


The time remaining until the end of the process is indicated on the display with flashing “ : “. By pushing **OK** for 10 seconds and releasing the button, you get back to the temperature display.

---

### Operation Shaking

For shaking, the mains switch on the back of the system must be switched on.

With the scaled slider, the shaking frequency can be set to a value between 0 and 1.500 rpm in steps of 250. For fine-tuning, in steps of 10, use the keys – and + . The currently entered rotational speed is indicated on the display for three seconds. The recommended min. rotational speed is 80 rpm, lower rotational speeds can also be set.

Please note that the Heating Thermo Shaker MHR 11 will shake immediately if the shaking frequency is set, independent of the start of the tempering.

If the rotational speed is to be displayed during a running process, decrease or increase the rotational speed with + and – by 10 rotations, the current rotational speed value is then displayed for three seconds.

## 3. Assembly of Accessories

### Anti-condensation plate BA 24 / 96

The anti-condensation plate BA 96 (800013000) is placed on the attached block, no assembly steps are required. For using the anti-condensation plate

BA 24 (800012900), remove the socket screw in the block with the enclosed screwdriver. Screw the enclosed threaded pin with isolation knob into the anti-condensation plate. Place the anti-condensation plate on the block and screw the threaded pin in deeper so that it takes hold of the bolt of the system and with that tightly attaches block and anti-condensation to the system.

## 4. Troubleshooting

### **The display remains blank**

Please check if the mains switch on the back side is switched on. If it is switched on but voltage is present at the socket, check the micro-fuse and replace if required (IEC 127-2/III, 250 V, 2 A time-lag). This fuse - and a spare fuse - are located in the IEC-bushing (in which the connection cable is plugged). The fuse box can be pulled out with a screwdriver.

### **There are excessive temperature fluctuations**

Check the seat of the changeable block by pulling it up. If it is too loose or if the block exhibits roughness, e.g. contamination, the heat / coolness is not transferred correctly.

### **Information on the capacitive touch display**

The touch display reacts on finger pressure, even when wearing thin latex gloves, but not on stylos. Due to its glass surface, the front panel is insensitive to dirt, chemicals and mechanical damage. Please avoid scratches in the coating since they can cause malfunctions. Please take note of the following cleaning instructions.

## 5. Maintenance

### Cleaning

Regularly clean the housing and the changeable blocks of the Thermo Shaker and BlockThermostates.

#### Precautions for avoiding electric shock

Electronic devices can cause electric shocks in case of an operating error. Never try to repair electric parts. Never open the housing.

- Switch off the system and disconnect it from the power supply before starting with cleaning or disinfection works.
- Never let get liquids inside the housing (ventilations slit).
- Do not perform spray disinfection.
- Do only connect the system with the power supply if it is completely dry.

The repair service may only be performed by staff authorized and trained by the manufacturer. A modification of the system is not permitted.

#### Caution when handling aggressive chemicals

Do not use aggressive chemicals like e.g. strong and weak bases, strong acids, formaldehyde, acetone, halogenated hydrocarbons or phenol for cleaning the system and its accessories.

- In case of contamination with aggressive chemicals, clean the system with a neutral detergent immediately.
- Use neither corrosive detergents nor aggressive solvents or abrasive polishing agents.

## **Cleaning**

1. Disconnect the system from the power supply before you start cleaning.
2. Clean all outer parts of the system with a mild detergent and a lint-free cloth.
3. Wipe off the detergent with distilled water.
4. Dry all cleaned parts.

## **Disinfection**

1. Disconnect the system from the power supply before you start disinfecting.
2. Let the system cool down.
3. Clean the system as described above.
4. Select a disinfection method compliant to the applicable local legal regulations and directives.
5. Wipe off all outer parts of the system with the disinfectant and a lint-free cloth.

## **Service**

If a technical problem arises, please contact your local dealer or the DITABIS service department. The DITABIS contact information can be found at [www.ditabis.com](http://www.ditabis.com). If required, the system will be returned for repair – please take note of our service guidelines, which can be found at [www.ditabis.com](http://www.ditabis.com) as well.

## **Decontamination before shipment**

If you send the system to the authorised technical service for repair or to your distributor for disposal, decontaminate all parts you want to send. Document the decontamination in a Decontamination Certificate (incl. serial number) and include it with the shipment. Please detach exchangeable blocks before shipment.

## 6. Technical Data

<b>Technical Data</b>	<b>MHR 11</b>
Temperature-working range	Amb. +3°C to +130°C
Temperature-adjustable range	0°C to +137°C
Accuracy / resolution	+/- 0.1°C / 0.1°C
Max heating time	11.5°C / min.
Shaking frequency	80 – 1,500 rpm
Orbit	3 mm round
Dimensions (without block) W x D x H	220 x 330 x 109 mm
Block - capacity	1 changeable block
Weight (without block)	6.5 kg
Power input	200 W
Electr. supply	230 V, 50 Hz (115 V / 60 Hz available)
Fuse	2 AT
Protection class	IP21
<b>Environmental conditions</b>	
Ambient air temperature in service (non-condensing)	3° – 50°C
Relative humidity	max. 80%

<b>Technical Data</b>	<b>TH 21</b>
Temperature-working range	Amb. +3°C to +130°C
Temperature-adjustable range	0°C to +137°C
Accuracy / resolution	+/- 0.1°C / 0.1°C
Max heating time	9.5°C / min
Dimensions (without block) W x D x H	220 x 330 x 109 mm
Block - capacity	2 changeable blocks
Weight (without block)	3.5 kg
Power input	200 W

Electr. supply 230 V / 50 Hz  
(115 V / 60 Hz available)

Fuse 2 AT

Protection class IP21

### Environmental conditions

Ambient air temperature in service (non-condensing) 3° – 50°C

Relative humidity max. 80%

## 7. Article description blocks & accessories

Blocks	Art. No.	Dimensions of the Containers
For micro test tubes	800010800 800010900 800011000 800011100	BM 02 for 96 x 0,2 ml conical & 8.container stripes BM 05 for 38 x 0,5 ml conical BM 15 for 24 x 1,5 ml conical BM 20 for 24 x 2,0 ml cylindrical
For sample tubes	800011500 800011600 800011700 800011800 800014200 800016100 800015200 800016000 800011900	BP 10 for 24 x D=10,3 mm, 46 mm deep, round bottom, lid BP 12 for 24 x D=12 mm, 20 mm deep, flat bottom BP 15 for 24 x D=15 mm, 30 mm deep, flat bottom BP 16 for 24 x D=16,5 mm, 46 mm deep, round bottom, lid BP 17.0 for 24 x D=17 mm, 30 mm deep, flat bottom BP 19.2 for 24 x D=18,7 mm, 25 mm deep, flat bottom BP 23.0 for 12 x D=22,5 mm, 56 mm deep, flat bottom BP 25.5 for 12 x D=23,0 mm, 56 mm deep, flat bottom BP 28 for 8 x D=28 mm, 40 mm deep, flat bottom
For centrifuge tubes	800012200 800012300	BZ 15 "Falcon" tube 14 x 15 ml, with isolation lid BZ 50 for "Falcon" tube 6 x 50 ml, with isolation lid
For PCR plates	800010400 800010300	BC 96 for 96-Well "V" bottom BC 84 for 384-Well
For micro plates	800012000 800010600	BV 96 for 96 x round or "V" bottom BF 96 for flat bottom
For deep-well plates	800010500	BD 96 for 96-Deep-Well plates


For other applications	800011200	BN 10 for 36 rectangular cuvettes, outer diameter 12.5 mm
	800012100	
	800014100	SO 10.4 for 24 x D=10,4 mm, 180 mm deep
	800012600	SO 12.0 for 24 x 12ml tubes
	800012700	SO 20.5 for 12 x 20ml tubes
For customized containers	800010700	BM 00 without boring, height 25 mm
	800011300	BO 37 without boring, height 37 mm
	800011400	BO 50 without boring, height 50 mm
For slides	800012400	BY 12 for 12 slides for hybridisation

## Accessories

Art. No.	Description
800012800	BI 01 Isolation lid
800013800	BT 01 Tempering tub for 1 block, transparent
800013900	BT 02 Tempering tub for 2 blocks, transparent
800012900	BA 24 Anti-condensation plate for 1 block (micro tubes)
800013000	BA 96 Anti-condensation plate for 1 block (PCR/ test plates)
800012500	BH 10 Heating top for the connection with cooling units
800013100	BR 05 Rack with holder f. tubes 0.5 ml
800014300	BR 15 Rack with holder f. tubes 1.5/ 2.0ml

## 8. Explanations

	Caution! Risk of electric shock!
	Caution!
	Caution! Hot surface!
	Important information
	Earth conductor
	Fuse

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