## Thermo Shakers and BlockThermostats with Smart Control

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# Content

<b>1. Basic Data</b>	<b>4</b>
Safety	4
Warranty	5
Disclaimer of warranty	5
Important information	5
Environmental conditions	6
Temperature accuracy	6
<b>2. First Steps</b>	<b>7</b>
Setup	7
<b>3. Operation</b>	<b>7</b>
Parameter setting	7
Programming function	13
<b>4. Assembly of Accessories</b>	<b>19</b>
Anti-Condensation Plate BA 24 / 96	19
Data transfer	20
<b>5. Troubleshooting</b>	<b>23</b>
<b>6. Maintenance</b>	<b>24</b>
Cleaning	24
Service	25
<b>7. Technical data</b>	<b>26</b>
<b>8. Article description blocks &amp; accessories</b>	<b>30</b>
<b>9. Explanations</b>	<b>31</b>

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## 1. Basic Data

The HLC by DITABIS Thermo Shakers and Block Thermostats are used for tempering and mixing of solutions in closed reaction tubes and plates. Please pay attention that the target temperature of the samples should not exceed the boiling point of the solution, which should be heated. Please consider this when setting the target temperature via the software of the Thermo Shaker or Block Thermostat. If there are any uncertainties regarding the liquids, which should be heated or cooled, please contact DITABIS. The HLC by DITABIS devices are meant for indoor usage only. Please use only HLC by DITABIS accessories or accessories recommended by HLC by DITABIS.

### Safety

The systems comply with the standards and directives mentioned in the applicable CE declaration.

#### **Please take note of the following safety measures:**

- » Do only connect the system to an earthed mains power socket of 230 V, 50 Hz.
- » Do not place hot blocks on inflammable or not heat-resistant surfaces.
- » Make sure to only use containers in the block which are suitable for the desired temperature range.
- » If liquid gets into the system, immediately pull out the plug and contact our service department to ensure complete safety.
- » Please do not touch or transport hot systems.
- » Contact with highly flammable fluids must be avoided.
- » The ventilation of the device must not be blocked. To ensure the proper ventilation at all times please consider a minimum distance of 10 centimetres on each side of the device.
- » Damaged mains cable should only be replaced with equal ones.

## Warranty

All functions of the systems have been tested thoroughly. After that, the system and the accessories were packed carefully in perfect condition. If, however, any damages or defects are detected during installation or setup, please contact your local dealer or the DITABIS service department. The DITABIS contact information can be found at [www.ditabis.de](http://www.ditabis.de) (see also our Terms and Conditions at [www.ditabis.com](http://www.ditabis.com)). If used under normal laboratory conditions according to the Operating Manual, DITABIS grants a warranty of one year, starting with the day of shipment.

## Disclaimer of warranty

In case of unspecified use of the device, the manufacturer assumes no responsibility.

## Important information

This manual complies with the standards and directives mentioned in the EN DIN 61010-1 Norm. IQ, OQ and PQ documents are available on request.

The systems are provided with a high-precision self-optimising temperature controller. Due to the self-optimisation it might occur that the set temperature is exceeded in individual cases.

Please note, that by heating your samples dangerous gases may be emitted. In this case, the device must be used with an extractor hood.

If the ThermoMixers are operated with a high shaking frequency, the vibrations of the system might transmit to the underground.

Please do consider this when placing the system.

Contact with highly flammable liquids must be avoided. The exchangeable block may be hot.

The ventilation of the device must not be blocked. To ensure the proper ventilation at all times please consider a minimum distance of 10 centimetres on each side of the device.

## **Environmental conditions**

The ideal area of operation of the devices is achieved at 80% relative air humidity at most and lays between +3°C - +50°C ambient air temperature. This range should not be underrun or exceeded.

## **Temperature accuracy**

Every HLC by DITABIS system is calibrated with a gauged high-precision temperature measuring device to temperature accuracy and reproducibility. This calibration is performed at a room temperature of +20,0°C at 10 different temperature calibration points and oil as a reference liquid. HLC by DITABIS reaches a very high accuracy by measuring the block temperature directly inside the block. It can, however, not be avoided that the room temperature affects the sample temperature in case of a large temperature difference between room and block temperature. The anti-condensation plate reduces this effect. We recommend using the anti-condensate plate BA 24 / BA 96 or the tempering tub BT 01 / BT 02. When setting a process time, please note that the sample temperature is always reached a little later than the displayed block temperature.

## 2. First Steps

### Setup

Plug the mains cable in the IEC socket on the back of the system and connect it with the mains power 230 V, 50 Hz via a Schuko socket. Please make sure that you are using earthed mains power sockets.

### Changeable blocks

At delivery, the changeable blocks are packed separately. Place the blocks on the tempering plates and take care of the proper alignment of the central threaded bolt and the two small corner bolts ensuring that the block lays accurately on the tempering plate. Then attach the block with the enclosed socket screw tightly at the bottom of the tempering plate to get the block fixed for shaking and to obtain even better temperature accuracy. Check the tight attachment of the block by pulling it up

Never start the process without the block being tightly attached. Do always use two blocks for two-block systems, even if only one of them is used. This ensures safety as well as a smooth shaking operation.

Switch on the system with the mains switch on the back of the system. The start menu will be displayed. The green LED turns on, on the right sight of the touch panel, as soon as a process is running.

## 3. Operation

### Parameter setting

---


#### Start menu

This menu is displayed after the power is switched on, after abortion or after the normal completion of a process. The last set values are retained even after the switch-off of the system.

---


25,6 °C	SP: 25,6 °C
0 rpm	<b>Start</b>
	SP: 0 rpm
	100% Short-Mix
	Program
time: not started	SP: endless B0

## Setup menu

Push the **left and right arrow keys at the same time**, then the setup menu is displayed. With **+** and **–** the values can be selected for setting.

The selected box flashes, with **ok** it will be permanently set for changing.

With the keys **+** and **–**, the input or the values can be changed. If an input is confirmed with **ok**, the box will highlighted flashing again. With confirming the button **Start menu**, the start menu is displayed again.

The following values can be adapted:

Language:	EN
Tempering unit:	°C
Block number:	B0
confirm process end:	Yes
rotation direct.:	left
goto:	Page2
HLC system Set-up	Startmenu

### Page 1:

Language: English, German, French, Spanish

Tempering unit: °C or °F

Block number (10-point calibration)

Confirm end:

**Yes:** at the process end, a signal is audible until the **ok** key is pushed. Until then, the process runs with the previous data.

**No:** the process ends without confirmation.

Contrast:	70%
Brightness:	90%
Display:	norm
Goto	Page 1
HLC system Set-Up	Startmenu

Rotation direction: right, left

## Page 2:

Contrast: Mentioned in %

Brightness: Mentioned in %

Display colour: Normal - inverse

## 10-point-calibration

25,6 °C	SP: 25,6 °C
0 rpm	<b>Start</b>
	SP: 0 rpm
	100% Short-Mix
	Program
time:	SP: endless
not started	<b>B0</b>

By default the standard calibration begins with the block number "B0". 10 temperature bases are underlying this calibration. The calibration is optimized for all DITABIS standard blocks. If you are using customized blocks you will be able to do individual calibrations in addition. For the calibration you need a gauged high-precision temperature measuring device.

-16,0:	-18,0	76,5:	74,5
2,5:	0,5	95,0:	93,0
21,0:	19,0	113,5:	110,5
39,5:	37,5	132,0:	130,0
58,0:	<b>56,0</b>	150,0:	148,0
B1			
HLC Block Set-Up		Startmenu	

Select block number in the set up menu and confirm with **ok**. You are now able to calibrate and save different individual calibrations, from standard block B1 up to 4. Therefore please select B1 to B4. Confirm again with **ok**. Now you will get to the calibration menu. In the first and third column you can see the temperature base. In the second and fourth column you can see the variable set values.

By pushing the **+** and **-** keys the entries and values can be changed. The selected field is blinking; by pushing **ok** the field is selected for changing.

## Example:

You have adjusted the device for 95,0°C in the start menu and are measuring only 85,0°C

-16,0:	-16,0	76,5:	74,2
2,5:	2,5	95,0:	<b>85,0</b>
21,0:	19,0	113,5:	109,4
39,5:	37,5	132,0:	126,6
58,0:	60,0	150,0:	145,0
B1			
HLC Block Set-Up		Startmenu	

with a calibrated thermometer inside your sample. Now you have to change the temperature set value in column four about 10°C.

By doing so you can calibrate every single base individually.

## Temperature set value

<b>25,6 °C</b>	SP: <b>25,6 °C</b>
<b>0 rpm</b>	<b>Start</b>
	SP: 0 rpm
	100% Short-Mix
	Program
time:	SP: endless
not started	<b>B0</b>

Navigate with **+** and **-** to the temperature input box and confirm with **ok**. The individual numbers can be selected with **arrow left** and **right**, their value can be changed with **+** and **-**. The input is confirmed with **ok** and the cursor automatically jumps to the **Start** button.

A set value change can also be made in the same way during a running process. After confirming the input with **ok**, the new set value is effective immediately. For the hundreds, the following values can be selected:  
 0 for temperatures up to +99,9°C  
 1 for temperatures from +100,0°C  
 - for temperatures below 0,0°C

## Important information

For cooling systems, the minimum temperature may not be smaller than the difference from room temperature mentioned below.  
 MKR 13 / TK 23: 16°C below r.t.  
 MKR 23: 11°C below r.t.

## Mix set value

Navigate with **+** and **-** to the shaking frequency input box and confirm with **ok**. The individual numbers can be selected with **arrow left** and **right**, its value can be changed with **+** and **-**. The input is confirmed with **ok** and the cursor

25,6 °C	SP: 37,0 °C
0 rpm	Start
	SP: 0 rpm
	100% Short-Mix
	Program
time: not started	SP: endless B0

automatically jumps to the button **Start**. The 1. Digit number cannot be selected and changed. A set value change can also be made in the same way during a running process. After confirming the input with **ok**, the new set value is effective immediately.

### Duration of the process (time setting)

25,6 °C	SP: 37,0 °C
0 rpm	Start
	SP: 0 rpm
	100% Short-Mix
	Program
time: not started	SP: <b>endless</b> B0

Navigate with **+** and **-** to the time input box and confirm with **ok**. The individual numbers can be selected with **arrow left** and **right**, its value can be changed with **+** and **-**. The input is confirmed with **ok** and the cursor automatically jumps to the button **Start**. For an endless process, select **00:00**.

### 100% Short-Mix (vortexing)

25,6 °C	SP: 37,0 °C
0 rpm	Start
	SP: 0 rpm
	<b>100% Short-Mix</b>
	Program
time: not started	SP: endless B0

Navigate with **+** and **-** to the button **100% Short-Mix** and confirm with **ok**. The system shakes with maximum mixing speed regardless of a process is started or not. During the whole shaking process **ok** must be pushed continuously..


### Select program mode

25,6 °C	SP: 25,6 °C
0 rpm	Start
	SP: 0 rpm
	100% Short-Mix
	<b>Program</b>
time: not started	SP: endless B0

Navigate with **+** and **-** to the button **Program** and confirm with **ok**. For detailed information on programming see page 10.

---


## Process without time setting

Only set a set temperature and a set shaking frequency but set the time on continuous **00:00**. Now, the process time is displayed continuously. It only starts to run if the temperature set value is reached. A set value change can also be made as described above during a running process. Changes are retained also after the completion of the process and are displayed in the start menu as current set values. For terminating the process, push **Stop**.


---

## Process with time setting

If a process time is set, the remaining time is displayed in the menu. The time only starts to run if the temperature set value is reached. This does also apply for subsequent temperature changes. A set value change can also be made as described above during a running process.

The process ends without a signal or continuous with a signal until it is terminated manually. The relevant modus can be set in the setup menu.

### For changes during the process

If set value changes have been made, you can select to save these changes and to have them displayed in the start menu as new set values.

### Stopping the process before the expiration of the time


If the process is stopped or changed before the expiration of the time, you can select whether the process really should be termi-

---

nated or if it should continue. Depending on the selection, the start menu or the selection to continue the process from the stop (time runs from stop) or to make a restart (time runs from zero) is displayed.

---

## Process end

Setting in setup:

### **Confirm end: Yes**

If a process time is preset, the expiration of the time is indicated by a flashing LED and a signal, the process continues. You can terminate the process with **ok**.

### **Confirm end: No**

If a process time is preset, only the information that the process is completed will be displayed after the expiration of the time. The device turns off automatically.

---

## Programming function


### Program menu

123456789	Startmenu
free: 12/30	Delete Start
temp(°C): 42,0	prog-no: 1
risingtime: 6	cycles: 1
speed: 1100	step: 4
pause(min): 0,5	hours: 01
mix(s) 1	minutes: 25

Navigate with **+** and **-** to the button **Programs** and confirm with **ok**. A total of 30 different program steps can be preset, distributed on max. 9 different programs.

Confirm the highlighted flashing program number if you want to change it. Confirm the progr. no. with **ok**. With **+** and **-**, the program numbers 1-9 can be selected. If you want to edit the displayed program, select the different adjustable parameters with **+** and **-**. For editing, confirm with **ok** and change the value. For final confirmation confirm with **ok** again.

---

## Setting of the individual paramater

<b>1</b> 2 3 4 5 6 7 8 9		Startmenu	
free: 12/30		Delete	Start
temp (°C): 42,0		prog-no:	<b>1</b>
risingtime: 6		cycles:	1
speed: 1100		step:	4
pause (min): 0,5		hours:	01
mix (s) 1		minutes:	25

By selecting "Programs" in the start menu, you will have a total number of 30 different (recurrent) program steps available. These can be selected in up to 9 different programs. You have the choice of f. ex. either 9 programs with 3 program runnings each or f. ex. 1 program with 30 program runnings. For normal operation the choice will be somewhere in between.

"**prog-no.**" is blinking and shows 1 as standard. By pressing **ok** any program number. between 1 and 9 can be selected using buttons **+** and **-**. Confirm your entry by pressing **ok** and the next field "cycles" will be blinking.

**Cycles** = Number of repetitions of the program, selection 1-9 possible without having to confirm "Start" again.

If no changes are to be made, select the next field by using the **-** (down) button. By pressing **ok** you can change the number of cycles from 1 to 9 using the **+** and **-** buttons. Confirm your entry by pressing **ok** and the next field "step" will be blinking.

**Step** = Section within a program. Selection 1-9 possible.

If no changes are to be made, select the next field by using the **-** (down) button.

---

For changing the program step displayed, press **ok**. A value between 1 and 30 can be selected using buttons **+** and **−**. Confirm your entry by pressing **ok** and the next field "hours" will be blinking. If in the field "minutes" no entries are to be made, select the next field by using the **−** (down) button.

**Hours** = Duration of the section in hours.

If no changes at "hours" are to be made, select the next field by using the **−** (down) button. For any changes press **ok**; the ones column is highlighted and can be changed with buttons **+** and **−**. For changing the tens column use button for marking and **+** and **−** for changing. Confirm your entry by pressing **ok** and the next field "minutes" will be blinking.

**Minutes** = Duration of the section in minutes.

For entries press **ok**; the tens column is highlighted and can be changed by using buttons **+** and **−**. For changing the ones column use button for marking and **+** and **−** for changing. Confirm your entry by pressing **ok** and the next field "temp °C" will be blinking.

### **Time "endless"**

If a program step (normally the last one) shall be running without a time setting until the program will be completed manually, i.e. endless, the hours and minutes must be set to 00.

**Temperature** = Set temperature of the section.


The temperature set value is 37°C as standard. If no changes are to be made, select the next field by using the **−** (down) button. To change this value press **ok**. Every single digit can be changed. The cursor is placed at the tens column of the set value. This value can be changed with buttons **+** and **−**. To change another digit in this data field, it must be marked by buttons **<** and **>**. The highlighted value will be changed with button **+** or **−**. After setting the value it must be confirmed with **ok**; the field "rising time" will be blinking.

**Rising time** = Time setting for tempering the sample to set temperature. This setting is only required if the sample is to be tempered slower than in presttings of the device. If the sample is to be tempered as fast as possible, set the value to 0. With "rising time" the heating or cooling speed compared with the normal speed

---

can be reduced. You have the possibility to determine a certain period (= rising time; with a max. of 99 minutes to choose) in which a set temperature shall be reached in the actual program step, if a slower than normal tempering is desired. As shown in the illustration beside, a temperature of 60,0°C in step 2 shall be reached within 20 minutes, starting from a temperature of 37,0°C in step 1. If no changes are to be made in "rising time" (fastest possible heating or cooling time), select the next field by using the – (down) button. For changing the standard value 0 press **ok**; the ones column is highlighted and can be changed with buttons **+** and **–**. For changing the tens column use button for marking and **+** and **–** for changing. Confirm your entry by pressing **ok** and the next field "speed" will be blinking.

#### Example of a temperature program:

**Rotational speed** = Shaking frequency with which the blocks are shaken. The last number (1. digit) can not be selected and changed.

For "speed" 0 is set as standard. If no changes are to be made (tempering without shaking), select the next field by using the – (down) button. To change this value press **ok**; every single digit can be changed. The cursor is placed at the hundreds column of the set value. This value can be changed with buttons **+** and **–**. To change another digit in this data field, it must be marked by buttons **<** and **>**. The highlighted value will be changed with button **+** or **–**. The last digit (ones column) cannot be chosen and changed. After setting the value it must be confirmed with **ok**; the field "pause" is blinking.

---


**Interval shaking** = Interval shaking is an alternative when the probe shall not be subject to permanent shaking. A resting time of up to 9.9 minutes can be set in "pause". The following short shaking time of up to 9 seconds is set in "mix sec". This is a reasonable choice for the last program running of a complete program with time setting "endless" during off working hours.

**Pause [min]** = Resting time of up to 9.9 minutes. If no intervals are to be set, enter the value 0.0 here.

Duration of pause is 0.0 as standard. If no changes are to be made, select the next field "step" by using the – (down) button. To change this value press **ok**; the ones column is highlighted and can be changed by using buttons **+** and **–**. For changing the 1/10 digit press button **>**; the digit is highlighted and can be changed with buttons **+** and **–**. Confirm your entry by pressing **ok** and the next field "mix sec" will be blinking.

**Mix [s]** = Short shaking time up to 9 seconds after resting time.

For "mix sec" 1 is set as standard. If no changes are to be made, select the next field "step" by using the – (down) button. For changing the shaking time press **ok** and change time by using buttons **+** and **–**. Confirm your entry by pressing **ok** and the next field "step" will be blinking.


---

### Entry of further program steps

If another program step shall be added to a chosen program, press **ok** at the blinking field "step" and increase the number of the last program step by 1. Further entries as already described.

### To change the program data entry

To change the chosen program, enter the program step number to be changed. Select the field to be changed by using buttons **+** and **–**. Press **ok** for changing. Changes are made as already described

and must be confirmed with **ok**. Select a corresponding field by using buttons **+** and **-** or **< / >**.

### Termination of program data entry

To complete a chosen program, mark either "Start" (to start the program immediately) - screen 6 will appear, or "Startmenu" (to save the program only) - screen 1 will appear. You can also select "prog-no." (for entry of another program). When "Delete" is selected, the **complete** program with all program steps as displayed is deleted

### Starting a program

1 2 3 4 5 6 7 8 9	Startmenu
free: 12/30	Delete <b>Start</b>
temp (°C): 42,0	prog-no: 1
rising time: 6	cycles: 1
speed: 1100	step: 4
pause (min): 0,5	hours: 01
mix (s) 1	minutes: 25

Select the relevant prog. no. and confirm the **Start** button in the program menu.

### Changing parameters during the program process

<b>25,6 °C</b>	SP: 25,6 °C
1100 rpm	<b>Stop</b>
	SP: 1100 rpm
	100% Short-Mix
	Program
time: 01:30	SP: 02:00
process running	P1-Z1-S1 B0

Change of temperature and rotational speed as well as performing a Short-Mix are possible during a running program.

The set time is the time of the running step, the remaining time is the time of the complete process. In the lower right part of the display, program no., cycle and step are indicated.

---

## Termination/ interruption of the program:

26,6 °C	SP: 26,6 °C
1100 rpm	<b>Stop process</b>
	SP: 1100 rpm
	100% Short-Mix
	cont. proc.
time: 01:30	SP: 02:00
process running	P1-Z1-S1 B0

If you push **Stop**, the options **stop process** and **continue process** are displayed.

26,6 °C	SP: 26,6 °C
0 rpm	<b>Start menu</b>
	SP: 1100 rpm
stopped and finished	
B0	

If you select **stop process**, the **Start menu** or **Save changes** is displayed if parameters have been changed during the process.

26,6 °C	SP: 26,6 °C
0 rpm	<b>Cont. from Stop</b>
	SP: 1100 rpm
	Proc. re-start
time rem.: 00:10	SP: 02:00
stopped at	P1-Z1-S1 B0

If you select **Continue with process** the selection **Continue from stop** or **Process Restart** is displayed (program starts from cycle 1, step 1).

---

## 4. Assembly of accessories

### Anti-Condensation Plate BA 24 / 96

The anti-condensation plate BA 96 (800013000) is placed on the attached block, no assembly steps are required. For using the anti-condensation plate BA 24 (800012900), remove the socket screw in the block with the enclosed

screwdriver. Screw the enclosed threaded pin with isolation knob into the anti-condensation plate. Place the anti-condensation plate on the block and screw the threaded pin in deeper so that it takes hold of the bolt of the system and with that tightly attaches block and anti-condensation to the system.

## **Data transfer**

A USB 1.1. interface for the communication with a PC is provided as a standard. Connect the system via the USB-port on the left side with a USB-cable to your PC. If this is not available at your computer, install the suitable USB-driver. You will find it under [www.ftdichip.com/Drivers/VCP.htm](http://www.ftdichip.com/Drivers/VCP.htm). These USB-drivers generate a new virtual COM-interface in your PC. The chip in the Smart Control is called FT232B.

### **Access process data via USB**

Start a terminal program (e.g. Hyper terminal) on your PC. This program you will find under Windows in "Programs" / "Accessories" / "Communication". The window "Hyper Terminal" is opened. Proceed systematically, please note that for the connection setting for "Bits per second" 115200 is selected. Select at "Transfer": "Record Text", enter the location and define the file name as .txt Start the Hyper Terminal and the process. Every minute, the following data are recorded in text format and separated by comma in the form shown on the right.

Since the data are separated by commas, you can easily create an Excel file and generate one or more curves from that. If you open the text file you selected for saving, you can see the complete process.

START	00:12,+61.5,0300,01,01,02,00:08,02,9,0.1
Setpoints	00:13,+64.7,0000,01,01,02,00:09,01,9,0.1
xx:xx,+37.0,0200,01,01,01,00	STOP
:01,00,9,0.2	
time	
hh:mm,tmp,rpm,prg,cycle,st	
ep,steptime hh:mm,rem.	
risetime_m,mix_s,pause_m	
00:00,+26.9,0000,01,01,01,0	
0:00,00,9,0.2	
00:01,+33.7,0220,01,01,01,0	
0:00,00,9,0.2	
00:02,+36.7,0000,01,01,01,0	
0:00,00,9,0.2	
00:03,+36.9,0000,01,01,01,0	
0:00,00,9,0.2	
START	
Setpoints	
xx:xx,+70.0,0300,01,01,02,00	
:12,10,9,0.1	
time	
hh:mm,tmp,rpm,prg,cycle,st	
ep,steptime hh:mm,rem.	
risetime_m,mix_s,pause_m	
00:04,+36.9,0000,01,01,02,0	
0:00,10,9,0.1	
00:05,+39.1,0000,01,01,02,0	
0:01,09,9,0.1	
00:06,+42.1,0000,01,01,02,0	
0:02,08,9,0.1	
00:07,+45.2,0300,01,01,02,0	
0:03,07,9,0.100:08,+48.4,030	
0,01,01,02,00:04,06,9,0.1	
00:09,+51.7,0000,01,01,02,0	
0:05,05,9,0.1	
00:10,+54.9,0000,01,01,02,0	
0:06,04,9,0.1	
00:11,+58.2,0300,01,01,02,0	
0:07,03,9,0.1	

## **Enter process data via USB**

For the process control of a HLC system with Smart Control via a USB-interface, four commands are possible and required:

### **Start**

Starts a process without program with the parameters indicated on the display. Time settings, if applicable, are deleted and overwritten by continuous since time settings must be programmed on the computer.

### **Stop**

Stops all processes and returns to the start menu, running programs / processes are aborted.

### **t=0370**

Changes the set temperature to a new value.

The temperature must always be entered with four digits.

t=0370 sets 37.0°C, t=0050 sets 5.0°C, t=-060 sets -6.0°C.

t= is the identification for a temperature value. The associated numerical value reflects the value in tenth of a degree.

The first digit after the = can be a 0, a 1 or a hyphen. Invalid temperature values are - like for the input at the system - replaced by the next valid value.

### **r=020**

Changes the set rotational speed to a new value.

The rotational speed must always be entered with three digits.

r=120 sets 1,200 1/min, r=045 sets 450 1/min.

r= is the identification for a rotational speed value. The associated numerical value reflects the value in 10 1/min.

Invalid rotational speed values are - like for the input at the system - replaced by the next valid value.

### **r=000**

Stops the rotational speed, e.g. if shaking should be stopped but tempering be continued.

Other inputs like e.g. r=0 are interpreted as rotational speed and cause the motor to rotate with the minimum rotational speed.

### **Every command must be confirmed with "CR"**

In case of unknown commands and/or in case of spelling errors, nothing is done. The inputs in the setup menu like colour, contrast and brightness of the display, language, rotation direction etc. can only be done at the operating panel of the system itself.

## **5. Troubleshooting**

### **The display remains blank**

Please check if the mains switch on the back side is switched on.

If it is switched on but voltage is present at the socket, check the micro-fuse and replace if required. (IEC 127-2/III, 250 V, 2 A time-lag).

This fuse - and a spare fuse - are located in the IEC-bushing (in which the connection cable is plugged). The fuse box can be pulled out with a screwdriver.

### **The system does not cool, heat or shake as set**

Check whether the display indicates the correct system type name at power-on. If not, please contact your local dealer or the HLC by DITABIS service department.

### **There are excessive temperature fluctuations**

Check the seat of the changeable block by pulling it up. If it is too loose or if the block exhibits roughness, e.g. contamination, the heat / coolness is not transferred correctly.

### **Information on the capacitive touch display**

The touch display reacts on finger pressure, even when wearing thin latex gloves, but not on stylos.

Due to its glass surface, the front panel is insensitive to dirt, chemicals and mechanical damage. Please avoid scratches in the coating since they can cause malfunctions. Please take note of the following cleaning instructions.

## **6. Maintenance**

### **Cleaning**

Regularly clean the housing and the changeable blocks of the Thermo Shakers and BlockThermostates.

#### **Precautions for avoiding electric shocks**

Electronic devices can cause electric shocks in case of an operating error. Never try to repair electric parts. Never open the housing.

- Switch off the system and disconnect it from the power supply before starting with cleaning or disinfection works.
- Never let get liquids inside the housing (ventilations slit).
- Do not perform spray disinfection.
- Do only connect the system with the power supply if it is completely dry.

The repair service may only be performed by staff authorized and trained by the manufacturer. A modification of the system is

not permitted.

### **Caution when handling aggressive chemicals**

Do not use aggressive chemicals like e.g. strong and weak bases, strong acids, formaldehyde, acetone, halogenated hydrocarbons or phenol for cleaning the system and its accessories.

- In case of contamination with aggressive chemicals, clean the system with a neutral detergent immediately.
- Use neither corrosive detergents nor aggressive solvents or abrasive polishing agents.

### **Cleaning**

1. Please disconnect the system from the power supply before you start cleaning.
2. Please clean all outer parts of the system with a mild detergent and a lint-free cloth.
3. Please wipe off the detergent with Aqua dest..
4. Please dry all cleaned parts.

### **Disinfection**

1. Please disconnect the system from the power supply before you start disinfecting.
2. Let the system cool down.
3. Please clean the system as described above.
4. Please select a disinfection method compliant to the applicable local legal regulations and directives.
5. Please wipe off all outer parts of the system with the disinfectant and a lint-free cloth.

### **Service**

If a technical problem arises, please contact your local dealer or the DITABIS service department. The DITABIS contact information can be found at [www.ditabis.com](http://www.ditabis.com). If required, the system will be returned for repair – please take note of our service guidelines, which can be found at [www.ditabis.com](http://www.ditabis.com) as well.

### Decontamination before shipment

If you send the system to the authorised technical service for repair or to your distributor for disposal, decontaminate all parts you want to send. Document the decontamination in a Decontamination Certificate (incl. serial number) and include it with the shipment.

## 7. Technical data

Technical Data	MKR 13
Temperature-working range	RT -16°C to +100°C
Temperature-adjustable range	-10°C to +105°C
Accuracy / resolution	+/- 0.1°C / 0.1°C
Max. heating time	6.0°C / min
Max. cooling time	12.0°C / min
Shaking frequency	200 – 1,500 rpm
Orbit	3 mm round
Dimensions (without block) W x D x H	220 x 330 x 144 mm
Block - capacity	1 changeable block
Weight (without block)	9.0 kg
Power input	130 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
Fuse	2AT
Protection class	IP21
Environmental conditions	
Ambient air temperature	3° – 50°C
in service (non-condensing)	
Relative humidity	max. 80%
Technical Data	MKR 23
Temperature-working range	RT -11°C to +70°C
Temperature-adjustable range	0°C to +80°C

Accuracy / resolution	+/- 0.3°C / 0.1°C
Max. heating time	3.5°C / min
Max. cooling time	6.5°C / min
Shaking frequency	200 – 1,200 rpm
Orbit	3 mm round
Dimensions (without block) W x D x H	220 x 330 x 144 mm
Block - capacity	2 changeable blocks
Weight (without block)	9.5 kg
Power input	200 W
Electr. supply	230 V / 50 Hz (115 V / 60 hz available)
Fuse	2AT
Protection class	IP21

### **Environmental conditions**

Ambient air temperature	3° – 50°C
in service (non-condensing)	
Relative humidity	max. 80%

<b>Technical Data</b>	<b>MHR 13</b>
Temperature-working range	RT +3°C to +130°C
Temperature-adjustable range	0°C to +137°C
Accuracy / resolution	+/- 0.1°C / 0.1°C
Max. heating time	11.5°C / min
Shaking frequency	200 – 1,500 rpm
Orbit	3 mm round
Dimensions (without block) W x D x H	220 x 330 x 109 mm
Block - capacity	1 changeable block
Weight (without block)	6.5 kg
Power input	200 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
Fuse	2AT

Protection class	IP21
------------------	------

**Environmental conditions**

Ambient air temperature	3° – 50°C
-------------------------	-----------

in service (non-condensing)

Relative humidity	max. 80%
-------------------	----------

**Technical Data****MHR 23**

Temperature-working range	RT +3°C to +130°C
---------------------------	-------------------

Temperature-adjustable range	0°C to +137°C
------------------------------	---------------

Accuracy / resolution	+/- 0.1°C / 0.1°C
-----------------------	-------------------

Max. heating time	9.5°C / min
-------------------	-------------

Shaking frequency	200 – 1,500 rpm
-------------------	-----------------

Orbit	3 mm round
-------	------------

Dimensions (without block) W x D x H	220 x 330 x 109 mm
--------------------------------------	--------------------

Block - capacity	2 changeable blocks
------------------	---------------------

Weight (without block)	7.0 kg
------------------------	--------

Power input	350 W
-------------	-------

Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
----------------	--

Fuse	2AT
------	-----

Protection class	IP21
------------------	------

**Environmental conditions**

Ambient air temperature	3° – 50°C
-------------------------	-----------

in service (non-condensing)

Relative humidity	max. 80%
-------------------	----------

**Technical Data****MHL 23**

Temperature-working range	RT +3°C to +130°C
---------------------------	-------------------

Temperature-adjustable range	0°C to +137°C
------------------------------	---------------

Accuracy / resolution	+/- 0.1°C / 0.1°C
Max. heating time	9.5°C / min
Shaking frequency	200 – 1,300 rpm
Orbit	3 mm linear
Dimensions (without block) W x D x H	220 x 330 x 109 mm
Block - capacity	2 changeable blocks
Weight (without block)	7.0 kg
Power input	350 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
Fuse	2AT
Protection class	IP21
<b>Environmental conditions</b>	
Ambient air temperature	3° – 50°C
in service (non-condensing)	
Relative humidity	max. 80%

Technical Data	TK 23
Temperature-working range	RT -16°C to +90°C
Temperature-adjustable range	-10°C to +105°C
Accuracy / resolution	+/- 0.3°C / 0.1°C
Max. heating time up	4.0°C / min
Max. cooling time	7.0°C / min
Dimensions (without block) W x D x H	220 x 330 x 144 mm
Block - capacity	2 changeable blocks
Weight (without block)	7.5 kg
Power input	130 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
Fuse	2AT
Protection class	IP21

#### **Environmental conditions**

Ambient air temperature in service (non-condensing)	3° – 50°C
Relative humidity	max. 80%

## 8. Article description blocks & accessories

Blocks	Art. No.	Dimensions of the Containers
For micro test tubes	800010800 800010900 800011000 800011100	BM 02 for 96 x 0.2 ml conical & 8-container stripes BM 05 for 38 x 0.5 ml conical BM 15 for 24 x 1.5 ml conical BM 20 for 24 x 2.0 ml cylindrical
For sample tubes	800011500 800011600 800011700 800011800 800014200 800016100 800015200 800016000 800011900	BP 10 for 24 x D=10.3 mm, 46 mm deep, round bottom, lid BP 12 for 24 x D=12 mm, 20 mm deep, flat bottom BP 15 for 24 x D=15 mm, 30 mm deep, flat bottom BP 16 for 24 x D=16.5 mm, 46 mm deep, round bottom, lid BP 17.0 for 24 x D=17 mm, 30 mm deep, flat bottom BP 19.2 for 24 x D=18.7 mm, 25 mm deep, flat bottom BP 23.0 for 12 x D=22.5 mm, 56 mm deep, flat bottom BP 25.5 for 12 x D=23.0 mm, 56 mm deep, flat bottom BP 28 for 8 x D=28 mm, 40 mm deep, flat bottom
For centri- fuge tubes	800012200 800012300	BZ 15 "Falcon" tubes 14 x 15 ml, with isolation lid BZ 50 for "Falcon" tubes 6 x 50 ml, with isolation lid
For PCR plates	800010400 800010300	BC 96 for 96-well "V"-bottom BC 84 for 384-well
For micro plates	800012000 800010600	BV 96 for 96 x round or "V"-bottom BF 96 for flat bottom
For deep- well plates	800010500	BD 96 for 96-deep-well-plates
For other applications	800011200 800012100  800014100	BN 10 for 36 rectangular cuvettes, outer diameter 12.5 mm BW 01 as tub for deep-well, stacked microtiter plates and other con- tainers, incl. isolation lid and unloading device SO 10.4 for 24 x D=10,4 mm, 180 mm deep

	800012600	SO 12.0 for 24 x 12ml test tubes
	800012700	SO 20.5 for 12 x 20ml test tubes
For custom- ized containers	800010700	BM 00 without drilling, height 25 mm
	800011300	BO 37 without drillings, height 37 mm
	800011400	BO 50 without drillings, height 50 mm
For Slides	800012400	BY 12 for 12 slides for hybridisation

## Accessories

Art. No.	Description
800012800	BI 01 Isolation lid
800013800	BT 01 Tempering tub for 1 block, transparent
800013900	BT 02 Tempering tub for 2 blocks, transparent
800012900	BA 24 Anti-condensation plate for 1 block (micro test tubes)
800013000	BA 96 Anti-condensation-plate for 1 block (PCR / test plates)
800013100	BR 05 Rack with holder f. tubes 0.5 ml
800014300	BR 15 Rack with holder f. tubes 1.5 / 2,0ml

## 9. Explanations

	Caution! Risk of electric shock!
	Caution!

	<p>Caution! Hot surface!</p>
	<p>Important information</p>
	<p>Earth conductor</p>
	<p>Fuse</p>

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