User manual for:

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Important safety instructions

1. To obtain full use of the cabinet, we recommend reading this instruction manual.
2. It is the user's responsibility to operate the appliance in accordance with the instructions given.
3. Contact your dealer immediately in case of any malfunctions.
4. Place the cabinet in a dry and ventilated place.
5. Keep the cabinet away from strongly heat-emitting sources and do not expose it to direct sunlight.
6. Always keep in mind that all electrical devices are sources of potential danger.
7. Do not store inflammable material such as thinner, gasoline etc. in the cabinet.
8. We declare that no asbestos nor any CFC are used in the construction.
9. The oil in the compressor does not contain PCB.

Unpacking and installation

Remove the wooden pallet and the packing. Stainless steel surfaces are supplied with a protection foil, which must be removed before installation

Electrical connecting

The cabinet operates on 230 V/50 Hz.

The wall socket should be easily accessible.

All earthing requirements stipulated by the local electricity authorities must be observed. The cabinet plug and wall socket should then give correct earthing. If in doubt, contact your local supplier or authorized electrician.

The flexible cord fitted to this appliance has three cores for use with a 3-pin 13-Amp or 3-pin 15-Amp plug. If a B.S. 1363 (13-Amp) fused plug is used, it should be fitted with a 13-Amp fuse.

The wires in this mains lead are coloured in accordance with the following code:
Green/Yellow: Earth, Blue: Neutral, Brown: Live.

The main electrical connections must be done by skilled electricians.

Start-up of the cabinet

Before use, we recommend that the cabinet is cleaned, see the section on maintenance and cleaning.

Important !

If the cabinet has been vertically placed during transport, please wait 2 hours before starting up the cabinet.

Thermostat

Display set temperature:

Press this key and the display will show the set temperature. Press the key again to return to normal reading.

Set new temperature:

Press this key for 3 sec. and the display shows the set temperature.

Press this key to increase the set temperature.

Press this key to lower the set temperature.

Press this key to save the new settings. The display will flash with the new value and will then return to normal reading.

Changing of parameters:

+

Press these keys continuously for 5 seconds to access the list of parameters. The display will show the number of the parameter and after 2 seconds the value.

Press this key to increase the value.

Press this key to lower the value.

Press this key to save the new settings. The display will flash with the new value and will then return to normal reading.

See table of parameters from page 11.

Alarm codes:

P1 Flashing in the display: indicates that the cabinet sensor is defective.
The cabinet will strive to keep the set temperature until it has been repaired.

HA Flashing in the display: indicates that the temperature is too high.

LA Flashing in the display: indicates that the temperature is too low.

Defrosting


We recommend defrosting the freezer every three month or if the layer of ice exceeds 5 mm.

Do never use metal or other sharp objects to scrape the internal sides as it can damage the freezer

1. Disconnect the freezer
2. Empty the freezer
3. Leave the freezer with the lid open
4. Use the plastic scraper to remove the ice from the inner liner
5. Remove the drain plug and open the drain on the back side of the freezer and place a tray under the drain
6. Clean the freezer with mild detergent and wipe thoroughly
7. Connect the freezer and it is ready for use.

Lock

The freezer is supplied with a lock in the handle

Maintenance and cleaning

Switch the freezer off at the socket.

The cabinet must be periodically cleaned. Clean the external and internal surfaces of the cabinet with a light soap solution and subsequently wipe dry. Stainless steel surfaces can be maintained using steel oil.

Do NOT use cleansers containing chlorine or other harsh cleansers, as these can damage the surfaces and the internal cooling system.

Clean the condenser and the compressor compartment using a vacuum cleaner and a stiff brush.

Do NOT hose the compressor compartment with water as this can cause short circuits and damage on the electrical parts.

Service

The cooling system is a hermetically sealed system and does not require supervision, only cleaning.

If the cabinet fails to cool, check if the reason is a power cut.

If you cannot locate the reason to the failure of the cabinet, please contact your supplier. Please inform model and serial number of the cabinet. You can find this information on the rating label which is placed inside the cabinet in the top right hand side.

Disposal

Disposal of the cabinet must take place in an environmentally correct way. Please note existing regulation on disposal. There may be special requirements and conditions which must be observed.

XR20C (SE-45 Series)

PAR.	DESCRIPTION	RANGE	DEFAULT	U.M.
Set	Set point	LS-US	-45	°C/°F
Hy	diFFerential. Relay compressor tripping differential. The compressor stops on reaching the Setpoint value, and restarts at temperature equal to the Setpoint plus the value of the differential.	0.1...25.0	2,0	°C/°F
LS	Lower SEt. Minimum possible setpoint value.	-55.0...HSE	-45	°C/°F
US	Higher SEt. Maximum possible setpoint value.	LSE...302	-5,0	°C/°F
Ot	CAlibration 1. Calibration 1. Positive or negative temperature value added to the value read by probe 1, based on "ot" parameter settings.	-9.9...9.9	0	°C/°F
P2P	Second probe presence: n=not present, y=present	n/y	N	flag
oE	CAlibration 2. Calibration 2. Positive or negative temperature value added to the value read by probe 2, based on "P2" parameter settings.	-9.9...9.9	0	°C/°F
odS	delay (at) On compressor. Delay time in activating the compressor relay after switch-on of instrument.	0...99	2	min
AC	Anti short cycle delay	0...50	5	Min
CON	On time (compressor). Compressor activation time in the event of faulty probe.	0...99	10	min
COF	OFF time (compressor). Compressor in disabled state time in the event of a faulty probe.	0...99	10	min
CF	Measurement units: °C=Celsius, °F=Fahrenheit	°C/°F	°C	°C/°F
rES	Resolution (only for °C): dE=decimal between -9.9 and 9,9; in=integer	dE-in	dE	flag
Lod	Default display: P1=thermostat probe, P2=evaporator probe, SP=Set point	P1,P2,SP	P1	flag
td	defrost type. Type of defrosting. El = electric defrost. in = reverse cycle defrost (hot gas)	EL/in	In	flag
dtE	Defrost termination temp.	-50...50	6	°C/°F
idF	defrost interval time. Interval between the start of two successive defrosting operations.	0...99	0	hours
MdF	Max length for defrost	0...99	0	min
dFd	Display during defrost: rt=real temp. it=start defrost temp., St=Set point, dF=label Df	Rt/it/SP/dF	St	flag
dAd	Max display delay after defrost	0...255	1	Min
dAF	Defrost delay after fast freezing	0...23	0	hours
ALc	Temperature alarm configuration: rE=relative to SET, Ab=absolute	rE/Ab	rE	Flag
ALU	Higher ALarm. Maximum temperature alarm. Temperature value. which if exceeded in an upward direction triggers the activation of the alarm signal.	AL...99	10	°C/°F
ALL	Lower ALarm. Minimum temperature alarm. Temperature value. which if exceeded in a downward direction, triggers the activation of the alarm signal.	-55...AU	-1	°C/°F
ALd	Temperature alarm delay	0...99	20	min
dAO	Exclusion of temperature alarm at startup.	0...23	1	hour
PbC	Kind of probe	Ptc/ntc	Ntc	flag