English

Deutsch

Français

Italiano

Español


Operating Instructions **FiveEasyPlus™** Conductivity meter FEP30

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METTLER TOLEDO

Operating Instructions
FiveEasyPlus™

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Español

Table of contents

1	Introduction	3
2	Safety measures	4
3	Installation	6
3.1	Unpacking	6
3.2	Installing the Electrode Arm (optional)	6
4	Display and Key controls	7
5	Operating the Instrument	9
5.1	Calibration	9
5.1.1	Selecting a standard	9
5.1.2	Performing a Calibration	9
5.2	Sample Measurement	9
5.2.1	Performing a Conductivity Measurement	9
5.2.2	Performing a TDS / salinity measurement	10
5.3	Setup	10
5.3.1	Set temperature correction coefficient	10
5.3.2	Set reference temperature	10
5.3.3	Set temperature unit	10
5.3.4	Set TDS factor	10
5.4	Using the Memory	10
5.4.1	Storing a Reading	10
5.4.2	Recalling from Memory	11
5.4.3	Clearing the Memory	11
5.5	Printing	11
5.5.1	Connection and Configuration	11
5.5.2	Printing a Measurement Results	11
5.5.3	Printing Calibration Results	12
5.5.4	Printing from Memory	12
5.6	Self-Diagnosis	12
5.7	Error Messages	12
6	Maintenance	14
6.1	Meter Maintenance	14
6.2	Disposal	14
7	Accessories	15
8	Specifications	16
9	Appendix	17
9.1	Conductivity standards	17
9.2	Examples of temperature coefficients (alpha-values)	17
9.3	Practical salinity scale (UNESCO 1978)	17
9.4	Conductivity to TDS conversion factors	18
9.5	Error Limits	18

1 Introduction

Thank you for purchasing this high quality METTLER TOLEDO meter. Ease of use, reliable results and ergonomic design – these are our goals for products.

The FiveEasy meters have an excellent price/performance ratio and many useful accessories are already included in the standard delivery of the instruments. Other accessories such as electrodes, solutions or electrode arms are of course also available from METTLER TOLEDO.

You will find much helpful functionality in these meters: Automatic endpoint format which freezes the measurement automatically when the value is stable and predefined conductivity standards are only few of many exiting features.

2 Safety measures

Measures for your protection

Risk of explosion

- Never work in an environment subject to explosion hazards! The housing of the instrument is not gas tight (explosion hazard due to spark formation, corrosion caused by the ingress of gases).

Risk of corrosion

- When using chemicals and solvents, comply with the instructions of the producer and the general lab safety rules!

Measures for your operational safety

Caution

- Never unscrew the two halves of the housing!
- Have the meter serviced only by METTLER TOLEDO Service!
- Any spillage should be wiped off immediately! Some solvents might cause corrosion of the housing.
- Avoid the following environmental influences:
 - Powerful vibrations
 - Direct sunlight
 - Atmospheric humidity greater than 80%
 - Corrosive gas atmosphere
 - Temperatures below 5 °C and above 40 °C
 - Powerful electric or magnetic fields

FCC Rules

This device complies with Part 15 of the FCC Rules and Radio Interference Requirements of the Canadian Department of Communications. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.


This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

3 Installation

3.1 Unpacking


Carefully unpack the meter. Keep the calibration certificate in a safe place.

Insert the right adapter clip into the power adapter slot:


3.2 Installing the Electrode Arm (optional)

The electrode arm can be fastened to FiveEasyPlus™ on the left or the right side.

- 1 Remove the electrode stand rubber cover on that side of the bench meter where you prefer the stand. Keep the cover in a safe place.
- 2 Insert the electrode stand firmly into the opening as shown.
- 3 Fix the opening of the arm upon the stand from above, adjust it to a moderate height.

4 Display and Key controls

- 1 Calibration icon - calibration in progress
- 2 Measurement icon - measurement or calibration is running
- 3 Setup icon - instrument is in the setup mode
- 4 TDS factor in the setup mode
- 5 Conductivity / TDS / Salinity reading / Cell constant in calibration process
- 6 Endpoint stability $\sqrt{\quad}$, Auto endpoint \sqrt{A}
- 7 Temperature during measurement or reference temperature in setup mode
- 8 Standard  / Memory number **MR** / Error index **Err**

	 Press & release	 Press & hold for 3 seconds
	<ul style="list-style-type: none"> • Start or endpoint measurement • Confirm setting, store entered value 	<ul style="list-style-type: none"> • Turn auto endpoint on \sqrt{A} / off $\sqrt{\text{off}}$
	<ul style="list-style-type: none"> • Start calibration 	<ul style="list-style-type: none"> • Review the latest calibration data
	<ul style="list-style-type: none"> • Meter on • Back to measurement screen 	<ul style="list-style-type: none"> • Meter off
	<ul style="list-style-type: none"> • Store current reading to memory • Increase value during setting • Scroll up through the memory 	<ul style="list-style-type: none"> • Recall stored data • Print current data set of stored data
	<ul style="list-style-type: none"> • Switch between conductivity, TDS and salinity measuring modes • Decrease value during setting • Scroll down through the memory 	<ul style="list-style-type: none"> • Enter Setup mode
	<ul style="list-style-type: none"> • Start self-diagnosis 	

5 Operating the Instrument

5.1 Calibration

5.1.1 Selecting a standard

When using the Five conductivity meter, you have to select a standard for calibration.

- 1 Press and hold the **Mode/Setup** key until the setup icon appears on the display and the current standard blinks.
- 2 Use **▲** or **▼** to select another standard and press **Read** to confirm your selection.
- 3 Press **Exit** to leave the setup mode.

The three predefined standards are

- 84 $\mu\text{S}/\text{cm}$
- 1413 $\mu\text{S}/\text{cm}$
- 12.88 mS/cm

Tables for automatic temperature compensation are programmed in the meter for each standard (see also "Appendix (page 17)").

5.1.2 Performing a Calibration

- 1 Place the electrode in the defined calibration standard and press **Cal**.
 - ⇒ The calibration icon and the measurement icon appear on the display. The meter endpoints according to the preselected endpoint mode after the signal has stabilized or after pressing **Read**. The standard value at the measurement temperature is displayed and stored; the measurement icon disappears from the display.
- 2 The cell constant is shown on the display for 3 seconds. To reject the calibration, press **Exit**, before the meter returns automatically after 3 seconds to the measurement screen.

Note

To ensure the most accurate conductivity readings, you should verify your cell constant with a standard solution once a day and recalibrate if necessary. Use always fresh standards.

5.2 Sample Measurement

5.2.1 Performing a Conductivity Measurement

- ▶ Press **Mode** to switch between conductivity, TDS and salinity measurement mode.
- 1 Place the conductivity sensor in the sample and press **Read** to start the measurement:
 - ⇒ The measurement icon appears on the display and the decimal point blinks. The display shows the conductivity of the sample.
 - 2 The automatic endpoint **A** is the default setting of the meter. When the signal has stabilized, the display freezes automatically, \sqrt{A} appears and the measurement icon disappears.
 - or -
 - To manually endpoint a measurement, press **Read**. The display freezes and $\sqrt{\quad}$ appears.

Note

- By pressing and holding **Read**, you can switch between the auto and manual endpoint modes.
- Stability criterion for conductivity measurement – The sensor input signal of the meter may not change by more than 0.4% from the measured average conductivity of the probe in 6 seconds.

5.2.2 Performing a TDS / salinity measurement

- ▶ Press **Mode** to switch between conductivity, TDS and salinity measurement mode.
- To perform a TDS / salinity measurement, follow the same procedure as for a conductivity measurement, see "Performing a Conductivity Measurement (page 9)".

5.3 Setup

5.3.1 Set temperature correction coefficient

- 1 Press **Setup** the key until the setup icon appears on the display and the current standard blinks.
- 2 Press **Read** to ignore standard setting; the meter changes to the α -coefficient setting automatically.
- 3 Use ▲ or ▼ to increase or decrease the temperature correction value.
- 4 Press **Read** to confirm your setting.
- 5 Continue with reference temperature setting or press **Exit** to go back to the measurement screen.

5.3.2 Set reference temperature

- ▶ After confirming the temperature correction coefficient, the reference temperature appears on the display. See "Set temperature correction coefficient (page 10)".
- 1 Use ▲ or ▼ to toggle between 25 °C and 20 °C.
- 2 Continue with temperature unit setting or press **Exit** to go back to the measurement screen.

5.3.3 Set temperature unit

- ▶ After selecting the reference temperature, the current temperature unit blinks. See "Set reference temperature (page 10)".
- 1 Use ▲ or ▼ to toggle between °C and °F
- 2 Press **Read** to confirm your selection.
- 3 Continue with TDS factor setting or press **Exit** to go back to the measurement screen.

5.3.4 Set TDS factor

- ▶ After selecting the temperature unit, the current TDS factor blinks.
- 1 Use ▲ or ▼ to increase or decrease the value.
- 2 Press **Read** to confirm your setting.
- ⇒ The meter will exit automatically to the measurement screen.

5.4 Using the Memory

5.4.1 Storing a Reading

- FiveEasyPlus™ meter can store up to 99 endpointed results.
- Press **STO** when the measurement has endpointed.

⇒ **M01** indicates that one result has been stored.

Note

If you press **STO** when **M99** is displayed, **FUL** indicates the memory is full. To store further data you will have to clear the memory (see "Clearing the Memory (page 11)").

5.4.2 Recalling from Memory

- 1 Press and hold **RCL** to recall the stored values from memory when the current measurement has endpointed.
- 2 Press **▲** or **▼** to scroll through the stored results. **R01** to **R99** indicates which result is being displayed.
- 3 Press **Read** to exit.

5.4.3 Clearing the Memory

- 1 Continue pressing **▲** or **▼** to scroll through the stored results until **MRCL** appears.
- 2 Then press **Read**, **CLr** blinks.
- 3 Press **Read** again to confirm the deletion or press **Exit** to return to the measurement mode without deleting the data.

5.5 Printing

5.5.1 Connection and Configuration

A printer can be connected to the RS-232 interface of FEP30. The use of the **RS-P25, RS-P26, RS-P28** printer is recommended, because they recognize the FEP30 and adjust the correct parameters automatically.

In case a different printer is used, the following parameters need to be set:

Baud rate:	1200 bps
Data bit:	8 bit
Parity bit:	None
Stop bit:	1 bit

5.5.2 Printing a Measurement Results

If a printer is connected to the FEP30, a print-out is automatically generated after each endpointed measurement or calibration.

The format for the print-out following a conductivity measurement is:

EP, Value, Unit, Temp.
 AE 85.1µS/cm 28.4°C ATC
 Ref.25°C,T.Comp.2.00%/°C

The details for the second line are:

```

AE 85.1µS/cm 25.0°C ATC
|           |           |----- ATC or MTC
|           |           |----- Temperature value
|           |----- conductivity value
|----- AE for Auto-End-Point,
|----- ME for Manual-End-Point
  
```

The details for the third line are:

Ref.25°C, T.Comp.2.00%/°C
| ----- Coefficient value
----- Reference Temperature

The format for the print-out following a TDS measurement is:

EP, Value, Unit, Temp.
ME 42.0mg/L 25.0°C MTC
Ref.25°C, T.Comp. 2.00%/°C

The format for a salinity measurement is:

EP, Value, Unit, Temp.
ME 0.04psu, 25.0°C MTC

The print-out in case of an error message is:

EP, Value, Unit, Temp.
Error2

5.5.3 Printing Calibration Results

The print-out for a successful calibration is:

Standard: 84µS/cm
Temperature: 25°C
Cell Const.: 1.176E-1/cm
ATC/MTC: MTC

5.5.4 Printing from Memory

When scrolling through the memory (see "Recalling from Memory (page 1 1)") you can print the entry that is currently viewed by pressing and holding the **STO/RCL** key for 2 s.

The printout format is following:

Recall 05:
EP, Value, Unit, Temp.
AE 85.1µS/cm 28.4°C ATC
Ref.25°C, T.Comp. 2.00%/°C

5.6 Self-Diagnosis

- 1 Press and hold **Read** and **Cal** simultaneously until the meter displays the full screen.
 - ⇒ Each icon blinks one after the other. This way you may check whether all icons are correctly shown. The next step is to check that the keys are functioning correctly. This requires user interaction.
- 2 When **b** blinks, five icons are displayed. Press the five keys in any order.
 - ⇒ Each time you press a key an icon disappears from the screen.
- 3 Continue to press the other keys until all the icons have disappeared.
 - ⇒ When the self-diagnosis has been completed successfully, **PAS** appears. If self-diagnosis fails, error message **Err 1** appears (see "Error Messages").

Note

You have to finish pressing all five keys within two minutes, otherwise **Err 1** appears and you will have to repeat the procedure.

5.7 Error Messages

Error 0	Memory access error	Reset to factory settings.
---------	---------------------	----------------------------

Error 1	Self-diagnosis failed	Repeat the self-diagnosis procedure and make sure that you finish pressing all five keys within two minutes.
Error 2	Measured values out of range	Check if the electrode is properly connected and placed in a sample solution.
Error 3	Measured standard temperature out of range (5...35 °C)	Keep the standard temperature within the range for calibration.
Error 4	Measuring temperature out of range (0...100 °C)	Check if the electrode is properly connected and keep the sample temperature within the range.
Error 9	The current data set has already been stored once	A measurement can only be stored once. Perform a new measurement to store a new data set.

6 Maintenance

6.1 Meter Maintenance

Never unscrew the two halves of the housing!

The Five series instruments do not require any maintenance other than occasional wipe with a damp cloth.

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is attacked by some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). Any spillage should be immediately wiped off.

6.2 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

7 Accessories

Item	Order No.
Power supply	51302950
Electrode arm	51302951
Rubber covers (for electrode arm hole)	51302952
Shorting plug BNC	51302859
84 $\mu\text{S}/\text{cm}$ standard solution, 500mL	51302153
1413 $\mu\text{S}/\text{cm}$ standard solution sachets, 30 x 20 mL	51302049
1413 $\mu\text{S}/\text{cm}$ standard solution, 250 mL	51350096
12.88 mS/cm standard solution sachets, 30 x 20 mL	51302050
12.88 mS/cm standard solution, 250 mL	51350098
LE703, conductivity sensor, 10 $\mu\text{S}/\text{cm}$ - 500mS/cm	51340335
LE 740, conductivity sensor, 0.01 $\mu\text{S}/\text{cm}$ - 500 $\mu\text{S}/\text{cm}$	51340336

8 Specifications

FiveEasyPlus™ conductivity meter FEP30		
Measurement range	Conductivity	0.0 µS/cm... 199.9 mS/cm
	TDS	0.1 mg/L ... 199.9 g/L
	Salinity	0.00 ... 19.99 psu
	Temperature	0 °C ... 100 °C
Resolution	Conductivity	Auto range
		0.10 µS/cm ... 19.99 µS/cm
		20.0 µS/cm ... 199.9 µS/cm
		200 µS/cm ... 1999 µS/cm
		2.00 mS/cm ... 19.99 mS/cm
	20.0 mS/cm ... 199.9 mS/cm	
	TDS	Auto range, same as conductivity
Salinity	0.00 ... 19.99 psu	
Temperature	0.1 °C	
Limits of error	Conductivity	± 0.5 % of measured value
	TDS	± 0.5 % of measured value
	Salinity	± 0.5 % of measured value
	Temperature	± 0.3 °C
Number of memory	99	
Power requirements	9 V DC, 0.4 W ± 10% For use with CSA certified (or equivalent approved) power source, which must have a limited circuit output.	
Size/weight	200 x 175 x 52 mm / 0.6 kg	
Display	Liquid crystal	
Signal Input	Mini-Din	
IP rating	For indoor use only	
Ambient conditions	Temperature	5...40 °C
	Relative humidity	5 %...80 % (non-condensing)
	Height	up to 2000 m
	Installation category	II
	Pollution degree	2
Materials	Housing	ABS/PC reinforced
	Window	Membrane
	Keypad	Membrane

9 Appendix

9.1 Conductivity standards

T (°C)	84 µS/cm	1413 µS/cm	12.88 mS/cm
5	53 µS/cm	896 µS/cm	8.22 mS/cm
10	60 µS/cm	1020 µS/cm	9.33 mS/cm
15	68 µS/cm	1147 µS/cm	10.48 mS/cm
20	76 µS/cm	1278 µS/cm	11.67 mS/cm
25	84 µS/cm	1413 µS/cm	12.88 mS/cm
30	92 µS/cm	1552 µS/cm	14.12 mS/cm
35	101 µS/cm	1667 µS/cm	15.39 mS/cm

9.2 Examples of temperature coefficients (alpha-values)

Substance at 25 °C	Concentration [%]	Temperature coefficient alpha [%/°C]
HCl	10	1.56
KCl	10	1.88
CH ₃ COOH	10	1.69
NaCl	10	2.14
H ₂ SO ₄	10	1.28
HF	1.5	7.20

α-coefficients of conductivity standards for a calculation to reference temperature 25 °C

Standard	Measurement temp.: 15 °C	Measurement temp.: 20 °C	Measurement temp.: 30 °C	Measurement temp.: 35 °C
84 µS/cm	1.95	1.95	1.95	2.01
1413 µS/cm	1.94	1.94	1.94	1.99
12.88 mS/cm	1.90	1.89	1.91	1.95

9.3 Practical salinity scale (UNESCO 1978)

The salinity is calculated according to the official definition of UNESCO 1978. Therefore the salinity S_{psu} of a sample in psu (practical salinity unit) at standard atmospheric pressure is calculated as follows:

$$S = \sum_{j=0}^5 a_j R_T^{j/2} - \frac{(T-15)}{1+k(T-15)} \sum_{j=0}^5 b_j R_T^{j/2}$$

a ₀ = 0.0080	b ₀ = 0.0005	k = 0.00162
a ₁ = -0.1692	b ₁ = -0.0056	
a ₂ = 25.3851	b ₂ = -0.0066	
a ₃ = 14.0941	b ₃ = -0.0375	
a ₄ = -7.0261	b ₄ = 0.0636	
a ₅ = 2.7081	b ₅ = -0.0144	

$$R_T = \frac{R_{\text{Sample}}(T)}{R_{\text{KCl}}(T)}$$

(32.4356 g KCl per 1000 g of solution)

9.4 Conductivity to TDS conversion factors

Conductivity at 25 °C	TDS KCl		TDS NaCl	
	ppm value	factor	ppm value	factor
84 µS/cm	40.38	0.5048	38.04	0.4755
447 µS/cm	225.6	0.5047	215.5	0.4822
1413 µS/cm	744.7	0.527	702.1	0.4969
1500 µS/cm	757.1	0.5047	737.1	0.4914
8974 µS/cm	5101	0.5685	4487	0.5000
12.880 µS/cm	7447	0.5782	7230	0.5613
15.000 µS/cm	8759	0.5839	8532	0.5688
80 mS/cm	52.168	0.6521	48.384	0.6048

9.5 Error Limits

Message	Description	Range not accepted
Err 2	Measured values out of range	C: > 199.9 mS/cm TDS: < 0.1 mg/L or > 199.9 g/L SAL: > 19.99 psu
Err 3	Calibration standard temperature out of range	T: < 5 °C or > 35 °C
Err 4	Temperature out of range	T: < 0 °C or > 100 °C

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pH, conductivity & DO meters
EL20, EL2, FG4 (FEPx0, FEx0, ELx0, FGx, ELx -Series)
For additional types, see page type code
For optional equipment, see page accessories

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

2006/95/EC Low voltage (LVD)
2004/108/EC Electromagnetic compatibility (EMC)

and that the standards have been applied.
und die Normen zur Anwendung gelangten.

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Chris Radloff
General Manager

Rolf Truttmann
Head SBU pH Lab

References of standards for this declaration of conformity, or parts thereof:
Harmonized standards of Europe and Switzerland:

Safety standards:

IEC/EN61010-1:2001
IEC/EN61010-1:2010

EMC standards:

IEC61326-1:2005 / EN61326-1:2006 (class B)
IEC61326-1:2005 / EN61326-1:2006 (Basic requirements)

Metrological standards:

--

IP standards:

IEC/EN60529:1991+A1:2000

Standards for Canada, USA and Australia:

CAN/CSA-C22.2 No. 61010-1-04
UL Std. No. 61010-1 (2nd Edition)

AS/NZS CISPR 11, AS/NZS 61000.4.3

EC - DECLARATION OF CONFORMITY

EG-Konformitätsklärung
KD-Nr.: 51302969 D

Doku-Nr.: 20070006

Type code Typenschlüssel	
other types of same construction: andere Typen/Modelle mit der gleichen Konstruktion:	
EL20, EL2, FG4	Tested types
Portable equipment FGx ELx	Portable equipment Series name for FiveGo Series name for Education Line (Design is different)
Bench equipment FEx0 ELx0	Bench equipment Series name for FiveEasy Series name for Education Line (Design is different)
Bench equipment FEPx0	Bench equipment, software and with RS-232 port Series name for FiveEasy Plus
x	Type of meter: 2 = ph, 3 = conductivity, 4 = DO (Sensor and software are different)
Where x in the model designation may be any number 0 to 9, denoting SELV/ELVEL secondary circuits or minor mechanical differences.	
Remarks Bemerkungen:	
Bench equipment: (0.6Kg / 200 x 175 x 52mm) Rating: 9Vdc, 0.4W	
Portable equipment: (0.18Kg / 169 x 82 x 36mm) Rating: 6Vdc, 5mA powered with batteries: 4 x AAA / LR3 1.5V or 4 x NiMH 1.2V rechargeable	

EC - DECLARATION OF CONFORMITY

EG-Konformitätserklärung
KD-Nr.: 51302969 D

Doku-Nr.: 20070006

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Bench equipment:

FiveEasy™ Power supply, FiveEasy™ electrode arm
Rubber covers (for electrode arm hole), Shorting plug BNC

Portable equipment:

FiveGo™ electrode clip, Rubber clip covers (for electrode clip slot)
Education Line Battery cover, Wrist strap
FiveGo™ carry bag, IP54 electrode caps
Sealing kit (O-rings for battery case)

pH 4.01 buffer sachets, 30 x 20mL, pH 4.01 buffer solution, 6 x 250mL
pH 7.00 buffer sachets, 30 x 20mL, pH 7.00 buffer solution, 6 x 250mL
pH 9.21 buffer sachets, 30 x 20mL, pH 9.21 buffer solution, 6 x 250mL
pH 10.01 buffer sachets, 30 x 20mL
pH 10.01 buffer solution, 6 x 250mL
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)


LE407, combination pH sensor, plastic, gel electrolyte
LE408, combination pH sensor, plastic, liquid electrolyte, refillable
LE409, combination pH sensor, glass, refillable
LE420, special pH sensor with PTFE junction, refillable
LE427, spear tip design, low maintenance
LE438, 3-in-1 pH sensor, plastic, ATC
LE501, Redox-sensor, glass
LE510, Redox-sensor, plastic
ATC probe, temperature sensor

84 µS/cm standard solution, 500mL
1413 µS/cm standard solution sachets, 30 x 20 mL
1413 µS/cm standard solution, 250 mL
12.88 mS/cm standard solution sachets, 30 x 20 mL
12.88 mS/cm standard solution, 250 mL
LE703, conductivity sensor
LE 740, conductivity senso

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

www.mt.com/phlab

For more information

Mettler-Toledo AG, Analytical

CH-8603 Schwerzenbach, Switzerland

Tel. +41 (0)44 806 77 11

Fax +41 (0)44 806 73 50

www.mt.com

Subject to technical changes.

© Mettler-Toledo AG 06/2012

51710770D

* 5 1 7 1 0 7 7 0 *