
OPTIKA

M I C R O S C O P E S
I T A L Y

Ver. 4.1.0

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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1.0 DESCRIPTION	page 4
2.0 INTRODUCTION	page 7
3.0 UNPACKING AND ASSEMBLY	page 7
4.0 USING THE MICROSCOPE	page 8
5.0 MAINTENANCE	page 10
6.0 ELECTRICS	page 11
7.0 RECOVERY AND RECYCLING	page 11

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

DANGER

The instructions on this symbol to avoid possible severe personal injuries.

WARNING

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

WARNING

Possibility of electric shock.

HOT!

Attention: high temperature surfaces. Avoid direct contact.



NOTE

Technical notes or usage tips.

1.0 DESCRIPTION


1.0 DESCRIPTION

B-500 ERGO MODEL

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily classroom and laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

3.0 UNPACKING AND ASSEMBLY

- 3.1** The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.
- 3.2** The objectives are located in individual protective vials. Remove the objectives from the vials and insert them into the microscope nosepiece in the order from the lowest magnification to the highest, in a clockwise direction from the rear.
- 3.3** Insert the eyepieces into the eye tubes. Place the observation head onto the top of the arm and tighten the lock-screw. If you have a trinocular microscope and want to use the photo port, insert your adaptor in the hole on the top of the observation head and fixate it with the screw on the side.

Connect the mains plug into the socket at the base

WARNING

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

WARNING

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

WARNING

AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

4.0 USING THE MICROSCOPE

4.1 Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

4.2 Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

4.3 Lamp settings

The microscope is fitted with a 20w illuminator (3w high-efficiency white led source for b-500 ergo model) . Insert the plug of the cable into the power socket and turn on the switch on the back of the main body. Turn the brightness adjustment knob to a brightness suitable for observation.

4.4 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. After finishing the interpupillary adjustment, turn the dioptric adjustment ring on the left eye tube to the "0" position at the graduation line on the eyetube.

4.5 Focus and diopter adjustment

Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus with an objective with a low magnification, and then lock the focus-stop knob. Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left dioptric adjustment ring to a sharp and clear image also with the other eye. Turn the tension-adjust-knob to get a suitable tension for the focus system. When the image is in focus, turn the nosepiece to choose the objective you need.

4.6 Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

4.7 Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

4.8 Field diaphragm

Adjust the field diaphragm to align the system and to control stray light.

4.9 Additional filters

The chromatic correction filter (blue) and the monochromator filter (green) can be inserted for definition of specimen parts. The use of the green filter is recommended for phase contrast observations.

4.10 Using the polarizer (optional)

The polarizing set consists of two filters; the analyser and the polarizer. The analyser is placed in the analyser holder and the polarizer on top of the illuminator.

4.11 Video capturing (optional)

Two kinds of observation heads, binocular (standard or ergonomic) or trinocular, are available for the b-500 series. Trinocular version has a dedicated photo-port, while binocular version can receive a camera in either eyepiece tube, via a photo/video adaptor. Before taking a picture or filming video, if necessary, obscure with a dark cloth both the camera viewfinder and the eyepieces and pull out the light path selector lever. Please refer to the adaptor and camera manuals for further details.

4.12 Using the phase contrast set (optional)

The phase contrast set has 4 objectives for phase contrast and one condenser, with 4 rings of phase, and one position for bright field. To use the phase contrast, remove the abbe condenser of the microscope and replace it with the phase contrast condenser.

If, during use, the phase contrast diaphragm rings (anuli) are dislodged, they can be realigned rotating the alignment screws placed on two sides of each ring, using the enclosed tools. In particular:

Select an objective corresponding to the ring to be realigned, and focus on a slide under bright field illumination.

The selector of the condenser has 5 selectable positions:

- | | |
|-----|---|
| 0 | → To use when it is necessary to be used the instrument in bright field |
| 10 | → To use with 10x phase contrast objective |
| 20 | → To use with 20x phase contrast objective |
| 40 | → To use with 40x phase contrast objective |
| 100 | → To use with 100x phase contrast objective |



Insert the phase ring to be realigned.

Remove an eyepiece, and replace it with the centring telescope.

Focus the centring telescope, without touching the microscope focus, in order to focus the phase ring and the phase plate simultaneously. In particular, the phase plate is visible as a dark ring, and the phase ring as a bright ring divided into sectors.

Turning the two phase ring centring screws, position the phase ring (bright ring) so that it is completely contained in the phase plate profile (dark ring). (See figure.)

Remove the telescope, and replace the eyepiece.

4.13 Dark field condenser (optional)

For dark field observation, the Abbe condenser in standard equipment must be replaced with a special condenser for dark field, showed in the picture.

5.0 MAINTENANCE

5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing).
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

5.3 How to change the lamp

- Before changing the lamp, you must pull the plug out of the electrical socket. To avoid being burnt, wait until the lamp has cooled down.
- Hold the main body and tilt the microscope until the bottom base can be seen. Loosen the door-lock screw of the small door on the bottom base to open the lamp housing. The lamp and lamp holder come out with the door.
- Pull out the old lamp. Hold the new lamp without touching it with your bare hands (e.g. using its packaging, or a cloth) and insert it into the holder. Replace/close the lamp holder.

B-500 ERGO: the high-efficiency white LED source in B-500 ERGO is designed to last many years without any performance decrease. Never try to replace it with another kind of LED source. In case of malfunctioning contact OPTIKA service.

5.4 If you need to send the microscope to Optika for maintenance, please use the original packaging.


Mains Power: 90-240 V, 50/60 Hz

Light source: 20 W halogen lamp
3W high-efficiency LED (for B-500 ERGO B-500PL)

7.0 RECOVERY AND RECYCLING

The appliance reports the symbol:

This symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials it contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, Optika Microscopes, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor you bought it from, or to other Optika Microscopes distributors. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96 / EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the “crossed bin” symbol indicated above.

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