
Heating Mantles

HM50C HM100C HM250C
HM500C and HM1000C

Instructions for use
Manuel d'utilisation
Istruzioni per l'uso
Instrucciones de funcionamiento
Bedienungsanweisung

Bibby The logo consists of the word "Scientific" in a bold, sans-serif font, with a stylized circular graphic element to its left.

Version 1.0

HEATING MANTLE

HM50C, HM100C, HM250C, HM500C and
HM1000C

Instructions for use

Introduction

Thank you for purchasing this piece of Stuart equipment. To get the best performance from the equipment, and for your own safety, please read these instructions carefully before use. Before discarding the packaging check that all parts are present and correct.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Stuart, the protection provided may be impaired.

Electrical Installation

THIS INSTRUMENT MUST BE EARTHED

Before connection please ensure that the line supply corresponds to that shown on the rating plate.

Model	Supply requirements	Power
HM50C	230V, 50/60Hz	80W
HM100C	230V, 50/60Hz	80W
HM250C	230V, 50/60Hz	170W
HM500C	230V, 50/60Hz	220W
HM1000C	230V, 50/60Hz	320W

All models are supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other. Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. This involves cutting off the moulded plug, preparing the cable and connecting to the re-wireable plug in accordance with its instructions.

**THIS OPERATION SHOULD ONLY BE
UNDERTAKEN BY A QUALIFIED ELECTRICIAN**

NOTE: Refer to the equipment's rating plate to ensure that the plug and fusing are suitable for the voltage and wattage stated.

The wires in the mains cable are coloured as follows:

BROWN	- LIVE
BLUE	- NEUTRAL
GREEN/YELLOW	- EARTH

The instruments are fitted with an IEC socket at the rear of the instrument for connection of the mains lead. The appropriate mains lead should be connected BEFORE connection to the mains supply.

Should the mains lead need replacement, a cable of 1mm² of harmonised code H05W-F connected to an IEC320 plug should be used.

**IF IN DOUBT CONSULT A QUALIFIED
ELECTRICIAN**

Safety Advice Before Use


- ❖ Never move or carry the unit while it is in use or connected to the mains supply.
- ❖ Do not touch the heating element or any glass vessel whilst in use.
- ❖ Only use flasks corresponding to the correct capacity of your heating mantle.
- ❖ Do not use the heating mantle to heat inflammable liquids.
- ❖ Do not immerse the unit in water or fluids.
- ❖ Always fill the glassware/vessel away from the heating mantle and ensure there is no liquid on the outside of the vessel before placing it in the mantle. If spillage does occur, disconnect from the mains supply and follow the instructions as detailed in Maintenance & Servicing.
- ❖ Do not cover the mantle while in use.
- ❖ Do not block or obstruct the ventilation slots.
- ❖ Do not leave the equipment switched on without a charged flask.

General Description


The Stuart heating mantles have been specifically designed for heating fluids in round bottomed flasks. Heating control is provided by a built-in solid state simmerstat (50, 100, 250 and 500ml mantles) or an energy regulator (1000ml mantle). The maximum element temperature is 450°C. The mantles are housed in an aluminium case with stove paint finish to give good chemical resistance. The units are provided with ventilation slots in the base and around the rim to allow convection cooling. All mantles also have a support rod positioning clamp at the side of the unit.

Controls

Front view

Rear view

Operation

1. Place the unit on a clean, dry, non-combustible, solid work surface with at least 300mm suitable clearance all around from other equipment.
2. With the mains electricity supply switched off, connect the mains lead with IEC plug to the mains input IEC socket. Connect the plug to the mains supply.
3. Place a charged, clean, dry glass vessel of the size indicated on the mantle rating plate label into the unit. Wherever possible, the glass vessel should be supported within the mantle by means of a support rod and clamp.
4. Switch on the mains electricity supply. Adjust the heating control knob to the required setting. Setting 1 is low heat up; setting 10 is maximum heat up.
NOTE: The mains power on' indicator will illuminate. The amber heating element on' indicator will illuminate/pulsate when the heaters are in operation.
5. When the process is complete, switch the heating control knob to the off position. Disconnect the mains electricity supply.
6. Remove the charged vessel. Handle hot charged vessels with care.

Maintenance & Servicing

⚠ WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a clean cloth. DO NOT USE SOLVENTS FOR CLEANING ANY PART OF THIS EQUIPMENT.

Fuse replacement

The mains fuse holder is located at the rear of the unit within the mains input IEC socket. Refer to the table below for the correct fuse type and rating for your model.

Model	Total heating power (W)	Fuse rating (A)
HM50C	75	0.5
HM100C	100	0.5
HM250C	200	1.25
HM500C	280	2.5
HM1000C	380	2.5

Spillage

In the event of spillage or glassware fracture, do not touch the mantle. Disconnect the product from the mains electricity supply and allow the product to cool. Wearing suitable hand protection, carefully remove any pieces of broken glassware. Wipe off all excess liquid from the mantle and surrounding area using an absorbent soft cloth. Drain off any residual fluid retained in the mantle.

In the case of excessive spillage/ flask fracture, invert the mantle and allow it to drain for minimum of one hour. Then proceed with the following drying out procedure. Place the complete mantle, correct way up, in a heated oven at 50°C for a minimum period of 40 hours. Before further use, the mantle must be subjected to electrical safety testing by competent service personnel.

Servicing

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

For a comprehensive list of parts required by service engineers conducting internal repairs, please contact the Service Department, quoting both the model number and serial number:

E-mail: service@bibby-scientific.com

Tel: +44 (0)1785 810475

Fax: +44 (0)1785 810471

Only spare parts supplied or specified by Stuart or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument. If in any doubt, please contact the Service Department of Bibby Scientific Ltd. or the point of sale.

For any other technical enquiries, please contact the Technical Support Department:

E-mail: stuarthelp@bibby-scientific.com

Tel: +44 (0)1785 810433.

Spares and Accessories

Description	Catalogue Number
Mains cable (UK)	AZ9165
Mains cable (EU)	AZ6747
F500mA fuse pack (pack of 10)	AZ9038
F1.25A fuse pack (pack of 10)	AZ9039
F2.5A fuse pack (pack of 10)	AZ9040

Warranty

Stuart warrants this equipment to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim, Stuart will replace any defective component or replace the unit free of charge. This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.


Bibby Scientific Ltd.
Beacon Road,
Stone, Staffordshire
ST15 0SA
United Kingdom
Tel: +44 (0) 1785 812121
Fax: +44 (0) 1785 813748
E-mail: equipment@bibby-scientific.com
www.stuart-equipment.com

Technical Specification

Description	Specification
Flask capacity:	
HM50C	50ml
HM100C	100ml
HM250C	250ml
HM500C	500ml
HM1000C	1000ml
Heating element construction	Thermal insulated element wire stitched into a cartridge construction
Maximum element temperature	450°C. Nominal maximum
Case construction	Aluminium
Thermal insulation	Rockwool mineral insulation
Mains input supply voltage	230V, 50/60Hz
IP rating	IPX0

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
These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these

standards and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Department of Bibby Scientific Ltd.

Declaration of Conformity

Laboratory Heating Mantle, Model HM50C, HM100C, HM250C, HM500C, HM1000C


These products comply with the requirements of the EU Directives listed below:

2004/108/EC EMC Directive.
2006/95/EC Low voltage Directive (LVD)

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326-1:2006 (Electrical Equipment for Measurement, Control and Laboratory use).
EN 61010-1: 2001
(Safety Requirements Electrical Equipment for Measurement, Control and Laboratory use)
EN 61010-2-010: 2003 (Particular Requirements for Laboratory Equipment for Heating of Materials).
EN 61010-2-051 : 2003 (Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for laboratory equipment for mixing and stirring)
FCC 47CFR Part 15B Class B

CE mark affixed '12

Signed:  (Mr S. Marriott)

Date: 1/2/12

Authority: Technical Director
Bibby Scientific Ltd

Bibby Scientific Ltd - Stone - Staffs - ST15 0SA - UK
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

INSPECTION REPORT

MODEL HM50C HM250C HM1000C
HM100C HM500C


ELECTRICAL SAFETY

- | | |
|---------------------|-------------------------------------|
| 1. Earth continuity | <input checked="" type="checkbox"/> |
| 2. Insulation | <input checked="" type="checkbox"/> |
| 3. Flash test | <input checked="" type="checkbox"/> |

FUNCTIONAL

- | | |
|----------------------|-------------------------------------|
| 1. Indicators | <input checked="" type="checkbox"/> |
| 2. Heating function | <input checked="" type="checkbox"/> |
| 3. Visual acceptance | <input checked="" type="checkbox"/> |

QUALITY CONTROL INSPECTOR

UK
Bibby Scientific Ltd.
Beacon Road, Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 810405
e-mail: sales@bibby-scientific.com
www.bibby-scientific.com

France
Bibby Scientific France SAS
ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bsf@bibby-scientific.fr
www.bibby-scientific.com

North & South America
Bibby Scientific US Inc.
7/a Techne Inc.
3 Terri Lane, Suite 10
Burlington, NJ 08016 USA
Toll Free (in NA): 800-225-9243
Tel: +1 609 589 2560
Fax: +1 609-589-2571
e-mail: labproducts@techneusa.com
www.techneusa.com

Middle East
Bibby Scientific Middle East Ltd.
PO Box 27842, Engomi 2433
Nicosia
Cyprus
Tel: +357 22 660 423
Fax: +357 22 660 424
e-mail: sales@bibby-scientificme.com

Italy
Bibby Scientific Italia Srl
Via Alcide de Gasperi 56
20077 Riozzo di Cerro al Lambro
Milano Italia
Tel: +39 (02) 98230679
Fax: +39 (02) 98230211
e-mail: marketing@bibby-scientific.it
www.bibby-scientific.it

Asia
Bibby Scientific - Singapore
Prudential Tower, Level 26
30 Cecil Street
Singapore 049712
Tel: +65 6631 2976
Fax: +44 (0)1785 810405
e-mail: info@bibby-scientific.com
www.bibby-scientific.com