

# Data Sheet (EN)

Translation of the german original

**WELCH**  
by Gardner Denver





Designation, Model	Order no.
<b>Diaphragm pump chemical resistant in an explosion-proof Design</b> Model:  <b>MPC 301 Zp Ex</b> <b>II 2G c IIB T4 X</b> <b>(10°C ≤ Ta ≤ 40°C)</b>	<b>4000481-04</b>

Figure	Description
 <p style="text-align: center;"><i>Inert gas connections</i></p>	<p>The two phase diaphragm pump consists of the pump casing, the drive unit and the drive motor.</p> <p>The pump casing contains the drive unit and two pump heads. Both pump head contains a diaphragm and the work valves. The pump heads are arranged in a horizontally-opposed pattern.</p> <p>The pump heads are driven via an eccentric shaft with a connecting rod.</p>

Technical Data		
Parameter	Unit	MPC 301 Zp Ex
<b>Ex-Marking</b>	-	 II 2G c IIB T4 X (10°C ≤ Ta ≤ 40°C)
<b>Type Examination Certificate no.</b>	<i>more information</i>	IBEXU04ATEXB017 X
<b>Pumping speed</b>	m <sup>3</sup> / h	2.3
DIN 28432 (at speed of 1500 rpm)	l / min	38
<b>Ultimate pressure</b> (at speed of 1500 rpm)	mbar	< 8
<b>Max. Inlet pressure</b>	bar	1
<b>Max. Outlet pressure</b> (absolute pressure)	bar	1.5 *)
<b>Inert flushing of the drive</b>	l / h	20 ±10 % **)
<b>Intake-/ Exhaust pressure</b>	-	Hose nozzle DN8 for hose inside diameter 8 mm
<b>Operating temperature</b>	°C	+ 10 to + 40
<b>Max. Operating gas temperature</b> (measured on the inlet of the vacuum pump)	°C	+ 60 *)
<b>Bearing</b>		maintenance-free
<b>Reference surface sound pressure level</b> DIN EN ISO 2151	dB (A)	< 44
<b>Voltage / Frequency</b>	V, Hz	400, 50
<b>Three-phase a.c. motor – model</b> (without motor protection switch, switch and cable)	-	 CD 71 L – 4 II 2G EEx de IIC T4 PTB 99 ATEX 1051
<b>Power</b>	W	180
<b>Operating mode</b>	-	S 1
<b>Type of protection</b> DIN EN 60529	-	IP 55
<b>Weight</b>	kg	22.9
<b>Dimensions (W/D/H)</b> (without connections)	mm	240 / 300 / 260


\*) If this value is exceeded, the information about the temperature class inside and outside is no longer applicable.


\*\*) The inert gas flushing must be monitored. When exceeding or falling below the tolerances, the pump must be switched off.

# Data Sheet (EN)

Translation of the german original

**WELCH**  
by Gardner Denver

Intake pressure- / Pumping speed - Diagram	Materials of the medium-affecting Pump parts	
	<b>Connection heads / Pump heads with Insert</b>	PTFE carbon-fibre reinforced
	<b>O-Rings</b>	EPDM
	<b>Screw fittings</b>	PVDF
	<b>Valves</b>	PEEK
	<b>Diaphragms</b>	Elastomer + PTFE-Layer
<b>Manifolds</b>	PP	

Circuitry of the Pump heads	Principle of Operation	
	<p>The two pump heads are connected in parallel, two further pumps in series.</p>	<p>Motor, eccentric shaft and connecting rod set the diaphragms in stroke movement. This changes the size of the space between the diaphragms and pump head (pump chamber). Increasing the size of the pump chamber opens the inlet valve while the outlet valve is closed (intake process). Decreasing the size of the pump chamber ejects the gas through the outlet valve. The valves are actuated by the gas being pumped.</p>

## Application

**The Diaphragm pump in an explosion-proof Design is intended to:**

- Aspirating, pumping and compressing neutral and aggressive gases and vapors.
- Generating a vacuum down to an ultimate pressure < 2 mbar.
- Aspirating, pumping and compressing an explosive atmosphere – comprising air and combustible gases, vapors and mists in any mix ratio – from zone 1 areas at risk of explosion (device category 2 according to ATEX).
- Installation and operation of the diaphragm pumps in zone 1 areas at risk of explosion (device category 2 according to ATEX).
- The permissible temperature classes and explosion groups of the atmospheres – comprising air and combustible gases, vapors and mists in any mix ratio – that are to be pumped or which are ambient are determined according to the specifications of standard EN 13 463-1. The diaphragm pumps mentioned here have been assigned to T4 IIB.

**Substances that tend to decompose spontaneously, such as acetylene C<sub>2</sub>H<sub>2</sub>, carbon disulfide CS<sub>2</sub> and explosives, lie outside the scope of the application of Atex-Directive 2014/34/EU.**



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax: 01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**Use the above details to contact us if this literature doesn't answer all your questions.**

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

