


Guide to Operations

***innova*[™] 43/43R**

Incubator Shakers

MANUAL NO: M1320-0050
Revision E
May 28, 2008



NEW BRUNSWICK SCIENTIFIC CO., INC.

BOX 4005 • 44 TALMADGE ROAD • EDISON, NJ 08818-4005

Telephone: 1-732-287-1200 • 1-800-631-5417

Fax: 732-287-4222 • Telex: 4753012 NBSCO

Internet: <http://www.nbsc.com> • E-mail: bioinfo@nbsc.com

INTERNATIONAL OFFICES:

BELGIUM

New Brunswick Scientific NV-SA
Stationsstraat 180/4
3110 Rotselaar
België/Belgique
Tel: 32 (0)16 56 28 31
Fax: 32 (0)16 57 27 53
E-mail: sales@nbsnv-sa.be

CHINA

New Brunswick Scientific Co., Inc.
A903, 904 Yin Hai Building
No. 250, Cao Xi Road
Shanghai 200235, P.R. China
Tel: 86 21 6484 5955 or 5966
Fax: 86 21 6484 5933
E-mail: nbschc@online.sh.cn

FRANCE

New Brunswick Scientific SARL
12-14, Rond Point des Champs Elysées
75008 Paris
France
Tel: 33 (0)1 5353 1511
Fax: 33 (0)1 5353 1557
E-mail: sales@nbsnarl.fr

GERMANY

New Brunswick Scientific GmbH
In Der Au 14
D-72622 Nürtingen
Deutschland
Tel: 49 (0)7022 932490
Fax: 49 (0)7022 32486
E-mail: sales@nbsgmbh.de

THE NETHERLANDS

New Brunswick Scientific BV
Kerkenbos 1101, 6546 BC Nijmegen
P.O Box 6826, 6503 GH Nijmegen
Nederland
Tel: 31 (0)24 3717 600
Fax: 31 (0)24 3717 640
E-mail: sales@nbsbv.nl

UNITED KINGDOM

New Brunswick Scientific (UK) Ltd.
17 Alban Park
St. Albans, Herts. AL4 0JJ
United Kingdom
Tel: 44 (0)1727 853855 or 0800 581331
Fax: 44 (0)1727 835666
E-mail: bioinfo@nbsuk.co.uk
Web: www.nbsuk.co.uk

**CAUTION!**

This equipment *must* be operated as described in this manual. If operational guidelines are not followed, equipment damage and personal injury *can* occur.

Please read the entire User's Guide before attempting to use this unit. Do not use this equipment in a hazardous atmosphere or with hazardous materials for which the equipment was not designed.

New Brunswick Scientific Co., Inc. (NBS) is not responsible for any damage to this equipment that may result from the use of an accessory not manufactured by NBS.

Copyright Notice

New Brunswick Scientific Company, Inc.
Box 4005 44 Talmadge Road
Edison, New Jersey 08818-4005

© Copyright 2007 New Brunswick Scientific Co., Inc.
All Rights Reserved.

Reproduction, adaptation, or translation without prior written permission from New Brunswick Scientific is prohibited.

Disclaimer Notice

New Brunswick Scientific Co., Inc. reserves the right to change information in this document without notice. Updates to information in this document reflect our commitment to continuing product development and improvement.

Manual Conventions



NOTE:

Notes contain essential information that deserves special attention.



CAUTION!

Caution messages appear before procedures which, if caution is not observed, could result in damage to the equipment.



WARNING!

Warning messages alert you to specific procedures or practices which, if not followed correctly, could result in serious personal injury.

Bold

Text in boldface type emphasizes key words or phrases.



This particular *Warning* message, whether found in the manual or on the unit, means **HOT SURFACE**—and therefore represents a potential danger to touch.



CRUSH WARNING!

Crush Warning messages alert you to specific procedures or practices regarding heavy objects which, if not followed correctly, could result in serious personal injury .



Biohazard-related messages.



WARRANTY

Every Instrument manufactured by the
New Brunswick Scientific Co., Inc. is warranted
to be free from defects in material and workmanship.

This apparatus is warranted for three (3) years against faulty components and assembly. During the first two (2) years of this warranty period, New Brunswick Scientific Co., Inc. (NBS) will also cover labor charges associated with the repair of this instrument. Our obligation under this warranty is limited to repairing or replacing the instrument or part thereof which shall, within three (3) years after the date of shipment, prove to be defective after our examination.

This warranty does not extend to any NBS products that have been subjected to misuse, neglect, accident or improper installation or application; nor shall it extend to products that have been repaired or altered outside the NBS factory without prior authorization from New Brunswick Scientific Co., Inc.

TABLE OF CONTENTS

1	OVERVIEW.....	11
2	INSPECTION & UNPACKING OF EQUIPMENT	13
2.1	INSPECTION OF BOXES	13
2.2	PACKING LIST VERIFICATION	13
2.3	UNPACKING OF EQUIPMENT	13
2.4	OUT OF BOX CONCERNS	13
2.5	WARRANTY REGISTRATION	13
3	PREPARING THE LOCATION.....	14
3.1	PHYSICAL LOCATION	14
3.2	ENVIRONMENT	14
3.3	ELECTRICAL REQUIREMENTS	14
3.4	SPACE REQUIREMENTS	15
3.5	NARROW DOORWAYS	15
4	INSTALLATION	17
4.1	LEVELING THE SHAKER	17
5	FEATURES	19
5.1	CONTROLS	19
5.2	LCD DISPLAY	20
5.3	CHANGING SCREENS	20
5.4	DISPLAY ICONS	21
5.5	ALARMS	22
5.6	LID.....	23
5.7	SPILL PAN/WATER RESERVOIR	23
5.8	SOFTWARE INTERFACES.....	24
5.9	INTERIOR LIGHT(S)	24
5.10	HEATER	24
5.11	REFRIGERATION (43R ONLY).....	24
5.12	SERVICE ACCESSIBILITY	25
6	GETTING STARTED	26
6.1	PLATFORM ASSEMBLIES	26
6.2	INSTALLATION OF PLATFORM	26
6.3	INSTALLATION OF CLAMPS	26
6.4	ELECTRICAL CONNECTIONS	27
6.5	OPTIONAL GASSING MANIFOLD KIT	28
6.6	OPTIONAL GERMICIDAL UV LAMP	30
6.7	OPTIONAL PHOTOSYNTHETIC LAMPS	30
6.8	OPTIONAL HUMIDITY MONITOR	31
6.9	FILLING THE WATER RESERVOIR	33

6.10	DRAINING THE WATER RESERVOIR.....	34
7	OPERATION	35
7.1	SAFETY PRECAUTIONS	35
7.2	OPENING THE LID	35
7.3	STARTING THE SHAKER.....	36
7.4	USING THE LCD SCREENS	36
7.4.1	<i>Display Screen</i>	36
7.4.2	<i>Summary Screen</i>	39
7.4.3	<i>Setup Screen</i>	40
7.4.4	<i>Lamps Screen</i>	41
7.4.5	<i>RS232 Screen</i>	43
7.4.6	<i>Calibrate Screen</i>	44
7.4.7	<i>Programs Screen</i>	45
7.5	PROGRAMMING THE SHAKER	45
7.5.1	<i>Timer Only</i>	45
7.5.2	<i>Programmed Steps</i>	46
7.5.3	<i>Creating a Program</i>	46
7.5.4	<i>Editing a Program</i>	49
7.5.5	<i>Running a Program</i>	50
7.6	MUTING THE AUDIBLE ALARM	50
7.7	TEMPERATURE OFFSET CALIBRATION	50
7.7.1	<i>Calculating the Offset Value</i>	51
7.7.2	<i>Setting the Offset</i>	51
7.8	USING CALSPEED	52
7.9	POWER INTERRUPTION.....	52
8	MAINTENANCE	54
8.1	ROUTINE MAINTENANCE	54
8.2	CLEANING EXTERNAL & INTERNAL SURFACES.....	54
8.3	BIOHAZARD DECONTAMINATION.....	54
9	SERVICE & ACCESSORIES	56
9.1	TROUBLESHOOTING	56
9.2	PRODUCT RETURNS.....	58
9.3	OPENING THE SERVICE COMPARTMENT	59
9.4	FUSE REPLACEMENT	59
9.5	REPLACING CABINET LIGHTS.....	60
9.6	REPLACING OPTIONAL PHOTOSYNTHETIC LAMPS.....	61
9.7	ADJUSTING LID TENSION	62
9.8	REPLACING OPTIONAL GERMICIDAL LAMP.....	63
9.9	BELT REPLACEMENT OR ADJUSTMENT	64
9.10	REPLACEMENT PARTS.....	67
9.11	ACCESSORIES	67
9.11.1	<i>Platforms</i>	67
9.11.2	<i>Flask Clamps for Platforms</i>	68

9.11.3	<i>Replacement Clamp Hardware Kits</i>	69
9.11.4	<i>Test Tube Racks & Other Accessories</i>	69
9.11.5	<i>Optional Gas Manifold Kit</i>	70
9.11.6	<i>Optional Photosynthetic Lamps</i>	70
9.11.7	<i>Optional UV Germicidal Lamp</i>	70
9.11.8	<i>Optional Humidity Monitor</i>	70
9.11.9	<i>Optional Utility Basket</i>	70
9.11.10	<i>Optional Remote Alarm</i>	71
10	DRAWINGS, GRAPHS & TABLES	72
10.1	SCHEMATICS	72
10.2	LOAD & SPEED GRAPHS	73
10.3	LIST OF DRAWINGS	75
10.4	LIST OF TABLES	76
11	SPECIFICATIONS	77
11.1	CERTIFICATIONS	78
11.1.1	<i>CE Declaration of Conformity</i>	79
12	APPENDIX: REMOTE PROGRAMMING	80
12.1	SETTING UP WITH HYPERTERMINAL	80
12.2	OVERVIEW OF COMMAND SETS	81
12.3	INDEX TO COMMAND CODES	82
12.4	SET COMMANDS	83
12.5	PROFILE CONTROL COMMANDS	83
12.6	REPORT REQUEST COMMANDS	85
12.7	SET/GET DATE & TIME	85
13	APPENDIX: PRODUCT RETURNS	87
13.1	RETURN PROCEDURE	87
13.2	RETURN AUTHORIZATION AND DECONTAMINATION CERTIFICATE	87
14	INDEX	89

1 OVERVIEW

The Innova 43/43R Incubator Shakers are large-capacity orbital shakers that utilize an eccentric counterbalanced drive mechanism. They provide horizontal plane rotary motion in either a 1-inch (2.54 cm) or a 2-inch (5 cm) diameter circular orbit, depending on the model. A Proportional/Integral (PI) microprocessor controller with instantaneous digital feedback controls the speed over the entire speed range.

The Innova 43R provides temperature control from 20° below ambient (as low as 4°C) to 80°C, and the Innova 43 from 5° above ambient to 80°C. Naturally, both these ranges depend on relative humidity and other ambient factors, as well as the options installed in the unit. Ambient temperature is measured at one meter from the exterior of the unit.

Erlenmeyer flasks (up to 6 liters in size) as well as a wide variety of tubes and plates can be accommodated using the New Brunswick Scientific shaker accessories described in Section 9.11. These are easily accessed on the platform.

To facilitate moving the unit, each shaker is equipped with casters that roll freely when the feet are raised.

The Innova 43/43R may be operated in the following ways:

- **Continuously:** at a set speed and temperature, until user intervention.
- **In a timed mode:** run at a set speed, time and temperature for a period of up to 99.9 hours, after which the shaker automatically shuts off.
- **Via the shaker's programmable controller:** run through multiple temperature and speed changes for an extended period of time.
- **Via computer through an RS-232 interface.**

See Section 7 for more details on these various modes of operation.

For safe operation, the Innova 43/43R shakers are designed with a safety switch that automatically stops the shaker mechanism when the lid is opened.

The Innova 43/43R is equipped with visual and/or audible alarms that alert the user to the following conditions:

- The end of a timed run
- Deviations from speed setpoint
- Deviations from temperature setpoint
- Power failure
- Lid open

To accommodate customer needs, a wide variety of platforms can be used with the Innova 43/43R:

- Universal platforms are the most flexible, providing hole patterns for flask clamps, test tube racks and other accessories.
- Dedicated platforms are supplied with flask clamps attached; they are designed solely and expressly for this purpose.
- Test tube racks, microplate holders, and test tube rack holders are also available (a universal platform is needed for all test tube racks and holders).

For further information on these accessories, see Section 9.11.

We also recommend that you consult the load and speed graphs provided in Section 10.2.

2 INSPECTION & UNPACKING OF EQUIPMENT

2.1 *Inspection of Boxes*

After you receive your order from New Brunswick Scientific, inspect the boxes carefully for any damage that may have occurred during shipping. Report any damage immediately to the carrier and to your local NBS Customer Service Department.

2.2 *Packing List Verification*

Verify against your NBS packing list that you have received all of the correct materials.

2.3 *Unpacking of Equipment*



CRUSH WARNING!

**Do not attempt to lift the Innova 43/43R by hand.
Always use a lifter or other suitable equipment when raising or handling the unit.**

To unpack your Innova 43/43R, you will need the following tools:

- Claw hammer
- Scissors to cut nylon strapping
- Tool to remove 3-inch (7.6 cm) metal staples

Save all packing materials and this user's guide/operating manual.

2.4 *Out of Box Concerns*

If any part of your order was damaged during shipping, is missing pieces, or fails to operate properly, please fill out the *Customer Satisfaction Form 6300*, included with your equipment, and return it by fax.

2.5 *Warranty Registration*

Please complete and return your warranty card or register electronically at our Website:
www.nbsc.com

3 PREPARING THE LOCATION

3.1 *Physical Location*

It is essential that the instrument be situated in an area where there is sufficient space for the shaker to clear walls and potential obstructions during operation. The surface on which the unit is placed must be smooth, level, and able to support the shaker under full load operating conditions.



CRUSH WARNING!

Do not attempt to lift the Innova 43/43R by hand.

Always use a lifter or other suitable equipment when raising or handling the unit.

NOTE:

The Innova 43/43R has casters and leveling feet. When you need to move the unit across the floor, make sure the feet are raised enough so that the casters can roll freely without scraping the floor with the feet. *At all other times, make sure the feet are down, to avoid any unintended rolling.*

3.2 *Environment*

The shaker is designed to operate optimally in the following ambient conditions:

- 10° to 35°C
- 20 to 80% Relative Humidity (non-condensing)

3.3 *Electrical Requirements*

The Innova 43/43R can be equipped to run on:

- 100 Volts, 50 or 60 Hz, 1500 VA maximum
- 120 Volts, 60 Hz, 1500 VA maximum
- 230 Volts, 50 Hz, 1500 VA maximum


Check your shaker's Electrical Specification Plate (located on the back of the unit) to be sure the electrical requirements of your unit match the output of your electric supply. If they do not match, contact your NBS representative.

3.4 Space Requirements

Allow at least four inches (10 cm) around the shaker for proper ventilation and for access to Power Switch and RS-232 port accessibility on the right side. Also allow enough room above the shaker for the lid to be fully open (*see Figure 1*).

Be sure to keep the power plug and power outlet easily accessible to facilitate unplugging the unit, as needed.

Figure 1: Innova 43/43R Front View



NOTE:

Be sure to allow at least four inches (10.5 cm) around shaker for ventilation, access to power cord (rear panel), and access to power switch and RS-232 port (right side).

3.5 Narrow Doorways

If you need to move the Innova 43/43R through a doorway that is less than 32.5 inches (82.5 cm) wide, remove the front window cover to reduce the depth of the unit to 30.3 inches (77 cm):

1. Carefully remove the five screws on the underside of the plastic window cover, retaining them for reuse.

2. *Being very careful not to hit the control knob or button*, slide the window cover upward to remove it from the spring clips on the front of the unit.

When the shaker is in the desired location, reinstall the window cover by reversing the procedure.

4 INSTALLATION



CRUSH WARNING!

Do not attempt to lift the Innova 43/43R by hand.
Always use a lifter or other suitable equipment when raising or handling the unit.

4.1 Leveling the Shaker


NOTE:

The Innova 43/43R has casters and leveling feet. When you need to move the unit across the floor, make sure the feet are raised enough so that the casters can roll freely without scraping the floor with the feet. *At all other times, make sure the feet are down, to avoid any unintended rolling.*

Make sure that the shaker is placed on a level surface, in its intended location. Lower all four adjustable feet as described below until they are solidly on the surface. If the shaker is not level, adjust the feet as needed to achieve leveling:

1. Immobilize the top lock nut against the unit with one wrench **whenever you adjust the foot**, to keep the threaded stud from falling out (*see Figure 2*).

Figure 2: Adjustable Foot (as shipped)




2. With a second wrench against the flats of the threaded stud, just above the foot. Rotate clockwise to lower the foot or counter-clockwise to raise the foot.

3. Place a level on the top of the unit. If necessary, make further adjustments by repeating all steps until the unit is level.
4. After installing a platform (*see Section 6*), fully load the shaker and do a test run at normal speed (*see Section 7*). Make additional leveling adjustments if necessary.

5 FEATURES

5.1 Controls

Figure 3: Front Panel (detail)



- **START/STOP SWITCH** This switch is used to start or stop the shaker. It will also activate the timer when a timed run is desired. If the unit is stopped and restarted, the timer automatically returns to the beginning of a run.
- **CONTROL KNOB** This knob is multifunctional. It is used to change screens, and to select and change operating conditions.
- **RS-232 PORT** See Section 5.8 for details.
- **POWER SWITCH** This rocker switch is a circuit breaker that turns power on and off to the entire Innova 43/43R.


 **NOTE:**

In addition to the power switch, the power cord is also used to conduct power or to break the power circuit to the shaker. Whenever power to the shaker may be a hazard (during cleaning, maintenance or service work), be sure to disconnect the power cord from the electrical outlet.

5.2 LCD Display

When the unit is powered up, using the Power Switch located on the right side of the machine (see Figure 3), the initial display screen will remain for a short time while the system boots up. Then the main screen, called **DISP** for Display, appears (see Figure 4). This screen will indicate the same parameters that were in effect when the power was turned off.

Figure 4: Display Screen



- For an explanation of the icons on the display screen, see Section 5.4.
- Turning the Control Knob will highlight functions and/or values that can be changed.
- Whenever you turn the Control Knob, the chamber light will go on (for details, see Section 5.9).

For more information on working in the Display Screen, see Section 7.4.1.

5.3 Changing Screens

You can change screens displayed by highlighting the screen name field in the lower left corner, pressing the Control Knob in until it clicks, rotating the Knob left or right (which also makes clicking sounds) to the desired screen and clicking the Control Knob in again. Table 1 describes the various screens:

Table 1: Screens

Screen Name	Meaning	Features/Modes
DISP	Display	Shows two user-selectable parameters* and actual values.
SUMM	Summary	Shows all parameters*, setpoints and actual values.
SET	Set-Up	Set day of week, set time, enable or mute alarm, lock or unlock operating parameters
LAMP	Lamps	Internal chamber light: ON (always on); OFF (always off); AUTO (default mode), light goes on and stays on when lid is open, shuts off 15 seconds after lid is closed, and goes on for 15 seconds when Control Knob is moved. <input type="checkbox"/> Photosynthetic lights (GRO): ON, OFF, NONE <input type="checkbox"/> UV Light (UV): ON, OFF, NONE
COMM	Communication (RS-232)	SET: set baud rate OFF disables RS-232 MONITOR: PC commands shaker to read setpoints and actual values on a schedule determined by PC software. Parameters are unlocked and can be changed by program or manually. SLAVE: PC controls shaker and logs data. TALK: Shaker sends setpoint and actual data to PC at one-minute intervals.
CAL	Calibration	Allows user to enter a temperature offset. Self-calibrates the speed sensor.
PROG	Program	Allows user to set up 1- 4 programs, each with 1-15 steps.

* see Table 2, Display Screen Parameters Optional ¹Not installed

Table 2: Display Screen Parameters

Parameter Name	Meaning
RPM	Shaking speed, in revolutions/minute
°C	Chamber temperature, in degrees Celsius
HRS	Programmed time remaining, in hours
<input type="checkbox"/> %RH	Relative Humidity, in percent.
<input type="checkbox"/> UV	Status of Ultraviolet germicidal lamp
<input type="checkbox"/> GRO	Status of Photosynthetic growth lamps

Optional

5.4 Display Icons

There are six icons that help identify operating conditions. Five of these icons are located at the bottom of the display, and they are visible, when applicable, in any screen you are viewing, except the Program subscreens. The sixth appears, when applicable, next to °C whenever the temperature is onscreen.

Table 3: Display Icons

<i>Icon</i>	<i>Explanation</i>
	This icon appears when audible alarms are muted.
	This icon appears when the possibility to make manual or programmed parameter changes is disabled (locked). This is controlled by settings on the SET screen. See Section 7.4.3
	This icon appears when the shaker lid is open.
	This icon appears when the Heater is on. See Section 5.10.
	This icon appears when a user-defined Program is running. See Section 7.5.5.
	This icon appears to the right of °C if the Temperature Offset feature is being used. See Section 7.7.

5.5 Alarms

If an alarm condition exists, the field in the lower right corner will alternate the Day and Time with characters indicating the nature of the alarm condition, accompanied by an audible alarm (unless muted):

Table 4: Alarms

<i>Indication</i>	<i>Description</i>
TEMP	The temperature deviates more than $\pm 1^{\circ}\text{C}$ from setpoint after achieving control temperature range. After lid is opened, alarm will be disabled for 5 minutes while chamber recovers to setpoint.
SPEED	The speed deviates more than ± 5 RPM from setpoint after achieving operating speed setpoint. After lid is opened, alarm will be disabled for 5 minutes while chamber recovers to setpoint.
POWER	Indicates unit is powering up (both at normal power-up and after power interruption); will flash until the Control Knob is moved.
HRS	Indicates when timed run is completed.

5.6 Lid

Access to the chamber of the Innova 43/43R is through the top-opening hinged lid. To open the lid, firmly press the pedal or lift up on the front handle (*see Figure 1*). The lid stands on its hinges until you choose to close it.

When the lid opens, the following will happen:

- Heater turns off
- Shaker stops
- Interior light goes on and, when it is in **AUTO** mode, will remain on for 15 seconds after the lid is closed
- UV germicidal lamp (if so equipped) turns off

5.7 Spill Pan/Water Reservoir

The Innova 43/43R is equipped with a spill pan and spill cover to protect the drive mechanism in the case of accidental spills and/or broken glassware. This pan can also be used as a water reservoir to humidify the chamber and to reduce evaporation. An optional factory-installed humidity monitor is also available.

The reservoir can be drained through the quick-connect valve on the front of the unit.

5.8 Software Interfaces

The RS-232 port is next to the Power Switch on the right side of the control panel (*see Figure 3*). It can be used to interface a computer to the shaker for control of operating conditions or data logging applications (*see Section 12*).

The customer is responsible for securing the proper driver to interface with the RS-232.

5.9 Interior Light(s)

When the **LAMP** screen is in its default **AUTO** mode, the interior (“chamber”) light is activated for 15 seconds whenever you turn the Control Knob. It will automatically shut off after 15 seconds of Control Knob inactivity.

The chamber light will also go on when the lid is open.

In addition, you can set the chamber light to be continuously **ON** or **OFF** by selecting either mode in the **LAMP** screen.

There are two additional light options for refrigerated units only: interior photosynthetic growth lamps (*see Section 6.7*) and a germicidal UV lamp located outside the chamber but in the airflow path (*see Section 6.6*).

5.10 Heater

The chamber temperature is sensed by a 1000 ohm platinum RTD. A 650W heater is controlled using pulse width modulation on a 2.5-second duty cycle. This cycle time is fast enough to prevent noticeable changes in air temperature due to the cycling.

Whenever the heater is on, the Heater On icon  will appear in the display. The heater automatically stops running when the lid is opened.

NOTE:

Operating the shaker for sustained periods at temperatures above 60°C may result in a shorter fan lifespan.

5.11 Refrigeration (43R Only)

The refrigeration system in the Innova 43R is a variable-capacity system carefully designed with self-checks to maintain the setpoint, to balance pressure within the system, and to prevent freezing on the evaporator surface.

When the shaker is powered up, there is a four-minute time delay prior to compressor start-up.

5.12 Service Accessibility

In the unlikely event that your Innova 43/43R should need service, all electronic boards, refrigeration, heating components are easily accessed, **by an authorized service technician**, behind the right side panel of the shaker.

6 GETTING STARTED


6.1 Platform Assemblies

The Innova 43/43R can be used with a variety of NBS platforms that will accept a wide range of clamps for flasks, test tubes, etc. A platform, which is **required for operation**, is a separate item, not included with the shaker assembly. Refer to Section 9.11 for details on available platforms and platform accessories.

6.2 Installation of Platform

Prior to use, a spill cover and platform must be installed on the unit. The shaker is shipped with four Allen head platform screws installed in the bearing housing (see Figure 5 below, which also shows the spill cover that you must install).

Figure 5: Installing Platform & Spill Cover



6.3 Installation of Clamps

Flask clamps purchased for use with universal platforms (see Section 9.11.2) require installation. Clamps are installed by securing the base of the clamp to the platform with the correct type and number of screws. All clamps are shipped complete with hardware.

Clamps for 2-, 2.8-, 4- and 6-liter flasks are shipped with an additional girdle to keep the flasks in place (see Figure 6). To install these double girdle clamps:

1. Place clamp on platform, secure in place with correct type of screws: #S2116-3051, 10-24 x 5/16-inch flat Phillips head, provided with the clamps (*see Figure 6a*).
2. Place the loose girdle around the upper portion of clamp body so that it is held in place by the legs of the clamp.
3. Insert the flask into the clamp, and push the girdle down so the rubber tubes are in contact with the platform and the flask.

Figure 6: Double Girdle Clamp Installation





Figure 6a: Clamp Fastener



6.4 Electrical Connections

Before making electrical connections, verify that the power source voltage matches the voltage on the electrical specification plate and that the **ON/OFF** switch is in the **OFF** position. The electrical specification plate is located on the rear panel of the unit near the power connector.

Connect the power cord to the power connector, then connect the other end to a suitable, grounded receptacle. Make sure there is enough clearance to disconnect the plug whenever necessary.

Figure 7: Rear Panel

6.5 *Optional Gassing Manifold Kit*

This option is factory-installed. The manifold delivers gas into the chamber via as many as 12 ports (see *Figures 8 & 9 on the following page*). The manifold can be adapted to the desired tubing configuration by adding or subtracting ports or by temporarily clamping off unused tubes. You may elect to use splitters (barbed Y-connectors) after the manifold to increase the number of flasks you can serve.

You determine the appropriate gas flow rate using a pressure regulator (which you supply) on the gas supply.

 **NOTE:**

Gas supply should be regulated and it must never exceed 15 PSIG.

50 feet of 1/16-inch (ID) sterilizable silicone tubing is supplied with the kit. Filters may be needed: 0.22 μ syringe filters (which you supply) can be fitted to the individual manifold ports to maintain a sterile barrier.

Figure 8: Gas Monitor (optional)





Figure 9: Gas Fitting (with optional Gassing Manifold)



SAFETY PRECAUTIONS!

- Never use the gas manifold with flammable gases.
- Never exceed 15 PSIG inlet pressure to the manifold.
- Make sure free lengths of tubing are secured by the hook and loop straps to keep them from contact with flasks or other moving parts.

6.6 Optional Germicidal UV Lamp




WARNING!

Never try to operate the UV Germicidal Lamp while the shaker lid is open.

This option, a germicidal ultraviolet lamp, is placed inside the Service Compartment (see *Figure 10*), outside the chamber to help reduce the risk of contamination. The lamp is identified on the display screen as **UV**.

Figure 10: Germicidal UV Lamps (optional)




The UV germicidal lamp is factory-installed and available on refrigerated units only. Operation instructions are provided with this option (see Section 9.11.7).

6.7 Optional Photosynthetic Lamps


This factory-installed option, available on refrigerated units only, provides nine photosynthetic growth lamps inside the chamber (see *Figure 11 on the following page*). They can be turned on and off manually, by the easily set programmable timer, or by computer. These lamps are identified on the display screen as **GRO**. The replacement bulb part number is P0300-0221.

The recommended operating temperature when using this option is 15 - 37°C; the maximum operating temperature is 70°C.

Figure 11: Photosynthetic Lamps (optional)

6.8 *Optional Humidity Monitor*

This optional factory-installed accessory allows you to monitor relative humidity levels in the chamber throughout your run.

Figure 12a: Humidity Monitor (optional)

When the humidity sensor is present in the chamber, **maximum shaker temperature is automatically limited to 60°C.**

Specifications:

Humidity Operating Range:	0 to 100% RH
Temperature Operating Range:	4 to 60°C
Accuracy:	See <i>Operating Range & Error Graph</i> , below
Hysteresis:	$\pm 1.5\%$ RH
Recovery Time after Saturation:	10 seconds
Long Term Stability:	0.5% RH per year

Figure 12b: RH Sensor Operating Range & Error Graph

Operating Range and Humidity Error




Figure 12c: RH Sensor Connection





Figure 12d: Humidity Sensor Mounting



CAUTION!

Operating the Innova 43/43R at temperatures above 60°C can result in permanent damage to the sensor.

Water can be added to the spill pan in order to elevate the humidity level in the chamber. See Section 6.9 for instructions.

6.9 Filling the Water Reservoir

If you choose to use the spill pan as a water reservoir to reduce evaporation and to raise the humidity level in the chamber:

1. Open the lid and temporarily remove the platform.
2. Make sure that the drain check valve is closed.



NOTE:

As you add water, do not allow the water to splash or flow into the recessed middle of the pan, where the bearing housing assembly is mounted. Pour water very slowly into the shallow area beyond the edge of the spill cover to protect the bearing housing.

3. Accessing the pan/reservoir from the left, the right or in front of the spill cover, **slowly** fill the reservoir with no more than 3 liters of distilled water. A long, narrow watering can or a flexible hose will make it easier to access the pan while protecting the bearing housing from accidental overflow.

At a 37°C setpoint, the chamber loses approximately 500 mL/hr from the pan. At a 25°C setpoint, and placed in a 25°C room, the chamber achieves a relative humidity equilibrium that is approximately 15% above the ambient humidity.

6.10 *Draining the Water Reservoir*

To drain water from the water reservoir/spill pan:

1. Attach the quick-connect drain fitting, direct it to a container or drain, and allow the water to gravity drain.
2. When the reservoir is empty, detach the fitting.

The spill pan/water reservoir drain is located in the front, on the left, under the humidity tray.

7 OPERATION

7.1 Safety Precautions

Before operating the shaker, verify that anyone involved with its operation has been instructed in both general safety practices for laboratories and specific safety practices for this apparatus.



The user is also responsible for following local guidelines for handling hazardous waste and biohazardous materials that may be generated from the use of this equipment.



If service should be required on a unit that is going to be returned to an NBS facility, it must be completely decontaminated and cleaned prior to its return, and a Returned Material Safety Sheet must be filled out to certify that you have complied. See Section 13, *Product Returns*.



It is the responsibility of the user to carry out appropriate decontamination procedures if hazardous material is spilled on or inside the equipment. Before using any cleaning or decontamination method other than those suggested by the manufacturer, users should check with New Brunswick Scientific that the proposed method would not damage the equipment.

This equipment is **not** “explosion-proof” and should never be used with flammable substances or used to grow organisms that produce flammable by-products.



CAUTION!

To prevent damage to the shaker and its contents, never run the shaker without a platform.

7.2 Opening the Lid

Open the lid by firmly pressing the foot pedal or lifting up on the front handle. The lid will stand open on its hinges until you swing it down to the closed position.

7.3 Starting the Shaker

To initially start the shaker, close the lid and turn the power switch (located on the righthand side of the control panel) to the **ON** position. The display will come on (first showing only **New Brunswick Scientific**, then briefly displaying the model number, **43** or **43R**, and the stroke, **1 inch** or **2 inch**, and then quickly moving into the Display screen), and the audible alarm will sound. To mute it, see Section 7.6.

When the shaker begins to operate, the LCD display will track the speed as it accelerates to the last entered setpoint. The shaking action may be started or stopped by pressing the Start/Stop button on the front panel.

 **NOTE:**

The shaker will not operate if the lid is open. This is indicated by the “lid open” symbol appearing in the bottom line of the display (see Figure 4 or Table 3).

7.4 Using the LCD Screens

7.4.1 Display Screen

When you turn the power on, this is the first screen to appear after the company title screen. The default display parameters are temperature (**°C**) and shaking speed (**RPM**).

You can change the displayed parameters. **To replace a parameter:**

1. Using the Control Knob, highlight the parameter that you wish to replace. For this example, we will replace **RPM** (see *Figure 13*).

 **NOTE:**


If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

 **NOTE:**

“UV” (UV Germicidal lamp) and “GRO” (photosynthetic growth lamps) will appear in this screen with the word “NONE” if your shaker is not equipped with these optional features.

Figure 13: Changing Display Parameter

2. Click the Control Knob in. **RPM** will flash.
3. Turn the Knob until the desired parameter appears in the highlighted field. For this example, we will select **HRS**.
4. Click the Knob in, to set and save the parameter (*see Figure 14*).

Figure 14: Changed Display Parameter

You can also use this screen to verify a setpoint, even though the values displayed here are actual (current) values. **To view a setpoint:**


1. Use the Control Knob to highlight the value (in this example, we will view the temperature setpoint, so we will highlight the current **°C**, which is **23**).
2. Click the Knob in to display the current setpoint, which will flash.

At this point you can modify the setpoint or click the Knob in again to return to the normal display, which will be the actual temperature.

To modify a setpoint in this screen:

1. Use the Control Knob to highlight the current value (we will continue to use the temperature as our example, so we will select **23**).
2. Click the Knob in to display the current setpoint (in this example, **38.5**—see *Figure 15*), which will flash.

Figure 15: Changing Setpoint



3. Turn or spin the Knob to reset the setpoint (in this example, turn the Knob left to decrease the setpoint to **37.0**).

 **NOTE:**

If you turn the Control Knob slowly, one click left or right will change the setpoint by an increment of one tenth of a degree Celsius (0.1°C). If you spin the Knob, the value will change by larger increments.

4. Click the Knob in to set and save this new setpoint.

 **NOTE:**

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

5. The display will automatically return to the actual value.

To move out of this screen and into the next:

1. Use the Control Knob to highlight **DISP**, then click the Knob in. **DISP** begins to flash.
2. Turn the Knob to the right until the next screen, Summary (**SUMM**), appears. If you turn too far and enter another screen, just turn the Knob back to the left to recapture the **SUMM** screen.
3. Click the Knob in to select the screen and to work in it.

7.4.2 Summary Screen

In this screen (see *Figure 16*), you can see both the current **ACTUAL** readings and the **SET**points for shaking speed (**RPM**), chamber temperature (**°C**), elapsed time in a programmed run (**HRS**) and, if you are using the optional Humidity Monitor (see *Section 9.11.8*), the percentage of relative humidity (**%RH**).

Figure 16: Summary Screen

PARAM	ACTUAL	SET
RPM	100	100
°C	45.1	45.0
HRS	0.0	0.0
%RH	50.0	N/A
SUMM	Th 16:18	

 **NOTE:**

The current day (Su, Mo, Tu, We, Th, Fr or Sa) and time always remains visible in the lower righthand corner.

The only elements you can modify in this screen are setpoints. **To change setpoints in this screen:**

1. Turn the Knob until the desired setpoint is highlighted, then click the Knob in. The setpoint will begin to flash.
2. Turn the Knob to the right to increase the number, or to the left to decrease it. One click left or right will increase the setpoint by an increment of one (one whole unit or one tenth unit, depending on the parameter). Move the Knob more rapidly (you can spin it) to change the value by larger increments.
3. Click the Knob in to set and save the new value.

 **NOTE:**

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

4. Repeat the above steps to change any or all of the other setpoints.

To move out of this screen and into the next:

1. Use the Control Knob to highlight **SUMM**, then click the Knob in. **SUMM** begins to flash.
2. Turn the Knob to the right until the next screen, Setup (**SET**), appears. If you turn too far and enter another screen, just turn the Knob back to the left to recapture the **SET** screen.
3. Click the Knob in to select the screen and to work in it.

7.4.3 Setup Screen

Here you can set the day of the week and the time (on a 24-hour clock). This screen also allows you to lock all of your settings from further changes, and to mute or enable the audible alarm.

Figure 17: Setup Screen

SETUP	
Day :	Thu
Hour/Min :	16:19
Lock :	Off
Mute :	Off
SET	Th 16:19

To change the day:

1. Turn the knob to highlight the day (**Thu** in the sample screen above), then click inward once. The day will flash.
2. Turn the knob left or right to select the day of choice: Sun, Mon, Tue, Wed, Thu, Fri or Sat.
3. Click the knob in to set and save your choice.

 **NOTE:**

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

To change the time (Hour/Min):

1. Turn the knob to highlight the time (**16:19** in the sample screen above), then click inward once. The time will flash.
2. Turn the knob left or right to change the time. Left moves backward, right moves forward in time. One click right or left changes by one minute; spin the knob to move more rapidly.
3. Click the knob once inward to set and save your choice.

To lock the settings:

1. Turn the knob to highlight **Lock**, then click inward once. The current status (**Off** in the sample screen above) will flash.
2. Turn the knob in either direction; the only other choice is **On**. Click once inward to select and save **On**, or continue turning to return to **Off**.
3. When you set **Lock** to **On**, a padlock icon (*see Section 5.4*) will appear at the bottom of the screen. This icon will remain on display through all main display screens until you turn the locking function off.

To mute the audible alarm:

1. Turn the knob to highlight **Mute**, then click inward once. The current status (**Off** in the sample screen above) will flash.
2. Turn the knob in either direction; the only other choice is **On**. Click once inward to select and save **On**, or continue turning to return to **Off**.
3. When you set **Mute** to **On**, the  icon will appear at the bottom of the screen. This icon will remain on display through all screens until you turn the muting function off.

To move out of this screen and into the next:

1. Use the Control Knob to highlight **SET**, then click the Knob in. **SET** begins to flash.
2. Turn the Knob to the right until the next screen, Lamps (**LAMP**), appears. If you turn too far and enter another screen, just turn the Knob back to the left to recapture the **LAMP** screen.
3. Click the Knob in to select the screen and to work in it.

7.4.4 Lamps Screen

In this screen (*see Figure 18*), you can turn the chamber light (**Chamber**), the optional UV germicidal (“decontamination”) lamp (**UV Decont**) and the optional photosynthetic growth lamps (**Growth**) on and off.

On means the light is always on, and **Off** means the light is always off, unless you add additional programming.

There is an additional mode for the chamber light: **Auto**. In **Auto** mode, the light will go on every time you activate the control knob or open the lid. This is the default mode.

 **NOTE:**

The Lamps screen, shown below, will always indicate Chamber light mode. If the shaker is not equipped with the optional UV germicidal lamp and/or photosynthetic growth lamps, UV Decont and/or Growth will say None.

Figure 18: Lamps Screen

<u>LAMPS</u>	
Chamber :	Auto
UV Decont :	Off
Growth :	Off
LAMP	Th 16:20

To change the mode setting for any of the lamps:

1. Turn the Control Knob to highlight the setting for the lamp of choice, then click the knob in. The current setting will flash (in the sample screen, we will use the **Chamber** light as an example).
2. Turn the Knob left or right until the desired mode setting appears (**Auto** in this example).
3. Click the Knob in to save the new setting.

 **NOTE:**

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

To move out of this screen and into the next:


1. Use the Control Knob to highlight **LAMP**, then click the Knob in. **LAMP** begins to flash.

2. Turn the Knob to the right until the next screen, RS232 (**COMM**), appears. If you turn too far and enter another screen, just turn the Knob back to the left to recapture the **COMM** screen.
3. Click the Knob in to select the screen and to work in it.

7.4.5 RS232 Screen

This screen (see *Figure 19*) is used only if you have connected a personal computer to the RS-232 port (see *Section 5.8*). Here you can select the RS-232 port's **Mode** and the **Baud Rate** appropriate to your PC.

Figure 19: RS232 Screen



To change the Communication Mode:

1. Turn the Control Knob to highlight the current setting (**Off** in the sample screen above), then click the Knob in. The current setting will flash.
2. Turn the Control Knob left or right until the desired mode setting appears (see *Table 5: Communication Mode, on the following page*).

Table 5: Communication Mode

Mode	Application
Off	The RS-232 port is not open to communication in either direction.
Slave	The shaker can be fully controlled from the computer.
Talk	The shaker sends current value reports to the computer once per minute.
Monit [Monitor]	The shaker responds only to "Report Requests".

3. Click the control knob inward once to save the new setting.

 **NOTE:**

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

To change the Baud Rate:

1. Turn the control knob to highlight the current setting (**19200** in the sample screen above), then click the knob inward once. The current setting will flash.
2. Turn the control knob left or right until the desired setting appears: 9600, 19200 or 38400. The setting you choose should match the baud rate of your computer.
3. Click the control knob inward once to save the new setting.


To move out of this screen and into the next:

1. Use the Control Knob to highlight **COMM**, then click the Knob in. **COMM** begins to flash.
2. Turn the Knob to the right until the next screen, Calibrate (**CAL**), appears. If you turn too far and enter another screen, just turn the Knob back to the left to recapture the **CAL** screen.
3. Click the Knob in to select the screen and to work in it.

7.4.6 Calibrate Screen

Use this screen (*see Figure 20 on the following page*) to create a temperature offset and to calibrate the shaking speed (*for details, see Sections 7.7, Temperature Offset Calibration, and 0, Using Calspeed*).


Figure 20: Calibrate Screen



7.4.7 Programs Screen

Use this screen (*see Figure 21*) to set as many as four operating programs for the shaker. Each program can have as many as 15 steps. See Section 7.5 below for complete details.

Figure 21: Programs Screen



7.5 Programming the Shaker

7.5.1 Timer Only

By setting an **HRS** setpoint in the **DISP** or **SUMM** screen, the shaker may be set to automatically stop after a preset time period of 0.1 to 99.9 hours.

If the time is set to 0.0, the shaker will operate continuously until either the lid is opened or the Start/Stop button is pushed.

7.5.2 Programmed Steps

The resident software for the Innova 43/43R can store up to four programs, each having as many as 15 steps. *Each step* can be programmed in one-minute increments, for total periods of one minute up to 99 hours 59 minutes each.

To enter the programming mode, use the Control Knob to select the **PROG** screen (see *Figure 21, repeated below*). At this point, you can **Run** a program, **Edit** a program, generate a **New** program or turn a program **Off**. **Off** is the default mode.

Figure 21: Programs Screen


PROGRAMS	
Program 1 :	Off
Program 2 :	Off
Program 3 :	Off
Program 4 :	Off
PROG	Th 16:23

7.5.3 Creating a Program

To write a **New** program:

1. Use the Control Knob to highlight the mode of Program 1 (in our sample screen, Figure 21, it is **Off**), then click the Knob in. The selected field begins to flash.
2. Turn the Knob until the field says **New**. Click the Knob in to select this mode. The screen for Program 1 - Step 1 will open (see *Figure 22*):

Figure 22: Program 1 – Step 1



NOTE:

If your shaker is not equipped with these optional features, “UV” and “GRO” will appear with the word “Off” in this screen but will not be programmable.

3. Turn the Knob to highlight the Time setting (00:00 in the sample screen above), then click the Knob in. The field will flash.
4. Turn the Knob until the desired running **duration** for this step (from 00:01, which means one minute, to 99:59) appears, then click the Knob in to save the setting. We will set Step 1 time to eight hours for this example (*see Figure 23*).
5. Turn the Knob to highlight the °C temperature setting (20.0 in the sample screen above), then click the Knob in. The field will flash.
6. To set the temperature desired (°C from 4.0 to 80.0) for the time period you have set, turn the Knob (left to decrease, right to increase). When the desired value appears, click the Knob in to save the setting. We will set Step 1 temperature to 37.0°C (*see Figure 23*).
7. Turn the Knob to highlight the **RPM** setting (**Off** in the sample screen above), then click the Knob in. The field will flash.
8. Turn the Knob to select the desired shaking speed (25 to 400 RPM) for this time period, then click the Knob in to save the setting. We will set Step 1 speed to 150 RPM (*see Figure 23*).

9. *If you do not have the optional **UV** lamp, skip this step.* If you have the UV lamp (default setting is **Off**) and you wish to have it on during the time period you have set, use the Control Knob to select the field, click the Knob to select the field (it will flash), turn the Knob to select **On**, then click the Knob in to save the setting. In our sample Step 1 (*see Figure 23*), however, **UV** remains **Off**.
10. *If you do not have the optional **GRO** lamps, skip this step.* If you have the photosynthetic lamps (default setting is **Off**) and you wish to have them on during the time period you have set, use the Control Knob to select the field, click the Knob to select the field (it will flash), turn the Knob to select **On**, then click the Knob in to save the setting. In our sample Step 1 (*see Figure 18*), however, **GRO** remains **Off**.

 **NOTE:**

DO NOT select “Save” yet!

11. To program Step 2 (*see sample screen in Figure 23*): use the Control Knob to highlight the arrows at the bottom left of the screen. Click the Knob to make the arrows flash, then turn the Knob to the right until Step 2 appears. Click the Knob in to work in this screen, and repeat Steps 3-10.

Figure 23: Program 1 – Step 2

PRG1 - STEP	01	02	
Time	08:00	03:00	Step number
°C	37.0	40.0	
RPM	150	225	
UV	Off	Off	
GRO	Off	Off	
	Save	Cancel	

Whatever step is shown on the right-hand side of the display is settable (in Figure 23, it is Step 2). To scroll between steps, select the arrows (bottom left), then turn the Control Knob clockwise or counterclockwise, then select the desired step.

 **NOTE:**

The time entered for each step is for that step only; it is not cumulative (i.e., the elapsed time from the start of the program).

12. Continue to program as many as 15 steps in the same manner. For our sample program, we have only three steps (*see Figure 24*).

Figure 24: Program 1 – Step 3

PRG1 - STEP	02	03	Step number
Time	03:00	12:00	
°C	40.0	4.0	
RPM	225	0	
UV	Off	Off	
GRO	Off	Off	
	Save	Cancel	

The three-step program we have illustrated is designed to start incubating cultures at 37°C and hold that temperature for eight hours, shaking at 150 RPM. After eight hours, Step 2 engages, increasing the temperature setpoint to 40°C for a temperature induction, holding that temperature for three hours and increasing the shaking speed to 225 RPM. After this interval, the temperature will be cooled to 4°C and held at that temperature for twelve hours; no shaking will occur during this period, as the speed is set to 0 RPM.

Please note that temperature reductions, even at temperatures above ambient (going from 40° to 30°C, for example), require refrigeration to be effective.

To save the entire program:

1. Once all your steps are set, use the Control Knob to highlight **Save** at the bottom of the screen, then click the Knob in. The field will flash.
2. Click the Knob in again to save the program. The display will read **Process Running – Saving Profile** for a few seconds, then return to the main Programs (**PROG**) screen.

If you wish, you can set Programs 2, 3 and 4 and save them in the same manner.


7.5.4 Editing a Program

Use the **Edit function** to open a program you have already created and saved, in order to modify your settings following the same procedures.

7.5.5 Running a Program

Use the **Run** function to turn a specific program on. Naturally, only one program can run at a time. When you change the mode to **Run**, the screen will show the **Run** icon:

Figure 25: Programs Screen – Running



To stop a program: You can abort a run any time by changing the program's function to **Off**.

7.6 Muting the Audible Alarm

The Innova 43/43R shakers have an audible alarm that is activated under predetermined conditions (see Table 4). It may be muted in the following way:

1. Turn the Control Knob until the **SET** screen is highlighted on the display. Click the Knob in to work in this screen.
2. Turn the Knob to highlight the **Mute** mode (the default setting is **Off**, which means the alarm is audible). *To mute the alarm*, click the Knob in. The field will flash.
3. Turn the Knob to change the setting to **On**, then click the Knob to save this selection.
4. To reactivate the audible alarm at any time, repeat steps 1-3, reversing **Off** and **On**.

7.7 Temperature Offset Calibration

The temperature probe and the temperature controller are calibrated together at the factory. The temperature probe measures the temperature of the air at the probe's location, near the return vent. The controller uses the probe input to adjust air temperature, up or down, to match the temperature setpoint.

Depending on various conditions within the chamber, such as flask placement and size, the heat produced by growing organisms, heat losses due to liquid evaporation from flasks, etc., the display temperature may differ from temperatures within the flasks themselves. You can calculate the correction value for this offset and program the shaker to display a corrected temperature.

7.7.1 Calculating the Offset Value

If you wish to have the temperature display (“Indicated Temperature”) match the temperature at a given point, or match the average of a series of points within the chamber (“Actual Temperature”), proceed as follows:

1. Let the unit equilibrate at or near the desired temperature, then record the Indicated Temperature.
2. Now record the Actual Temperature.
3. Calculate the temperature correction value using this formula:
Actual Temperature – Indicated Temperature = Temperature Offset Value.
4. To set the Temperature Calibration Offset, follow the procedure outlined in Section 7.7.2 below.

7.7.2 Setting the Offset


To set the temperature calibration **Offset**:

1. Use the Control Knob to enter the **CAL** screen (*see Figure 20, repeated below*).
2. Turn the Control Knob to highlight the current setting (**0.0** in the sample screen above), then click the Knob in. The current setting will flash.
3. Turn the Knob (left for negative settings or right for positive settings) to display the desired setting. As you turn the Knob, each click represents one tenth of a degree Celsius (0.1°C).
4. When you reach the desired setting, click the Control Knob in to save the new value.
5. When the Temperature Offset is set to any value other than 0, the  icon will appear next to °C in the **DISP** and **SUMM** screens.

NOTE:

If you highlight an item, change it, but do not save your selection, after a few seconds the screen will revert to its previous setting.

Figure 20: Calibrate Screen



7.8 Using Calspeed

The **Calspeed** function, which is set in the **CAL** screen (see Figure 20 above), is used to calibrate the speed of the shaking mechanism. Calibrated at the factory, speed does not need to be recalibrated until a major operating component (e.g., drive belt) is changed.

NOTE:

Prior to calibrating the speed, make sure that the platform is properly secured to the subplatform, and any flasks present are secured.

When the shaker is running, the **CAL** screen shows the indicated RPM value. If you wish to calibrate the speed, set the speed to a value that can be measured—a setpoint of 250 RPM works well. The use of a strobe is recommended for accuracy. If, after measuring the actual speed, you wish to adjust the indicated value:

1. Click the Control Knob in.
2. Set the new value.
3. Click the Knob in again to save the setting.
4. Turn the shaker **OFF**, wait a few seconds, then turn the shaker back **ON**.

7.9 Power Interruption

In the event of a power failure, the Innova 43/43R Shakers are equipped with an **automatic restart** function. The shaker's non-volatile memory retains all stored information.

If the shaker was in operation prior to the power interruption, the shaker will begin to operate at its last entered setpoints. The alarm/POWER display will flash, indicating that a power interruption has occurred. Turn the Control Knob in any direction to acknowledge the visual alarm. The flashing will stop.

8 MAINTENANCE



WARNING!

When cleaning the unit, always turn off the shaker and disconnect the power cord from the power supply.

8.1 Routine Maintenance

No routine maintenance schedule is required for the Innova 43 and 43R.

To ensure that your shaker retains its attractive appearance, an occasional cleaning, using a cloth with conventional household (non-abrasive) cleaner is recommended (*see Section 8.2 below for more details*).

We also suggest that the area around the shaker be vacuumed or swept to remove dust and other debris, ensuring proper air flow in and around the shaker.

8.2 Cleaning External & Internal Surfaces

The unit may be cleaned using a damp cloth or any standard household or laboratory cleaner to wipe down its outer surfaces. Do not use abrasive or corrosive compounds to clean this instrument, as they may damage the unit and void the warranty.

If Biohazard decontamination is required, see Section 8.3 below.

8.3 Biohazard Decontamination



It is the responsibility of the user to carry out appropriate decontamination procedures if hazardous material is spilled on or inside the equipment. Before using any cleaning or decontamination method other than those suggested by the manufacturer, users should check with New Brunswick Scientific that the proposed method would not damage the equipment.

Commercially available household bleach solutions, when diluted at a 1:10 ratio, are effective in routine decontamination of the instrument. The method for decontaminating a spill depends upon the nature of the spill.

Spills involving fresh cultures or samples known to have low concentrations of biomass should be flooded with decontamination solution and soaked for 5 min before cleanup. Spills involving samples with high concentrations of biomass, or involving organic matter, or occurring in areas warmer than room ambient temperature, should be exposed to decontamination solution for *at least one hour* before cleanup.

**WARNING!**

Personnel involved in the cleanup of any spill should wear gloves, safety glasses, and a laboratory coat or gown during the cleanup process. Respiratory protection should be considered for spills where aerosolization is suspected.

9 SERVICE & ACCESSORIES

9.1 Troubleshooting

If any problems occur with your shaker, do not attempt to perform any service on the unit other than specified in this manual. Unauthorized servicing may void the warranty. Please contact your local NBS Customer Service Department

In any correspondence with NBS, please refer to the model number and serial number of your unit. This information is on the electrical specification plate, located on the rear panel of the unit, above the power connector.

There are some problems, however, that you can investigate and correct yourself. Refer to the following Troubleshooting Guide:

Symptom(s)	Probable Cause(s) & Solution(s)
Shaker does not run.	No power; display is not on; power cord is not plugged in and/or power switch is off: plug in power cord (to working electric outlet), and turn on power switch.
	Lid is open—look for Open Lid icon on display: close the lid firmly.
	Lid is closed but Open Lid icon is on display: lid switch is broken; call for service.
	On/Off switch is not working: call for service.
	Fuse(s) burned out: check and replace as needed.
	Shaking speed has been set to Zero by program running (look for Run icon in display) or by computer interface: reset shaking speed.
	Defective main board: call for service
	Defective display controller board: call for service.
	Jammed shaking mechanism: call for service
	Defective motor: call for service
	Drive belt out of alignment or worn: call for service.
Shaker runs slowly and/or no speed indication.	Fuse(s) burned out: check and replace as needed.
	Incorrect speed calibration: recalibrate shaking speed.
	Defective main board: call for service.
	Defective motor: call for service.
	Drive belt out of alignment or worn: call for service.

...continued...

Symptom(s)	Probable Cause(s) & Solution(s)
Shaker does not run at set speed.	Shaker is running in Program mode (look for Run icon in display).
	Shaker speed has been changed by RS-232 command/ computer interface.
	Shaker is overloaded: remove some contents & balance load.
	Defective motor: call for service.
	Drive belt out of alignment or worn: call for service.
Operating noise	Load out of balance: unload all contents from platform, then reload.
	Loose component(s) in Platform, subplatform, and/or drive assembly: call for service.
Incubator does not reach set temperature.	Shaker running in Program mode: Run icon in display.
	Temperature setpoint changed by RS-232/computer interface command.
	Heater fuse blown: replace.
	Compressor fuse blown: replace.
	Compressor over-pressure switch activated: call for service.
	Chamber fans not working: call for service
	Ambient temperature too high or too low: cool or heat the room as needed.
	Lid is not completely closed (even though Open Lid icon may not be on display): open it and reclose it firmly.
	Line voltage is too low.
	Frequency on line voltage if set incorrectly: reset.
	Defective heater: call for service.
Incorrect temperature indication.	Temperature Offset has been programmed (look for [adjacent to C on display).
	Defective RTD assembly: call for service.
	Defective main board: call for service.
Incorrect relative humidity %	Defective humidity probe or broken contact: call for service.

...continued...

Symptom(s)	Probable Cause(s) & Solution(s)
Photosynthetic GRO lamp does not work.	Shaker is running in Program mode (look for Run icon in display), and program calls for GRO lamp to be off.
	GRO lamp mode has been changed by RS-232 command/computer interface.
	Bulb is burned out; replace (<i>see Section 9.6</i>).
	Fuse is burned out: replace (<i>see Section 9.4</i>)
	Check ballast voltage.
	Check cable connections.
	Check switch settings.
UV germicidal lamp does not work.	GRO lamp mode has been changed by RS-232 command/computer interface.
	Shaker is running in Program mode (look for Run icon in display), and program calls for UV lamp to be off.
	UV lamp is burned out; replace (<i>see Section 9.8</i>)
	Fuse is burned out: replace (<i>see Section 9.4</i>)
	Check ballast voltage.
	Check cable connections.
	Check switch settings.
UV lamp mode has been changed by RS-232 command/computer interface.	

9.2 Product Returns

Should you need to return your Innova 43/43R to NBS for any reason, first contact Customer Service to obtain a Returned Material Authorization (RMA) number. This number must appear on the outside of the shipping container, otherwise NBS Receiving will refuse to accept the shipment.

In addition, you must also certify that the instrument being returned has been thoroughly cleaned and decontaminated. **A form for this purpose is provided in Section 13** of this manual; it can also be downloaded from our website (www.nbsc.com). A copy of this completed Return Authorization and Decontamination Certificate must be attached to the outside of the container, with a second copy packed inside with the instrument.

9.3 Opening the Service Compartment



WARNING!

Before opening the Service Compartment, always turn off the shaker and disconnect the power cord from the power supply.

The Service Compartment (*see Figure 26*) contains the shaker's electronic and temperature control components. The user should access the Service Compartment only to replace a fuse—*see Section 9.4 for Fuse Replacement instructions*). All other service intervention are to be performed by **authorized service technicians only**.

9.4 Fuse Replacement



WARNING!

When replacing fuses, always turn off the shaker and disconnect the power cord from the power supply.

There are four fuses that can be replaced by the user. These are located inside the Service Compartment. To access the fuses (*see Figure 26*):


1. Turn the power off and unplug the shaker.
2. Remove the four screws and washers that hold the right side service panel grille in place. Set them aside for reuse.
3. Fuses are numbered 1 through 4 (*see Table 6*); access the fuse you wish to replace by using a coin or a blade screwdriver to turn and release the spring-loaded cap.
4. Replace the fuse with a new one of the same type and rating (*see Table 6 below*):

Table 6: Fuses

Fuse No.	Function	Type	Rating
F1	Heater	Slo Blo®	6.3 A
F2	☐GRO Lamps (photosynthetic)		3.15 A
F3	☐UV Lamp (germicidal)		3.15 A
F4	Compressor		10 A

☐ optional equipment

Figure 26: Fuse Replacement



9.5 Replacing Cabinet Lights



WARNING!

Before installing or removing the quartz bulbs:


- (1) Make sure power is disconnected from the unit,
- (2) Wear plastic gloves to protect the bulb from skin oils and to protect your hands just in case the bulb should break,
- (3) Before removal, make sure the bulb has fully cooled.



NOTE:

Cabinet lights are connected in series. When one burns out, both lights will be off. It is therefore prudent to replace both bulbs at the same time.

Be sure to follow all the **WARNINGS** and the **NOTE** above as well as the procedure outlined below when you replace the cabinet bulbs with reference to Figure 27 below:

Figure 27: Replacing Cabinet Lights

1. Make sure power is disconnected from the shaker.
2. Open the lid.
3. Wearing plastic gloves, remove the first used bulb, then carefully insert the new bulb's prongs fully into the connector base.
4. Repeat step 3 for the second bulb.

9.6 Replacing Optional Photosynthetic Lamps




WARNING!

Before removing or installing the photosynthetic "GRO" lamp, make sure power is disconnected from the unit.

To replace the photosynthetic GRO lamp bulbs, with reference to Figure 28:

1. Make sure power is disconnected from the shaker.
2. Open the lid.
3. Remove the two thumbscrews from the diffuser holder.

Figure 28: Replacing Photosynthetic Lamps

4. Remove the diffuser by placing a finger underneath and rotating the diffuser in the slot at the back of the photosynthetic light assembly. All 9 lamps should now be accessible to you.
5. To remove the lamp for replacement, rotate the lamp 90°, taking care to support the lamp before it drops down.
6. To replace the lamp, orient the pins of the new lamp vertically. Simultaneously, gently press both sides of the lamp into the slots.
7. Rotate the new lamp 90°. The lamp should catch in the holders on both ends. Make sure both ends are secure.
8. Put the diffuser back in place, reinstalling the diffuser holder with its two thumbscrews.

9.7 Adjusting Lid Tension

In the event the lid does not stay in the fully open position, you may want to adjust the lid tension: tighten the 5/16-inch nuts on the hinges to increase spring tension.



WARNING!

The following procedures are provided for your information only. Do not attempt to perform these service interventions yourself unless you are an authorized service technician.

9.8 Replacing Optional Germicidal Lamp



WARNING!


Before installing or removing the germicidal “UV” lamp, the authorized technician will:

- (1) Make sure power is completely disconnected from the unit,
- (2) Wear plastic gloves to protect the lamp from skin oils and to protect your hands just in case the lamp should break,
- (3) Before removal, make sure the lamp has fully cooled.

To replace the UV lamp, your service technician will follow all the **WARNINGS** above and the procedure and **WARNINGS** outlined below, with reference to Figure 29:

1. Make sure power is disconnected from the shaker.
2. Open the Service Compartment by removing the side panel.
3. Free the heater cable harness and the temperature probe from their tie wraps and Permagem[®].

Figure 29: Germicidal (UV) Lamp Replacement



NOTE that wires, tie wraps and Permagem[®] are not shown above.


4. Remove the evaporator cover from the Service Compartment, paying careful attention to the drain tubes.
5. Wearing plastic gloves, remove the germicidal lamp by pressing it toward the spring-loaded tombstone. The lamp has a single pin at each end.
6. Still wearing plastic gloves, install a new lamp.
7. Pull the heater cable harness and temperature probe back through the holes in the Service Compartment.

8. Replace the fan cover, making sure the heater cable harness wires and temperature probe have been pulled through. Adjust the temperature probe position to ensure that the probe tip can protrude 3 – 3.25 inches (7.6 – 8.2 cm) from the fan cover face.
9. Secure the heater cable harness and temperature probe with tie wraps, and add Permagum[®] to fill the hole.
10. Reinstall the side panel that closes off the Service Compartment.
11. Use the plastic fitting to hold the temperature probe, ensuring that its tip is 3 – 3.25 inches (7.6 – 8.2 cm) from the fan cover face.

9.9 Belt Replacement or Adjustment


To gain access to the drive belt, your service technician will follow these steps with reference to Figures 30 & 31 below:

1. Make sure that the power is switched off and the shaker is unplugged.
2. Open the lid and remove the platform.
3. Using a 7/32-inch Allen key (hex wrench), remove the four bolts that hold the subplatform in place, and remove the subplatform.
4. Disconnect the motor wire by unplugging its connector.
5. With a wrench, remove the nuts that hold the bearing housing to the bottom weldment.
6. Position the counterweight so it is closest to you, and tie it in place with tie wraps or strong cord, **in order to prevent it from dropping against your fingers or causing you to lose your grip.**
7. Remove the bearing housing with care: it weighs 47 lbs. / 21.3 kg.
8. Turn the bearing housing over to access the belt, motor and pulley.

Figure 30: Drive Assembly

To replace the drive belt, follow these steps with reference to Figures 30 & 31:

1. Loosen the four hex nuts that hold the motor assembly in place and remove the old belt.
2. Guide the belt onto the two pulley grooves.
3. Tighten the motor nuts.
4. Check the belt adjustment by applying finger pressure to the belt midway between the two pulleys. The belt should deflect approximately 3/8-inch (9-10 mm).
5. If the deflection is not right, loosen the nuts, readjust the belt, then retighten the motor nuts.

Figure 31: Belt Replacement & Adjustment**If an adjustment is required:**

1. Loosen the four hex nuts on the motor assembly.
2. Move the motor assembly until the belt is tight.
3. Tighten the hex nuts and recheck the belt tension by exerting pressure on the belt. The belt should deflect approximately 3/8-inch (9.5 mm).

Reassembly:

1. In reverse order, reinstall the bearing housing. Once it is securely in place, remove the tie wraps/cord.
2. Secure the wiring to ensure no interference during rotation.
3. Reinstall the subplatform and platform.
4. Recalibrate speed (*see Section 0*).

9.10 Replacement Parts

In Table 7 below you will find commonly needed replacement parts for your Innova 43/43R shaker:

Table 7: Replacement Parts

Part Description	Quantity	NBS Part No.
Spare Parts Kit (includes 1 motor drive poly-V belt, two each hinge assemblies, halogen bulbs, and 10A, 6.3A & 3.15A fuses)	1	M1320-6000
Spare Parts Kit (includes 1 motor drive poly-V belt, and two each 10A, 6.3A & 3.15A fuses)	1	M1282-6000
Flask Clamp screws, 10-24 x 5/16-inch flat Phillips head (+)	25	S2116-3051
Test Tube Rack screws, 10-24 x 1/2-inch pan head, and washers (#10 plain and #10 lock washers)	10 each	M1289-0800
Motor Drive Poly-V Belt, 28 inches long	1	P0700-7070
Fuse, 10A	1	P0380-3162
Fuse, 6.3A	1	P0380-3630
Fuse, 3.15A	1	P0380-3460
Chamber Light Bulbs	2	P0300-0470
UV Germicidal Bulb (if so equipped)	1	P0300-0351
Photosynthetic Bulb (if so equipped)	1	P0300-0221
UV Germicidal Light Ballast (if so equipped)	1	P0420-5113
Photosynthetic Light Ballast (if so equipped)	1	P0420-5192

9.11 Accessories

When ordering accessories, you may be asked to provide the model number and serial number of your shaker. This information is on the electrical specification plate, located on the rear panel of the unit. The serial number is also labeled in the lower right corner of the front panel, below the chamber lid seal.

9.11.1 Platforms

Table 8 lists the various solid 18- x 30-inch (45.7 x 76 cm) platforms available for your Innova 43/43R. Consult Section 10.2 for load and speed graphs.

Table 8: Available Platforms

Description	NBS Part No.
Universal Platform	M1250-9920
50 mL Erlenmeyer flask Dedicated Platform	M1191-9908
125 mL Erlenmeyer flask Dedicated Platform	M1191-9909
250 mL Erlenmeyer flask Dedicated Platform	M1191-9910

...continued...

Description	NBS Part No.
500 mL Erlenmeyer flask Dedicated Platform	M1191-9911
1L Erlenmeyer flask Dedicated Platform	AG-1
2L Erlenmeyer flask Dedicated Platform	AG-2
2.8L Fernbach flask Dedicated Platform	AG-28
4L Erlenmeyer flask Dedicated Platform	AG-4
6L Erlenmeyer flask Dedicated Platform	AG-6

To add your choice, following is a list of flask capacities for the Universal Platform and for the dedicated platforms designed solely for one flask size:

Table 9: Platform Flask Capacities

Flask Type	Universal Platform Capacity*	Dedicated Platform Capacity*
10 ml Erlenmeyer Flasks	183	N/A
25 ml Erlenmeyer Flasks	92	N/A
50 ml Erlenmeyer Flasks	92	108
125 ml Erlenmeyer Flasks	39	60
250 ml Erlenmeyer Flasks	30	40
500 ml Erlenmeyer Flasks	18	24
1 L Erlenmeyer Flasks	12	15
2 L Erlenmeyer Flasks	8	12
2.8 L Erlenmeyer Flasks	6	6
4 L Erlenmeyer Flasks	6	6
5 L Erlenmeyer Flasks	6	6
6 L Erlenmeyer Flasks	4	4

*Clamps are sold separately.

9.11.2 Flask Clamps for Platforms

The following clamps are available for the matching Dedicated Platform; all may be used on the Universal Platform:

Table 10: Flask Clamps

Clamp Size	Part Number
10 ml Erlenmeyer Flask	ACE-10S
25 ml Erlenmeyer Flask	M1190-9004
50 ml Erlenmeyer Flask	M1190-9000
125 ml Erlenmeyer Flask	M1190-9001
250 ml Erlenmeyer Flask	M1190-9002
500 ml Erlenmeyer Flask	M1190-9003
1 L Erlenmeyer Flask	ACE-1000S
2 L Erlenmeyer Flask	ACE-2000S

...continued...

Clamp Size	Part Number
2.8 L Fernbach Flask	ACFE-2800S
4 L Erlenmeyer Flask	ACE-4000S
5 L Erlenmeyer Flask	ACE-5000S
6 L Erlenmeyer Flask	ACE-6000S

9.11.3 Replacement Clamp Hardware Kits

NBS flask clamps come complete with mounting screws. Additional screws are available separately in packs of 25 (*see Section 9.10*).

9.11.4 Test Tube Racks & Other Accessories

Table 11: Racks & Trays

Accessory Description		NBS Part Number	Platform Capacity
Adjustable angle Test Tube Rack for tubes 8 – 11 mm diameter	80 tube capacity	M1289-0100	7
	60 tube capacity	M1289-0010	9
	48 tube capacity	M1289-0001	9
Adjustable angle Test Tube Rack for tubes 12 - 15 mm diameter	60 tube capacity	M1289-0200	7
	44 tube capacity	M1289-0020	9
	34 tube capacity	M1289-0002	9
Adjustable angle Test Tube Rack for tubes 15 –18 mm diameter	42 tube capacity	M1289-0300	7
	31 tube capacity	M1289-0030	9
	24 tube capacity	M1289-0003	9
Adjustable angle Test Tube Rack for tubes 18 – 21 mm diameter	30 tube capacity	M1289-0400	7
	23 tube capacity	M1289-0040	9
	18 tube capacity	M1289-0004	9
Adjustable angle Test Tube Rack for tubes 22 – 26 mm diameter	22 tube capacity	M1289-0500	7
	16 tube capacity	M1289-0050	9
	13 tube capacity	M1289-0005	9
Adjustable angle Test Tube Rack for tubes 26 - 30 mm diameter	20 tube capacity	M1289-0600	7
	16 tube capacity	M1289-0060	9
	12 tube capacity	M1289-0006	9
Microplate holder rack (stacked)	3 deep well or 9 standard	M1289-0700	16
Microplate holder rack (single layer)	5 deep well or standard	TTR-221	4
Angled Test Tube Rack Holder* for user-supplied test tube racks that are 4-5 in. (10-13 mm) wide and up to 15 in. (38 mm) long.		TTR-210	4
Angled Test Tube Rack Spacer Bar* for use with TTR-210 to accommodate test tubes racks that are less than 5 in. (13 mm) wide.		TTR-215	NA

* Universal Platform Required

9.11.5 Optional Gas Manifold Kit

Table 12: Gas Manifold Kit (Factory-Installed)

<i>Description</i>	<i>NBS Part Number</i>
Gas Manifold Kit	M1320-0500

9.11.6 Optional Photosynthetic Lamps

Table 13: Photosynthetic Growth Lamps (Factory-Installed)

<i>Description</i>	<i>NBS Part Number</i>
Photosynthetic Growth Lamp Assembly (9 lamps)	M1320-0300

9.11.7 Optional UV Germicidal Lamp

Table 14: Germicidal Lamp (Factory-Installed)

<i>Description</i>	<i>NBS Part Number</i>
UV Germicidal Lamp Assembly	M1320-0400

9.11.8 Optional Humidity Monitor

Table 15: Humidity Monitor (Factory-Installed)


<i>Description</i>	<i>NBS Part Number</i>
Humidity Monitor (for use with spill pan as water reservoir)	M1320-0600

9.11.9 Optional Utility Basket

Table 16: Utility Basket

<i>Description</i>	<i>NBS Part Number</i>
Utility Basket	M1320-0700

Figure 32: Utility Basket (optional)




9.11.10 Optional Remote Alarm

The Innova 43/43R can be equipped with a factory-installed remote alarm component (NBS part number M1320-8029). When it is hooked up to your relay and receiving equipment, this device will send notification of an alarm condition to the remote location you choose.


10 DRAWINGS, GRAPHS & TABLES

10.1 Schematics

Figure 33: Control Schematics




NOTE:
 For a larger version of this drawing, contact your sales or service representative.

Figure 34: Refrigeration Schematic


10.2 Load & Speed Graphs

The following graphs explain the maximum recommended shaking speed according to the load you put on your shaker's platform and the diameter of the shaker's orbit (1-inch or 2-inch).

**Figure 35: Load & Speed for 1-Inch Orbit Shakers
(20% Full)**



**Figure 36: Load & Speed for 2-Inch Orbit Shakers
(20% Full)**



10.3 List of Drawings

Figure	Description	Page
1	Innova 43/43R Front View	15
2	Adjustable Foot	17
3	Front Panel (Detail)	19
4	Control Panel (Right Side)	20
5	Installing Platform & Spill Cover	26
6	Double Girdle Clamp Installation	27
6a	Clamp Fastener	27
7	Rear Panel	28
8	Gas Manifold (optional)	29
9	Gas Fitting (with optional Gas Manifold)	29
10	Germicidal UV Lamps (optional)	30
11	Photosynthetic Lamps (optional)	31
12a	Humidity Monitor (optional)	31
12b	RH Sensor Operating Range & Error Graph	32
12c	RH Sensor Connection	32
12d	RH Sensor Mounting	33
13	Changing Display Parameter	37
14	Changed Display Parameter	37
15	Changing Setpoint	38
16	Summary Screen	39
17	Setup Screen	40
18	Lamps Screen	42
19	RS232 Screen	43
20	Calibrate Screen	45, 52
21	Programs Screen	45, 46
22	Program 1– Step 1	47
23	Program 1– Step 2	48
24	Program 1– Step 3	49
25	Programs Screen – Running	50
26	Fuse Replacement	60
27	Replacing Cabinet Lights	61
28	Replacing Photosynthetic Lamps	61
29	Germicidal (UV) Lamp Replacement	63
30	Drive Assembly	65
31	Belt Replacement & Adjustment	66
32	Utility Basket (optional)	71
33	Control Schematic	72
34	Refrigeration Schematic	73
35	Load & Speed for 1-Inch Orbit Shakers	74
36	Load & Speed for 2-Inch Orbit Shakers	74

10.4 List of Tables

Table	Description	Page
1	Screens	21
2	Display Screen Parameters	21
3	Display Icons	22
4	Alarms	23
5	Communication Mode	44
6	Fuses	59
7	Replacement Parts	67
8	Available Platforms	67
9	Platform Flask Capacities	68
10	Flask Clamps	68
11	Racks & Trays	69
12	Gas Manifold Kit	70
13	Photosynthetic Growth Lamps	70
14	Germidical [UV] Lamp	70
15	Humidity Monitor	70
16	Utility Basket	70
17	COM Parameters for HyperTerminal	80

<i>Innova 43 & 43R Incubator Shaker</i>	
CHAMBER DIMENSIONS	24 inches deep X 34.5 inches wide X 19 inches high from top of platform (61 cm deep X 87.6 cm wide X 48 cm high from top of platform)
CHAMBER LID	Foot-activated lid opening.
PLATFORM	Aluminum, 18" X 30" (46 X 76 cm). Select universal or dedicated styles
WEIGHT	
43	460 lbs (209 kg)
43R	510 lbs (231 kg)

 **NOTE:**

Use of baffled flasks will significantly reduce maximum speed for any shaker. We may be able to improve this maximum speed by using an alternative motor pulley. Contact your NBS representative for more information.


11.1 Certifications

The Innova 43 and 43R have been tested to ETL standards, to comply with UL and CAN/CSA electrical safety standards (see “ETL Regulatory Standards” in the specifications table).

As attested in the *CE Declaration of Conformity* reproduced on the following page, the Innova 43 and 43R also conform to the appropriate CE standards (see also “CE Regulatory Standards” in the specifications table).

11.1.1

CE Declaration of Conformity

 New Brunswick Scientific	
CE	CE
DECLARATION OF CONFORMITY	
New Brunswick Scientific, Hereby declares that the product(s) listed below conform to the European Union directive and standards identified in this declaration.	
<u>Product(s)</u>	
Innova 44/44R and 43/43R	
<u>EU Directive(s)</u>	
Low Voltage (73/23/EEC/93/68/EEC) Electromagnetic Compatibility (89/336/EEC/93/68/EEC)	
<u>Standard(s)</u>	
IEC61010-1: 2001 (2 nd Edition) EN61010-1: 2001 (2 nd Edition) EN61010-2-010 (1992) +A1: 1996 EN61010-2-081 (2002) EN61326: 1997+A1: 1998 EN61000-3-2: 1995 EN61000-3-3: 1995	
The conformity assessment procedures were performed at the following Testing Lab.:	
Intertek Testing Services, 40 Commerce Way, Totowa, NJ 07512	
The technical documentation relevant to the above equipment will be held at:	
New Brunswick Scientific Company PO Box 4005 44 Talmadge Road Edison, New Jersey 08818-4005 U.S.A Tel. (732) 287-1200 Fax. (732) 287-4222	
 _____ Lee Eppstein VP of Science & Technology	 _____ Date

12 APPENDIX: REMOTE PROGRAMMING

12.1 Setting Up with HyperTerminal

You can use the Microsoft® HyperTerminal program on your personal computer to communicate with your shaker:

1. Before you set up HyperTerminal, be sure to check the **Mode** and **Baud Rate** parameter settings (see Section 7.4.5).
2. While you are in the RS232 screen (see Figure 15), set the **Mode** parameter to **Slave**, and select the appropriate **Baud Rate**. The default setting, as you see in your screen, is 19200; your other choices are 9600 or 38400.
3. Connect the RS-232 cable between your computer and the shaker's RS-232 port (see Figure 4).
4. Start HyperTerminal, and follow the steps for Setting a New Connection. When prompted, select the COM port you will be using, then set the COM parameters (see Table 17 below):

Table 17: COM Parameters for HyperTerminal

Baud Rate	Enter the setting on your shaker (see Step 2 above)
Data Bits	8
Parity	None
Stop Bits	1
Flow Control	None

5. To be sure that you have established a connection with your shaker, turn the shaker **OFF**, open the HyperTerminal connection, then turn the shaker back **ON**. Once a connection has been made, the shaker will send an **OK** to HyperTerminal.
6. You may wish to send a trial command to the shaker (see *Sample Test below*). We suggest that you enable *Caps Lock* on your computer, because all commands must be sent as capital letters.

Sample Test:

1. To request the current setpoint, type **RS** at the terminal. When you type **R**, you should see the **R** on the computer display. This is an echo character, sent back by the Innova 43/43R. When you type the **S**, you should see the echo **S** displayed at the computer.
2. Use the **←Enter** key to send the carriage return control character. This will move the cursor back over the first character.
3. **RS** is a Request command, so your shaker will send back the data corresponding to the command.

12.2 Overview of Command Sets

Using an optional RS-232 interface, the user has the capability to control various functions of the Innova 43/43R shakers by computer. A detailed command set follows this introduction.

The serial interface uses an echo character protocol that is implemented to indicate that a valid character has been received by the shaker. For every character sent by the user computer, the Innova 43/43R will return a matching character. This echo character should be read by the computer and compared to the character that was sent.

NOTE:

When you are sending a command, do not allow the pause between any two characters of your message to last longer than 10 seconds. If it takes any longer, the shaker's serial port will reset itself.

In the following example, the **Master** (your computer) will command (the **C** in **CS**) the shaker (**Slave**) to set its speed to **150 RPM**:

Master sends:	Innova 43/43R responds:
CS<Space>150<CR>	CS<Space>150<CR>

In the next example, the **Master** (your computer) will request that the shaker (**Slave**) report (the **R** in **RV**) the actual parameters:

Master sends:	Innova 43/43R responds:
RV<CR>	First the echo: RV<CR>
	Then the report: <Arg1><Tab><Arg2><Tab><Arg3>...<Arg7><CR><LF>

NOTE:

On power-up, even if the shaker is not set in Slave mode at the time, the shaker will always send these four characters: OK<CR><LF>. Before starting any communications, be sure to clear the computer's receive buffer.

This is the command set format:

<Command><Space><Arg1><Space><Argn...><Terminator>

Most Commands are defined as two capital characters using standard ASCII format (see the command sets provided below).

If a command has additional arguments, the ASCII space control character separates each argument. The command is terminated using the ASCII carriage return control character; an optional line feed control character may also be included but must follow the carriage return character.

When a command is sent that returns data, the characters returned are in ASCII format. Each additional argument returned will be separated by the ASCII space or tab control character; refer to the selected command for the actual format. Data returned is terminated with a carriage return line feed control character sequence.

Section 12.3 provides an Index to the 2-character command codes.

Sections 12.4-12.7 provide command sets for your use for Set Commands, Profile Control Commands, Report Request Commands and Set/Get Date & Time.

12.3 *Index to Command Codes*

Command Code	Meaning
CL	Command Grow Lamp
CS	Command Speed
CT	Command Temperature
CU	Command UV Lamp
PC	Clear Profile
PM	Profile Mode
PR	Read Profile Step in Memory
PS	Profile Stop/Start
PW	Write Profile Step to Memory
RI	Report Software Version
RP	Report Parameter List
RS	Report Setpoint Values
RV	Report Measured Values
=D	Set Date and Time
?D	Get Date and Time

12.4 Set Commands

Code	Meaning	Format
CL	Command Grow Lamp	CL <Space><Arg1><CR> where Arg1 = Set Grow Lamp state (0 = off, 1 = on)
CS	Command Speed	CS <Space><Arg1><CR> where Arg1 = Agitation setpoint (####)
CT	Command Temperature	CT <Space><Arg1><CR> where Arg1 = Temperature setpoint (###.#)
CU	Command UV Lamp	CU <Space><Arg1><CR> where Arg1 = Set UV Lamp state (0 = off, 1 = on)

12.5 Profile Control Commands

Code	Meaning	Format
PC	Clear Profile (clears profile step)	PC <Space><Arg1><Space><Arg2><CR> where Arg1 = Profile number (range 1 to 4) and Arg2 =Step number (range 1 to 15) NOTE: If Arg2 is not included, all 15 steps of the selected profile will be cleared.
PM	Profile Mode (returns status of current profile run/hold mode)	PM <CR> There are no parameters to this command. Return data format: <Arg1><Space><Arg2><Space><Arg3><CR> where Arg1 =Run/hold status (0 = end, 1 = run), Arg2 =Profile number and Arg3 =Step number NOTE: If a profile has not been selected to run, a 0 will be returned for Profile number and Step number.

...continued...

Code	Meaning	Format
PR	Read Profile step in memory	<p>PR<Space><Arg1><Space><Arg2><CR> where Arg1= Profile number (range 1 to 4) and Arg2=Step number (range 1 to 15)</p> <p>Return data format: <Arg1><Space><Arg2><Space><Arg3><Space>...<CR><LF> where Arg1=Profile number (range 1 to 4), Arg2=Step number, Arg3=Temperature setpoint, Arg4=Agitation setpoint, Arg5=(Reserved, reads 0.0), Arg6=Step hours, Arg7=Step minutes, Arg8=UV lamp, and Arg9=Grow lamp</p>
PS	Profile Stop/ Start (stops or starts a program profile)	<p>PS<CR> With no other parameters, any profile currently running will be stopped.</p> <p>PS<Space><Arg1><Space><Arg2><CR> where Arg1=Profile number (range 1 to 4), and Arg2=Step number (range 1 to 15).</p> <p>NOTE: If Arg2 is not included, the profile will start in Step 1.</p>
PW	Write Profile step to memory	<p>PW<Space><Arg1><Space><Arg2><Space><Arg3><Space>...<CR> where Arg1=Profile number (range 1 to 4), Arg2=Step number (range 1 to 15), Arg3=Temperature setpoint, Arg4=Agitation setpoint, Arg5=(Reserved); Enter 0.0, Arg6=Step hours (range 0 to 99), Arg7=Step minutes (range 0 to 59), Arg8=UV lamp (0 = off, 1 = on), and Arg9=Grow lamp (0 = off, 1 = on)</p>

For example, the following **PW** command:

PW<Space>**25.0**<Space>**150**<Space>**0.0**<Space>**2**<Space>**30**<Space>**0**
 <Space>**1**<CR>

will set Program 1, Step 1 to the following conditions:

- Temperature setpoint 25.0°
- Agitation setpoint 150 rpm
- Step Time 2 hours 30 minutes
- UV Lamp OFF
- GRO Lamp ON

12.6 Report Request Commands

Code	Meaning	Format
RI	Report Software Version	RI <CR> Return data format: <Arg1><Tab><Arg2><CR><LF> where Arg1 =Display module software version, and Arg2 =Control module software version
RP	Report Parameter List (ASCII text header)	RP <CR> Return data format: <Arg1><Tab><Arg2><Tab><Arg3>...<Arg6><CR><LF> where Arg1 ="Shaker Speed", Arg2 ="Temperature", Arg3 ="% Relative Humidity", Arg4 ="Reserved", Arg5 ="Grow Lamp Status", and Arg6 ="UV Lamp Status"
RV	Report Measured Values	RV <CR> RS <CR> Return data format: <Arg1><Tab><Arg2><Tab><Arg3>...<Arg6><CR><LF>
RS	Report Setpoint Values	where Arg1 = Shaker Speed, Arg2 =Temperature, Arg3 =% Relative Humidity, Arg4 =(Reserved, reads 0.0), Arg5 =Grow Lamp Status, and Arg6 =UV Lamp Status

12.7 Set/Get Date & Time

Code	Meaning	Format
=D	Set Date & Time	=D <Space><Arg1><Space><Arg2><Space><Arg3><Space>...<Arg7><CR>
?D	Get Date & Time	?D <CR> Return data format: <Arg1><Tab><Arg2><Tab><Arg3><Tab>...<Arg7><CR><LF> where Arg1 =Hours (range 0-23), Arg2 =Minutes (range 0-59), Arg3 =Seconds (range 0-59), Arg4 =Year (range 00-99), Arg5 =Month (range 01-12), Arg6 =Date (range 1-31), Arg7 =Day (range 1-7; 1=Monday, 2=Tuesday...7=Sunday)

For example, this **=D** command:

```
=D<Space><01><Space><30><Space><00><Space><04><Space><07>  
<Space><04><Space><07><CR>
```

will set the real date and time on your shaker to:

1:30:00 am, 2004, July 4, Sunday.

13 APPENDIX: PRODUCT RETURNS

13.1 Return Procedure


As explained in Section 9.2, should you need to return your Innova 43/43R to NBS for any reason, first contact Customer Service to obtain a Returned Material Authorization (RMA) number. This number must appear on the outside of the shipping container, otherwise NBS Receiving will refuse to accept the shipment.

In addition, you must also certify that the instrument being returned has been thoroughly cleaned and decontaminated. A form for this purpose is provided on the following page; you can photocopy it and fill it out by hand. It can also be downloaded from our website (www.nbsc.com), if you prefer to fill it out electronically.

A copy of the completed Return Authorization and Decontamination Certificate must be attached to the outside of the container, with a second copy packed inside with the instrument.

13.2 Return Authorization and Decontamination Certificate

A sample form for you to copy and fill out is provided on the following page.



New Brunswick Scientific Return Authorization and Decontamination Certificate

Contact New Brunswick Scientific for an RMA number prior to returning any equipment, then complete this form and attach it to the outside container of the equipment being returned to our facility. In addition, please enclose a completed, duplicate copy of this form with the returned item.

Returned Material Authorization (RMA) Number _____

Equipment being returned: Model Number _____ Serial Number _____

Reason for its return:

This equipment (check all that apply):

New Product

Never used

Biohazards

Not used

Used, but decontaminated with

Hazardous Chemicals

Not used

Used, but decontaminated with

Radioactive Materials

Not used

Used, but decontaminated with

I certify that the equipment described above has been thoroughly cleaned and decontaminated of all chemical, biological and radioactive contaminants and also certify that the returned unit is safe for unprotected human contact.

By: _____
Signature Print name

Title: _____ Date: _____

Company: _____

Address: _____

Phone: _____ Fax: _____ email: _____

Form 2847

14 INDEX

A

Accessories, 67
 Racks & Trays, 69
Actual Temperature, 51
Adjusting Lid Tension, 62
Alarms, 11, 22, 77
Arrows Icon, 21
ASCII Format, 81
ASCII Text Headers, 85
Audible Alarm
 Disabling, 50
Audible Alarm Muted Icon, 20
Audible Alarms, 22
Automatic Restart, 52, 77

B

Baud Rate
 Setting the, 43
Belt Replacement, 64
Biohazards, 54

C

Calibrate Screen, 44
Calibrating RPM, 52
Calspeed, 52
CAUTION
 Symbol for, 5
Chamber Light
 Auto, 42
 On/Off, 41
Changing a Setpoint, 39
Changing Baud Rate, 44
Changing Communication Mode, 43
Changing Screens, 38
Changing the Day, 40
Changing the Time, 41
Clamp Hardware Kits, 69
Clearance, 15

COM Parameters, 80
COMM Screen, 43
Command Codes, 81
Command Format, 81
Communication Mode, 44
Continuous Operation, 11
Control Knob, 20
Copyright Notice, 5
Creating a Program, 46
Customer Satisfaction Form 6300, 13

D

Day, 40
Declaration of Conformity, 79
Decontamination, 54
Decontamination Certificate, 87
Decontamination Certificate, 88
Dimensions
 Chamber, 77
 Overall, 77
Disclaimer Notice, 5
Display Icons, 21
Display Parameters, 20
Display Screen, 36
 Changing Parameters, 36
 Modifying Setpoints, 37
Display Screens
 Changing the, 20
Door
 Opening the, 35
Door Open Icon, 20, 21, 36
Drawings
 List of, 75
Drive, 77
Drive Belt
 Removing the, 64
 Replacing the, 64

E

Edit Program, 49

Electrical Requirements, 14, 77
Environment, 14

F

Firmware
 Operation via, 11
Flask Capacities, 68
Flask Clamps, 68
 Installation of, 26
Front View, 15
Fuses, 59

G

Gas Manifold Kit, 28, 70
Germicidal Lamp, 70
Glide-Up Door, 23
Growth
 On/Off, 41
Growth Lamps, 70

H

Heater, 24, 77
Heater On Icon, 20, 21
HRS Alarm, 23
Humidity Monitor, 23, 70
HyperTerminal, 80

I

Icons, 21
Incomplete Order, 13
Index of Drawings, 75
Index of Tables, 76
Indicated Temperature, 51
Inspection
 Of Boxes, 13
Installing Flask Clamps, 26
Interior Light, 24
International Offices
 List of, 3

L

LAMP, 24
Lamps Screen, 41
LCD Display, 20, 77

Leveling
 One Shaker, No Base, 17
Lid Tension
 Adjusting the, 62
Load & Speed Graphs, 73
Locked Icon, 21
Locking Settings, 41
Loudspeaker Icon, 21, 41

M

Manual Conventions, 5
Missing Parts, 13
Modifying a Setpoint, 37, 39
Monitor Mode, 44
Motor, 77
Moving Screen to Screen, 38
Mute, 41, 50
Mute Icon, 20, 21
Muting the Alarm, 41

N

Narrow Doorway, 15
New Program, 46
NOTE
 Symbol for, 5

O

Operation Modes, 11
Overview, 11

P

Packing List
 Verification of, 13
Padlock Icon, 21, 41
Parameter Locked Icon, 20
Photosynthetic Lamps, 24, 30, 31, 70
Platform
 Available Types, 67
Platform Assemblies, 26
Platform Screw Locations, 26
Platforms, 12
POWER Alarm, 23
Power Failure, 52
Power Interruption, 52
Problems, 56

Product Returns, 58, 87
Program Running Icon, 21
Programming, 46
Programs Screen, 45

R

Refrigeration, 24
Refrigeration Schematic, 73
Remote Alarm Option, 71
Remote Programming, 81
Replacement Parts, 67
Replacing UV Light, 63
Report Request Commands, 85
Resetting the Day, 40
Resetting the Time, 41
Return Authorization, 87
Return Authorization, 88
Returned Material Authorization Number,
58, 87
Returned Material Safety Sheet, 58
RMA, 58, 87, 88
RS-232, 24, 77, 81
RS-232 Interface
Operation via, 11
RS232 Screen, 43
Run Icon, 20, 21, 50

S

Safety, 77
Safety Precautions, 35
Safety Switch, 11
Sample Program Screen, 48
Save Program, 49
Screen Parameters, 21
Service Accessibility, 25
Set Commands, 83
Setpoint Retention, 77
Setting the Baud Rate, 43
Setup Screen, 40
Changing the Day, 40
Changing the Time, 41
Locking Settings, 41
Muting the Alarm, 41
Shaking Speed, 77
Shipping Damage, 13

Slave Mode, 44
Software Interface, 24
Space Requirements, 15
Specifications, 77
SPEED Alarm, 23
Spill Pan, 23
Filling the, 33
Spills, 54
Summary Screen, 39
Changing Setpoints, 39
Symptoms, 56

T

Tables
List of, 76
Talk Mode, 44
TEMP Alarm, 23
Temperature Offset, 45, 50, 52
Setting the, 51
Temperature Offset Calibration, 51
Temperature Offset Value, 51
Temperature Range, 77
Test Tube Racks, 69
Time, 41
Timed Functions, 45
Timed Run, 11
Troubleshooting, 56

U

Unpacking
Tools for, 13
Utility Basket, 70
UV Decont
On/Off, 41
UV Lamp, 24, 70

V

Visible Alarms, 22

W

WARNING
Symbol for, 5
Warranty, 6, 7
Warranty Registration, 13
Water Reservoir, 23

Filling the, 33

Wavy Arrows Icon, 21