Meter for pH/Redox/Ion/DO/Conductivity

www.store.mt.com/search?stq=30997001

English

1 Information for Download

Scan the QR code to be directed to www.mt.com, where you will find the reference manual, user manual, and declaration of conformity.

2 Safety Information

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the power plug for damage and replace them if damaged.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts.

Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

3 Intended Use

This instrument is designed to measure pH, redox potential, conductivity, ions concentration, and dissolved oxygen (DO). It has been developed for indoor operation, and any other type of use is considered unintended. The meter complies with the IP52 standard. Never operate it in an environment subject to explosion hazards.

4 Technical Data

The instrument is operated using an AC adapter. The AC adapter is suitable for all supply line voltages ranging from 100 to 240 V AC $\pm 10\%$ and 50-60 Hz. Instrument should be operated under ambient conditions with a temperature range of 5 °C to 40 °C, and a relative humidity between 5% and 80% (non-condensing). The installation must adhere to category II, with a

pollution degree of 2, and can operate at altitudes up to 5000 meters above sea level. Please avoid powerful vibrations, direct exposure to sunlight, environments containing corrosive gases, and strong electric or magnetic fields to ensure proper operation and longevity.

5 Installation

Installing the Instrument

After unboxing all the components, connect the terminal to the meter using the high speed data (HSD) terminal cable. Press the metal strip on the connector and forcefully push the connector into the socket on both the terminal and the meter side. A click sound can be heard when a proper connection is made. Then, connect the measurement module, if there is any optional module required, and sensors. Then, plug the power supply cable into the power socket of the meter and the country specific plug on the adapter into your power outlet. Afterwards, switch ON the instrument by pressing the power button.

Note

For more information, please refer Chapter 4: Installation in the reference manual by scanning a QR code provided in the document.

Connecting a Sensor

The DES sensor is compatible with digital measurement module (DMM), while non-DES sensors must be connected to analog measurement module (AMM).

Sensor	Module
DES Sensor	DMM
Non-DES Sensor	AMM

For detailed connection instructions, please refer the user manual for the sensors.

6 Maintenance

Instrument Maintenance

Maintenance of the instrument is crucial for ensuring its accuracy and longevity. Never open the housing of the meter by yourself. The instrument requires occasional wiping with a dry and clean cloth. The housing is made of acrylonitrile butadiene styrene (ABS) and polycarbonate (PC) and is sensitive to solvents such

as toluene, xylene, and methyl ethyl ketone (MEK). Use water with mild detergent and wipe it with a clean and dry cloth. Wipe off any spills immediately.

Sensor Maintenance

Regular maintenance is crucial for extending the lifespan of sensors. Always ensure that the pH, Ion, Redox sensors are protected with the wetting cap filled with InLab® Storage Solution when not in use. The conductivity and DO sensors should be kept dry. Additionally, it is important to change the complete reference electrolyte regularly to prevent crystallization. Don't allow bubbles to form inside the sensor, particularly near the junction.

7 Transport and Storage

Store and transport the instrument within the temperature range -25 °C to 70 °C. Keep the instrument dry. If the instrument is not used for long time then keep the instrument in the carrying case to protect it from damages.

8 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this equipment may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this equipment in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this equipment. Should this equipment be passed on to other parties, the content of this directive must also be passed on to the other party.

9 Declaration of Conformity

The declaration of conformity, which details compliance with applicable standards and regulations, can be easily accessed by scanning the QR code provided in this document.

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© 05/2025 METTLER TOLEDO. All rights reserved.
Mettler-Toledo GmbH