

# Densito – the hand-held density pro

Whether you check the purity of water or the quality of salad oil: Densito always shows the result within seconds – directly in the measurement units you want to use. Using the oscillating tube method combined with precise temperature measurement and an elegant user interface, Densito gives you reliable measurement results, easily – just immerse the sampling tube, pull the trigger and read the final result!

PortableLab™ – Lab power in your hands! Densito is delivered in a special case that contains everything you need for successful measurements: Pipettes for sample handling, cleaning towels, density standards, flasks for samples and cleaning solutions, operating instructions and a CD-ROM with data transmission software, a tutorial to get started and more.



PortableLab™



**Density:**  
Absolute and relative

The density of the sample can be displayed in g/cm<sup>3</sup>, lb/gal (US) and lb/gal (IP). The specific gravity is calculated using the density of water at the measurement temperature or at any other temperature (e.g. 20 °C or 4 °C).



**Electroplating and photo:**  
Automatic temperature compensation

Temperature compensation is the key to fast and reliable control of electroplating or photographic solutions. With Densito, the procedure is very simple – just choose the sample type and perform the measurement. It is an easy matter to check different solutions one after another, quickly.



**Alcohol:**  
Replaces 4 sets of hydrometers

Densito replaces 4 complete sets of hydrometers (w/w%, v/v%, US proof, UK proof) and gives you accurate results in record time!




**Batteries:**  
Effortless acid test  
  
Densito performs battery acid tests without need of a calculator or lookup table. It displays the sulphuric acid concentration directly in weight%.



**Food:**  
The sugar expert  
  
Brix%, °Baumé or °Plato: Densito determines whichever you want.



**Petroleum products:**  
API compliant measurements  
  
The built-in API tables (product groups A, B and D) make Densito the ideal choice for quality control of crude oil, gasoline, petroleum products and lubricants.



**Custom applications:**  
Define your own calculations

If you need concentration determinations such as heavy Baumé grades, light Baumé grades, Twadell grades or milk grades, you can teach Densito the necessary calculation formula.



# Functional design

## Temperature compensation

In order to perform accurate density measurements, the result must be corrected for sample temperature. Densito not only measures the sample temperature very accurately, it also lets you choose a temperature correction coefficient before each measurement. This makes it easy to quickly measure many different types of samples.

## Plain language display

The backlit dot matrix display shows the results in large digits and offers a plain language and easy-to-navigate user interface.

## User guidance

Thanks to the pictogram-labeled keys, Densito is easy to learn and use.

## Controlled sampling

With Densito, the sampling speed can be adjusted to the task at hand: Slow for vis-  
cous samples to avoid the formation of  
air bubbles or very fast for efficient rins-  
ing of the cell.

## Clearly visible measurement cell

The most common reasons for bad results are  
air bubbles or impurities in the measurement  
cell. Densito's cell is clearly visible at all times.  
Potential problems are apparent at once.

## External sampling

Sampling can be difficult with  
samples of very high viscosity  
or samples that degas easily. In  
cases like this, Densito lets you  
connect an external syringe.

## Save and transfer your data whenever you want

Densito saves up to 1100 results  
including sample identification, meas-  
urement unit, temperature correction  
coefficient, instrument identification,  
date and time. You can transfer the  
data to a PC and printer any time  
using the infrared interface. The PC  
software to do this comes with the  
instrument.

## PortableLab™

### Refracto 30GS

portable  
high-end  
refractometer

### Refracto 30PX

portable  
refractometer

### Densito 30PX

portable  
density meter

### X-mate<sup>Pro</sup>

portable  
pH/  
multiparameter-  
instrument

# Densito 30PX

## Density meter



## PortableLab™

### Specifications Densito 30PX

**Measurement principle:** Density measurement using the oscillating tube method • **Measurement range:** 0 to 2 g/cm³, Resolution: 0.0001 g/cm³ • **Accuracy:** ± 0.001 g/cm³ • **Measurement units:** Density, specific gravity, temperature compensated density, temperature compensated specific gravity, Brix%, alcohol (w/w%, v/v%, US proof and UK proof), °Baumé, °Plato, API (tables A, B and D), % sulphuric acid (w/w%), user defined units • **Temperature:** measurement range: 0 – 40 °C, resolution 0.1 °C, display: °C or °F • **Ambient temperature range:** 5 °C – 35 °C • **Temperature compensation:** Automatic (Brix%, alcohol, °Plato, API, % sulphuric acid) or by using a user-defined temperature correction coefficient. Up to 10 temperature compensation coefficients can be stored in the instrument. • **Calibration:** With dry air or pure water (supplied) or other density standards • **Data memory:** Capacity for up to 1100 results (measured value, sample identifica-  
tion, temperature compensation coefficient, date and time) • **Display:** Backlit LCD matrix • **Interface:** Infrared Interface for data transfer to printer and PC, IrDA or RS232C protocol • **Weight:** approx. 360 g • **Batteries:** 2 x LR3 1.5 V batteries. Type AAA, approx. 90 hours battery life • **Materials:** Measurement cell: Borosilicate glass. Housing: PBT and PET; Sampling pump: PP and PTFE; Materials with sample contact: PTFE, PPS, borosilicate glass and PP.

AT	Mettler-Toledo Ges.m.b.H., A-1100 Wien Tel. +43-1-604 19 80, Fax +43-1-604 28 80
AU	Mettler-Toledo Ltd., Port Melbourne, 3207 Victoria Tel. +61-1300 659 791, Fax +61-3-9645 3935
BE	N.V. Mettler-Toledo S.A., B-1932 Zaventem Tel. +32-2-334 02 11, Fax +32-2-334 03 34
CA	Mettler-Toledo, Inc., Mississauga, Ontario L4V 1J2 Tel. +1-800-638 8537, Fax +1-905 612 8094
CH	Mettler-Toledo (Schweiz) AG, CH-8606 Greifensee Tel. +41-1-944 45 45, Fax +41-1-944 45 10
CN	Mettler-Toledo (Shanghai) Ltd., Shanghai 200233 Tel. +86-21-6485 0435, Fax +86-21-6485 3351
CZ	Mettler-Toledo, s.r.o., CZ-10000 Praha 10 Tel. +420-2-72 123 150, Fax +420-2-72 123 170
DE	Mettler-Toledo GmbH, D-35353 Giessen 11 Tel. +49-6-41 50 70, Fax +49-6-41 507 319
DK	Mettler-Toledo A/S, DK-2600 Glostrup Tel. +45-4327 0800, Fax +45-4327 0828
ES	Mettler-Toledo S.A.E., E-08908 L'Hospitalet del Llobregat Tel. +34-93-223 76 00, Fax +34-93-223 76 01
FR	Mettler-Toledo s.a., F-78222 Viroflay Cedex Tel. +33-1-3097 1717, Fax +33-1-3097 1773
HK	Mettler-Toledo HK Ltd., Cheun Sha Won Kowloon Tel. +852-2744 1221, Fax +852-2744 6878
HR	Mettler-Toledo, d.o.o., HR-10000 Zagreb Tel. +385-1-292 06 33, Fax +385-1-295 81 40
HU	Mettler-Toledo, Kft, H-1139 Budapest Tel. +36-1-288 40 40, Fax +36-1-288 40 50
IN	Mettler-Toledo India Pvt. Ltd. 400 072 Mumbai Tel. +91-22-285 708 08, Fax +91-22-285 750 71
IT	Mettler-Toledo S.p.A., I-20026 Novate Milanese Tel. +39-2-333 321, Fax +39-2-356 29 73
KZ	Mettler-Toledo C.A. 480009 Almaty Tel. +7-3272-980 834, Fax +7-3272-980 835
MY	Mettler-Toledo (M) Sdn. Bhd. 40150 Shah Alam Tel. +60-3-7845 5773, Fax +60-3-7845 8773
MX	Mettler-Toledo S.A. de C.V., México C.P. 06430 Tel. +52-55-554 75700, Fax +52-55 554 12366
NL	Mettler-Toledo B.V., NL-4004 HA Tiel Tel. +31-344 63 8363, Fax +31-344 63 8390
PL	Mettler-Toledo, Sp. z o.o., PL-02-924 Warszawa Tel. +48-22-651 92 32, Fax +48-22-651 71 72
RU	Mettler-Toledo C.I.S., 101000 Moscow Tel. +7-095-921 92 11, Fax +7-095-921 63 53
SE	Mettler-Toledo AB, S-12008 Stockholm Tel. +46-8-702 50 00, Fax +46-8-642 45 62
SG	Mettler-Toledo (S) Pte. Ltd., 139959 Singapore Tel. +65-689 000 11, Fax +65-689 000 12
SK	Mettler-Toledo s.r.o, SK-82 104 Bratislava Tel. +421-43 42 74 96, Fax +421-43 33 7 190
SI	Mettler-Toledo, d.o.o., SI-1236 Trzin Tel. +386-1-562 18 01, Fax +386-1-562 17 89
TH	Mettler-Toledo (Thailand), Huy Kwang Bangkok 10310 Tel. +66-2-723 0300, Fax +66-2-719 6479
UK	Mettler-Toledo Ltd., Leicester, LE4 1AW Tel. +44-116-235 70 70, Fax +44-116-236 63 99
US	Mettler-Toledo, Inc., Columbus, OH 43240-2020 Tel. +1-614-438 4511, Fax +1-614-438 4900
L.A.	Mettler-Toledo, Latin America, Miami, FL 33122-1528 Tel. +1-305-351 2442, Fax +1-305-351 2425

Subject to technical changes.  
© 06/2003 Mettler-Toledo GmbH  
Printed in Switzerland 51724230A

METTLER TOLEDO

PortableLab™





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.