CLD-01

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

DESCRIPTION

CLD-01 is a double flexible-arm illuminator, with white LED source. It offers low power consumption, low operating temperatures, and a colour temperature similar to sunlight.

SPECIFICATIONS

Light source:	High power 14W LED
Chromaticity coordinates (x,y):	0.320, 0.320 (pure white)
Illuminance:	70000 lux (at 10 cm distance)
Power supply:	External power supply provided. Input: 110/240 Vac 50/60 Hz; output: 6 Vdc 5A
Accessories:	Provided with focusable lens on each arm. Also available: CL-17 : lens pair with built-in linear polarizer

INSTRUCTIONS

Connect the power supply (included) to the jack located on the back of the instrument.

Bend the flexible arms, aiming the light in the desired direction.


The front lenses can be moved in order to optimally focus the output light on the sample.

The intensity of the light can be adjusted by using the dedicated knob.

Clean with a soft dry cloth. Do not use solvents.

WARNING: Do not stare directly into the output lenses.

Art.13 Dlsg 25 july 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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