

SELECT FUSION OPERATOR'S MANUAL

ISSUE – F

17.03.08



INTRODUCTION

The *SELECT FUSION* Range.

Provides dual quality laboratory water for analytical and Life science applications direct from a potable water supply. The units have been designed to provide both Grade I+ and Grade II water quality according to **BS EN ISO 3696 “Water For Analytical Laboratory Use”**

There are 3 models in the range: -

SELECT FUSION 40/80/160

Producing 6, 12 and 24 litres/hr*

- *Outputs based on 25 °C, 60 psi (4Bar) feedwater.

Reverse Osmosis provides the principle water treatment technology. Pre-filtered mains water is supplied to the reverse osmosis module. Purified water passes through the membrane as permeate, the rejected impurities are flushed to waste via the Drain (concentrate) stream. The permeate is further polished via a primary Puripac PP8 cartridge and fed to the internal 20 litre storage tank.

The quality of water available for dispense from the tank will be at Grade II level. To achieve Grade 1+ quality, the stored water is recirculated in a loop around a second polishing cartridge and is also irradiated with a dual wavelength (185/254nm) Ultra violet lamp to photo-oxidise any residual Total Organic carbon and deactivate any bacteria present.

To maintain purity the high quality dispense valve is fitted with a sterile 0.2um filter to prevent ingress of air borne particles and bacteria from contaminating the dispensed water.

The system status is constantly monitored for: -

- Water quality - MΩ-cm
- Temperature - °C
- Flowrate - L/min
- Total Organic Carbon - ppb as C

These parameters can be viewed via the graphic display.



In all cases where this symbol is used, documentation needs to be consulted in order to find out the nature of the potential HAZARD and any actions that have to be taken.

INTENDED USE



The unit is not for use in explosive atmospheres.

The unit is for indoor use only and is not to be washed down.

The unit is not classed as a Medical Device.

The unit should only be fed only from a potable drinking water supply as covered by the EC Directive “**Relating to the quality of water intended for human consumption**”. Refer to Section 1.3.1 for details of feedwater conditions

The unit should only powered from a single-phase electrical supply. Refer to Section 1.2.1 for details of required electrical supply.

CONTENTS

SECTION 1 – SPECIFICATION

- 1.1 Standard Features
- 1.2 Electrical Specifications/Connections
- 1.3 Feedwater Specification
- 1.4 Treated Water Specification
- 1.5 Water Services Connections
- 1.6 Weights and Dimensions
- 1.7 Environmental
- 1.8 Standards Applied
- 1.9 Caution and Warning Statements
- 1.10 General Safety Information

SECTION 2 – PROCESS DESCRIPTION

- 2.1 Operational Overview
- 2.2 Disinfection

SECTION 3 – INSTALLATION AND COMMISSIONING

- 3.1 Environment
- 3.2 Unpacking
- 3.3 Installation Kit
- 3.4 Wall Mounting Instructions (Optional)
- 3.5 Raw Water Pre-filter Assembly
- 3.6 Cartridge Packs
- 3.7 Preliminary Checks
- 3.8 Installation
- 3.9 Cartridge Pack Location
- 3.10 0.2 um Point of Use Filter Location.
- 3.11 Handling the Unit.

SECTION 4 – OPERATIONAL PROCEDURES – QUICK REFERENCE

- 4.1 Operator Interface
- 4.2 Starting the Unit for the first time.
- 4.3 Starting during normal operation
- 4.4 Stopping the unit during normal operation.
- 4.5 Stopping the unit in an emergency
- 4.6 Dispensing water from the unit
- 4.7 Standby Operation

SECTION 5 – OPERATING PROCEDURES – DETAILED REFERENCE

- 5.1 Control Philosophy
- 5.2 Menu Structure
- 5.3 Clean Routine

- 5.4 Dispense Options
- 5.5 User Settings Menu
- 5.6 Processing Screen Information
- 5.7 Safety Features
- 5.8 Error and Alarm Messages

SECTION 6 – DISINFECTION

- 6.1 Procedure

SECTION 7 – MAINTENANCE

- 7.1 Replacing the pre-filter cartridge
- 7.2 Replacing the Puripac PP-8 cartridge
- 7.3 Replacing the NCP media Pack
- 7.4 Replacing 0.2um Point of Use Filter
- 7.5 Replacing the Ultra Violet Lamp
- 7.6 Replacing Fuses/s
- 7.7 Replacing main PCB Lithium Battery
- 7.8 Cleaning External Surfaces
- 7.9 Preventative Maintenance Guide.

SECTION 8 – TROUBLE SHOOTING

- 8.1 General Fault Conditions
- 8.2 Alarm Messages

SECTION 9 – CONSUMABLES AND SPARES

- 9.1 Consumables
- 9.2 Recommended Spares
- 9.3 Serviceable Spares
- 9.4 Disposal of Consumables
- 9.5 Accessories
- 9.6 Consumable Shelf life

SECTION 10 – CE DECLARATION

SECTION 11 – CONTACT US

- Appendix A** Pushfit Water Connectors
- WEEE Declaration
- Material Safety Data Sheets:
- Process Flow diagram

HISTORY

Issue	Date of Change	Description of Change	Approved By
A	18.03.05	Manual Release	<i>M. R. Bosley</i>
B	08.06.05	Text changes	<i>M. R. Bosley</i>
C	09.06.06	CE D of C Appendix E 7.6.06	<i>M. R. Bosley</i>
D	06.11.06	Contents page add 7.10, new Flow diagrams 40/80/160 & 40BP/80BP/160BP	<i>M. R. Bosley</i>
E	05.09.07	CE D of C Appendix F WEEE Declaration	<i>R. King</i>
F	07.03.08	7.9 & 9.1 Consumables - Part No R041146 amended to R011146 for External CI Carbon Filter. Addition of Section 4.7 Standby Operations	<i>J. Hallett</i>

SECTION



SPECIFICATION

This section provides details of the **SELECT FUSION** Unit electrical, mechanical and process specification as well as relevant standards complied with and general safety information. The section is divided into the following sub-sections.

SECTION CONTENTS

- 1.1 Standard Features
- 1.2 Electrical Specifications / Connections
- 1.3 Feedwater Specifications
- 1.4 Treated Water Specification
- 1.5 Water Services Connections
- 1.6 Weights & Dimensions
- 1.7 Environmental
- 1.8 Standards Applied
- 1.9 Caution & Warning Statements
- 1.10 General Safety Information

1.0 **SPECIFICATION**

1.1 **STANDARD FEATURES**

- Grade I+ and Grade II water quality available direct from mains water supply.
- Flowrate, T.O.C., Temperature and Water quality monitoring.
- Wall or bench mountable.
- System status alarms.
- Automatic cleaning/sanitisation procedure.
- Easy change disposable cartridge packs.
- Optional boost pump variant for users with low mains pressure.
- Multiple water take offs.
- Microprocessor controlled with graphical user interface.

1.2 **ELECTRICAL SPECIFICATIONS / CONNECTIONS**

1.2.1 **Mains supply**

Electrical Supply	KW Rating	Current Draw (Amps)	External Fuse Rating
Single phase 100-240v ±10%, 50 Hz and Earth	0.1	1.0	5 Amps (conforming to BS1362)

1.2.2 **Fuse Rating / Type**

IEC module fuse type: - (20x5mm), T5AH250V, conforming to IEC 127

Number per unit = 2



The mains supply is double pole/live & neutral fused.

1.2.3 **Serial Port Connection Details (see Fig-13)**

Application: Uploading of new software

Connection to external printer or PC

Connector Type: RJ45 to 9 pin D-type socket. Contact Purite for supply of Serial Port lead.

1.2.4 Alarm Port connection Details (see Fig-13)

Application: Volt free alarm output e.g. connection to BMS system

Connector type: Din 3-pin Plug (max rating 34Vdc / 24ac, 1 amp).
Max. Lead length: 10 mtrs. Contact Purite for supply of Alarm lead.

1.3 FEEDWATER SPECIFICATION.

1.3.1 Feedwater Flowrates

Model	Max Feed Flowrate L/hr	Drain Flowrate L/hr	Output @ 10°C (Ltr/hr)	Output @ 25 °C (Ltr/hr)
Fusion 40	66	60	3.6	6
Fusion 80	72	60	7.2	12
Fusion 160	84	60	14.4	24

1.3.2 Feedwater Quality

The unit has been designed to only operate on a Potable water supply conforming to current **EC Directive “Relating to the quality of water intended for human consumption”**, but with the following additions.

Pre Filtration	Filtered to 5 Micron
Free Chlorine	< 0.1 ppm as Cl ₂
Total Dissolved Solids (Max)	1000 ppm
Temperature	1-40°C (33.8 – 103.9° F)

1.3.3 Feedwater Pressures

Feedwater Pressure (Max)	(90 psi)
Feedwater Pressure (Min)	(30 psi)*
Design Pressure	(60 psi)

- At 30 psi rated outputs specified in 1.3.1 will be reduced by up to 50% when compared to design pressure.

1.4 TREATED WATER SPECIFICATION

Parameter	High Purity Dispense	Tank
BS EN ISO 3696	Grade 1+	Grade II
ASTM D1193-99e1	Type-1	Type III
Resistivity @ 25 °C	15-18.2MΩ-cm	1-15 MΩ-cm
PH	6-8	6-8
TOC (ppb as C)	<30	<50
Bacteria	<1 cfu/ml	<99%rejection
Particles	0.2 um Free	N/a
Organics	<0.001AU@254nm	N/a
Residual Solids (ppm)	<0.1	<0.5

1.5 WATER SERVICES CONNECTIONS

<u>Description</u>	<u>Connection</u>
Drain	8mm Pushfit
Mains Water In	8mm Pushfit
Overflow	12mm Pushfit

1.6 WEIGHTS & DIMENSIONS

Unit	Weight kg	Height (mm)	Width (mm)	Depth (mm)
Fusion 40/80/160	22 (Dry) 42 (Working)	630	440	548
Fusion 40/80/160 with boost pump	24 (Dry) 44 (Working)	630	440	548

1.7 ENVIRONMENTAL

Room storage and operating temperature range	5 to 40°C (41-104°F)
Relative Humidity	30 to 80%
Max Altitude	2000m
Transport and Storage temperature (limited by RO membranes)	-5 to 85°C (with frost protection liquid – 40 to 85°C)
RFI/EMI Radiation	The 'EMC' environment must be within the limits to which the unit has been tested, see section 1.7. Care must be taken not to have sources of RFI/EMI, which are liable to cause electromagnetic disturbance to the unit. If the unit is affected by such disturbances, the sources should be suppressed or relocated.

1.8 STANDARDS APPLIED

EMC	<p>BS EN 61326:1998/IEC 61326-1 :1997;Class A Electrical equipment, for measurement, control and laboratory use EMC requirements.</p> <p>BSEN 61000-3-2:1995, Incorporating Amendments 1 & 2. Mains Harmonic Emissions.</p> <p>BS EN 61000-3-2:1995, Incorporating Amendment 1 Mains Flicker Emissions.</p>
LVD	<p>BS EN 61010-1:1993 , Incorporating Amendment 1 Safety requirements for electrical equipment, for measurement control and laboratory use.</p>

1.9 CAUTION AND WARNING STATEMENTS

- These instructions provide information to ensure safe and continuing operation of the equipment and that safe working practices can be adopted as required. The manual should be read and understood before the equipment is placed into service.
- **Purite Limited** reserves the right to make engineering refinements to the equipment that may not be described herein. Any questions that cannot be answered specifically by these instructions should be addressed to **Purite** or their agents for response.
- **Purite** will not accept any responsibility for any equipment supplied or the actions of such equipment or associated system when the customer has made a modification that is considered by Purite to question the integrity of the original design philosophy.
- If the unit's performance becomes impaired and any remedial work appears to be outside the scope of this manual, then seek advice from Purite's **Service Department, Tel. +44 (0) 1844 211555**. Quoting the unit's serial number.
- The unit must not be dismantled unless carried out by Purite Service Department personnel or authorised trained personnel. On no account must the unit be connected to the electrical supply with the top control cover removed.
- Always refer to the Safety Data Sheets, in Appendix A, before handling any of the recommended cleaning disinfection adaptors or consumable cartridges.
- There is the potential for sensitive equipment/devices located in close proximity to the Fusion unit to be affected by electromagnetic or other interference generated from the unit. If affected by interference the relevant equipment/device should be relocated.
- The use of mobile phones in close proximity to the Fusion unit should be avoided where possible.

- The 'Caution' symbol is used throughout this manual to highlight where particular care must be taken to ensure the safety of the operator, and the protection provided by the equipment is, not impaired.

1.10 GENERAL SAFETY INFORMATION

- **Explanation of symbols and references**



Danger This symbol refers to any immediate dangers that may threaten the safety and life of persons. Failure to observe these notices will have severe consequences on health and safety, including life-threatening injuries.



Warning This symbol refers to a possible danger that threatens the safety and life of persons.



Caution This symbol refers to a possibly hazardous situation. Failure to observe these references may result in minor injuries and/or damage to property.

i

This symbol points out important information for working with the system in a proper manner. Failure to observe these references may result in malfunctions in the system or impact on the environment.

- **Additional safety requirements**

Country-specific requirements standards and regulations must be observed.

- **Usage in accordance with intended purpose.**

The **SELECT FUSION** units are used to purify, potable water. The units must only be operated with water supplied in accordance with the quality described in **Section 1.3** and operated in accordance with the parameters specified in **Section 1.7**. The units must not be operated unless in proper working order. Any malfunctions must be rectified immediately.

- **Operating staff**

Only persons who have read and understood these Operation Instructions should be permitted to operate the unit. When operating the units, it is particularly important to observe the safety information strictly.

- **Residual dangers**



Danger

Electrical Shock

Do not touch electrical components with wet hands.

Before performing tasks on parts of electrical system, isolate the system from the electrical power supply.

Mechanical force

Some parts of the system could be under pressure of up to 60 psi. Always release the pressure from the unit before repairs and maintenance tasks are carried out.

- **Bringing the system to a stop in the event of an emergency**

- Turn off the electrical supply and/or remove mains plug.
- Shut off the water supply.

After remedying the fault:

- Open the water supply.
- Turn on the electrical supply.
- Restart/operate the unit via the main front panel.

- **Safety information for maintenance tasks**

The operator must take care to ensure that authorised and qualified professionals who have been sufficiently informed for the task at hand by thoroughly studying the Operating Instructions perform all maintenance, inspection and assembly tasks. Professionally trained staff must properly perform these tasks.

The system must be shut down and protected from being placed in operation again unintentionally before all repair and maintenance tasks have been completed. It is essential to observe the procedure described in these Operating Instructions for shutting down the system.

Before beginning tasks on the electrical equipment of the system, a check must confirm that power has been disconnected from the corresponding section of the system. In addition, the system must be secured to prevent it from being turned on again unintentionally.

- **Disposing of system parts and operating materials**

When they need to be discarded, consumables must be disposed of according to local requirements. Refer to Appendix A – WEEE Declaration for disposal of electrical and electronic equipment supplied by Purite.

- **Unauthorised conversion and manufacturing replacement parts**

Conversion or modification of the system is only permitted with the approval of the manufacturer. The same applies to making changes in the programming for the control system. Original replacement parts and accessories authorised by the manufacturer enhance safety. Use of other parts will void the warranty.

- **Warranty claims and liability**

This product corresponds to the state of the art and was designed and manufactured in accordance with applicable rules of the technology, after which it was subjected to a quality control process. If there should nevertheless be any grounds for complaint, please direct requests for replacement to the manufacturer of this product in accordance with the general terms and conditions of sale and delivery.

Manufacturers Name: Purite Limited
Bandet Way
Thame
Oxon OX9 3SJ
England

Telephone: +44 (0) 1844 217141

Fax: +44 (0) 1844 218098

SECTION

2

PROCESS DESCRIPTION

This section provides general information regarding the process of purification utilised by the ***SELECT FUSION***. The section is divided into the following sub-sections.

SECTION CONTENTS

- 2.1 Operational Overview
- 2.2 Disinfection

SELECT FUSION



Fig-1
(Front view)



Fig-2
(Side view)



Fig-3
(Rear view)

2.0 PROCESS DESCRIPTION

2.1 OPERATIONAL OVERVIEW

All external hydraulic and electrical connections are positioned at the rear of the unit, to protect them from accidental damage.

Potable feedwater enters the unit via an external 5 micron, carbon impregnated, prefilter; protecting the inlet solenoid valve from particulates and the reverse osmosis membrane from possible chlorine degradation. The filter removes up to 95% of all particles > 5 micron and will dechlorinate to < 0.1 ppm on 0.5 ppm feedwater supply.

The Pre-filtered, water enters the reverse osmosis module/s. Purified water passes through the membrane and is termed; the permeate (85-95% of salts removed). The water containing the rejected salts passes across the membrane and out of the RO module to drain and is termed, the concentrate.

The concentrate flowrate is controlled by a fixed orifice restrictor located in the Drain line and is set at approximately 60l/hr.

The, permeate is further purified by a Puripac-8 cartridge containing mixed ion-exchange resin and activated carbon. The deionised water then feeds directly into the integral 20 litre storage tank. Water can be dispensed from this tank via the front mounted bib tap.

The quality of water drawn from the tank will have a quality conforming to Grade II standard.

To produce Grade 1+quality, water is drawn from this tank by a recirculation pump at a rate of up to 1 litre/min and further polished via an NCP media pack containing high purity semiconductor grade nuclear grade resin and activated carbon. The polished water is then further treated by a dual wavelength Ultra violet light that photo-oxidises any dissolved organics, thus reducing the TOC content of the water and deactivates any living bacteria present.

The high purity water can be dispensed from the system via the automatic dispense valve which is again protected from external contamination by a sterile 0.2 um, point of use filter.

The microprocessor control system constantly monitors the system performance and water qualities. If at any time the preset parameters are exceeded the system will respond with various alarm messages prompting actions. Under certain serious alarm conditions the unit will automatically shutdown to prevent injury and to protect the unit from damage. Refer to **Section 5.3** for details of all operating alarm conditions and safety features.

2.2 DISINFECTION

Regular disinfection of the unit is recommended to maximise the life of the RO membrane/s and ensure constant performance. If the unit is regularly used disinfection would be recommended every 6 months.

The process of disinfection is semi-automatic and requires only the substitution of the media pacs with a cleaning adaptor cartridge, which contains the required disinfecting chemical agent.

The graphical display will give prompts at every stage of the disinfection routine to guide

Refer to **Section 6.0** for details on how to carry out disinfection on the **SELECT FUSION** Unit.

SECTION

3

INSTALLATION

This section provides the recommended method of installation for the ***SELECT FUSION***. It is divided into the following sub-sections.

SECTION CONTENTS

- 3.1 Environment
- 3.2 Unpacking
- 3.3 Installation Kit
- 3.4 Wall Mounting Instructions (Optional)
- 3.5 Raw water Pre-filter Assembly
- 3.6 Cartridge Packs
- 3.7 Preliminary Checks
- 3.8 Installation
- 3.9 Cartridge Pack Location
- 3.10 0.2µm Point of use Filter Location
- 3.11 Handling the Unit.

3.0 INSTALLATION

3.1 ENVIRONMENT

The unit should be installed in a clean and dry environment as specified in **Section 1.6**. The unit can either be wall or bench mounted. Refer to **Section 3.4** for details on how to install the wall mounting kit.

3.2 UNPACKING

Remove all packaging materials and ensure the following items have been provided.

- Select Fusion Unit
- Installation kit
- Select Wall mount kit (Optional)
- Raw water pre-filter housing and wall mount bracket, c/w filter element.
- PP-8 Puripac cartridge
- NCP media pack
- User Manual

3.3 INSTALLATION KIT

The installation kit comprises of the following items:-

- **3 off security clips**
- **1 off 12- 8mm stem reducer**
- **2 off 1/2" hose tail fittings**
- **1 off 1/2"- 1/4" reducing bush**
- **1 off 12mm stem elbow**
- **3 off 8mm stem elbow**
- **1 off 8mm-1/4" straight adaptor**
- **4 mtrs 1/2" blue braided hose**
- **7 mtrs 8mm natural, nylon tubing**
- **3 mtrs 12mm natural, nylon tubing**
- **1 roll of PTFE sealing tape**



Fig-4 Installation Kit

3.4 WALL MOUNTING INSTRUCTIONS (OPTIONAL)

Check **Section 1.6** for working weights of the Select Fusion units before attempting to wall mount. It is vital that the supporting wall is sound and is of load bearing, brick or concrete construction and of adequate thickness to secure to and carry weights of potentially up to, **44Kgs.**

Wall fixings must be of good quality and of the largest size to match the 10mm diameter holes provided in the wall mount brackets and of adequate length to ensure a strong retention. All fixing points must be used and be sound and tightened before attaching the unit. If unsure contact Purite Customer Service Department for advice.


3.4.1 Fitting

A Select Fusion Unit can be heavy and difficult to handle at installation; therefore, it is advisable to have additional personnel and/or suitable equipment to assist with the installation of this unit.

Note:

When required, Select units are specifically built to fit a wall mounting bracket, therefore, if at a later date it is required to remove the unit to stand on a bench top, it **MUST** be fitted with 2 off Foot Brackets. Please contact the Purite Service Department for further information.)

FIG-5 WALL MOUNT BRACKET



3.4.2 Installation

- 3.4.2.1** Fix the wall-mounting bracket to a suitable wall with reference to the Safety First section; ensuring the bracket is upright and square. (Tip: Make the central fixing first then when square and level, mark through the other holes using bracket as a template.).
- 3.4.2.2** Remove plastic cover on left hand side of unit by loosening screws at front edge, ease away and remove backwards to disengage from supporting lugs at rear.
- 3.4.2.3** With the cover removed, place the unit onto the bracket and align the holes on the rear left hand side chassis panel with the threaded holes located on the bracket as indicated in (Fig- 5). (**Tip: connect up all water service lines before fixing the unit to the bracket**).
- 3.4.2.4** Use two M8x20 Hex Set Screws with washers supplied to retain the unit. Tighten only to finger tight at this stage (Figs 6 & 7).




Fig-6 Rear Fixings




Fig-7 Top Fixing

- 3.4.2.5** Now open the front door of the unit to locate the two fixing holes (one each side), which align with holes in wall mounting bracket in base of unit.
- 3.4.2.6** Bolt through both with two M8x50 Hex Set Screws, washers and nuts supplied (Figs 8 & 9 for closer detail of fixing points).




Fig-8 Lower mounting holes




Fig-9 LHS lower mounting point

3.4.2.7 Finally tighten all securing fasteners and refit side cover.

3.5 Raw water Pre-filter Assembly

As pre-treatment, the Select Fusion comes supplied with an external pre-filter, which should be connected, between the mains water supply and the unit.

The filter is designed to remove particulates and chlorine from the incoming supply, which if left untreated, could irreversibly damage the reverse osmosis membrane/s.

The filter should under normal running conditions be replaced every 3 months. (Refer to Section 7.1 for details)

The assembly comes complete with filter housing and one filter element. The Housing can be wall mounted by using the “L” shaped bracket.

To fit the filter, simply unscrew the bowl in a clockwise direction. Place the filter into the bowl ensuring that it locates over the spigot at the bottom.

When refitting the make sure the black O-ring seal is present on the top of the bowl. If the filter leaks check that the O-ring has not moved out of position.

Fig-10 Pre-filter Assembly

Pre-filter housing, complete with Pressure gauge, filter element, wall mount bracket and fixing screws.



3.6 Cartridge Packs

The Select Fusion is supplied with two consumable media packs.

Fig-11 Rear view of Cartridges



Fig-12 Select Fusion media packs



3.7 PRELIMINARY CHECKS

There are three main items to check before installing the unit:

- ❖ **The electrical supply**
- ❖ **The Potable feedwater supply**
- ❖ **Drainage**

3.7.1 Electrical Supply



A 100-240V, 50 Hz, Single Phase, earthed, supply is required and should be provided via a 3-pin socket or a switched fused spur fitted with a 5-amp fuse. The socket or switched spur must be easily accessible, to provide a suitable means of electrical isolation.

Connection Detail

1. Mains Inlet socket
2. Level control
3. Serial Port
4. Alarm Output
5. Fuse holder




Fig-13 Electrical connections

3.7.2 Internal Level Control.

This is an internal connection which links the level switches in the integral tank with the main control assembly. **Do Not unplug this connector.**

3.7.3 Serial Port.

The Serial Port connection can be used to link directly to a PC or Laptop for the downloading of stored operational data or for uploading of new software code.

The uploading of new software must only be carried out by Purite Engineers or Purite Approved Personnel.

3.7.4 Alarm Output

This volt free output can be used to connect to a Building Management System (BMS) to provide a remote common alarm warning. Refer to **Section 1.2.4** for details of connector type. An Alarm lead can be supplied, contact Purite Service Dept. for details.

3.7.5 Potable Feedwater Supply

A potable feed water supply should be provided, terminated with a suitable isolation valve. A mains pressure exceeding 90psi should be reduced, using a suitable pressure regulator valve, to 60 psi (Design Pressure), Refer to **Section 9.4** for details of Pressure regulator. The minimum pressure the unit will operate on is, *30 psi. Refer to **Section 1.3.2** for details of feedwater quality requirements.

NB * At 30 psi the rated design output of the unit will be reduced by as much as 50%

3.7.6 Drainage

A suitable, unrestricted, drain is required, capable of handling up to 60 l/hr. Refer to **Section 1.3.1 Feedwater Flowrates**, for Drain flows.

3.8 INSTALLATION.

Services connection panel details



Fig-14 Rear view of unit showing all services connections




Fig-15 Water services connections

3.8.1 Connection Details

The **Select Fusion** unit is fitted with 8mm push fit connections on the rear of the unit for the “**Mains Water In**” and “**Drain**”. The “**Overflow**” connection is a 12mm push fit connection.

Remove red transit plugs from, “**Mains Water in**”, “**Drain**” and “**Overflow**”. To remove transit plug push in black collet and withdraw red transit plug (see appendix-A)


Mains Water In (Via pre-filter)

Fig-16 Mains Water In Connection



The “**Mains Water In**” port should be connected directly to the outlet of the 10” blue pre-filter housing using the ½” blue tubing adapted down to the 8mm tubing supplied.

Fig-17 Pre-filter Connection To Unit



Using 1/2" blue reinforced tubing, fit 5/16 - 3/8" hose tail to one end and secure with security clip. Push the plain end into an 8mm stem elbow and push the stem end into the "**Mains Water In**" connection on the rear of the unit. Using the adaptors supplied in the installation kit connect up to mains water supply, ideally via an isolating valve.

3.8.2 Drain

Fig-18 Drain Connection



Using an 8mm stem elbow push the plain stem part into the "**Drain**" port. Then into the elbow, push in one end of a length of 8mm tubing. Run the tubing to a suitable drain, within 2 meters of the unit ensuring the tubing is unrestricted and does run more than 1 meter above the unit

3.8.3 Overflow

Fig-19. Overflow Connection detail



Using a 12mm Stem elbow supplied, push the plain stem end into the "**Overflow**" connection on the rear of the unit. To the female end of the elbow, push in a length of the 12mm tubing. Run the tubing to a suitable drain, within 2 meters of the unit, ensuring the tubing is unrestricted.

3.9 Cartridge Packs Location



Fig-20 Cartridge Pac positions

The Fusion unit has two media packs fitted. A PP-8 cartridge which is positioned in the left hand compartment and an NCP media pack which is fitted to the right hand compartment.

NB: Ensure that the two black sealing plugs have been removed from the top ports on the cartridge. See Fig 11.

For details of how to fit new replacement cartridge packs Refer to, **Section 7.2 & 7.3** For re-ordering details refer to **Section 9.1 Consumables and Spares**

3.10 0.2um Point of Use Filter



Fig-21 Point of Use Filter

Before starting the unit ensure that the Point of use filter has been fitted to the dispense point

For details on how to install the filter refer to **Section 7.4** in the Maintenance section of this manual.

For re-ordering details refer to **Section 9.1**, in the Consumables and Spares Section.

3.11 HANDLING THE UNIT

The following points should be considered when handling the unit:

- The dry weight of the unit can be as much as 24kgs. If the unit is to be moved over short distances, two people should be employed.
- If the unit is to be transported over some distance by foot then a trolley or other suitable device should be used.
- Never pick the unit up by the side covers; they are not designed to be load bearing. Always support the weight of the unit by holding it by its base.
- Always ensure the tank is fully drained down and the two cartridge media packs have been removed before transporting.

SECTION

4

OPERATING PROCEDURES QUICK REFERENCE

This section provides a quick reference for operating the ***SELECT FUSION***. For more detailed instructions and full procedures refer to Section-5. This section is divided into the following sub-sections.

SECTION CONTENTS

- 4.1 Operator interface
- 4.2 Starting the unit for the first time.
- 4.3 Starting during normal operation.
- 4.4 Stopping the unit during normal operation
- 4.5 Stopping the unit in an emergency
- 4.6 Dispensing water from the unit
- 4.7 Standby Operation

4.1 OPERATOR INTERFACE




Fig-22 Keypad display

4.1.1 Keypad




Press Keypad Button	Function	Button No.
	<ul style="list-style-type: none"> - Main Stop/Start button. - General Function button 	①
	<ul style="list-style-type: none"> - Menu select Button - General Function button 	②
	<ul style="list-style-type: none"> - Point of use dispense button 	③
	<ul style="list-style-type: none"> - Up / Down buttons, moves cursor to allows selection of different parameters and to increase/decrease values, within a menu screen 	④ & ⑤





4.2 STARTING THE UNIT FOR THE FIRST TIME



Before starting please ensure that the unit has been installed and set up as detailed in Section 3.0. Failing to set the unit up in the correct manor may cause irreversible damage.

- Ensure you have completed the following installation checks.
- Mains Water In, Drain and Overflow connections are correct as detailed in Section 3.8. **(Water supply off)**
- Power supply is connected as per Section 3.7.1 **(Switched Off)**
- Pre-filter is fitted in the external filter housing, as per section 3.5
- Both Cartridge Packs are fitted, as per Section 3.9
- 0.2 um Point of use Filter is fitted, as per 3.10
- Unit is securely mounted.
- Operators manual has been read and understood.

Display	Actions	Comments	Step
	<ul style="list-style-type: none"> • Switch power on to Unit. • Turn on water supply. • Check for leaks. 	When power is first switched on to unit Select Fusion screen displayed for a few seconds, displaying software version.	1
	<ul style="list-style-type: none"> • Press buttons ④ & ⑤ to set values. • Press button ② to change from hours to mins. • Once set select EXIT by pressing button ① 	Follow on screen prompts to set time.	2
	<ul style="list-style-type: none"> • Press buttons ④ & ⑤ to set values. • Press button ② to change from days to mths to yr. • Once set select EXIT by pressing button ① 	Follow on screen prompts to set date.	3

Display	Actions	Comments	Step
 <p>NCP Pac replacement Due: 00/00/00 EXIT RESET</p>	<ul style="list-style-type: none"> Press button ② to RESET, NCP pac replacement date 6 mths in future. Then select YES to confirm. 	<p>Follow on screen prompts to set future NCP pac replacement date.</p>	4
 <p>PP8 replacement due: 00/00/00 EXIT RESET</p>	<ul style="list-style-type: none"> Press button ② to RESET, PP8 pac replacement date 6 mths in the future. Then select YES to confirm. 	<ul style="list-style-type: none"> Follow on screen prompts to set future PP8 pac replacement date. On selecting RESET display reverts to “POWER ON screen. 	5
 <p>POWER ON 14:25:12 02/03/05 START MENU</p>	<ul style="list-style-type: none"> To start the unit select START. For details of menu functions refer to section 5 of this manual. 	<ul style="list-style-type: none"> On selecting START screen will change to PROCESSING screen. After a few seconds the screen will display the process mimic. The integral tank will start to fill after approx. 10-15 mins. Ignore all alarms at this stage. 	6
 <p>Pre Qual 11 MΩ PROCESS MIMIC</p>	<ul style="list-style-type: none"> Pressing buttons ④ & ⑤ will display:- Out Qual Pre Qual U Qual Temp Flowrate TOC 	<ul style="list-style-type: none"> By pressing buttons ① & ② will return the display to the PROCESSING screen. Refer to section 5 for details of displayed parameters. 	7

Display	Actions	Comments	Step
<div style="border: 1px solid black; background-color: #90EE90; padding: 5px; width: fit-content;"> <p>PROCESSING</p> <p>Pre Qual 11 MΩ</p> <p>STOP MENU</p> </div>	<p>Press button 2 to select MENU features or press button 1 to STOP the unit.</p>	<p>See section 5 for details of menu features.</p>	<p>8</p>

4.3 STARTING THE UNIT DURING NORMAL OPERATION.

- Setting of the time and date is only required when the unit is first started from new.
- To start the unit once commissioned follow instructions from step 6 above.

4.4 STOPPING THE UNIT DURING NORMAL OPERATION

- To stop the unit during normal operation simply press button, **1**, twice.

4.5 STOPPING THE UNIT IN AN EMERGENCY.

- If you need to stop the unit in an emergency as a result of a severe leak or fire, then switch the unit off at the power supply and turn off the water supply. **(Do not put yourself at risk)**

4.6 DISPENSING WATER FROM THE UNIT

- The **SELECT FUSION** unit produces two grades of purified water. Grade II water is stored in the integral 20 litre tank.
- To dispense water from the tank simply press the lever on the white bib tap located at the front of the unit as indicated in **Fig-23** below.



Fig-23 Bib Tap Operation

- To dispense the higher quality Grade I+ water simply press button **5** on the keypad. **See Fig- 22**
- The automatic dispense valve will only operate if there is sufficient water in the integral to tank to run the internal recirculation pump.
- This dispense point has three modes of operation, **Latched, Hold** and **volume**.
- Refer to **Section 5.4** for details on all three dispensing modes and how to set them up.

4.7 STANDBY OPERATIONS

4.7.1 Manual Standby Operation

When dispensing of purified water is not required. The unit can be manually put into a standby condition (recirculation of purified water stops), and put back into processing via the menu options within the “Processing” screen. When in a standby condition the recirculation circuit can be quick started by pressing the dispense button. During standby, there will no recirculation of the purified water for 50 minutes. After 50 minutes the purified water will be recirculated for 10 minutes. This cycle will continue while the unit remains in a standby condition i.e. no demand for dispensed purified water.

4.7.2 Intelligent Standby Operation

The unit will automatically go into standby condition (recirculation of purified water stops) when no purified water has been dispensed for 30 minutes. The unit can be taken out of standby by simply pressing the dispense button, this action will start the recirculation circuit providing purified water for dispensing. During standby, there will be no recirculation of the purified water for 50 minutes. After 50 minutes the purified water will be recirculated for 10 minutes. This cycle will continue while the unit remains in a standby condition i.e. no demand for dispensed purified water.

SECTION

5

OPERATING PROCEDURES- DETAILED REFERENCE

This section describes the full procedures for all aspects of the ***SELECT FUSION***. It provides more details to that contained in the Quick Reference Section. This Section is divided into the following subsections:

SECTION CONTENTS


- 5.1 Control Philosophy
- 5.2 Menu structure
- 5.3 Clean Routine
- 5.4 Dispense Options
- 5.5 User Settings Menu
- 5.6 Processing Screen Information
- 5.7 Safety Features
- 5.8 Error and Alarm Messages


5.1 CONTROL PHILOSOPHY

The operation of the **SELECT FUSION** is described as follows.

- 5.1.1 The feedwater enters the unit via an external 5-micron carbon impregnated prefilter. The prefilter removes up to 95 % of all particles >5micron and dechlorinates to <0.1 ppm free Chlorine.
- 5.1.2 The filtered and dechlorinated water now enters the unit via a 150micron stainless steel mesh type screen filter to protect the following inlet solenoid from any large particulates.
- 5.1.3 From the inlet valve the water pressure is boosted by an internal boost pump to 60 psi, the optimum operating pressure for the Reverse Osmosis (RO) membrane.
- 5.1.4 Water now enters the (RO) module. The water that passes through the membrane is termed, Permeate. The water containing the rejected salts, up to 98% of the original incoming salts is directed to drain at a rate of up to 60 l/hr.
- 5.1.5 After leaving the (RO) membrane the Permeate passes through the Puripac PP8 cartridge containing mixed bed ion-exchange resin and activated carbon.
- 5.1.6 The quality of the treated water is monitored and displayed as “**Pre Qual**” and has a quality of between 1-10 MΩ-cm. The water now fills the integral tank, the level being controlled by the two upper level sensors, which are linked to the operation of the inlet solenoid.
- 5.1.7 When the level of water in the treated water reaches the lower of the 2 upper level sensors the internal recirculation pump will start. Water is drawn by the pump and passes it firstly through the NCP media pack containing semiconductor, mixed bed nuclear grade resin and activated carbon. The quality, flowrate, TOC and Temperature of the final purified water being monitored and displayed on the user interface.
- 5.1.8 The recirculation pump is protected from running dry by a low level sensor located in the bottom of the tank.
- 5.1.9 From the NCP media pack the highly polished water at a quality of between 15-18.2 MΩ-cm, displayed as “**Out Qual**,” flows through the automatic 3-way dispense valve.
- 5.1.10 Water can be dispensed via this valve, 3 options of dispense are available, Hold, Latched and Volume.
- 5.1.11 To provide sterile water at point of use the valve is fitted with a 0.2-micron absolute rated capsule filter with protective filling bell.
- 5.1.12 Water from the integral tank can be dispensed via the white bib tap.
- 5.1.13 To maintain quality the water is recirculated around this loop at a rate of up to 1 l/min.
- 5.1.14 The Graphical screen displays all alarms, quality, and status readings along with a mimic of the process.


Refer to the process flow schematics in the Appendix.





5.3 CLEAN ROUTINE


- To access the “**Clean Routine**”, from the “Power On” screen, select “**Menu**”, to display the screen below. Then step down and select “**Clean Routine**”.



- On entering the Clean routine you are asked to confirm, “**Are you Sure**”, as once you have selected “Yes” the unit will be locked in to the clean cycle for approx 2hours.
- On selecting “**Yes**” you will be asked to select between a **Loop Clean or Full Clean**”. A Loop clean will only disinfect the integral tank and polishing loop. A full clean will disinfect in addition to the loop clean, the (RO) modules.
- For full instructions on how to carryout the disinfection cycles please refer to **Section 6 Cleaning and Disinfection**.

5.4 DISPENSE OPTIONS.

- To access “**Dispense Options**”, from the “Power On” screen select “**Menu**”, to display the screen below. Then step down and select “**Dispense Options**”.



5.4.1 On selecting “**Dispense Options**” you will be presented with the following choices: -

Dispense Method
Dispense Volume
Dispense Print

- On selecting “**Dispense method**” you will be presented with 3 choices of dispense: -

Manual Hold
Manual Latched
Volume

- If “**Manual Hold**” is selected on pressing the Dispense button water will be dispensed, on releasing the dispense button water flow will stop.
- If “**Manual Latched**” is selected on pressing the dispense button water will be dispensed, on releasing the button water will continue to flow. To stop dispense the button has to be pressed for a second time.
- If “**Volume**” is selected you will then need to set the volume to be dispensed by pressing EXIT and then selecting “**Dispense Volume**”. By using the ▲ ▼ arrow buttons program in the volume required.

From 0 – 1.00 litres the volume increases by 10mls
From 1.1 litres to 10.00 litres the volume increases by 100 mls
From 10.01 litres the volume increase by 1.0 litres up to a maximum dispense volume of 20 litres.

5.4.2 On selecting “**Dispense Print**” you can select this function to be Enabled or Disabled. If Enabled every time water is dispensed from the unit and if the unit is connected to a printer, a “Validation Report” will be printed out. (See Contents of a typical Validation Report print out below):

VALIDATION REPORT	
Time	:11:42:30
Date	:21/01/05
Unit ID Number	:1
Unit Serial Number	:26745
Outlet Quality	:18.2M
Uncompensated Quality	:29.8M
TOC	:5 ppb
Temperature	:21.0 degC
NCP pack Replacement Due	
PP8 pack Replacement Due	

5.5 USER SETTINGS MENU

- To select the User settings menu press, **“MENU”** from the “Power On” screen to display the screen below. Then step down and select **“Settings Menu”**.



- The table below shows all parameters within the User Settings Menu, which can be accessed via the “Power On” screen. To change any of the settings or values follow the on screen prompts and by using the ▲ ▼ buttons on the keypad.

Parameter	Range/Options	Factory Default Value
Select Language	English	English
Set Time	24 hr clock	Current Date
Set Date	DD/MM/YY	Current GMT
Audio Alarm Enable	Enable/Disable	Disabled
Set Quality Alarm	1.0MΩ – 18.2 MΩ	1.0MΩ
Set Temp Alarm	35 °C - 20 °C	35 °C
NCP Pack Date	6 mths Range	00/00/00
PP8 Pack Date	6 mths Range	00/00/00
UV Hours Date	8000 hrs	0 hrs
Data Log Enable	Enable / Disable	Disable
Dump Log Data	N/a	N/a
Set Unit ID	01-99	01
Auto Restart	Enable / Disable	Disabled

5.5.1 Select Language

The Select Fusion has French/English language available.

5.5.2 Set Time

The Select unit incorporates a real time, 24 hr clock. The clock will have to be adjusted for BST. The clock will be factory pre set to the correct GMT. For overseas users of the equipment the time settings may have to be altered as necessary to match individual time zones.

5.5.3 Set Date

The correct date will be pre-set into the unit, check before use and adjust if necessary. Leap year calendar changes will have to be adjusted for.

5.5.4 Audio Alarm Enable

The Unit incorporates a buzzer alarm, which will sound should any of the alarms be activated, see **Section 5.7** for list of alarms. The alarm can be permanently disabled or can be muted each time by pressing button **2** and selecting “**Mute Alarm**” from the menu list.

5.5.5 Set Quality Alarm

This alarm refers to the quality at which the NCP media pac can be pre-set to indicate replacement. The set point is factory pre set to alarm at 1M Ω .

NB: The capacity quoted for the NCP is to 1M Ω . If the set point is raised then there will be a reduction in cartridge capacity.

5.5.6 Set Temp Alarm

If the treated water temperature exceeds the pre-set temperature valve, a “**Temperature High**” alarm message will be displayed. The default temperature setting for the unit is set at 35°C and can be accessed via the “**User Settings**” menu **Section 5.5**.

5.5.7 NCP Pack Date

Following installation of a new NCP Pack the date of replacement must be reset. A date of 6 months in the future is set on pressing the **RESET** button.

5.5.8 PP8 Pack Date

Following installation of a new PP8 Puripac the date of replacement must be reset. A date of 6 months in the future is set on pressing the **RESET** button.

5.5.9 UV Hours Date

Following replacement of a new UV lamp the future date of the replacement lamp must be reset. On entering this screen (see below) on pressing "RESET" the totalled hours run will reset to "0" hrs and the date will change to 12 mths in the future. The Hours run figure will rise according to actual hours run and number of lamp on/off's. This also applies to the forecasted date.



5.5.10 Data Log Enable/ Dump Log Data

The internal microprocessor constantly records system data and stores it. This data can be downloaded to a compatible PC. For further information on this feature contact Purite Limited.

5.5.11 Set Unit ID

If there is more than one unit on site it may be useful to give them all different I.D. numbers. This feature allows you to identify up to 99 units.

5.5.12 Auto Restart

Can be set to "Enabled" or "Disabled". If Enabled is selected the unit will automatically restart following an interruption in the power supply and return to the start of "Processing" or Process Standby (depending on what mode the unit was in prior to the loss of power).

5.6 “PROCESSING” SCREEN INFORMATION



By using the ▲▼ arrow buttons you can select any of the **SELECT FUSION’s**, 6 performance indicators. These will be displayed at the top of the screen. The indicators available are: -

Indicators	Units	
“Flowrate”	L/min	Display of recirculation/dispense flowrate.
“Pre Qual”	MΩ	Water quality feeding Integral Tank
“Out Qual”	MΩ	Quality available at point of dispense
“U Qual”	MΩ	As “Out Qual” but with no compensation for temperature (Pharmaceutical Users)
“Temp”	°C	Dispense water temperature
“TOC”*	ppb	Total Organic Carbon indicator for dispense water. Range: <5ppb <10 ppb <30 ppb

* The message “**TOC Wait**” will be displayed if “TOC” is selected within 2 minutes of starting the unit. This is to allow the TOC monitor to complete its start up routine.

5.7 SAFETY FEATURES

The **SELECT FUSION** has a number of safety features pre programmed into its microprocessor control system designed as self protection should a fault occur.

See Table below for details of all safety features: -

Displayed Alarm Message	Safety feature
“Low Pressure”	If the unit is fitted with a boost pump and the boosted pressure falls to < 5 psi for more than 3 seconds the unit will shut down. If the “Audio Alarm” has been enabled in “User Settings” it will also sound. This feature will protect the boost pump from dry running.
“High Pressure”	If the unit is fitted with a boost pump and the boosted pressure exceeds 90 psi the unit will shut down. If the “Audio Alarm” has been enabled in “User Settings” it will also sound.
“Fit PP8 (Right) Pack”	If the left hand side PP8 Puripac has not been fully engaged or works loose during normal operation the unit will shut down, switching off the feedwater supply and recirculation pump, so preventing flooding. If the “Audio Alarm” has been enabled in “User Settings” it will also sound.
“Fit NCP (Left) Pack”	If the right hand side NCP media pack has not been fully engaged or works loose during normal operation the unit will shut down, switching off the feedwater supply and recirculation pump, so preventing flooding. If the “Audio Alarm” has been enabled in “User Settings” it will also sound.
“Fit Both Packs”	If both pack have not been fully engaged or work loose during normal operation the unit will shut down, switching off the feedwater supply and recirculation pump, so preventing flooding. If the “Audio Alarm” has been enabled in “User Settings” it will also sound
“Tank Low Level”	The internal recirculation pump is protected from running dry by a low level sensor located at the bottom of the 20 litre tank. At low level the screen mimic for the tank will show empty, the corresponding “Tank Low Level Alarm” will be raised and the recirculation pump will be switched off. The pump will restart when the level reaches the intermediate level switch.

For details of how to resolve and clear the above, refer to **Section 8 Trouble shooting.**


5.8 ERROR MESSAGES AND ALARMS

In conjunction with the safety features the unit will also display many “Error & Alarm” messages giving warnings of failed sensors, prompts for changing consumable items, which have exceeded their pre set exchange dates and indication that the water quality is outside desired limits. Details of these can be found in the table below.

Displayed Alarm/Error Message	Secondary Message	Reason for Alarm
“Temp Sensor Error”	Value Displayed “==.=”	Temp Sensor lead being disconnected or faulty or water temp is <1 or >100 Deg C
“Line Cell Error”	Value Displayed “==.=”	Outlet line cell sensor lead disconnected or faulty or final water quality > 60 MΩ
“Keypad Error”	-	Keypad lead disconnected or faulty
“Time Date Not Set”	-	On setting up unit for the first time these must be set.
“Change UV Lamp”	UV mimic will flash	Faulty or disconnected UV sensor lead or lamp has failed or life expectancy of lamp has been exceeded.
“Temperature High”	-	Recirculated water temperature has exceeded pre set value.
“Poor Quality”	-	Outlet Quality has exceeded pre set quality limits.
“NCP Pack Rep Due”		Forecasted 6 mths pack life has been exceeded
“PP8 Pack Rep Due”	PP8 mimic will flash	Forecasted 6 mths pack life has been exceeded
“Clock Values Invalid Check PCB Battery”	-	Internal main PCB battery voltage fallen below min value.

For further details of how to diagnose and clear alarms refer to, **Section 5 User Settings, Section 8 Troubleshooting and Section 7 Maintenance.**

SECTION



DISINFECTION

This section provides details on how to carry out a Disinfection on the ***SELECT FUSION*** Unit. The section is divided into the following sub-sections.

SECTION CONTENTS

- 6.1 Procedure

6.0 **DISINFECTION**

- It is recommended that your **SELECT FUSION** unit is disinfected typically twice per year to ensure it maintains consistent performance.
- The unit can perform two different disinfection routines.

- (i) **Loop Clean**
- (ii) **Full Clean**

- If a “Loop Clean” is selected then only the polishing loop will be disinfected including the integral tank.
- If a “Full Clean” is selected then the unit will carryout a “Loop Clean” plus disinfection of the Reverse Osmosis membranes.
- The disinfecting agent is contained within a “Disinfection pac Adaptor” (Refer to **Section 9.2** for ordering details). The chlorinated agent is in tablet form.
- For a “Loop Clean” only one adaptor will be required, for a “Full Clean” two adaptor pacs will be required.
- To perform a disinfection refer to **Section 6.1**.
- The clean cycles generally last for about 60-90 minutes.
- Please read the material safety data sheet, which can be found in the Appendix, for the Disinfection Pac adaptor before handling.
- If any alarm messages are displayed during the clean then refer to **Sections 5.7 & 5.8** for details, or **Section 8** for corrective actions.

6.1 **PROCEDURE**

- 6.1.1 Switch the unit on and select “**Menu**”.
- 6.1.2 From the displayed Menu listing select “**Clean Routine**”.
- 6.1.3 You will be asked to confirm, “**Are you sure?**”
- 6.1.4 Select “**Yes**”. Once selected the unit will be locked into the routine until it is completed.
- 6.1.5 You will then be asked to select either “**Loop Clean**” or “**Full Clean**”; select your desired choice.
- 6.1.6 Now simply follow the on screen prompts which will indicate into which position the disinfection pac adaptor/s must be placed.
- 6.1.7 Once the pacs have been fitted the process is fully automatic.
- 6.1.8 At the end of the routine the pacs can be removed and disposed off.
- 6.1.9 If a new Puripac PP8 cartridge and NCP media pac is to be fitted follow instructions detailed in **Sections 7.2 & 7.3**
- 6.1.10 If a new 0.2um Point of Use filter is to be fitted, then follow the instructions detailed in Section 7.4.
- 6.1.11 For guidance of when to disinfect the unit and replacement of all consumables, refer to the Preventative Maintenance Guide in Section 7.9

SECTION



MAINTENANCE

This section describes the recommended procedures for replacing consumables on the **SELECT FUSION** unit. This section is divided into the following sub-sections.

SECTION CONTENTS

- 7.1 Replacing the Pre-filter cartridge
- 7.2 Replacing the Puripac PP-8 cartridge.
- 7.3 Replacing the NCP cartridge pack
- 7.4 Replacing the 0.2 um point of use filter.
- 7.5 Replacing the Ultra violet Lamp.
- 7.6 Replacing Fuse/s
- 7.7 Replacing PCB Lithium battery.
- 7.8 Cleaning of External Surfaces.
- 7.9 Preventative Maintenance Guide

7.0 MAINTENANCE

This section provides details of how to carry out routine maintenance tasks on your **SELECT FUSION** unit. The tasks only relate to those involving the changing of consumable items.

For tasks that fall outside the changing of consumables, always refer to the supplier of the equipment for advice and help.

Section 11 “Contacting Us”, provides contact details should you need to contact Purite Limited.

Safety Information For Maintenance Task

The operator must take care to ensure that only authorised personnel who have been sufficiently informed for the task in hand, by thoroughly studying the Operating Instructions, perform all maintenance tasks.

The unit must be shut down and protected from being placed in operation again unintentionally before any maintenance tasks have been completed. It is essential to observe the procedure described in these Operating Instructions for shutting the unit down.

Before beginning any task on the electrical equipment in the unit, a check must confirm that the power has been disconnected and isolated from the unit. In addition, the unit must be secured to prevent it from being turned on again unintentionally.

If at anytime you are unsure please contact Purite or your supplier to seek assistance using the contact details in Section 11.

7.1 REPLACING THE PRE-FILTER CARTRIDGE.

Fig-26



- Switch of the Unit.
- Isolate mains water supply to Filter Housing.
- Relieve pressure in housing by pressing Red Button on top of Bowl as indicated. **(See Fig-26)**
- Unscrew bowl clockwise. **(NB. Bowl will be full of water)**
- Empty bowl and remove/dispose filter element.
- Place the new filter (R011146) into the bowl, and ensure it locates over the bottom spigot.
- Refit the bowl and tighten to hand tight.
- If the bowl leaks check that the black O-ring seal is in place.

7.2 REPLACING THE PURIPAC PP-8 CARTRIDGE



Fig-27



Fig-28

- Switch off the unit and isolate the mains water supply.
- Open the front door cover; the **Puripac PP-8 cartridge** is located to the left.
- Press down the **Blue** lever on top of the holder and gently pull away the exhausted cartridge. **(See Fig-27)**. There may be some drops of water lost from the cartridge.
- The cartridge can now be disposed of, refer to **Section 9.4** for advice on disposal.
- To refit a new pack simply remove it from its packaging. Be sure to remove the two black plugs used to seal the two ports, **(see Fig-11)**.
- Then slide the pack back into the holder. The locking lever will “Click” indicating the pack is secure.
- Switch on the electrical supply and turn on the water. Select “**Start**” from the menu screen. The unit should now ask you to confirm that a new pack has been fitted followed by a request to reset the next exchange date, just follow the simple on screen commands.
- **NB** If the new pack has not located correctly, the sensors located in the pack holder will set off an alarm and raise a corresponding alarm message on the screen, i.e., “**Fit PP8 (Left) Pack**”. Refer to **Section 8 Troubleshooting** should this occur.

7.3 REPLACING THE NCP CARTRIDGE PACK

- Simply follow the same procedures as laid down for the **Puripac PP_8 cartridge** in 7.2 above.
- The NCP cartridge pack is located to the right of the Puripac cartridge and can be released by depressing and holding down the, **White**, lever.
- **NB** If the new pack has not located correctly; the sensors located in the pack holder will set off an alarm and raise a corresponding alarm message on the screen, i.e., “**Fit NCP (Right) Pack**”. Push the cartridge fully into the retainer. Refer to **Section 8 Troubleshooting** should this not cure the alarm.

7.4 REPLACING 0.2 µm POINT OF USE FILTER.

Fig-29



- Fitting of a new filter can be carried out without turning the unit off.
 - Simply hold the body of the capsule and un-screw clockwise.
 - Once unscrewed dispose of the filter and remove the new one from its packaging. Refer to Section 9.4 for advice on disposal.
 - Wind about 4 turns of PTFE tape around the thread on the filter.
-
- Screw the new filter anti-clockwise back into the dispense valve. (NB: Be careful not to cross thread the filter)
 - The bleed nozzle should point out to the right hand side of the unit.
 - To bleed the filter of air, open the bleed screw a few turns, remove the lower bell-housing cap and dispense approx 500 mls of water, tapping the filter gently.
 - When all of the air has gone close the bleed screw.

7.5 REPLACING THE ULTRA VIOLET LAMP

NB: The replacing of the Ultra violet (UV) lamp should only be carried out by trained personnel. Please ensure the instructions have been read and fully understood before attempting to replace the lamp. The lamp should only require replacing once per year (8000 hrs). The unit will alarm, automatically, warning that the lamp should be replaced.

Fig-30



- Switch off the unit and isolate from electrical supply.
- Open the front door cover to reveal the side cover fixing screws. **(see Fig-30)**



Fig-31

- Unscrew the side cover using a Phillips or cross head screwdriver. **(See Fig-31)**
- Pull the cover backwards to release it from its rear hangers.
- Store to one side carefully.



Fig-32

- With the side cover removed the UV assembly chamber can be accessed. **(See Fig-32)**



Fig-33

- Unlock the black bayonet cap on top of the chamber.
- Gently pull out the lamp about 2", using the lead. **(See Fig-33)**
- By holding both the lamp and its lead disconnect them.
- Now fully withdraw the lamp.
- (NB: as the lamp contains Mercury (Hg) dispose of it accordingly).

- Remove the new lamp from its packaging. **(DO NOT HOLD THE LAMP BY THE GLASS TUBE, IF YOU ACCIDENTALLY TOUCH THE GLASS WIPE OFF ANY MARKS WITH ALCOHOL, ISOPROPYL ALCOHOL OR METHYLATED SPIRITS)**
- Lower the lamp back into the chamber, reconnecting the lead and then refit the bayonet cap.
- Replace the side cover and switch the power back onto the unit.
- To reset the UV lamp's future replacement date, refer to **Section 5.5.9 "UV Hours Date"**.

7.6 REPLACING FUSE/S

- Isolate the unit from the electrical supply and unplug the mains lead from the rear of the unit.
- With the aid of a small flat bladed screwdriver pries out the fuse carrier from the IEC module. The two fuses will now be visible.
- Remove the fuse/s and insert replacement fuse/s of the correct rating. Refer to **Section 9.2**, into carrier and push back into IEC module.
- Plug the mains lead back into the unit, switch the power back on and restart the unit.
- If the fuse/s blow again, then contact Purite or your supplier. Refer to; **Section 11 Contacting Us**, for details.




Fig-34 Fuse holder

7.7 REPLACING PCB LITHIUM BATTERY



It is recommended that the replacement of the battery is only carried by Purite approved Engineers or by Purite personnel. The Battery has a 5-year life, but it is recommended that it should be changed every 3 years.

- Switch off the unit and turn off the water supply, isolate from the electrical supply and unplug the mains lead from the rear of the unit.
- Open the front door to gain access to the two Top cover retaining button head socket screws.
- Unscrew using a 3mmAF Allen key.
- Lift off the lid.
- The Lithium battery is located at the top right hand edge of the main circuit board. **See Fig-35.**
- To remove the battery gently pries up the spring retainer and slide out he battery. **See Fig-36.** Dispose of the battery according to local regulations.
- Refit the new battery ensuring the +ve side is facing upwards



Incorrect fitting of the battery could cause irreversible damage to the main PCB. Ensure +ve side of battery is facing upwards. Always use recommended battery Part No. R083085.

- Replace top cover, switch the electrical power back on and turn on the water supply.

Fig-35





Fig-36



7.8 CLEANING OF EXTERNAL SURFACES

- Use a clean damp cloth to wipe the exterior surface of the unit.
- DO NOT allow excessive liquid on to the key pad.
- Do not use any Ketone-based solvents on the covers or front display.
- Industrial methylated spirits or iso-propyl alcohol/water based cleaners/surface disinfectants can be used to remove more persistent marks/stains and to disinfect.

SECTION



TROUBLESHOOTING

This section provides a guide to first line checks that can be undertaken should you encounter a problem with your **SELECT FUSION**. It is divided into the following sub-sections:

SECTION CONTENTS

- 8.1 General Fault conditions
- 8.2 Alarm messages

8.0 TROUBLESHOOTING

8.1 GENERAL FAULT CONDITIONS

Symptom/s	Possible causes	Actions
Unit does not power up	<ol style="list-style-type: none"> 1. Fuse blown in mains lead 2. Fuse blown in IEC module 3. Electrical fault in Fusion 	<ol style="list-style-type: none"> 1. Check 5 amp fuse in lead 2. Check fuses in IEC module (refer to Sec. 7.6) 3. Contact supplier
Make up rate into integral tank declining	<ol style="list-style-type: none"> 1. Feedwater pressure low. 2. Feedwater temp change. 3. Consumable dates exceeded 4. (RO) membrane fouled 	<ol style="list-style-type: none"> 1 & 2. Check feed conditions. (Refer to Sec.1.3) 3. Check all replacement dates and replace where necessary. 4. Contact supplier with details.
"Pre Qual" quality poor	<ol style="list-style-type: none"> 1. Puripac PP8 exhausted 2. Change in feedwater supply. 3. Linecell calibration drifted. 	<ol style="list-style-type: none"> 1. Check PP8 2. Check feedwater conductivity. 3. Contact supplier to recalibrate linecell.
"Out Qual" poor	<ol style="list-style-type: none"> 1. NCP media pac exhausted. 2. Linecell calibration drifted. 	<ol style="list-style-type: none"> 1. Check NCP pac 2. Contact supplier to recalibrate linecell.

8.2 ALARM MESSAGES

CAUSE	DISPLAY	ACTION
Temp probe lead disconnected or faulty.	TEMP. SENSOR ERROR	Check Sensor lead. Replace if error persists.
Line cell lead disconnected or faulty.	LINE CELL ERROR	Check sensor connections. Replace if error persists.
Boost Pump variants only. Pump faulty or loss of feed water. Low pressure switch fault.	LOW PRESSURE	<ul style="list-style-type: none"> - Check Boost Pump operation. - Check Feedwater supply pressure. - Replace pre-filter. - Contact supplier
Boost Pump variants only. <ul style="list-style-type: none"> - Increase in Feed Water Pressure. - PP-8 Cartridge expired 	HIGH PRESSURE	<ul style="list-style-type: none"> - Check Feedwater Pressure must be <30psi. - Change PP8 Pack - Contact supplier if problem persists.

PP8 pack not fitted or not fully engaged.	FIT PP8 (LEFT) PACK	<ul style="list-style-type: none"> - Fit pack. (Refer to Sec 7.2) - Check pack is fully engaged and retainer pin locked.
Puripac not fitted or not fully engaged.	FIT NCP (RIGHT) PACK	<ul style="list-style-type: none"> - Fit NCP Pack. (Refer to Sec. 7.3) - Check Pack is fully engaged and retainer pin located.
Puripac PP8 or NCP pack both not fitted or both not fully engaged.	FIT BOTH PACKS	<ul style="list-style-type: none"> - Re-fit Pack/s. - Check Pack/s are fully engaged and retainer pin located.
On Start Up time and date may not be set. Main PCB battery needs replacing.	TIME/DATE NOT SET	<ul style="list-style-type: none"> - Select Set Time and Set Date function from settings menu and reset. - Replace Lithium PCB battering (Ref to Sec. 7.7)

Pre set lamp life expired	CHANGE UV LAMP	Replace UV lamp (refer to Sec. 7.5)
6 Month Operating Life exceeded.	NCP PACK REP DUE	<ul style="list-style-type: none"> - Change NCP Pack. (Refer to Sec 7.3) - Reset Pack Rep date. (Refer to Sec 5)
<ul style="list-style-type: none"> - UV lamp not alight. - Sensor fault. 	UV FAILURE	<ul style="list-style-type: none"> - Check lamp is alight. - Check UV lamp hrs run. - Contact Supplier.
Water in Process Tank below low level. Usage exceeding make up rate.	TANK LOW LEVEL	<ul style="list-style-type: none"> - Check feedwater supply and (RO) performance. - Check tank level (visual). - Close all Take Offs.

Link between Keypad and Main PCB Faulty	KEYPAD ERROR	<ul style="list-style-type: none"> - Contact Purite CSD.
Temperature of Purified water above set limit.	TEMPERATURE HIGH	<ul style="list-style-type: none"> - Check Temp Setting. (Refer to Sec 5) - Drain off water from system. - Check feedwater temp. - Change to Intermittent running (Refer to Sec. 5)
Outlet quality above set limit.	POOR QUALITY	<ul style="list-style-type: none"> - Check quality setting in setting menu.(Refer to Sec. 5) - Change NCP cartridge.
The Puripac PP8 pre-set replacement date has expired.	PP8 PACK REP DUE	<ul style="list-style-type: none"> - Replace PP8 pack (Refer to Sec. 7.2) - Reset expiry date (Refer to Sec. 5)

NOTE: If remedial actions do not resolve the problem, turn off the water supply, isolate the unit from the electrical supply and contact either your authorised supplier or call Purite Service Department on 01844 211555. DO NOT ATTEMPT ANY REPAIRS WITHOUT FIRST CONTACTING PURITE LIMITED OR AUTHORISED SUPPLIER.

SECTION

9

CONSUMABLES & SPARES

This section contains details of all recommended spare parts and consumables used on the **SELECT FUSION**. It is divided into the following sub-sections:

SECTION CONTENTS

- 9.1 Consumables
- 9.2 Recommended spares
- 9.3 Serviceable spares
- 9.4 Disposal of consumables
- 9.5 Accessories.
- 9.6 Consumable Shelf Life

9.0 CONSUMABLES AND SPARES

9.1 CONSUMABLES

Equipment Description	Re-order Code
PP-8 Select Puripac	L998521
NCP Cartridge Pack	L991540
0.2 micron Point of use Filter	R090015
External 10" Carbon Filter	R011146

9.2 RECOMMENDED SPARES

5 amp Fuses	R081414
Disinfection pac adaptor	L998549
185/254nm UV lamp	R081634
Lithium PCB Battery	R083085

9.3 SERVICEABLE SPARES

Fusion 40 Membrane Pack	M996001
Fusion 80 Membrane Pack	M996002
Fusion 160 Membrane Pack	M996003

9.4 DISPOSAL OF CONSUMABLES

All consumables should be disposed of according to local/national regulations. Refer to Appendix A – WEEE Declaration for disposal of electrical and electronic equipment supplied by Purite. Refer to Material Safety Data Sheets for further details.

9.5 ACCESSORIES

Wall mount Kit	L998400
Raw water Pressure Regulator	L991110
Mains Lead	R081404
Installation Kit	PM00901

9.6 CONSUMABLE SHELF LIFE*

Puripac PP-8	24 months
NCP media pack	24 months
0.2µm POU Filter	36 months
CI Carbon Filter	Indefinite

***Shelf life only guaranteed if consumables are not removed from their packaging and are stored in a clean dry environment at room temperature away from sources of heat or direct sunlight.**

SECTION

10

CE DECLARATION

This section contains a copy of the CE Declaration Certificate for the ***SELECT FUSION***.

SECTION CONTENTS

10.0 CE DECLARATION



Purite Limited,
Bandet Way, Thame
Oxon, OX9 3SJ, UK.
Tel: +44 (0)1844 217141
Fax: +44 (0)1844 218098
Web site: www.purite.com
Email: mail@purite.com
Company Registration
No. 1464412

CE DECLARATION OF CONFORMITY

We, the Manufacturer:

Purite Limited, Bandet Way, Thame, Oxon, OX9 3SJ, United Kingdom
Declare under sole responsibility that the SELECT/PRESTIGE product range covered by
this declaration, conforms with the essential protection requirements of the following
Directives:

Applied Council Directives:

2004/108/EC - Electromagnetic Compatibility Directive (EMC).

2006/95/EC - Low Voltage Directive (LVD).

(Standard Applied: EN 61010-1:2001-02, Safety requirements for electrical equipment for
measurement, control and laboratory use Part 1: General requirements).

The **SELECT/PRESTIGE** range covered by this declaration is as follows:

HP 40, 80, 160, 320, 640 with and without BOOST PUMP
BIO 40, 80, 160, 320, 640 with and without BOOST PUMP
PRESTIGE DESCALE PD50/50H/100/100H/150/250
FUSION 40,80,160 with and without BOOST PUMP
BECKMAN HP320
NEPTUNE ANALYTICAL, LIFE SCIENCE, ULTIMATE

Purite Limited's liability under this declaration is limited to that set forth in the current Purite
Limited Terms and Conditions of Sale.

Signature: 

Full Name: R S Keep
Position: Director
Date of issue: 05.08.07

Appendix F



SECTION

11

CONTACTING US

This section contains details of customer contact points at Purite Limited.

SECTION CONTENTS

11.0 Contacting Us


APPENDIX

APPENDIX A

- (i) Pushfit Water Connectors
- (ii) WEEE Declaration
- (iii) Material Safety Data Sheets
- (iv) Process Flow diagrams

How to make a connection

To make a connection, the tube is simply pushed in by hand; the unique patented John Guest collet locking system then holds the tube firmly in place without deforming it or restricting flow.





Users in the United Kingdom who wish to discard electrical and electronic equipment that was supplied by Purite should contact B2B Compliance on 01691 676 124.

Users in other European Union countries should contact their Producer, who will be the organisation in their country that supplied the product.

Information on disposal in other countries outside the EU

This symbol is only valid in the European Union. If you wish to discard this product, please contact your local authorities or dealer and ask for the correct method of disposal.

MATERIAL SAFETY DATA SHEET

1. IDENTIFICATION OF THE PREPARATION AND THE COMPANY

PRODUCT NAME: SELECT PURIPAC 4/8/16 CARTRIDGE

SUPPLIER: Purite Limited
Bandet Way
Thame Industrial Estate
Thame
Oxon
OX9 3SJ
Tel no: 01844 217141

2. COMPOSITION:

Ingredient Name

(Hazard classification based on data for mixture).

Sulfonated copolymer and styrene and divinylbenzene in the hydrogen form (CAS# 069011-20-7)

Trimethylamine functionalised copolymer of styrene and divinylbenzene in hydroxide form (CAS# 069011-18-3).

Steam Activated Carbon - Chemical Formula C - CAS No 7440-44-0 (or 64365-11-3).

UN Number 1362 - EINECS Number 231-153-3 (or 264-846-4).

3. HAZARDS IDENTIFICATION:

Irritating to eyes.

4. FIRST AID MEASURE:

Never give fluids or induce vomiting if patient is unconscious or is having convulsions.

Inhalation

No adverse effects anticipated by this route of exposure.

Skin Contact

Wash off in flowing water or shower.

Eye Contact

Irrigate with flowing water immediately and continuously for 15 minutes. Consult medical personnel.

Ingestion

No adverse effects anticipated by this route of exposure incidental to proper industrial handling.

5. FIRE FIGHTING MEASURES:

Extinguishing Media

Carbon dioxide. Dry chemical. Foam.

Hazardous combustion products

Nitrogen, sulphur and carbon oxidation products.

Protection of Fire-fighters

Wear positive-pressure self-contained breathing apparatus and protective fire fighting clothing (includes fire fighting helmet, coat, pants, boots and gloves).

Specific Fire or Explosive Hazards

Product is not combustible until moisture is removed and then resin starts to burn in flame at 230 deg. C.

6. ACCIDENTAL RELEASE MEASURES:

Personal precautions

Spills may cause very slippery surfaces. Wear adequate personal protective equipment, see section 8. Exposure controls/Personal Protection.

Environmental Precautions

Prevent contamination of surface and ground water.

Methods of cleaning up

Contain and sweep up the material. Collect in suitable and properly labelled containers. Recover if possible, or dispose of according to applicable regulations, see Section 13, DISPOSAL CONSIDERATIONS.

7. HANDLING AND STORAGE:

Handling

Static electricity can accumulate on dry beads. May cause very slippery surfaces if spilt.

Storage

Avoid dehydration.
Store between 1 and 50 deg. C.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION:

Exposure Guidelines

None established.

Engineering Controls

Good general ventilation should be sufficient.

Personal Protective Equipment

Respiratory Protection

No respiratory protection should be needed.

Protective Clothing

Use gloves, impervious to this material, when prolonged or frequently repeated contact could occur. For brief contact, no precautions other than clean body-covering clothing should be needed.

Eye/Face protection

Use Chemical goggles.

9. PHYSICAL AND CHEMICAL PROPERTIES:

Appearance	Solid beads
Colour	White to dark amber
Boiling point/range	N/A
Freezing point/range	N/A
Water Solubility	N/A
Vapour pressure	N/A
Specific gravity	1.04 - 1.4.
Flash point	N/A
Flammability	N/A
Auto-ignition temp.	> 500 deg. C
Flammability-LFL	N/A
Flammability-UFL	N/A

10. STABILITY AND REACTIVITY:

Chemical Stability

Stable under normal handling and storage conditions, see Section 7, Handling and Storage.

Conditions to avoid

The severity of the reaction with oxidising materials can vary from slight degradation to an explosive reaction. Oxidising agents such as nitric acid, attack organic ion exchange resins under certain conditions. Before using strong oxidising agents, consult sources knowledgeable in handling such materials.

Materials to avoid

Strong oxidising agents.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Ingestion

Single dose oral toxicity is believed to be low. No hazards anticipated from swallowing small amounts incidental to normal handling operations.

Skin

Skin absorption is unlikely due to physical properties.

Inhalation

Vapours are unlikely due to physical properties.

Irritation

Skin

Prolonged exposure may cause skin irritation. May cause more severe response if skin is abraded (scratched or cut).

Eyes

May cause severe eye irritation. May cause moderate corneal injury. Effects are likely to heal.

Other information

No specific data available, however, repeated exposures are not anticipated to cause significant adverse effects.

12. ECOLOGICAL INFORMATION

Material is expected to degrade only slowly in the environment. May change pH of waters. Possible adverse effects on aquatic organisms are expected to result primarily from physical/mechanical effects rather than toxicity.

13. DISPOSAL CONSIDERATIONS

Any disposal practice must be in compliance with all local and national laws and regulations. Consider contaminants when disposing of used resins.

14. TRANSPORT INFORMATION

Product is not classified for any mode of transportation.

15. REGULATORY INFORMATION:

EEC Classification and user label information

Classification according to UK chemicals (Hazard information and Packaging) regulations 1994 or CHIP.

Hazard Symbol:

Xi - Irritant

Risk Phrases:

Irritating to eyes (R36)

Safety Phrases: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. (S26). Wear eye/face protection. (S39).

Chemical Name: Not required.


16. OTHER INFORMATION

No other information.

17. DISCLAIMER:

The data given above is based upon information received by us and is believed to be correct as general guidance for our customers, who should nonetheless seek confirmation regarding the application of such information as well as safety and/or suitability in relation to the products supplied either alone or in combination with other products. Given that the above information is for general guidance only, users (customers and otherwise) assume all responsibility and liability of whatsoever nature arising from the handling and use of the product and the application of such information (whether or not the product is used alone or in combination with another product) and users must satisfy themselves by specific enquiries and confirmation as to application of the above information and/or the safety and/or suitability of the product and by reference to their own tests and knowledge.

Signature:



Date:

17.02.05

MATERIAL SAFETY DATA SHEET

1. IDENTIFICATION OF SUBSTANCE / PREPARATION AND COMPANY

Product Name DISINFECTION PAC ADAPTOR
Synonym (s)
Company Identification PURITE LIMITED
BANDET WAY
THAME
OXON OX9 3SJ
Telephone: 01844 217141

2. COMPOSITION / INFORMATION ON INGREDIENTS

INGREDIENT NAME	CAS NO	EEC NO
Sodium Dichloro Isocyanurate dihydrate	51580-86-0	220-767-7

Classification Xn.NR22-31-36/37-S0/S3

3. HAZARDS/RISKS IDENTIFICATION

Residual Hazardous Substances: Contact with acids liberates chlorine gas.
Health Risk: Harmful if swallowed.
Environmental Risk: Very toxic to aquatic organisms and may cause long-term adverse effects in the aquatic environment.

4. FIRST AID MEASURES

General information

Skin contact Irritation Drench the skin with plenty of water. Remove contaminated clothing and wash before reuse. If large areas of the skin is damaged or if irritation persists seek medical attention.
Eye contact Irritation Irrigate thoroughly with water for at least 10 minutes. Obtain medical attention.
Inhalation Irritation of breathing passages. Shortness of breath, wheezing, choking and chest pains. Remove the individual from exposure, rest and keep warm. In severe cases, or if recovery is not rapid or complete seek medical attention.
Ingestion Irritation or burns to gastrointestinal tract. Nausea, vomiting and abdominal pain. Rinse mouth with water. DO NOT INDUCE VOMITING. If patient is conscious give water to drink. If patient feels unwell seek medical attention.
Immediate Treatment Antidote Symptomatic treatment.

5. FIRE FIGHTING MEASURES

Suitable extinguishing media Self contained breathing apparatus.
Special protective equipment for fire-fighters In the event of fire, wear self-contained breathing apparatus and full protective gear.
Unsuitable Extinguishers: Use extinguisher suitable to cause of fire. Keep fire exposed containers cool with water spray.
Hazardous Combustion Products: Toxic fumes, chlorine, nitrogen trichloride.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions Wear appropriate PPE – see section 8
Environmental precautions Prevent entry into drains and water courses.
Measures for cleaning Sweep up and place in suitable labelled dedicated containers. Hold for waste disposal.

7. HANDLING & STORAGE

Advice on safe handling	Avoid contact with skin and eyes. Do not breathe dusts. Keep away from organic materials, reducing agents and other oxidizers.
Storage conditions	Temperatures range – cool. Humidity range – Dry Keep away from suitable storage media. Store in original closed containers. See section 10

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure limit(s)	Normal ventilation is sufficient to control exposure.
Eye protection	Facemask and safety glasses are recommended as well as eyewash facility.
Hand protection	Half facemask plus goggles or full-face mask, both fitted with chlorine gas and dust/mist filters. PVC gloves. Always wash thoroughly after handling chemicals.
Respiratory protection	
Hygiene measures:	
Skin Protection:	Overalls, boots and impermeable apron or suite.

9. PHYSICAL & CHEMICAL PROPERTIES

Colour	White tablets.
Boiling point/range	240 – 250°C
Odour	Sharp, chlorine/bleach like.
Flammability	Not flammable in normal conditions of handling and use.
Solubility in water	0.96 25% @ 30°C.
Relative density/pH	Oxidising agent – 6.0 – 6.4 (1% solution)

NOTE: The physical data presented above are typical values and should not be construed as a specification.

10. STABILITY & REACTIVITY

Hazardous reactions	Decomposes at 240°C. Contact with small amounts of water may result in exothermic reactions and liberation of toxic fumes.
Materials to avoid	Other oxidizers, nitrogen containing compounds, organics.
Hazardous decomposition products	Nitrogen trichloride, Chlorine, Oxides of nitrogen. Oxides of carbon.
Conditions to avoid	Damp.

11. TOXICOLOGICAL INFORMATION

Effects	Harmful if swallowed, causing irritation or burns to the gastrointestinal tract resulting in nausea, vomiting, diarrhoea, abdominal pain, ulceration and bleeding. Inhalation of dust will cause irritation of the respiratory system, inhalation of decomposition products may cause severe respiratory damage. Contact with eyes and skin will cause irritation or chemical burns.
LD₅₀ 735 mg/kg oral-rat	> 2 g/kg skin-rabbit.

12. ECOLOGICAL INFORMATION

Environmental Effects:

Mobility	Soluble in water.
Aquatic Toxicity	LD ₅₀ , 96h, Rainbow trout = 0.22 ppm: LD ₅₀ , Daphnia magna = 0.55 mg/l.

13. DISPOSAL CONSIDERATIONS

Removal of product	Via an authorised waste disposal contractor to an approved waste disposal site, observing all local and national regulations.
Container	As substance.

14. TRANSPORT INFORMATION

Classification for ROAD and Rail Transport (ADR/RID)	Not regulated (Not dangerous for transport).
Classification for SEA transport (IMO-IMDG)	Not regulated (Not dangerous for transport).
Classification for AIR transport (IATA/ICAO)	Not regulated (Not dangerous for transport).

Contains: Environmentally hazardous substance, solid, N.O.S. (troclesene, sodium dehydrate).
UN Number 3077, Packing group 111, ADR Class/Item 9, 12^o(C), HI No 90
Emergency Action code 2X.

Primary Hazard Other dangerous substance.

15. REGULATORY INFORMATION

Classification & labelling have been performed according to EU directives 67/548/EEC and 99/45/EC including amendments

Hazard symbol and Indication of danger SODIUM DICHLORO ISOCYANURATE DIHYDRATE.
HARMFUL AND DANGEROUS FOR THE ENVIRONMENT.

EU EINECS (EINECS)

RISK PHASES 22-31 – 36/37 – 50/53.
THIS PRODUCT IS NOT HAZARDOUS ACCORDING TO EEC DIRECTIVES 67/548/EEC AND 99/45/EC. INDUSTRIAL GRADE PRODUCTS ARE INTENDED FOR USE IN ANALYTICAL AND MEDICAL APPLICATIONS WITHOUT EXTENSIVE PURIFICATION.

US Toxic Substances Control Act (TSCA) THIS PRODUCT SATISFIES ALL THE REQUIREMENTS OF THE EUROPEAN INVENTORY OF EXISTING CHEMICAL SUBSTANCES (EINECS).
ALL COMPONENTS OF THIS PRODUCT ARE IN COMPLIANCE WITH THE INVENTORY LISTING REQUIREMENTS OF THE US TOXIC SUBSTANCES CONTROL ACT (TSCA) CHEMICAL SUBSTANCE INVENTORY.
EEC NUMBER 220-676-7

16. OTHER INFORMATION


Further information


The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process unless specified in the text.

Signature: 

Date: 17.02.05


OPERATING MANUAL FOR SELECT FUSION






MATERIAL		PURITE LTD.	
FINISH		BANDEY WAY, THAME INDUSTRIAL ESTATE, THAME, OXON. OX9 3SL ENGLAND.	
ALL BURRS AND SHARP CORNERS TO BE REMOVED		TEL : 01844 217141 FAX : 01844 216088	
TOLERANCES		TITLE	
0-500mm	± 0.5	FLOW DIAGRAM	
500 +	± 1.0	SELECT FUSION 80	
HOLES TO BE	± 0.1	SCALE	X
UNLESS OTHERWISE STATED		DRAWN	J.COLE. L998301
DIMENSIONS IN MILLIMETRES		APPROVED	
		DATE	5.4.04
		DRG NO.	W3891
		ISSUE	A
		ISSUE	A

UNCONTROLLED COPY
UPDATES WILL NOT BE MADE TO THIS COPY




MATERIAL		—	
FINISH		—	
ALL BURRS AND SHARP CORNERS TO BE REMOVED			
TOLERANCES			
0-50mm	± 0.5		
50-100mm	± 1.0		
HOLES TO BE	± 0.1		
UNLESS OTHERWISE STATED			
DIMENSIONS IN MILLIMETRES			
PURITE LTD. BANGET WAY, THAME INDUSTRIAL ESTATE, THAME, OXON. OX9 3SL, ENGLAND. TEL : 01844 271141 FAX : 01844 218068		TITLE FLOW DIAGRAM SELECT FUSION 160	
ISSUE	DATE	MODIFICATION	APP'D BY
A	13.6.08	TIME ISSUE 1	
SCALE		DRAWN	J.COLE
		APPROVED	
	DATE	ORG NO.	ISSUE
	5.4.04	W3892	A

UNCONTROLLED COPY
UPDATES WILL NOT BE MADE TO THIS COPY




MATERIAL		—	
FINISH		—	
ALL BURRS AND SHARP CORNERS TO BE REMOVED			
TOLERANCES			
0-500mm	± 0.5		
500 +	± 1.0		
HOLES TO BE ± 0.1			
UNLESS OTHERWISE STATED			
DIMENSIONS IN MILLIMETRES			
THIS DESIGN IS THE PROPERTY OF PURITE LTD. AND IS ONLY ALLOWED TO BE USED BY EXPRESS PERMISSION AND LICENCE FROM PURITE LTD.			
PURITE LTD. BARRET WAY, FRAME INDUSTRIAL ESTATE, FRAME, OXON. OX8 3SU, ENGLAND. TEL : 01844 217141 FAX : 01844 218086		TITLE FLOW DIAGRAM SELECT FUSION 40 BP	
ISSUE	DATE	MODIFICATION	APPRO BY
A	13.6.06	was issue 1	
SCALE	X	PART NO.	ORDER NO.
DRAWN	J.COOLE	L988303	
APPROVED		DRG NO.	ISSUE
DATE	5.4.04	W3909	A

UNCONTROLLED COPY
UPDATES WILL NOT BE MADE TO THIS COPY



MATERIAL		PURITE LTD.	
FINISH		BANDEY WAY, THAME INDUSTRIAL ESTATE, THAME, OXON. OX9 3SL. ENGLAND.	
CORNERS TO BE REMOVED		TEL : 01844 217141 FAX : 01844 218088	
TOLERANCES		TITLE	
0-500mm	± 0.5	FLOW DIAGRAM	
500 +	± 1.0	SELECT FUSION 80 BP	
HOLES TO BE ± 0.1		SCALE	X
UNLESS OTHERWISE STATED		DRAWN	J. COLE
DIMENSIONS IN MILLIMETRES		APPROVED	
		DATE	5.4.04
		ISSUE	A
		ORDER NO.	
		PART NO.	L98B304
		DRG NO.	W3910
		ISSUE	A

UNCONTROLLED COPY
UPDATES WILL NOT BE MADE TO THIS COPY



MATERIAL		PURITE LTD.	
FINISH		BANDLET WAY, THAME INDUSTRIAL ESTATE, THAME, OXON. OX9 3SL, ENGLAND.	
ALL BURRS AND SHARP CORNERS TO BE REMOVED		TEL : 01844 217141 FAX : 01844 216088	
TOLERANCES		TITLE	
0-500mm	± 0.5	SELECT FUSION 160 BP	
500-1000mm	± 1.0	SCALE	X
1000mm+	± 1.5	DRAWN	J.COOLE
UNLESS OTHERWISE STATED		APPROVED	
DIMENSIONS IN MILLIMETRES		DATE	5.4.04
		ISSUE	A
		DATE	13.6.06
		MODIFICATION	
		APPD BY	
		ISSUE	A
		DATE	13.6.06
		WAS ISSUE 1	
		PERMISSION AND LICENCE FROM PURITE LTD.	
		IT IS THE PROPERTY OF PURITE LTD. AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS WITHOUT THE WRITTEN	
		THIS DESIGN IS THE PROPERTY OF PURITE LTD. AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS WITHOUT THE WRITTEN	

UNCONTROLLED COPY
UPDATES WILL NOT BE MADE TO THIS COPY