

UV-VIS Spectrophotometer

# UV-1900



SHIMADZU



# Navigate Your Way

## Easy to Operate, Obtain Answers Easily and Rapidly

Easy-to-use user interface design

Ergonomic touch-screen display

## Advanced Regulatory Compliance

Validation functions enable checks in accordance with Pharmacopeia (JP, USP, and EP) to be performed easily

In combination with LabSolutions™ DB/CS, comply with FDA 21 CFR Part 11 and PIC/S GMP guidelines

## High Performance to Meet Diverse Needs

Performance at the highest level in its class, provides advanced function than UV-1800


Ultra-fast scan performance, capable of obtaining high-accuracy spectra in just a few seconds

# Easy to Operate, Obtain Answers Easily and Rapidly



## Easy-to-Use Interface Grasp the Current Status and Operating Procedures at a Glance

The UV-1900 on-screen user interface includes large, easy-to-see icons deployed on a black background, so the instrument settings are evident at a glance. In addition, the large, easy-to-see icons improve intuitive understanding, which enables users to quickly become familiar with the operations. Furthermore, the user interface is designed to minimize transitions between windows, so users do not get confused during the operations.






## Hardware Design Based on Ergonomics

The control panel design is based on ergonomics and positioned at the very best viewing angle for the user. Users can operate easily with any posture. In addition, a stylus pen is equipped next to the panel, so users can switch easily between operations using their fingers or the stylus pen.



## Navigation Tabs Improve Usability

In quantitation mode on the UV-1900, the stages of the entire measurement process and the current status are always shown on the display. As a result, users know immediately what to do in the next step.




# Advanced Regulatory Compliance


## Full Support for Pharmacopeia, GLP/GMP, FDA 21 CFR Part 11 and Other Regulations

### Instrument Validation Functions Compliant with JP, USP, and EP

This instrument can not only run checks for nine JIS items, but also those stipulated in the Japanese Pharmacopoeia (JP), United States Pharmacopeia (USP), and the European Pharmacopoeia (EP). Naturally, the hardware is also compliant with the specifications required by each Pharmacopeia. In addition, the check conditions can be saved. As a result, once the conditions are saved, checks can be performed easily just by calling them up as needed. Check results can also be saved.



Main Screen



Parameters Setting Screen



During Testing Screen

The Check Results Screen shows a table of test results for '[EP]Control of Wavelengths'. The table has columns for Test W. (nm), Meas. W. (nm), Diff. (nm), Criteria (nm), Result, and 1/4. The data is as follows:

Test W. (nm)	Meas. W. (nm)	Diff. (nm)	Criteria (nm)	Result	1/4
241.21	241.80	0.59	1.00	Pass	▲
287.30	287.80	0.50	1.00	Pass	▲
361.31	361.55	0.24	1.00	Pass	▲
536.63	536.95	0.32	3.00	Pass	▼


Check Results Screen

## Resolution of 1 nm, the Highest in its Class

In addition to achieving a resolution of 1 nm, the highest in its class, by using a monochromator with a Czerny-Turner mounting, the UV-1900 also features a compact, bright optical system. The instrument is more than capable of meeting the wavelength resolution required in the European Pharmacopoeia.

## Improved Security Functions

An external control security function has been added to provide more support for compliance with regulations. Three user authority levels, "Administrator", "Developer", and "Operator", can be set for instrument users.



## Support for FDA 21 CFR Part 11, PIC/S GMP Guidelines and Other Regulations and Guidelines

Ensuring the integrity of data (database management), including the user management, user authority management, and data audit trails required for compliance with FDA 21 CFR Part 11, PIC/S GMP guidelines, and other ER/ES regulations, is possible.


### LabSolutions DB UV-Vis or UVProbe / LabSolutions DB System

The system allows for data management and user management with a database. Compliant with ER/ES regulations, the system is optimally configured for customers using a PC.



### LabSolutions CS UV-Vis\* or UVProbe / LabSolutions CS System (Network System)

The system is optimally configured for customers who want to manage data on a server together with LC and GC data for ER/ES compliance.



\*: coming soon

# High Performance to Meet Diverse Needs




## Low Stray Light

Stray light is at 0.5 % max. (198 nm), twice as low as the performance level of the UV-1800. With this stray light reduction, accurate measurements are possible up to the vicinity of 2 Abs even in the ultraviolet region. In addition, high-concentration samples can be quantified accurately.

The figure on the right is a calibration curve for acetic acid, created with absorbance at 200 nm.


The correlation coefficient is 0.9997, and correct measured values are obtained even in the vicinity of 2 Abs. Linearity will be lost in the high absorbance region due to the stray light.



## Ultra-Fast Scan

Spectra can be acquired as fast as 29,000 nm/min. Ultra-fast scan is effective in tracking chemical reactions in a short time. In addition to the absorbance change at specified wavelengths, spectra can also be acquired in a short time with the UV-1900. Therefore, more detailed behavior can be investigated by observing spectra with the UV-1900.

The figures below show the analysis of particle agglomeration process when salts are added to silver nanoparticles. Measurements of the 300 to 700 nm region were performed in ultra-fast scan mode. In addition to the decrease of absorbance at 400 nm and the increase of absorbance at 480 nm, the temporal changes of spectra can also be observed.




## High Reproducibility and Repeatability Accuracy

The photometric repeatability accuracy is 0.0002 Abs max. (0.5 Abs and 1.0 Abs), an improvement of five times over the performance level of the UV-1800. With this high photometric repeatability accuracy, variance in the measurement results is suppressed, enabling more accurate quantitation and the detection of low-concentration samples.

The figure on the right is a calibration curve for caffeine, created with absorbance at 273 nm. The calibration curve has an Abs = 0.0528 Conc., the lower limit of quantitation determined from the standard deviation is 0.0051 mg/L at a point where it would be 0.051 mg/L<sup>Note</sup> for the UV-1900.

Note: One method of determining the lower limit of quantitation is to use ten times the standard deviation. This is an actual measured value and is not guaranteed.

No.	Absorbance of Blank Solution (273 nm)
1	-0.00001
2	0.00001
3	-0.00002
4	0.00002
5	0.00001
6	-0.00003
7	0.00001
8	-0.00004
9	0.00001
10	0.00005
Standard Deviation $\sigma$	0.000025





# A Diversity of Measurement Modes

## Photometric


Measures the photometric value at a single wavelength or multiple (up to eight) wavelengths.

Reset K [ 1.2000]

No.	Sample Name	Abs	K * Abs	7/10
003	SAMPLE3	0.2493	0.2992	
004	SAMPLE4	0.4547	0.5456	
005	SAMPLE5	0.4549	0.5459	
006	SAMPLE6	0.6746	0.8095	
007	SAMPLE7	0.6746	0.8095	
008	SAMPLE8	0.6743	0.8092	
009	SAMPLE9	0.8801	1.0561	
010	SAMPLE10	0.8801	1.0561	



## Spectrum

Measures a sample spectrum using wavelength scanning.




## Quantitation

Generates a calibration curve from the measurement of standards, and then calculates the concentrations of unknowns.



[Activity Calculation Range] 5.0 s - 40.0 s

No.	Sample Name	Init.(Obs)	ΔAbs Activity	2/2
001	SAMPLE1	0.9356	-0.6813	0.5723
002	SAMPLE2	1.0489	-0.6451	0.5419



Results 3/13

SAMPLE5	
A1( 260.0 ) =	0.3269
A2( 230.0 ) =	0.3094
Ab( 320.0 ) =	0.1501
Abs Ratio = 1.1099	
DNA Conc. = 6.1265	
Protein Conc. = 15.750	

## Kinetics

Measures absorbance changes as a function of time, and obtains the enzymatic activity values. The kinetics measurement method or the rate measurement method can be selected.

## Time Course

Measures changes over time in photometric values at a specified wavelength.


## Biomethod

Quantifies DNA or protein concentrations.

# Applications

## Foods

This is an example of the analysis of food dyes. By using ultra-fast scan mode, the time needed for routine spectral checks can be shortened. The 350 nm to 900 nm region can be measured at 1 nm intervals in approx. 4 seconds.



## Pharmaceuticals and Life Sciences

This is an example of the analysis of  $\lambda$ DNA. Trace quantities (on the order of a few  $\mu$ L) can be measured by combining the instrument with NanoStick and TrayCell™.


Using TrayCell, a calibration curve for 4  $\mu$ L of  $\lambda$ DNA was obtained correctly in the range between 25  $\mu$ L and 500  $\mu$ L.



TrayCell





NanoStick



## Chemistry

This is an example of the analysis of colored cellophane tape. Materials can be confirmed quantitatively by using LabSolutions UV-Vis and color measurement software.



## Optional Software

# Control with LabSolutions UV-Vis Software



The UV-1900 can be controlled using LabSolutions UV-Vis software. LabSolutions UV-Vis is a next-generation Shimadzu UV control software pursuing efficiency of analysis. The simple design layout enables even first-time users to perform operations easily. A new spectrum evaluation function automates the measurement, analysis, and printing of reports to significantly enhance the efficiency of routine analysis. In addition, it achieves the easy transfer of measurement data. Users can export the measurement data in text format and import it into other software for analysis with Excel®.

(A separate USB cable is required to connect with a computer.)


For more details, refer to LabSolutions UV-Vis brochure (C101-E147).

Note: LabSolutions UV-Vis is the latest optional software. UV-1900 is equipped with UVProbe software as standard.

## Simple Design


From the start, the software's user-friendliness allows users to perform operations with ease.

With extensive features, LabSolutions UV-Vis meets a wide range of users' expectations.




Simple Main Window

Clear and simple layout of the graphs and measurement results makes it easy to read.



Instrument Control Panel

The instrument control panel that brings together the measurement functions enables automatic measurement, analysis and reporting.



Easy-to-follow Configuration Window

Large icons make it easy for users to understand and operate.


## Four Measurement Modes

It permits four measurement modes: spectrum, quantitative, photometric, and time course.


Users can open multiple measurement modes at the same time, so that data analysis can be performed in one mode while collecting data in another mode.

## Report

Easily create report layouts.




Print from the measurement window with one click.



Report is printed.

# Spectra Evaluation Function

In addition to providing measurement and analysis results, judgment results are also provided.  
With this feature, LabSolutions UV-Vis enables users to maintain a product's quality.




Quality of the sample can be determined with comprehensive judgment at a glance.

Multiple evaluation criteria can be set.

Detail Settings Window  
Evaluation method can be selected from a wealth of choices.


## Easy Transfer of Measurement Data

Users want to export measurement data immediately in text format, and import for analysis in other software, such as Excel.




### To Analysis Software

Automatically generates a text file when the spectra data are saved. It can be immediately imported into other software.



### To Excel


Real-time transfer of the spectrum waveform to Excel during measurement.  
No need to create a CSV file each time.




### Matrix Output

Outputs multiple spectra to one text file.  
Easy to import data into multivariate analysis software.

Line up wavelengths vertically.



Line up wavelengths horizontally.



Select how the data is ordered



## Control with UVProbe Software

## Spectrum Module

### Photometric Module (Quantitation)

### Kinetics Module (Time Course Measurement)

Report Generator

In addition to peak detection, area calculation and other data processing functions, UVProbe is equipped with various functions including: security functions that limit each user limited to specific functions, a data history log function, and an instrument audit trail function.

[illegible]

**Main Table**

File	Activity	Initial Reading	ID	Multiple Mils	Factor	Correct
Raw-Cuts	0.1566	0.1377	0.0087	-2.1758	1.0000	1.0000
Raw-Cuts	0.3188	0.2250	0.0087	-4.4195	1.0000	1.0000
Raw-Cuts	0.2615	0.1943	0.0077	-3.6261	1.0000	1.0000
Raw-Cuts	0.2176	0.1583	0.0087	-3.5973	1.0000	1.0000
Raw-Cuts	0.1737	0.1401	0.0084	-2.4091	1.0000	1.0000
Raw-Cuts	0.4000	0.3545	0.0087	-5.5455	1.0000	1.0000

**Method** LaserSharp-Buck  $\frac{1}{V} \frac{K_m}{V_m} \left[ \frac{1}{S} \right] \frac{1}{V_m}$

File	EQ	Weight	1	EQS	to
Raw-Cuts	5.7100	1	0.1566	0.1750	6.3681
Raw-Cuts	13.3000	1	0.3188	0.2762	3.1370
Raw-Cuts	10.0000	1	0.2615	0.1906	3.6274
Raw-Cuts	8.0000	1	0.2176	0.1290	4.5847
Raw-Cuts	6.4000	1	0.1737	0.1498	5.7572
Raw-Cuts	20.0000	1	0.4000	0.3550	2.5500

**Description**

Description	Date/Time	User
Scan motor initialization - Passed	01/19/00 PM 4:43:02	Administrator
Wavelength origin 1 search - Passed	01/19/00 PM 4:44:02	Administrator
Wavelength origin 2 search - Passed	01/19/00 PM 4:44:02	Administrator
Wavelength origin 3 search - Passed	01/19/00 PM 4:44:02	Administrator
Wavelength origin 4 search - Passed	01/19/00 PM 4:44:02	Administrator
Wavelength origin 5 search - Passed	01/19/00 PM 4:44:02	Administrator
Baseline - approx. 3500 nm - 6000 nm at 500 nm	01/19/00 PM 4:44:02	Administrator

**500.000 nm 0.000 Abs.**

**Cell Blank** **Auto Zero** **Go To W** **Start** **Discard**

**Standard Curve**

Sample ID	Type	Ex	Conc	VL550	VL550	Result
1	Oye A	Standard	1.000	0.136	0.121	0.013
2	Oye B	Standard	2.250	0.149	0.134	0.016
3	Oye C	Standard	3.900	0.283	0.252	0.032
4	Oye D	Standard	4.750	0.440	0.388	0.061

**Sample Data - (Active)**


Sample ID	Type	Ex	Conc	VL550	VL550	Result
1	Sample A	Unknown	3.291	0.264	0.216	0.026
2	Sample B	Unknown	3.067	0.240	0.224	0.026
3	Sample C	Unknown	3.246	0.412	0.384	0.020
4	Sample D	Unknown	3.923	0.427	0.387	0.046
5	Sample E	Unknown	3.801	0.349	0.311	0.030
6	Sample F	Unknown	2.272	0.308	0.282	0.014
7	Sample G	Unknown	3.471	0.281	0.250	0.031
8	Sample H	Unknown	2.272	0.189	0.134	0.014
9	Sample I	Unknown	1.415	0.768	0.793	0.017
10	Sample J	Unknown	3.471	0.768	0.793	0.017

**Description**  
 Gene master initialization - Passed  
 Wave-length origin 1 search - Passed  
 Instrument History

The screenshot shows the 'Spectrum' application window. The title bar reads 'Spectrum - report'. The menu bar includes 'File', 'Edit', 'View', 'Format', 'Tools', 'Window', and 'Help'. The toolbar contains various icons for file operations, editing, and analysis. The main window displays a 'Spectrum Peak Pick Report' for 'Data Set: Storage 1 - RawData-020 - D:\Program Files\SHIMADZU\MSProber\Data\CermetMetWtl SalsynData.spc'. The report includes a plot area with a gray background and axes ranging from 0.00 to 1.00. Below the plot, there are two empty rectangular boxes labeled 'Notes/Parameters' and 'Data'. The status bar at the bottom indicates 'Scan 1 / 1', 'Run: 655.050 in', and 'MIM'.

**Spectrum Peak Pick Report** 01/26/01 06:06:01

Data Set: Storage 055829 PM - RawData - D:\Program Files\SHIMADZU\Analysis\DATA\HATCH\055829 PM.SPC

Abs.  nm.

Measurement Properties  
Wavelength Range (nm): 300.0 to 500.0  
Scan Speed: Fast  
Sweeping Interval: 0.2

Instrument Properties  
Measuring Mode: Abs.  
Slit Width: 2.0

No.	WV	Wavelength	Abs.	Description
01	335	335.0	0.23	
02	340	340.0	0.80	
03	345	345.0	0.50	
04	350	350.0	0.20	
05	355	355.0	0.20	
06	360	360.0	0.20	
07	365	365.0	1.40	
08	370	370.0	0.20	
09	375	375.0	0.20	
10	380	380.0	0.20	
11	385	385.0	0.20	
12	390	390.0	0.20	
13	395	395.0	0.20	
14	400	400.0	0.20	
15	405	405.0	0.20	
16	410	410.0	0.20	
17	415	415.0	0.20	
18	420	420.0	0.20	
19	425	425.0	0.20	
20	430	430.0	0.20	
21	435	435.0	0.20	
22	440	440.0	0.20	
23	445	445.0	0.20	
24	450	450.0	0.20	
25	455	455.0	0.20	
26	460	460.0	0.20	
27	465	465.0	0.20	
28	470	470.0	0.20	
29	475	475.0	0.20	
30	480	480.0	0.20	
31	485	485.0	0.20	
32	490	490.0	0.20	
33	495	495.0	0.20	
34	500	500.0	0.20	

Page 1/1



# Accessories

## Film Holder

(P/N 204-58909)

Used in transmittance measurement of thin samples such as films and filters. Holds thin samples, such as films and filters, for analysis.



## Long-Path Rectangular Cell Holder

(P/N 204-23118-01)

Holds two rectangular cells with an optical path length of 10, 20, 30, 50, 70, or 100 mm.



## Sipper Unit

Model	P/N	Standard Sample Volume
Sipper Unit 160L (Standard Sipper)	206-23790-51	2.0 mL
Sipper Unit 160T (Triple-Pass Sipper)	206-23790-52	1.5 mL
Sipper Unit 160C (Constant-Temperature Sipper)	206-23790-53	2.5 mL
Sipper Unit 160U (Supermicro Sipper)	206-23790-54	0.5 mL

Four types of sipper units with different flow cells are available. The stepping motor-driven peristaltic pump ensures reliable and smooth aspiration of sample solution.

(Direct driving is possible from the UV-1900 so no interface is required.)

## CPS-100 Cell Positioner, Thermoelectrically Temperature Controlled

(P/N 206-29500-\*\*)

This attachment permits measurement of up to six sample cells under constant-temperature conditions. Combination of this attachment and the Kinetics mode provides measurement of temperature-sensitive enzyme kinetics of one to six samples.

- Number of cells: 6 on the sample side (temperature-controlled)  
1 on the reference side (temperature not controlled)
- Temperature control range: 16°C to 60°C
- Temperature display accuracy (difference from the true value):  $\pm 0.5^{\circ}\text{C}$
- Temperature control precision (variation of temperature):  $\pm 0.1^{\circ}\text{C}$
- Ambient temperature: 15°C to 35°C

Note: Square cells (P/N 200-34442) are not included, please purchase separately.

A USB adapter CPS

(P/N 206-25234-91) is required.




## Test Tube Holder

(P/N 207-23510-41)

Holds test tube instead of sample compartment.

- Specifications:  
Outside diameter:  $\phi 15$  to 18 mm  
Height: 84 to 115 mm  
Note: Test tube is not included.





## Multi-Cell Sample Compartment

(P/N 206-69160-41)

Holds up to six 10-mm square cells on the sample side. No temperature control capability.

- Number of cells: 6 on the sample side  
1 on the reference side

Note: Square cells are not included, please purchase separately.



Note: The use of a Solenoid Valve (Fluoropolymer) (P/N 204-06599-01) and the SWA-2 Sample Waste Unit (206-23820-58) are recommended when strong acids, strong alkalis, or organic solvents are to be measured.

## TCC-100 Thermoelectrically Temperature Controlled Cell Holder

(P/N 206-29510-\*\*)

Uses Peltier effect for controlling the sample and reference temperature, so no thermostated bath or cooling water is required.

- Number of cells: One each on the sample and reference sides (temperature-controlled)
- Temperature control range: 7°C to 60°C
- Temperature display accuracy (difference from the true value):  $\pm 0.5^{\circ}\text{C}$
- Temperature control precision (variation of temperature):  $\pm 0.1^{\circ}\text{C}$

Note: Square cells (P/N 200-34442) are not included, please purchase separately.



## Typical Specifications

# UV-1900

## UV-VIS Spectrophotometer

The UV-1900 is a double-beam UV-Vis spectrophotometer using Shimadzu's original LO-RAY-LIGH™ diffraction grating technology. In addition to its high optical performance, the UV-1900 features high resolution, low stray light, high reproducibility, and an ultra-fast scan function. It also has an easy-to-use interface on a color touch-screen display. The UV-1900 is designed to meet the needs of both high performance and usability.



## Hardware Specifications

Item	Specification
Wavelength range	190 to 1,100 nm
Spectral bandwidth	1 nm (190 to 1,100 nm)
Wavelength display	0.1 nm increments
Wavelength setting	0.1 nm increments (1 nm increments when setting scanning range)
Wavelength accuracy	± 0.05 nm at D2 peak 656.1 nm, ± 0.3 nm for entire range
Wavelength repeatability	± 0.1 nm
Wavelength slew rate	About 14,500 nm/min
Wavelength scanning speed	3,000 to 2 nm/min 29,000 nm/min when survey scanning
Lamp interchange wavelength	Automatic interchange linked to wavelength. The interchange wavelength can be set freely in the range of 295 to 364 nm (0.1 nm increments).
Stray light	Less than 0.02% at 220 nm (NaI) Less than 0.01% at 340 nm (NaNO <sub>2</sub> ) Less than 0.5% at 198 nm (KCl)
Photometric system	Double beam optics
Photometric range	Absorbance: -4 to 4 Abs Transmittance: 0% to 400%
Photometric accuracy	± 0.002 Abs at 0.5 Abs ± 0.0025 Abs at 1.0 Abs ± 0.006 Abs at 2.0 Abs (measured using NIST930D/NIST1930 or equivalent.)
Photometric repeatability	Less than ± 0.0001 Abs at 0.5 Abs Less than ± 0.0001 Abs at 1 Abs Less than ± 0.0005 Abs at 2 Abs

The specifications shown here represent the average performance of the UV-1900. These specifications are typical values, not guaranteed values. The guaranteed specifications are listed in a separate publication.

Item	Specification
Baseline stability	Less than 0.0003 Abs/Hr (700 nm, one hour after light source turned ON)
Baseline flatness	Less than ± 0.0006 Abs (1,100 to 190 nm, one hour after light source turned ON)
Noise level	Less than 0.00003 Abs (700 nm)
Light source	20-W halogen lamp and deuterium lamp Built-in light source auto position adjustment
Monochromator	LO-RAY-LIGH grade blazed holographic grating in Czerny-Turner mounting
Detector	Silicon photodiode
Sample compartment	Internal dimensions: W110 × D250 × H115 mm Distance between light beams: 100 mm
Power requirements	AC 100, 120, 220, 230, 240 V, 50/60 Hz, 140 VA
Environmental requirements	Temperature: 15°C to 35°C Humidity: 30% to 80% (without condensation; 70% max. at 30°C or higher)
Dimensions	W450 × D501 × H244 mm
Weight	16.6 kg
Output device	USB memory (optional) Data files saved in text format or UVPC format. UVPC-format files can be read directly by UVProbe and LabSolutions™ UV-Vis.
PC compatibility	UVProbe software (standard) LabSolutions UV-Vis software (optional) External control possible via USB.
Display	24-bit color touch screen
Supported languages	Japanese, English, Chinese, Spanish (Mexico), Portuguese (Brazil).

## Software Specifications

Measurement mode	Specification
Photometric mode	Single-wavelength measurement 1. Photometric modes: T% or Abs 2. Quantitation using K-factor method 3. Data table storage and recall functions Multiple-wavelength measurement 4. Photometric modes: T% or ABS 5. Measurements at up to eight designated wavelengths (set in 0.1 mm increments) 6. Data calculation at up to four wavelengths (difference or ratio between two wavelengths, calculation between three wavelengths, etc. ) is possible.
Spectrum mode	1. Measurement modes: ABS, T%, E 2. Number of repeat scans: 1 to 99 3. Recording system: Selection between single spectrum and data overlay 4. Data storage and recall 5. Data processing: Peak/valley detection, arithmetic operations, differentiation, smoothing, area calculation, point picking, data reading at cursor-specified point
Quantitation mode	1. Measurement methods: 1-wavelength, 2-wavelength, 3-wavelength, and 1st to 4th derivative methods 2. Quantitation methods: Automatic concentration calculation using K-factor Automatic concentration calculation using single-point calibration curve Multi-point calibration curve method (1st to 3rd order regression curves) 3. Measurement parameters: Number of standards (2 to 10) Number of repeat measurements (1 to 10 times) to obtain a mean value for quantitation.
Kinetics mode	1. Measures absorbance changes as a function of time and calculates the enzymatic activity values. 2. Measurement time: 1 to 9,999 sec/min 3. Measurement methods: 1-wavelength, 2-wavelength, multi-cell, and rate measurements
Time scan mode	1. Measures changes in measured values as a function of time. 2. Measurement mode: ABS, T%, E 3. Measurement time: 1 to 9,999 sec/min 4. Data processing functions (same as spectrum mode)

Measurement mode	Specification
Biomethod mode	DNA/Protein Quantitation 1. Calculation of DNA/protein concentration and absorbance ratio DNA concentration = $K1 \times A1 - K2 \times A2$ Protein concentration = $K3 \times A2 - K4 \times A1$ 2. Factors and measurement wavelengths can be set freely. 3. Background correction is possible. Quantitation of proteins 1. Quantitation methods: Lowry method, BCA method, Biuret method, CBB method (Bradford method), UV method
Maintenance	1. Baseline correction 2. Lamp usage time display and reset. 3. Security settings Functions can be restricted according to the user. 4. Instrument validation functions: 1) Compatible with 9 JIS items Wavelength accuracy, wavelength repeatability, resolution, stray light, photometric accuracy, photometric repeatability, baseline flatness, baseline stability, noise level. 2) Semi-automatic validation Validation inspections conducted interactively while inserting and removing inspection jigs. 3) Fully automatic validation Automatic validation inspections from measurement to evaluation and printout. 4) Setting inspection parameters and pass/fail criteria Authority to make changes can be protected by password access. 5) Detailed printout of results. 6) Bulk printout of results. 7) Equipped with method in accordance with Pharmacopeia (JP, USP, EP).
Shared functions	1. Automatic setting of measurement mode after instrument initialization. It's possible to specify standby and parameter files in the parameter setting window for each measurement mode. 2. Selection of displayed number of decimal places Absorbance: 3 or 4 decimal places Transmittance: 1 or 2 decimal places 3. Number of files that can be saved (built-in memory) Measurement parameters: 100 files max. Tabular data: 15 files max. Curve data: 16 files max. Validation condition: 10 files max. Validation result: 3 files max. 4. Setting of integration time (for fixed-wavelength measurement) 5. PC control Spectrophotometer can be controlled by an external PC. This function is also used when performing operation with the standard UVProbe software and the optional LabSolutions UV-Vis software provided. *A USB cable is required.

LO-RAY-LIGH and LabSolutions are trademarks of Shimadzu Corporation.



Shimadzu Corporation

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2018

First Edition: April 2018



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.