
Homogenisers

SHM1 SMH2

Instructions for use
Manuel d'utilisation
Istruzioni per l'uso
Instrucciones de funcionamiento
Bedienungsanweisung

Homogenisers

SHM1 & SHM2

Instructions for use

Introduction

Thank you for purchasing this piece of Stuart equipment. To get the best performance from the equipment and for your own safety, please read these instructions carefully before use. Before discarding the packaging check that all parts are present and correct.

This equipment is designed to operate under the following conditions:

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuation not exceeding 10%
- ❖ Over-voltage category II IEC60364-4-443
- ❖ Pollution degree 2

If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Stuart, the protection provided may be impaired.

Safety Advice Before Use

- ❖ Do not allow the motor to be submerged in any liquid.
- ❖ Do not operate the product if it has been damaged in any way.
- ❖ Do not use attachments not recommended by Stuart.
- ❖ Do not attempt to stop movement by hand.
- ❖ Do not attempt to modify any part of the product.
- ❖ Keep this product away from heated surfaces.
- ❖ Do not run the generator probes in air, only in the presence of fluid.

General Description

The Stuart SHM1 and SHM2 tissue homogenisers are variable speed, hand-held or retort-stand mounted homogenisers. They combine a high-speed, high-torque motor with a choice of autoclavable rotor-stator generator probes. The model SHM1 tissue homogeniser is suitable for most laboratory applications such as tissue homogenisation, production of emulsions and suspensions, cell disruption and protein extraction. It is lightweight and can be hand-held or mounted to a retort stand for lengthier processing. The SHM2 is a powerful, versatile and durable homogeniser that can handle virtually any processing application, making it ideal for tougher samples and larger volumes.

The SHM1 is supplied with a post clamp assembly for mounting it to a retort stand. A retort stand assembly, catalogue number SHM/STAND, is available as an accessory for both models. Both models are also supplied with a tool kit for dismantling the rotor-stator generator probes for easy cleaning.

Theory of Operation

The rotor shaft part of the probe is coupled directly to the drive motor of the homogeniser. When attached, the rotor shaft can spin up to 35,000rpm; this makes up the rotor part of the rotor-stator generator probe. The tube/collar assembly which fits over the rotor shaft is attached to the motor housing but does not spin. This is the stator part of the rotor-stator generator probe.

As the rotor knife spins within the tube/collar assembly, it creates a pumping action, pulling the sample into the open end of the probe and forcing it out through the windows in the tube. This sets up a shearing action, reducing the particle size of the sample. The speed differential between the rapidly moving portion and the relatively stationary portion of the sample sets up a second force called cavitation which pulls the sample apart, further reducing the particle size.

Processing efficiency can be affected by:

- ❖ Amount of material processed vs. size and speed of the generator probe.
- ❖ Container geometry (round vessels encourage swirling, while fluted or cornered vessels disrupt flow patterns for more effective mixing/processing).
- ❖ Processing speed vs. optimal speed.
- ❖ Size and type of material and flow characteristics (material particles must be small enough to be carried into the generator head for optimal processing).

Preparation for Use

1. DO NOT modify the plug or cord that is provided. If the plug will not fit the outlet, have the proper outlet installed by a qualified electrician. Damaged or worn power cords should be repaired or replaced immediately by a qualified electrician.
2. If using a retort stand assembly, attach the unit to the stand and clamp firmly in position (see details below). Place the stand on a firm, level, non-slip surface ensuring that there is sufficient free space on all sides without coming into contact with anything during use.
3. Make sure the power switch is in the OFF position before connecting to the electricity supply.

Operation: SHM1

Stand-mounting the SHM1 (optional)
Please follow the steps given below:

1. Secure the post by twisting it into the base plate until securely tightened.
2. Loosen the locking knob and slide the post clamp assembly down over the end of the post until the clamp is at the desired height. Lock in place.
3. Slide the motor drive into the opening on the end of the clamp. Push down firmly until it locks into place.
4. To remove the SHM1, press the black button on the back of the opening of the post clamp assembly and lift the motor from the cradle.

Attaching the rotor-stator generator probes

The rotor-stator generator probes use a quarter-turn bayonet mount. Simply push the probe into the motor housing as far as possible, turn clockwise and release. Remove the protective cap from the tip of the probe prior to use. Removal of the probe is the opposite of installation.

If using plastic disposable rotor-generator probes, the adapter, catalogue number SHM/ADAPT, will be required. The adapter also uses a quarter-turn bayonet mount. Push the adapter into the motor housing as far as possible, turn clockwise and release. Next, attach the plastic probes as described above.

Processing the sample

NOTE: As a safety precaution, it is recommended that safety glasses are worn when using the homogenisers to protect the user from potential splashes.

1. Remove the protective cap from the tip of the probe.
2. Insert the probe to a depth of one third to one half of the liquid height.

NOTE: Heavy sediments may require deeper immersion and this can be optimised by observation during processing. The top hole of the probe should not be immersed in the sample as this could result in liquid being drawn into the lower motor bearing.

When using disposable plastic probes, the probes should not be immersed above the cone leading up to the motor drive during processing, since liquid may otherwise be drawn up the shaft and contaminate the coupling.

3. Ensure the speed control switch is in the far left-hand position (refer to Fig 1).


Fig 1

4. Switch the power switch to ON and slide the speed control button to the required processing speed.

For most samples, three to four short bursts at top speed should be sufficient for complete homogenisation of the sample.

5. When processing is complete, switch the power switch to OFF and remove the probe for cleaning. Replace the protective cap on the end of the probe when not in use.

Cleaning

For cleaning of the rotor-stator generator probes, please see the section below on disassembly of the probes.

The motor drive housing should be wiped off after use, especially when concentrated and potentially damaging liquids are used during processing. Never use solvents to clean the unit or accessories.

Operation: SHM2

Stand-mounting the SHM2 (optional)

Please follow the steps given below:

1. Secure the post by twisting it into the base plate until securely tightened.
2. Loosen the locking knob and slide the post clamp assembly down over the end of the post until the clamp is at the desired height. Lock in place.
3. Twist the cross rod into the back of the motor housing until tight.
4. Loosen the second knob on the post clamp assembly and insert the cross rod with the motor attached (approximately 25mm of the cross rod should protrude from the back of the post clamp).
5. Lock in place.

Attaching the rotor-stator generator probes

The rotor-stator generator probes use a quarter-turn bayonet mount. Simply push the probe into the motor housing as far as possible, turn clockwise and release. Remove the protective cap from the tip of the probe prior to use. Removal of the probe is the opposite of installation.

If using plastic disposable rotor-generator probes, the adapter, catalogue number SHM/ADAPT, will be required. The adapter also uses a quarter-turn bayonet mount. Push the adapter into the motor housing as far as possible, turn clockwise and release. Next, attach the plastic probes as described above.

Processing the sample

NOTE: As a safety precaution, it is recommended that safety glasses are worn when using the homogenisers to protect the user from potential splashes.

1. Remove the protective cap from the tip of the probe.
2. Insert the probe to a depth of one third to one half of the liquid height.

NOTE: Heavy sediments may require deeper immersion and this can be optimised by observation during processing. The top hole of the probe should not be immersed in the sample as this could result in liquid being drawn into the

lower motor bearing. When using disposable plastic probes, the probes should not be immersed above the cone leading up to the motor drive during processing, since liquid may otherwise be drawn up the shaft and contaminate the coupling.

3. Ensure the speed control dial is set to the lowest setting (refer to Fig 2).


Fig 2

4. Slide the power switch to ON and adjust the rotary speed control dial to the required processing speed. Speed settings are given in the table below. These are approximate and actual speed will vary depending on probe size and sample processed.

	Speed setting					
Approx.	1	2	3	4	5	6
rpm	4,000	5,000	7,500	10,000	13,500	25,000

For most samples, three to four short bursts at top speed should be sufficient for complete homogenisation of the sample.

5. When processing is complete, slide the power switch to OFF and remove the probe for cleaning. Replace the protective cap on the end of the probe when not in use.

Cleaning

For cleaning of the rotor-stator generator probes, please see section below on disassembly of the probes.


The motor drive housing should be wiped off after use, especially when concentrated and potentially damaging liquids are used during processing. Never use solvents to clean the unit or accessories.

Disassembly and Cleaning of Rotor-stator Generator Probes

Stainless steel probes: disassembly


The stainless steel probes are easily disassembled for cleaning. Refer to Fig 3 and follow the instructions given below.

1. Insert the proper size probe tool into the knife end of the probe to hold the knife (the tool is marked 5, 7, 10 or 20 which corresponds to the diameter of the probe in mm).

2. Insert the second probe tool over the driven end of the generator tube.


3. Holding both tools, carefully turn the upper probe tool anti-clockwise, unscrewing the shaft from the knife.


4. Remove the knife and pull the shaft out of the tube.

5. Remove the upper bearing from the shaft. 195mm and longer probes have mid-bearings which can be autoclaved inside the tube.

6. If cleaning a 5mm or 7mm probe, thread the appropriate bearing press onto the end of the shaft. Insert the shaft back into the generator tube and press out the lower bearing.

7. Remove the bearing press from the shaft before cleaning the probe parts.

NOTE: Inspect the bearing for wear or black particles. Clean any dirt from the bearing before reassembly. Replace the bearing if worn or damaged. New bearings must be run-in before initial use, placing the probe in clean water, as follows:

1. Run 4 minutes on lowest motor setting.
2. Run 5 minutes on medium motor setting.
3. Run 1 minute on highest motor setting.
4. Disassemble generator probe.
5. Clean bearings of run-in dust.
6. Reassemble generator probe.


Fig 3


Autoclaving the stainless steel probes

The generator probes may be sterilised by autoclaving up to 121°C. If the probe contains all Teflon bearings, then it can be autoclaved as an assembly.


Stainless steel probes: reassembly

1. Insert the lower bearing into the knife end of the generator tube.

2. For 5mm and 7mm probes, use the appropriate bearing press to set the lower bearing by using it to gently push the bearing into position in the tube. Remove the bearing press, leaving the bearing in place.


3. Slide the upper bearing onto the shaft so that the flanged side of the bearing is resting against the bottom of the shaft spacer.


4. Insert the shaft into the driven end of the generator tube.

5. Place the upper probe tool over the driven end of the generator tube. Holding it in place, invert the probe so that the knife end points up.

6. Insert the knife into the processing end of the generator tube and turn clockwise to begin threading it onto the shaft.

7. Using the other probe tool, select the correct size for the probe and insert it into the knife. Turn the rotor shaft clockwise until it stops. DO NOT TIGHTEN OR TORQUE THE KNIFE as this may cause it to bend or break.

8. The generator should spin freely, if not, disassemble and reassemble the probe to ensure the bearings are fully seated.
9. Run the generator probe in clean water for 5 minutes at slow speed, 5 minutes at medium speed then 1 minute at full speed. The knife is self-tightening and running the generator in water will tighten the knife onto the shaft to the proper torque.

Disposable plastic probes

Non-sterile disposable plastic probes are ideal for situations where molecular contamination between samples cannot be tolerated. The plastic probes are manufactured of polycarbonate and an amorphous thermoplastic polythermide and are designed to be disposed of after a single use. However they can be cleaned and autoclaved up to seven times for repeated use if required. The probes have good chemical resistance to weak acids, chlorides, hypochlorite solution (5%) and many other chemicals. The probes are not recommended for extended use with phenol or guanidine-based reagents and should be disposed of after a single use in these chemicals. If any cracking or brittleness of the probe is detected, it is recommended that it is discarded.

NOTE: The plastic disposable probes are not pre-sterilised. If this is a requirement please use one of the following methods of sterilisation:

- ❖ Autoclave up to 121°C
- ❖ Ethylene oxide
- ❖ Formaldehyde (g)
- ❖ Glutaraldehyde (2%)
- ❖ Gamma irradiation up to 5MRAD

UV irradiation is NOT recommended.

For biological disinfection of the disposable rotor-stator generator probes the following reagents can be used:

- ❖ Alcohol (70%)
- ❖ Hypochlorite (5%)
- ❖ Formaldehyde (40%)

Troubleshooting

Problem	Corrective action
The SHM1/SHM2 is plugged in and turned on but is not functioning.	Check power cord connectors. Check wall socket for power to the outlet.
The motor is turned ON and makes a "buzzing" sound but is not working.	The brushes may be worn or not installed correctly. Please contact the Service Department.
Motor unit operating speed declines, stalls intermittently, or stops completely.	The brushes may be worn or not installed correctly. Please contact the Service Department.
Teflon bearing wears quickly.	Fluid level may be too low in the tube. Immerse the probe deeper into the fluid.
Excessive splashing in sample tube.	Fluid level too low for tube size.

Maintenance, Servicing & Repair

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent. Do not use harsh or abrasive cleaning agents. Spillages or splashes should be cleaned up straightaway after isolating the unit from the mains electricity supply.

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

For a comprehensive list of parts required by service engineers conducting internal repairs, please contact the Service Department, quoting both the model number and serial number:

Tel: +44 (0)1785 810475

Fax: +44 (0)1785 810471

Only spare parts supplied or specified by Stuart or its agents should be used. Fitting of non-approved parts may affect the performance of the safety features designed into the instrument.

For any other technical enquiries, please contact the Technical Support Department:
E-mail: stuarthelp@bibby-scientific.com
Tel: +44 (0)1785 810433.

Accessories

Description	Catalogue Number
Retort stand assembly	SHM/STAND
Probe adaptor (for use with disposable plastic probes)	SHM/ADAPT
Generator probe 5mm x 75mm (flat bottom)	SHM/5
Generator probe 7mm x 115mm (saw tooth)	SHM/7
Generator probe 10mm x 115mm (saw tooth)	SHM/10
Generator probe 20mm x 195mm (saw tooth)	SHM/20
Generator probe, plastic disposable 7mm x 110mm (hard tissue)	SHM/DISP

Other generator probes may be available upon request. Please contact equipment@bibby-scientific.com

Warranty

Stuart warrants the SHM1 and SHM2 to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. The stainless steel rotor-stator generator probes are covered by a warranty period of one (1) year. In the event of a justified claim, Stuart will replace any defective component or replace the unit free of charge. This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorized personnel.

Bibby Scientific Ltd.
Beacon Road,
Stone, Staffordshire
ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
E-mail: equipment@bibby-scientific.com
www.stuart-equipment.com

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these

standards and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Bibby Scientific Ltd.

Declaration of Conformity

Homogenisers, Models, SHM1/UK, SHM1/EURO, SHM2/UK, SHM2/EURO

This product complies with the requirements of the EU Directives listed below:

2004/108/EC EMC Directive.
2006/95/EC Low voltage Directive (LVD)

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN 61326-1:2006 (Electrical Equipment for Measurement, Control and Laboratory use).
EN 61010-1: 2001 (Safety Requirements Electrical Equipment for Measurement, Control and Laboratory use)
EN61010-2-051: 2003 (Particular Requirements for Laboratory Equipment for Mixing and Stirring)

Compliance Certificates and Full Reports.
09-378A, 09-378B, 09-378C

The above certificates and reports, from a independent test house, are available upon request.

CE mark affixed '09.

Signed:  (Mr C. Warren)

Date: JUNE 2009

Authority: Technical Manager
Bibby Scientific Ltd

Bibby Scientific Ltd - Stone - Staffs - ST15 0SA - UK
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

INSPECTION REPORT

MODEL SHM1 ☐ SHM2 ☐


ELECTRICAL SAFETY

- | | |
|---------------------|-------------------------------------|
| 1. Earth continuity | ☒ |
| 2. Insulation | ☒ |
| 3. Flash test | ☒ |

FUNCTIONAL

- | | |
|----------------------|-------------------------------------|
| 1. Indicators | ☒ |
| 2. Speed control | ☒ |
| 3. Visual acceptance | ☒ |

QUALITY CONTROL INSPECTOR

Bibby Scientific France SAS
ZI du Rocher Vert - BP 79
77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bsf@bibby-scientific.fr
www.bibby-scientific.com

Bibby Scientific Italia Srl
Via Alcide de Gasperi 56
20077 Riozzo di Cerro al Lambro
Milano Italia
Tel: +39 (0)2 98230679
Fax: +39 (0)2 98230211
e-mail: marketing@bibby-scientific.it
www.bibby-scientific.it

Bibby Scientific Ltd
Beacon Road Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
e-mail: info@bibby-scientific.com
www.bibby-scientific.com

Bibby Scientific US Ltd
3 Terri Lane Suite 10
Burlington NJ 08016
USA
Tel: 800-225-9243
Fax: 609-589-2571
www.bibby-scientific.com

Bibby Scientific (Asia) Ltd
Room 607 Yen Sheng Centre
64 Hoi Yuen Road Kwun Tong
Kowloon Hong Kong
Tel: +852 3583 1581
Fax: +852 3583 1580
e-mail bibby@bibby-scientificasia.com
www.bibby-scientificasia.com