UV-C rays passive cabinet
for PCR samples preparation

aura PCR

OPERATING MANUAL

PART 1: USER

PART 2: SERVICE AND MAINTENANCE

INSTRUCTION TRANSLATED FROM THE ORIGINAL
80.021.261

Revision 00

February 2024

BioAir S.p.A.

Head offices: Via Lombardia, 12 - 27010 Siziano (PV) Italy

+39 0382 6672.1 - +39 0382 610924

www.bioair.it - info@bioair.it

CONTENTS

LABELLING	6
MEANING OF SYMBOLS	6
PURPOSE AND USE OF THE MANUAL	7
PART 1: USER	9
1. ASSESSMENT OF RISKS	11
2. DEFINITION	13
2.2 DEFINITION	13
3. INSTALLATION	14
3.1 CHOICE OF LOCATION	14
3.2 POSITIONING	14
3.3 CONNECTION TO THE MAINS	14
3.3.1 Mains socket	14
3.3.2 Mains supply cord	15
3.3.3 Mains appliance inlet and main fuses	15
3.4 FINAL PROCEDURES	15
4. TECHNICAL DATA	16
4.1 SPECIFICATIONS	16
4.2 USE ENVIRONMENTAL CONDITIONS	16
4.3 TRANSPORT AND STORAGE CONDITIONS	16
4.4 PESO E DIMENSIONI	16
4.5 MATERIALS	16
4.6 PERFORMANCES	17
5. OPERATING PRINCIPLE	18
6. CABINET DESCRIPTION	18
6.1 CONTROL BOX	18
6.1.1 Front view	18
6.1.2 Rear view	18
6.2 SAFETY UV-C DEVICE	19
6.3 UV-C LAMPS	19
6.3.1 PROTECTION AGAINST UV-C RADIATIONS	19
7. OPERATING INSTRUCTIONS	20
7.1 SWITCHING ON THE CABINET	20
7.2 SETTING OF TIMER	20
7.3 SETTING OF LIGHT SWITCH	20
7.4 ACTIVATION OF UV-C CYCLE	21
7.5 INTERRUPTION OF UV-C CYCLE	21
7.6 SWITCHING OFF THE CABINET	21
7.7 GENERAL GUIDELINES	22
7.7.1 Features of location	22
7.7.2 Working posture	22
7.7.3 Working procedure	22
8. CLEANSING AND ROUTINE MAINTENANCE	23
8.1 CLEANSING	23
8.1.1 Painted surfaces	23
8.1.2 Stainless steel surfaces	23
8.1.3 Glass surfaces	23
8.1.4 Base plate	23
8.1.5 UV and LED lamps	23

PART 2: SERVICE AND MAINTENANCE	24
9. INSTALLATION	26
9.1 <i>HANDLING AND UNPACKING</i>	26
9.2 <i>CHOICE OF LOCATION</i>	26
9.3 <i>POSITIONING</i>	26
9.4 <i>CONNECION TO THE MAINS</i>	27
9.5 <i>PERFORMANCE TEST</i>	27
9.5.1 Check of UV-C emissions efficiency:	27
9.5.2 Advised additional instruments:	27
10. TROUBLE SHOOTING	28
11. WIRING DIAGRAM	29

Dear Customer!

Thank you very much for purchasing a BioAir product. You have just acquired an high quality cabinet and we require you take a few minutes to read this operating and maintenance manual in order to become acquainted with your instrument and use it at the best of its performances.

We are confident that you have purchased an excellent product and that you will be completely satisfied with it. Should you have any enquiry, pls do not hesitate to contact your local Distributor or, directly us at BioAir, we will be delighted to assist you.

Your BioAir Team

If you require any further information or explanation of the information contained in this manual, and for any request relevant to:

- service
- spare parts
- technical data

contact:

WITHIN ITALY

SAT Italia

+ 39 0382 6672.1

+ 39 0382 610924

sat@bioair.it

ABROAD

contact your local distributor or

service@bioair.it

For any additional information, please contact the manufacturer

BioAir S.p.A.

Via Lombardia, 12
27010 Siziano (PV)

+ 39 0382 6672.1

+ 39 0382 610924


info@bioair.it

and check the manufacturer's label for the code and the serial number of the cabinet.

LABELLING

On every cabinet is placed a label as per the below example of *Pic. 1* with the main technical data of the cabinet.

- Logo
- Name of manufacturer
- Address of manufacturer
- CE Mark (for European Market only)
- Model
- Code
- Voltage, Frequency, Power
- Warning (3)
- Serial number
- Year of manufacturing
- Standard conformity (1)
- Competent Body Mark (2)
- Factory use code (4)

Pic. 1

Do not remove this label.

For any enquiry relevant to service or spares parts, please always refer to the labelling of the unit and report those data to us.

MEANING OF SYMBOLS

GENERAL WARNING.
CAREFULLY READ THE MANUAL BEFORE USING.

WARNING FOR UV-C RADIATIONS HAZARD

PROVISIONS SIGN

MAIN GROUND TERMINAL BLOCK

Separated collection for electric and electronic components.
According to 2012/19/EU DIRECTIVE.

DON'T USE CHLORINE PRODUCTS

PURPOSE AND USE OF THE MANUAL

Instructions translated from the original.

In case of a disagreement between the translation and the original version of this manual or a notice or disclaimer, the original Italian version will prevail.

This manual is intended to explain the correct use of this cabinet with front opening, equipped with UV-C light for the execution of all operations relevant to the Polymerase Chain Reactions (PCR) for DNA preparation in the safest conditions for your product.

The manual consists of two parts:

1. Part 1 is specifically designed to instruct the user.
2. Part 2 is specifically designed to train service personnel only. Instructions listed in this part are not designed to user.

NOTE: *service personnel who will operates on this cabinet should be however authorized.*

This Manual shall never be parted from the cabinet in order to be quickly available for use, service, maintenance and routine control.

The instructions given in this Manual shall be carefully read and understood before switching this equipment on, before any service operations and before you operate it, in order to become familiar with its operation and function, and to avoid misuse.

While every care has been taken in the production of this manual, it is not designed to act as your procedure manual or act as a substitute for validation of your own procedures.

Nothing contained in this manual can be considered as a warranty, either expressed or implied, including, not in a restrictive way, the suitability warranty for any special purpose.

Nothing contained in this manual can be interpreted as a modification or confirmation of the terms of any purchase contract.

ALL RIGHTS ARE RESERVED ACCORDING TO THE INTERNATIONAL COPYRIGHT CONVENTIONS,

the reproduction of any part of this manual, in any form, is forbidden without the prior written authorization of *BioAir S.p.A* (hereinafter referred to as *BioAir*).

BioAir reserves the right to modify or update the content at any time without prior notice.

aura PCR

PART 1: USER

1. ASSESSMENT OF RISKS

AURA PCR is a UV-C rays passive cabinet for PCR samples preparation.

NOTE 1: notwithstanding conformity with the standard referred, during operating of the cabinet several residual risks remain and were assessed in order to take adequate measures to prevent the following hazards arising from normal use and foreseeable misuse:

- General
- Misuse
- Electric shock
- UV-C radiations

The analysis were done as per EN 61010-1:2010, 7.3.3 and Annex J.

The residual risks have been assessed as per the following Tab. 1:

Tab. 1 – Risk Category

CATEGORY	ASSESSMENT	DESCRIPTION
1	Broadly acceptable.	This Category fulfils the requirement for tolerable risk.
2	As low as reasonably practicable.	This Category does not automatically fulfil the requirement for tolerable risk. If possible, reduce these risk to Category 1. <i>If not possible, then the instructions should contain a description of the risk so that the Responsible Body can take appropriate steps to protect the safety of operators.</i>
3	Intolerable.	This category contains risks that are not tolerable risks.

NOTE 2: no residual risk of this cabinet belongs to Category 3.

NOTE 3: read the following warnings before operating this cabinet.
Additional warnings shall be shown at the specific paragraphs.

NOTE 4: other risks associated with transportation, installation, electrical connection, installation of accessories, as well as maintenance and/or service in general, are also described in Part 2 in the relative paragraphs.

THE MANUFACTURER DECLINES ALL CRIMINAL AND CIVIL LIABILITY FOR DEFECTS DUE TO MISUSE, UNAUTHORISED USE, TAMPERING OR ANY GENERAL USE OF THE EQUIPMENT THAT IS NOT PERMITTED OR THAT DOES NOT COMPLY WITH THE OPERATING INSTRUCTIONS CONTAINED IN THIS MANUAL.

WARNING!

1.1. Misuse hazard. Risk category: 1

- **Misuse or a use which doesn't comply with the instructions of this manual compromises the safety of the cabinet.**

Taken measures to reduce the risk:

- Warning marking (caution) is provided on the cabinet.
- Adequate information for user.

PRECAUTIONS TO BE TAKEN:

- *The cabinet shall be used by authorized and adequately trained personnel only.*
- *Read carefully this manual before using the cabinet.*
- *Keep this manual always and easily available.*
- *DO NOT allow unauthorized people or general public to approach the cabinet.*

WARNING!

1.2. UV-C radiations hazard. Risk category: 1

- **The cabinet is equipped with UV-C light.**
- **Direct exposure to UV-C radiations may cause damage to the eyes (photoconjunctivitis and photokeratitis) and skin (erythema and skin cancer).**
- **UV-C radiations generate ozone which may be an inhalation hazard.**
- **UV radiation causes degradation of certain materials that may be used in cabinet construction. The breakdown of some materials may result in the emission of toxic vapours.**

Taken measures to reduce the risk:

- The walls are made of tempered safety glass which contains 98% of UV-C radiations.
- The closure of the front access contains inside the cabinet at least 98% of UV-C radiations.
- A microswitch allows the lighting up of the UV-C light only when the front access is completely closed.
- The supplied UV-C lamp complies with IEC/EN 62471.
- Warning marking (caution) is provided on protective glass.
- Adequate information for user.

PRECAUTIONS TO BE TAKEN:

- *Prevent access to the place where the cabinet is installed to unauthorized people or general public.*
- *Authorize the use of this cabinet only to trained and informed personnel.*
- *For long and continuous exposure is recommended an additional non-transparent shielding.*
- *Use UV-C lamps recommended by BioAir only.*

WARNING!

1.3. UV-C radiations hazard. Risk category: 1

- The front window is made of glass and supports an impact energy equal to 4 J.

The test was performed according to EN 61010-1, paragraph 8.2.2, with smooth steel sphere of mass 500g ± 25g in a vertical drop of 80 cm on a horizontal surface.

- Minimum required impact energy value = 5 J
Lower values are acceptable up to 1 J, if additional protective measures are taken.
- **Breakage of the front window or walls compromises the protection against UV-C radiation, with risk for the safety of the operators.**

Taken measures to reduce the risk:

- Warning marking (caution) is provided on protective glass.
- Adequate information for user.

PRECAUTIONS TO BE TAKEN:

- *Prevent access to the place where the cabinet is installed to unauthorized people or general public.*
 - *Authorize the use of this cabinet only to trained and informed personnel.*
- In case of breakage of the glass during UV-C cycle protect immediately eyes and skin and moreover:**
- *keep people away from the cabinet;*
 - *switch immediately off UV-C function (or do not run it);*
 - *shut down the cabinet and prevent reusing; and*
 - *call servicing organization.*

2. DEFINITION

2.1 DEFINITION

Aura PCR is a passive UV-C cabinet designed to protect the operations of production of recombinant DNA by risks of cross contamination due to the presence of DNA in aerosol in the air.

After handling of material, the DNA contained in any aerosol sprayed inside the cabinet can be neutralized by means of UV-C radiation.

The cabinet has been designed to protect the product handled, and DO NOT provides at all any protection to the operator and the environment from any potentially microbiologically hazardous product handled therein.

Therefore DO NOT handle potentially biohazard material inside this cabinet.

THE MANUFACTURER DECLINES ALL CRIMINAL AND CIVIL LIABILITY FOR DEFECTS DUE TO MISUSE, UNAUTHORIZED USE, TAMPERING OR ANY GENERAL USE OF THE EQUIPMENT THAT IS NOT PERMITTED OR THAT DOES NOT COMPLY WITH THE OPERATING INSTRUCTIONS CONTAINED IN THIS MANUAL.

3. INSTALLATION

WARNING!

- Unpacking, positioning, assembling and connecting to the mains shall be carried out solely by specialized and authorized staff, following the instructions listed in the PART 2 of this manual (Service and maintenance).

3.1 CHOICE OF LOCATION

WARNING!

- Not suitable place may compromise the performance of the cabinet.

- *Install the cabinet in places with no air conditioning grilles and aisles, which are not subject to air drafts and without frequent foot traffic.*

3.2 POSITIONING

WARNING!

- The cabinet should be fitted on a support stand (not supplied).
- Use of this accessory is under total responsibility of the customer.

- *Place the cabinet on a stand or bench suitable.*
- *This operation should be carried out only by qualified service personnel.*

3.3 CONNECTION TO THE MAINS

NOTE: *this operation should be carried out only by authorized personnel.*

The cabinet is designed to operate in the range defined by the table of specifications (see 4.1).

Check the mains supply to match with the one declared on the data plate of the cabinet.

3.3.1 Mains socket

A near-based (max 2 m) grounding mains socket shall be available to connect the cabinet.

Grounding must comply with local regulations.

Check that the mains voltage matches what shown on the manufacturer's label on the cabinet.

WARNING!

- The mains plug of the cord is the means for disconnection and it should be visible and easy-accessible.

- *Connect the cabinet to a dedicated power line, with a near-based grounding socket (max 2 m), according to the length of the power cord (no extension is allowed).*
- *Neither objects nor furniture should be hidden the plug in order to can quick disconnect the cabinet in case of emergency.*

3.3.2 Mains supply cord

The mains supply cord is a detachable flexible cord with IEC 60320 C13 female mains connector and CEE 7/7 mains plug (French-German) moulded at either end (Pic. 2). It complies with IEC 60227 or IEC 60245 and it is approved by a Competent Body.

WARNING!

- Use of unsuitable mains supply cord compromises the safety of the cabinet.
- Replacement of the mains supply cord (i.e. due to damaged cord or mains plug not compatible with mains socket), should be done with an approved one only, with same specifications as per Tab. 2.

Tab. 2 – Mains supply cord specifications

Type of cable:	H05VV-F
Numbers of wires:	3
Section (mm ²):	1
Max length (m):	3
Cable conformity:	IEC 60799
Mains connector conformity:	IEC 60320-1
Mains plug conformity:	IEC 60884-1
Approvals:	European

Pic. 2

3.3.3 Mains appliance inlet and main fuses

The appliance inlet (Pic. 3) is placed on the left rear of the control panel and consists to a IEC 60320 C14 male connector with a double fuse holder, EMI filter and 2 poles switch.

The main fuses are replaceable from the outside, located in the inlet appliance.

A label close to the appliance inlet reports the values of the fuses.

Fuses rating are also available in the table of specifications (see 4.1).

Pic. 3

3.4 FINAL PROCEDURES

Once installation is over, carry out the final procedures before operating the cabinet:

- check the main switch be turned off;
- carefully clean the cabinet inside and outside following instructions reported at relevant paragraph (see 8.1);

Should any problem arise, please contact your nearest BioAir distributor.

4. TECHNICAL DATA

4.1 SPECIFICATIONS	
Marks of conformity:	
Reference Standard:	IEC 61010-1:2010 / EN 61010-1:2010 IEC 61326-1:2012 / EN 61236-1:2013
Electrical insulating/protection class [IEC 61140]:	I
Mains supply voltage:	220-230 V~ 50/60 Hz
Main fuses:	F2A L, 250V
LED lamps (W):	2 x 7
CRI (Color Rendering Index):	80
Lamp colour temperature (K):	4000
Average life at 90% yield (h):	40000
UV-C lamps (W):	3x 15 T8
Required power line (W):	100
Absorbed power (W):	100
Sustained impact maximum energy of the glass [EN 61010-1, clause 8.2.2] (J):	4
Window glass UVC radiations retention (%):	98
4.2 USE ENVIRONMENTAL CONDITIONS	
Use:	indoor
Altitude (m):	up to 2000
Temperature (°C):	from 10 to 35
Maximum relative humidity (%):	80 for temperatures up to 31 °C, decreasing linearly to 50 at 40 °C
Max mains supply voltage fluctuations (%):	up to ±10
Transient overvoltage category:	II
Pollution degree:	2
4.3 TRANSPORT AND STORAGE CONDITIONS	
Ambient temperature (°C):	from -5 to 45
Relative humidity (%):	up to 90
Atmospheric pressure (mbar):	from 800 to 1060
4.4 PESO E DIMENSIONI	
Weight (kg):	41
Overall L x D x H (mm):	650 x 545 x 730
Maximum front aperture L x H (mm) :	555 x 430
Working space L x D x H (mm):	550 x 470 x 570
4.5 MATERIALS	
Upper structure:	cold rolled steel, stove enamel coated RAL 7035
Central structure, front panel and service tray:	tempered safety glass - 6 mm
Working surface:	PE

4.6 PERFORMANCES	
Illuminance [EN 12469] (lux):	>700
Sound level [EN ISO 3744] (dB[A]):	NA
Surface power density of single UV-C lamp at 1 m ($\mu\text{W}/\text{cm}^2$):	49,0
UV-C spectral peak (nm):	253,7
UV-C lamp average life (h):	8000


5. OPERATING PRINCIPLE

The cabinet is equipped with UV-C light, suitable to neutralize the DNA present in its own workspace.

The front opening allows the operator to perform PCR procedures reducing the possible dispersion of aerosols containing DNA in the outside environment.

6. CABINET DESCRIPTION

1. Control box
2. Appliance inlet/main switch/main fuses/ EMI filter
3. Service hinged tray
4. Data plate
5. Front panel
6. Working surface
7. Mains pilot light
8. UV-C cycle timer
9. Microswitch for UV-C enablement
10. Holes or extra utilities
11. Hook for front panel opening
12. Knob for front panel closure


Pic. 4

6.1 CONTROL BOX

6.1.1 Front view

On the right side of the control panel front there are the following controls:

1. electronic timer for UV-C cycle, also acts as power on indicator
2. switch/selector lighting (optional)


Pic. 5

6.1.2 Rear view

On the right side of the control panel back there is the appliance inlet which the following controls:

3. mains inlet connector
4. main fuses box
5. main switch


Pic. 6

6.2 SAFETY UV-C DEVICE

A microswitch (Pic. 7) located on the lower left side of the frontal aperture, allows the UV-C cycle to start on only when the front opening is completely closed.

Pic. 7

6.3 UV-C LAMPS

The cabinet contains 3 lamps UV-C suitable to reclaim their inner zone from DNA left over from previous work sessions.

The disinfection must be performed after any work session.

The efficiency of this method depends on right exposure time of the cabinet work zone to UV-C radiations.

WARNING! UV-C radiations hazard.

- See ASSESSMENT OF RISKS, warning and provision 0

NOTE: *the fact that UV lamps are producing light is not directly related to their actual neutralising power on DNA aerosols and microorganisms. Their actual efficiency may be tested with proper devices.*

On request it is possible to perform a test of efficiency of radiation emitted by the UV-C: contact your BioAir distributor.

6.3.1 PROTECTION AGAINST UV-C RADIATIONS

The protection supplied by the cabinet windows to the user is enough to guarantee a high level of safety during light start, programming and for occasional and short standing in the visible or reflected area of the UV-C lamp due to work needs.

WARNING!

We do not recommend long stays or workstations that can be reached by direct or reflected UVC light.

PRECAUTIONS TO BE TAKEN:

- *If it is essential to remain in the same room during the operation of the UVC lamp, it is the customer's responsibility to carry out a risk analysis to take the appropriate measures.*

7. OPERATING INSTRUCTIONS

7.1 SWITCHING ON THE CABINET

1. Turn on the cabinet by the main switch (Pic. 6) located on the right side of the back view of the control box (Pic. 4).

The timer (1) of Pic. 5 comes on, generates two series of beeps showing the time set to 0. The LED to the left of the time remains off.

The lighting comes on if the front window is open or the switch (2) of Pic. 5 is in AUTO position.

The cabinet is ready to use.

7.2 SETTING OF TIMER

KEY	DESCRIPTION
MIN	Increasing of time (minutes)
SEC	Increasing of time (seconds)
START/STOP	Start of cycle
	Stop of cycle
	Alarm mute
MIN+SEC	Reset of time

Pic. 8

Maximum time setting is 99 minutes and 59 seconds.

The LED next the time displayed turns on at the countdown start.

The set time remains stored until it is reset or turning off of the cabinet.

7.3 SETTING OF LIGHT SWITCH

At the center of the control panel is located a selector (Pic. 10) identified with the beside label (Pic. 9), which allows the turning on of the light as per the following scheme:

Pic. 9

	MAN	AUTO
MEANING	MANUAL TURNING ON OF THE LIGHT	AUTOMATIC TURNING ON OF THE LIGHT
DESCRIPTION OF THE FUNCTION	At the end of the UV-C cycle, the lighting will remain off until you open the front glass panel.	At the end of the UV-C cycle, the lighting comes on automatically.

Pic. 10


7.4 ACTIVATION OF UV-C CYCLE

WARNING! UV-C radiations hazard.

- See ASSESSMENT OF RISKS, warning and provision 0

1. After turning on the cabinet (7.1), open the front access and hang the knob to the hook on the roof of the control panel.
2. Rotate the inner shelf upwards in vertical position, so that it does not interfere with the radiation emitted by the UV-C lamps (Pic. 11).
3. Close the front access by locking on the base, with the lower glass panel knob.
4. Set on the timer (Pic. 8) the cycle time required by the keys MIN and SEC (see 7.2).

Pic. 11

 **NOTE:** normally an exposure between 5 and 30 minutes is enough to obtain the neutralization of DNA in aerosol. For more details, refer to the existing literature on this topic.

5. Press once the key START/STOP.

UV-C light immediately comes on and it will start the countdown that will end the cycle at the end of the set time.

If the switch MAN-AUTO is on AUTO position, the lighting will go off.

At the end of the cycle the timer generates a series of beeps for 30 seconds, mutable by pressing the key START/STOP.

NOTE 1: lighting is not switchable on during UV-C cycle.

NOTE 2: the opening up of the front access will switch off the lamps UV-C and switch on the lighting, but it does not stop the set cycle.

7.5 INTERRUPTION OF UV-C CYCLE

Should it be necessary stop the UV-C cycle, proceed as follows:

1. Press the START / STOP.

The countdown stops showing the remaining time of the cycle.

To continue with the remaining time press the key START / STOP again.

NOTE: to cancel the cycle, you have to reset the set time by pressing keys MIN and SEC and do if necessary a new time setting.

7.6 SWITCHING OFF THE CABINET

1. Switch off the cabinet by the main switch (Pic. 6) located on the right side of the back view of the control box (Pic. 4).

WARNING!

- The means of disconnection of the cabinet from the mains supply is only the plug of the mains supply cord.

7.7 GENERAL GUIDELINES

NOTE: all the involved people to the use of this cabinet shall be informed on relating operating procedures and restrictions and shall be obliged to carefully follow the user instructions of this manual.

In order to achieve better comfort for operators, *BioAir* suggests the following recommendations:

7.7.1 Features of location

- Chair should assure fully support of the body;
- If the chair has adjustable low back support, match the chair backrest to support the lower part of the spine;
- If different operators may use the cabinet, don't assume the settings are properly set for everybody;
- Vary the leg position frequently and avoid placing objects that limit the legroom under the cabinet.

7.7.2 Working posture

- Keep forearms, wrists and hands aligned in a straight neutral position: avoid bending or angling the wrist;
- Working intensely for long time poses risks: get up from the chair periodically and take brief walks.

7.7.3 Working procedure

- Operate the cabinet according to what prescribed by its definition (see **Errore. L'origine riferimento non è stata trovata.** and 2.1).
- DO NOT use the working area for storage of material.
- Request an efficiency test report of the cabinet after installation, components replacement, any repositioning or dislocation, and in any case at least once a year.
- Carry out the work according to a specific programme in which working steps, needed materials and accessories (to be put in the cabinet) are listed.
- Put inside the cabinet only essential materials and equipment, before operating the cabinet, so that nothing would pass through the front aperture during operation.
- Wash hands and arms always with a germicidal soap before and after work.
- Wear suitable clothing, long sleeves gowns with elastic cuffs, and rubber gloves.
- Operate in the innermost part of the working area.
- DO NOT operate the cabinet if there are damages to the glass parts.
- Restrict activities inside the room where the cabinet is located, not to cause drafts.
- Clean and disinfect the inner surfaces after each work session (see *par. 8*).

8. CLEANSING AND ROUTINE MAINTENANCE

WARNING!

- The cleaning and maintenance should be performed in the absence of power supply.
- *Disconnect the power supply by means of the power cord plug.*

8.1 CLEANSING

8.1.1 Painted surfaces

- Use common detergents, like soap solutions.
- *DO NOT USE solvents or other paint aggressive products.*

8.1.2 Stainless steel surfaces

- Specific products are recommended for normal cleaning.
- The disinfection of the inner surfaces of the working space shall be carried out by using disinfectant liquids, following the relating instructions to prevent misuse.
Decontaminate by swabbing with a disinfectant chlorine-free product.

WARNING!

- Chlorine products damage stainless steel surfaces.
- *DO NOT use products containing chlorine for stainless steel surfaces.*

8.1.3 Glass surfaces

- Disinfection of the glass window should be carried out with isopropyl alcohol.

WARNING! See ASSESSMENT OF RISKS, warnings and provisions 1.2 and 1.3.

8.1.4 Base plate

- Use specific products for polyethylene.

8.1.5 UV and LED lamps

- Use a woolen cloth moistened with water or an ammonia light solution.
- *DO NOT use oily products or waxes.*

PART 2: SERVICE AND MAINTENANCE

WARNING!

- All the operations described in this part of the manual should be performed only by qualified technicians, properly trained authorized.
- The operator is not allowed in any of these operations.

BioAir declines all the responsibility for damage resulting from non-compliance with these provisions.

INFORMATION FOR TRAINED AND AUTHORIZED SERVICE STAFF ONLY

9. INSTALLATION

WARNING!

- Unpacking, positioning and connection to the mains shall be carried out solely by specialized staff, following the instructions listed in this manual.

NOTE: the operations performed during installation should comply with the law provisions of the safety at work.

9.1 HANDLING AND UNPACKING

WARNING!

- Most breakages occur while these steps.

The cabinet is shipped on a dedicated pallet with protective wrapping film and outer cardboard.

To handle and lift the cabinet is necessary to use a hydraulic fork-lift truck (trans-pallet) with lifting forks up to a height of about 1 meter. The length of the forks should be at least 80 cm.

Remove the packing material only when the pallet is safely placed on the floor.

Disassemble the cardboard case and remove the polyethylene foil which covers the cabinet.

9.2 CHOICE OF LOCATION

WARNING!

- Not suitable installation place may compromise the performance of the cabinet.

The installation place should meet the following provisions:

- possibility of limiting access and shield the area;
- an available near-based (max 2 m) grounding mains socket;
- air-drafts-free ;
- no air conditioning louvers close to the cabinet; and
- no subject to foot traffic.

9.3 POSITIONING

WARNING!

- The cabinet is to be placed on a stand or a bench: check that they are adequate to support the weight of the cabinet.
- The choice of stand or bench to use is under total responsibility of the customer.

WARNING!

- Lifting must be carried out with the aid of a mechanical means or by at least 2 persons (the cabinet net weight is 41 kg).

WARNING! Risk of harm to body parts.

- Do not catch the cabinet neither by the top, nor by the glass structure.

The cabinet shall be levelled.

INFORMATION FOR TRAINED AND AUTHORIZED SERVICE STAFF ONLY

9.4 CONNECION TO THE MAINS

The cabinet is supplied with a detachable power supply cord.

It is recommended to connect the cabinet to a dedicated power line.

Mains supply cord, mains socket and main fuses specifications are available in the *PART 1* of this manual from 3.3.1 to 3.3.3.

See also 4.1, in the *PART 1*, for complete information about specifications.

Make sure that the voltage indicated on the cabinet tag is consistent with the mains voltage.

Grounding must comply with local regulations.

9.5 PERFORMANCE TEST

WARNING!

- Make performance tests with adequately (precision and accuracy) and freshly (max 1 year) calibrated instruments.
- Calibration certificate of the instruments for tests shall report readings in the range of the measurement to carry out.

9.5.1 Check of UV-C emissions efficiency:

NOTE: *this test is optional and it is performed only at the specific request of the customer.*

Instruments: luxmeter/radiometer with UV-C probe.

9.5.2 Advised additional instruments:

- digital multimeter capable of measuring true RMS AC voltage;

INFORMATION FOR TRAINED AND AUTHORIZED SERVICE STAFF ONLY

10. TROUBLE SHOOTING

WARNING! Electric shock hazard.

- The means of disconnection of the cabinet from the mains supply is only the plug of the mains supply cord.

ATTENZIONE! UV-C radiations hazard.

- See ASSESSMENT OF RISKS, warnings and provision 1.2 and 1.3
- During a technical service (eg in the case of actions for failure to UV system power-on or breaking of the glass) the taken protection on the cabinet can not be guaranteed.

- *Use a screen (or wear protective clothing and goggles/masks adequate) to protect skin and eyes from direct radiation of the UV-C lamps.*
- *Make sure there are no unprotected people around the cabinet.*

Trouble shooting table

Problem	Suggested action
1. The timer doesn't come on. No lamps turn on.	<ul style="list-style-type: none"> • Check the main switch at the mains appliance inlet X1 • Check the mains socket • Check the main fuses F10 e F20 at the mains appliance inlet X1
2. The timer comes on. Lighting lamp doesn't turn on. UV-C lamps are switchable on.	<ul style="list-style-type: none"> • Check the switch S1 MAN/AUTO • Check the LED lamps • Check the power supply for the LED lamps • Check the electronic timer TEMP1
3. The timer comes on. Lighting lamp is switchable on. UV-C lamp are not switchable on.	<ul style="list-style-type: none"> • Check the front window is well closed. • Check microswitch S2 • Check the UV-C lamps UV1-UV2-UV3 • Check the electronic ballasts R1.UV e R2.UV for UV-C lamps • Check the electronic timer TEMP1

INFORMATION FOR TRAINED AND AUTHORIZED SERVICE STAFF ONLY

11. WIRING DIAGRAM

The wiring diagram will be provided to the service staff adequately trained by *BioAir*.

Contact BioAir to request the necessary training and technical documentation.

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