

up to 60ml/min flow rate

MODEL: EV045

- Variable Speed
- Manual Control, Reversible
- Fast Prime Switch

Compact cased pump to provide variable speed/flow rate with fast prime switch. Accepts three different tube sizes.

Controls: Direction switch
Fast prime push button
Rotary speed control
Potentiometer

Pumphead: M045 Polycarbonate

Rotor: 2 roller Polycarbonate
3 roller Stainless Steel

Rollers: Nylatron

Motor: Permanent Magnet

Power Supply: 110/230v 50/60Hz 20w

Tube: Autoprene
Silicone
Viton
Tygon



Weight: 1.7 kgs

Dimensions: W205 x H110 x D195mm

Enclosure: IP30. Chemical resistant
Epoxy Polyester Coating

FLOW RATES ml/min	
TUBE	FLOW RATE (BASED ON 2 ROLLER ROTOR)
1.6mm	0 - 12 ml/min
3.2mm	0 - 36 ml/min
4.0mm	0 - 60 ml/min

WOLF LABORATORIES LIMITED

Tel: +44(0)1759 301142 Fax: +44(0)1759 301143

Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk