

OVERHEAD STIRRER

BOECO OVERHEAD STIRRER OSD-20

- ▶ With digital display
- ▶ Robust, slimline, ergonomic design
- ▶ With constant power-drive
- ▶ Two speed ranges for universal use from 60-2.000 rpm
- ▶ Push-through agitator shafts (only when stationary)



BOECO OVERHEAD
STIRRER OSD-20

Technical details

Stirring quantity (H ₂ O):	20 L.
Max.viscosity:	10.000 mPas
Motor rating input/output:	70 / 35 W
Output at stirring shaft:	26 W
Max. ON-time	100%
Max. torque at chuck (100rpm)	150 Ncm
Speed range I (50 Hz)	60 - 500 rpm.
Speed range II (50 Hz)	240 - 2.000 rpm
Speed display:	digital
Chuck range:	0,5 - 10 mm
Diameter:	13 mm
Length of extension arm:	160 mm
Protection class:	IP 20
(acc. to DIN EN 60529)	
Dimensions in mm (WxDxH):	88 x 212 x 294
Net weight:	approx 3,1 kg

The OSD-20 conforms to the CE mark and all relevant EC standards and DIN norms and has been manufactured under ISO 9001:2000

Code

Description


BOE 8359300	BOECO overhead stirrer OSD-20, 220-240V, 50/60 Hz,
-------------	---

Accessories

BOE 3160100	Plate stand R 1826
BOE 2657700	Boss head clamp R 182
BOE 0741000	Propeller stirrer, R 1342 Stirrer ø 50 mm, Shaft ø 8 mm, Shaft length: 350 mm



R 1826



R 182



R 1342