

Mini Roll-In Incubator



- "Flush-to-Floor" design accomodates 20 position Cell Production Roller Apparatus
- Optional shelves accommodate Bellco Magnetic Stirrers & Rockers
- Forced air circulation facilitates even temperature distribution in chamber
- Accessory outlets in chamber provide power for supplemental equipment
- Control panel features main power switch, convenience outlet power switch, safety thermostat, and microprocessor-based digital temperature controller
- Temperature Range: + 9 to 50°C
- Temperature Uniformity: $\pm 1.0^{\circ}\text{C}$
- Accepts optional Over-Temperature Alarm and Temperature Recorder Modules
- Left hand opening door (standard)
- Special models available with ports, solid and right hand opening doors
- Exterior Dimensions: 104cm x 81cm x 130cm (40x32x51")
- Chamber Dimensions: 89cm x 66cm x 102cm (21 cu.Ft.; 35x26x40")
- Shipping Weight: Approx. 350 pounds

Stock No	Power (V/Hz/W)
7728-06005	<u>100-1155/50-60-1500</u>
7728-65220	<u>230/50/1500</u>
Optional Shelves	
7728-40035	SST Shelf
7728-40050	Ultra Vented HD SST Shelf