

## Roll-In Incubator



- "Flush-to-Floor" design accommodates 45 position standard or 52 position programmable CPRA
- Optional shelves accommodate Bellco Magnetic Stirrers & Rockers
- Forced air circulation facilitates even temperature distribution in chamber
- Accessory outlets in chamber provide power for supplemental equipment
- Control panel features main power switch, convenience outlet power switch
- safety thermostat, and microprocessor-based digital temperature controller
- Temperature range: 9°C above ambient to 50°C ( $\pm 2^\circ\text{C}$ )
- Accepts optional Over-Temperature Alarm and Temperature Recorder Modules
- Left hand opening door (standard)
- Special models available with ports, solid and right hand opening doors
- Exterior Dimensions: 104cm x 81cm x 221cm (41x32x87")
- Chamber Dimensions: 89cm x 66cm x 193cm (40 cu. Ft.; 35x26x76")
- Shipping Weight: Approx. 700 pounds

| Stock No   | Power (V/Hz/A)   | Description        |
|------------|------------------|--------------------|
| 7728-09100 | <u>100/50/15</u> | Single Window Door |
| 7728-09005 | <u>115/60/15</u> | Single Window Door |
| 7728-95220 | <u>230/50/10</u> | Single Window Door |
| 7728-59005 | <u>115/60/15</u> | Double Window Door |
| 7728-59522 | <u>230/50/10</u> | Double Window Door |

### Accessories

|            |   |
|------------|---|
| 7728-40035 | Stainless Steel Shelf*                        |
| 7728-40050 | UltraVented Heavy Duty Stainless Steel Shelf* |