

Programmable Top Drive Roller Apparatus with and w/o Battery Backup



- Brushless motor & controller for smooth operation and easy maintenance
- Omron Controller
- Alarm jack and RS485 network port (RS232 and 4-20mA also available)
- 6-3/8" center to center accommodates up to 120mm diameter glass or plastic roller bottles
- Easy-to-remove roller bar assembly (tool provided)
- Front Locking Casters and 4 bumpers
- Speed Range: 0.08 to 7.0 rpm
- Battery level indicator flashes during battery mode (battery backup models only)
- Minimum of 18 hours battery backup depending upon load and rpm
- Rotation switch for CW and CCW operation
- Detachable cord /fused power inlet module
- All Position Drive standard on all models
- Optional audible/visual motion alarm (factory install only)
- Dimensions:
53 Position: 76"H x 32"W x 26"D (193 x 81 x 66cm)
86 Position: 76"H x 48"W x 26"D (193 x 122 x 66cm)

Stock No.	Description	Volts/Hz/Amp	Battery Back-up	Replaces
7400-10115	53 Position	115/60/3.0	no	7635-11000
7400-10230	53 Position	230/50/2.0	no	7635-12000
7400-20115	53 Position	115/60/3.0	yes	7635-11115
7400-20230	53 Position	230/50/2.0	yes	7635-11230
7420-10115	86 Position	115/60/3.0	no	7635-11000
7420-10230	86 Position	230/50/2.0	no	7635-12000
7420-20115	86 Position	115/60/3.0	yes	7635-11115
7420-20230	86 Position	230/50/2.0	yes	7635-11230

Roll-In Incubator



- "Flush-to-Floor" design accommodates 45 position standard or 52 position programmable CPRA
- Optional shelves accommodate Bellco Magnetic Stirrers & Rockers
- Forced air circulation facilitates even temperature distribution in chamber
- Accessory outlets in chamber provide power for supplemental equipment
- Control panel features main power switch, convenience outlet power switch
- safety thermostat, and microprocessor-based digital temperature controller
- Temperature range: 9°C above ambient to 50°C ($\pm 2^\circ\text{C}$)
- Accepts optional Over-Temperature Alarm and Temperature Recorder Modules
- Left hand opening door (standard)
- Special models available with ports, solid and right hand opening doors
- Exterior Dimensions: 104cm x 81cm x 221cm (41x32x87")
- Chamber Dimensions: 89cm x 66cm x 193cm (40 cu. Ft.; 35x26x76")
- Shipping Weight: Approx. 700 pounds

Stock No	Power (V/Hz/A)	Description
7728-09100	<u>100/50/15</u>	Single Window Door
7728-09005	<u>115/60/15</u>	Single Window Door
7728-95220	<u>230/50/10</u>	Single Window Door
7728-59005	<u>115/60/15</u>	Double Window Door
7728-59522	<u>230/50/10</u>	Double Window Door

Accessories

7728-40035	Stainless Steel Shelf*
7728-40050	UltraVented Heavy Duty Stainless Steel Shelf*