

## Modular Roller Apparatus



- Modular construction— add a deck at a time
- to 4.0 RPM standard: faster and slower models available upon request
- Analog speed indicator
- Locking speed control knob
- Accepts all Bellco Roller Bottles and 110 mm and 120 mm diameter, glass and plastic\* bottles up to 51 cm body length

Stock No.	Description	Voltage	Bottle Positions	Dimensions W x D x H
7630-75000	Base Only	<u>115</u>	5	75 x 63 x 37cm (30x24x14.5")
7630-75011	Base and 1 Deck	<u>115</u>	10	75 x 63 x 55cm (30x24x21.6")
7630-75012	Base and 2 Decks	<u>115</u>	15	75 x 63 x 73cm (30x24x28.7")
7630-75013	Base and 3 Decks	<u>115</u>	20	75 x 63 x 91cm (30x24x35.8")
7630-75014	Base and 4 Decks	<u>115</u>	25	75 x 63 x 109cm (30x24x43")
7630-75015	Base and 5 Decks	<u>115</u>	30	75 x 63 x 127cm (30x24x50")
7630-75016	Base and 6 Decks	<u>115</u>	35	75 x 63 x 145cm (30x24x57")
7630-75017	Base and 7 Decks	<u>115</u>	40	75 x 63 x 163cm (30x24x64")
7630-75018	Base and 8 Decks	<u>115</u>	45	75 x 63 x 180cm (30x24x71")
7630-75220	Base Only	<u>230</u>	5	75 x 63 x 37cm (30x24x14.5")
7630-75221	Base and 1 Deck	<u>230</u>	10	75 x 63 x 55cm (30x24x21.6")
7630-75222	Base and 2 Decks	<u>230</u>	15	75 x 63 x 73cm (30x24x28.7")
7630-75223	Base and 3 Decks	<u>230</u>	20	75 x 63 x 91cm (30x24x35.8")
7630-75224	Base and 4 Decks	<u>230</u>	25	75 x 63 x 109cm (30x24x43")
7630-75225	Base and 5 Decks	<u>230</u>	30	75 x 63 x 127cm (30x24x50")
7630-75226	Base and 6 Decks	<u>230</u>	35	75 x 63 x 145cm (30x24x57")
7630-75227	Base and 7 Decks	<u>230</u>	40	75 x 63 x 163cm (30x24x64")
7630-75228	Base and 8 Decks	<u>230</u>	45	75 x 63 x 180cm (30x24x71")
7630-75010	Additional Deck		5	75 x 63 x 18cm (30x24x7")

\* 7630-50555 All Position Drive is recommended for use with plastic bottles (1 per deck required)