

Bellco Roller Drum with Digital Display



- All purpose roller drum for blood, tissue culture, cells and virus applications
- Digital Display has 0.55" LED
- Locking speed knob
- Adjustable foot provides 0-7° drum tilt
- 3 speed ranges available: 0.2 to 3.0; 0.3 to 12.0; 0.5 to 50.0rpm
- Preset speed range from .5 to 50rpm
- Belt & Pulleys easily removed and reset when changes in speed range are necessary
- Four drums available to accomodate standard tube sizes
- Dimensions (base only): 12W x 13D x 17.5"H

Stock No.	Description	V/Hz/Amp	# Tubes	Tube O.D.
7736-10115	Drive Unit Only	115/60/2	~ ~ ~	~ ~ ~
7736-10230	Drive Unit Only	230/50/1	~ ~ ~	~ ~ ~
7736-20164	Drum Only	~ ~ ~	164	16
7736-20351	Drum Only	~ ~ ~	351	16
7736-20025	Drum Only	~ ~ ~	156	25
7736-20050	Drum Only	~ ~ ~	128	30