

## Adjustable Hanging Bar Spinner Flasks



The Bellco Adjustable Hanging Bar Spinner Flask Complete is designed for the efficient growth of suspension cultures and growth of cells on microcarrier beads. These autoclavable flasks are manufactured from Type 1 Class A Borosilicate glass with a flat bottom (no dimple), ideal for bar magnetic stirring. The PTFE impeller assembly rides on glass bearing for friction free operation. The 25 & 50ml have black phenolic white rubber lined screw caps. The 100 to 3L have Polypropylene caps & pouring rings. Clearly defined graduations are on each flask. The Adjustable Hanging Bar Spinner Flask complete consist of; Flask, Impeller Assembly, Impeller Glass Shaft & Cap Assembly.

### Special Notes

Flask is supplied with sidearm caps

Impeller Assembly consists of impeller magnet & magnet holder

Impeller Shaft consists of glass rod ( add 2-3" to extend through cap)

Cap assembly consists of screw cap. liner & stopper for 25ml & 50ml or compression fitting for sizes 100ml to 3L

For additional side arm screw caps. see caps

For Magnetic Stirrers, see Bench Top Equipment

For Stainless Steel Shafts, see Stainless Steel Assemblies

Note for 25ml & 50ml: Due to the smaller size of these vessels, probes and sensors are not available to fit the smaller caps and side arms

Make sure that the vessel is vented so that pressure will not build up and cause vessel to explode.

### Ordering Information

Stock No	Size (ml)	Center Cap (mm)	Sidearm Cap (mm)	O.D. (mm)	Hgt. (mm)
1967-00025	25*	38	18	38	120
1967-00050	50*	51	18	55	110
1967-00100	100	70	32	69	135
1967-00250	250	70	32	85	145
1967-00500	500	100	45	110	165
1967-01000	1000	100	45	135	240
1967-03000	3000	100	45	165	290