

Pipette Stands

When the pipette is not in use, it should be stored in an upright position in order to avoid contamination from dirty benches. Biohit offers stands for all of its pipettes.

The stylish Biohit Linear Stand is designed for all Biohit mechanical and electronic pipettes, and especially for mLINE, Proline Plus, Proline mechanical and ePET electronic pipettes. The stand is also compatible with a wide range of other manufacturers' pipettes.

The compact Carousel stands are available in different models, both for storage only and for charging.



Linear stand

eLINE

ePET

Cat. No.	Item
725600	Biohit Carousel Stand for 6 mechanical pipettes
725610	Biohit Pipette Holder for one Proline Plus / mLINE pipette
721000	Carousel Stand for 5 Proline mechanical pipettes
725620	Linear Stand for all Biohit pipette models

Cat. No.	Item
730981	eLINE Charging Stand for one pipette*
730991	eLINE Charging Carousel for 4 pipettes*
510001	Proline Charging Stand for one pipette*
511601	Proline Charging Carousel for 4 pipettes*
710999	ePET Stand for one pipette

* Includes Universal AC-adaptor (incl. Euro, US/Jpn, UK and China plugs)



Proline electronic

Proline mechanical

All mechanical pipettes

Reagent Vessel

The autoclavable and durable reagent vessel is made from polypropylene and is chemically resistant to all common reagents.

Cat. No.	Qty/Unit	Capacity
783500	16	120 ml

