

BIOHIT

Innovating for Health


eLINE®



A single innovation makes all the difference

eLINE Pro and Lite Dispensers hand you the most convenient and ergonomic approach to multiple dispensing. The automated liquid handling operation offers accurate and repetitive dispensing also with viscous, infectious and high vapour pressure liquids. The unique TipGuide selects and displays the most suitable tip size which saves time and eliminates the possibility of an incorrect tip selection. Additionally, a smart feature included in the whole eLINE product range is the unique electronic tip ejection, which reduces the risk of RSI (Repetitive Strain Injury).





Save your time, let the TipGuide make the tip choice for you

- No time consuming volume calculations
- Just set the desired volume and number of dispensings, the TipGuide selects and displays the suitable tip size

eLINE Pro and Lite Dispensers offer

- TipGuide – automatic guide for correct tip selection
- Outstanding ergonomics – dispensing with minimum force
- Unique electronic feather light tip ejection
- Comprehensive range of liquid handling protocols (eLINE Pro) with easy programming
- Contamination-free dispensing with positive displacement principle tips
- Volume range from 1 µl to 50 ml
- 5 speeds for aspiration and dispensing
- Charging indicator
- Direct charging option enabling dispensing while charging
- Compatibility with eLINE charging stand and carousel
- Compatibility with other dispenser tip brands

Operating modes	Pro	Lite
Multiple Dispensing with the Tip Guide (d)	x	x
Multiple Dispensing (d)	x	
Pipetting (P)	x	
Diluting (dd)	x	
Sequential Dispensing (Sd)	x	
Automatic Multiple Dispensing (Ad)	x	
Multi-Aspirating (SA)	x	
Custom mode (CST)	x	
- GL (select service intervals)		
- Sr (select the desired tip range)		
- SET (returns default settings)		
5 speeds for aspirating and dispensing	x	x

eLINE Dispenser Ordering details

Cat. No.	Item
730701	Biohit eLINE Pro Dispenser
730801	Biohit eLINE Lite Dispenser
730981	Biohit eLINE Charging Stand, with AC-adaptor
730991	eLINE Charging Carousel for 4 pipettes, with AC-adaptor
731001	eLINE replacement battery



New member of the eLINE family!

Biohit has several worldwide patents and patents pending for its products.