

## Ergomate – ergonomic pipetting system


When pipetting, ISO 8655 recommends touching off samples on the wall of receiving vessels, especially with small volumes. The new Biohit Ergomate removes the need to hold the plate at any specific angle and allows for safe dispensing, resulting in easier, safer, and more precise pipetting.



Ergonomic angle for proper positioning of wrist and hand for comfortable pipetting.

Micro plate grids for quick identification and sample management

Double-sided. Reverse side featuring numbering from 1-12.



When reversed plate holder creates an ideal aerosol barrier for use in hoods or on bench top surfaces

Micro Plate Grids

Plate holder can be reversed and used for flat surface pipetting or as an incubation station with lid

Ideal for safe transporting of plates



NEW! See page 80 for details

UNIQUELY  
BIOHIT!



Cat. No.	Item
723100	Ergomate Starter Kit, includes 3 Ergomate Plate holders and grids
723101	Ergomate Micro Plate Grids double-sided pack of 10

