

CLINICUT 60 CRYOSTAT

More than just a cryostat, the Clinicut 60 offers a systematic approach to the preparation, cutting and storage of fresh tissue specimens. Its unique specimen quick freezer is not simply a 'cold-spot' inside the chamber, but an insulated work station for tissue preparation and also enables specimens to be prepared in various shaped moulds.

With the addition of the compact, durable and reliable 3020 microtome, this cryostat's thoughtful design has also made the onerous task of regular cleaning so much easier.

Available as a manual or as a motorised rewind version, rapid frozen sections really are made easy with the Clinicut 60.



- *Small 'footprint' - takes up less space in the laboratory*
- *Thermostatically controlled chamber temperature to -35°C with fully programmable automatic defrost and digital display*
- *Excellent, rapid specimen quick freezer to -45°C with specimen moulds and heat dissipater*
- *Simple to set up - easy to use*
- *Accepts a wide range of knives and disposable blades*
- *Cuts excellent sections and incorporates the Easi-Set anti-roll plate*