

Captair® store AVP 804 (Up to 100 liters)



THE AIR OF THE LABORATORY IS PURIFIED

The air released by the filter contains 100 times less pollutants than the acceptable values defined by official authorities.



STANDARD EQUIPMENT :

- **Electronic "fan failure" alarm** for the continuous monitoring system of the ventilation.
- **Sampling port** for testing the saturation level of the filter.



VERY SILENT TURBINE FAN

Only 49 dBA



GOOD VISION OF THE STORAGE

Thanks to lockable transparent UV resistant door.



EASY TO SET-UP CLOSE TO THE USERS WORKSTATION

Does not require any ducting.




SPILL TRAY FILLED WITH AN ABSORBING MAT


To contain potential spill from broken containers.

2 Models adapted to your specific storage

captair® store 804 AVP
for the storage of acids
or bases or solvents



captair® store 804 AVPS
for the storage of acid-solvents
or acid-bases or bases-solvents




External Dimensions: H2150 x W800 x D505mm


Captair® store AVP 804 (Up to 100 liters)

2 Models for a greater storage capacity

Captair® store 804 AVPD
for the storage of acids-solvents
or acid-bases or acid-solvents



Captair® store 804 AVPSD
for the storage of acids or bases
or solvents in different compartments



External Dimensions: H2150 x W1600 x D505mm