

Ceti Star



Specifications : Star 24 ED	
Optical head	- binocular, inclined at 45° - dioptic adjustment on one eyepiece
Eyepiece	10x/20mm Widefield
Magnification	20x and 40x with turret changer
Working distance	75 mm
Focusing device	- rack and pinion with safety stop - fixed on vertical column
Illumination	- episcopic: fixed under the optical head - diascope: built-in the rectangular base - 2 individual On/Off switches
Stage	- frosted and black/white contrast plate - 2 stage clamps
Mains power	220-240V/50-60Hz (euroconnector)
Converter	built-in low voltage (12V-electronic)
Bulbs	- episcopic: Tungsten 12V/10W - diascope: halogen 12V/10W
Supplied with	2 rubber eyecaps, dust cover and instruction manual
Packing	styrofoam in cardboard box
Packing dimensions	34x18x35cm (WxDxH). Gross weight: 4,0kg