

## Cold Light Source – 150W



- Fiber optic illuminator
- Main power: 220V AC 50 Hz (euro connector))
- Light intensity: regulation with electronic control 0-100%
- Colour temperature: 3200°K
- Heat filter
- Bulb: 21V 150W halogen reflector bulb (500 burning hours)
- Cooling system: axial ventilator
- Includes focussing lenses
- Packing dimensions: 38x26x21 cm ((WxDxH). Gross Weight: 3,6 Kg

### Spare Parts for Cold light source

7060.2101 Bulb halogen reflector 21V 150W