

Ceti Polaris Polarimeter



Polaris 9100.000	
Rotating polarimeter 360° with 180° left and 180° right rotation	
Double vernier scale	accuracy 0.05°
Monochromatic light source of 589 nm	
Mains power	220-240V/50-60Hz, on/off switch
Reactor for Sodium bulb	
Bulb	sodium - 20Watt
Supplied with	1 Sodium bulb, 2 fuses, 3 test tubes (100 - 200 and 220mm), screwdriver and instruction manual
Packing	two styrofoam shelves in a cardboard box
Packing dimensions	60x20x40cm (WxDxH)
Gross weight	8,0kg