

## Coda Tower



Designed for filtering out the air and gas contaminants in the laboratories.

Four stage filter system for optimum absorption and oxidation of a wide variety of gases and particulate contaminants.

- Removes up to 99.97% of all contaminants
- Highly efficient four stage filter system
- Air flow creates an effective circulated air path
- Low noise level
- No draft
- Simple and economical
- Stable and easy to move on castor wheels

### **Specifications**

- Height : 198cm
- Depth : 53cm
- Width : 66cm
- Weight : 82Kg

Effective for rooms up to 27 square metres