

-40°C

LOW TEMPERATURE FREEZERS STAINLESS STEEL + WHITE GALVANIZED

LTFE

- Digital display
- Visual alarm
- Low energy consumption
- External contact for remote alarm
- Adjustable grid dividers
- Pull-out defrost drain for easy defrosting
- Lock and castors
- 100 % HCFC/CFC free
- White exterior

LTFS

- Digital display
- Low energy consumption
- Lock and castors
- Integrated handle
- 100 % HCFC/CFC free
- Stainless steel exterior



All stainless

LTFS 600



LCD display



LTFE 410



Digital display

STRONG, EASY TO CLEAN FREEZER WITH EXTRA TEMPERATURE STABILITY

LTFS SERIES: This -40 degrees series in stainless steel is ideal for all the users who prioritize easy cleaning and the right size. Standard features are lock, baskets and an elegant digital display without possibility of remote monitoring. A strong construction of the freezer provides long lifetime of the product as well as the 100 mm insulation ensures temperature stability.

LTFE SERIES: LTFE models are high performance freezers with lock, digital display and contact for remote monitoring. Flexible grid dividers can be configured to suit your individual requirements. As many of our other products, LTFE series has low energy consumption, a low noise level and compact design.



	GROSS CAPACITY (L)		NET CAPACITY (L)		TEMP. RANGE (°C)		AMBIENT TEMP. (°C)		WEIGHT (KG)		DIMENSIONS EXTERIOR (W×D×H)			DIMENSIONS INTERIOR (W×D×H)			REFRIGERANT		POWER SUPPLY (V)		FREQUENCY (HZ)		POWER CONSUMPTION (KWH/24H)		INSULATION (MM)	
LTFS 600	368	319	-12/-44	25	77	1575×750×960	1503×729×839	R404	230	50	4,50	100														
LTFS 700	429	372	-12/-44	25	82	1775×750×960	1703×729×839	R404	230	50	6,23	100														
LTFE 110	130	130	-10/-45	30	43	725×655×865	520×450×650	R404	230	50	3,84	100														
LTFE 210	226	226	-10/-45	30	54	1055×655×865	850×450×650	R404	230	50	5,03	100														
LTFE 310	299	299	-10/-45	30	64	1305×655×865	1100×450×650	R404	230	50	6,95	100														
LTFE 410	358	358	-10/-45	30	67	1505×655×865	1300×450×650	R404	230	50	7,54	100														
LTFE 510	416	416	-10/-45	30	77	1705×655×865	1500×450×650	R404	230	50	8,33	100														

We reserve the right to change specifications without notice. Subject to confirmation and availability.