

# Cooler for Bioruptor™

## Technical Specification



Ref. : BioAcc-Cool

Temperature Range	0 °C /80 °C
Temperature Sensor	Pt-100
Control System	Programmable Microprocessor
Temperature Set and Display Sensitivity	0.1 ° C
Temperature variation at 37 °C	<±0.4 °C
Temperature Fluctuation	±0.1 °C
Timer	100 hours +Hold Position
Pump capacity	5 lt./min.
External circulation	Optional
Tank volume,liters	14
Power consumption	1250 W
Power supply	230 V ,50 Hz.
Internal material	Stainless Steel
External material	Epoxy-Polyester painted stainless steel
Internal Dimensions (WxDxH)mm	240x300x200
External Dimensions (WxDxH)mm	445x340x590
Packing Dimensions (WxDxH)mm	560x410x670
Net / Packed weight kg.	26/28