

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

ESTA Storage freezers - F series



Model	F28	F38	F48	F58	F68
width	750	1040	1305	1535	1770
depth	730	730	730	730	730
height	840	840	840	840	840
litres	203	313	413	499	588
cu.ft.	7.2	10.4	14.6	17.6	20.8
oC	-20	-20	-20	-20	-20

Locking handle

Defrost drain

Replaceable gaskets

Climate Class 4 (30 deg c ambient)

External thermometer and thermostat

Hammered aluminum internal liner