

Incubator Hood TH 15

Constant temperatures; also for high vessels (e.g. 1000 ml Erlenmeyer flasks)
Can be combined with all compact shakers or used as independent unit

Basic Equipment: Coated metal housing with 3 windows, the front door opens upwards
Safety cut-off of the shaker if combined with a Bühler Compact Shaker
Standard temperature controller 5°C above ambient to +50°C



Accessories
① Shelf TH 15
Description see page 20

No. 0052 126

Temperature range: 5°C above ambient to +50°C
: to +60°C option
Temperature accuracy: + / - 1% of the set value
Inner dimensions: 490 x 410 x 310 mm (w x d x h)

Technical Data
Temperature sensor: PT 100
Air circulation: 240 m³/h
Electrical supply: 230 V or 115 V, 50 / 60 Hz
please indicate in case of order

External Dimensions (w x d x h): 510 x 450 x 455 mm
Weight: 21 kg

TH 15 No. 6161 000

Option
Extended temperature range up to + 60°C No. 0052 105
Installation only possible in our works
Description see page 20

Accessories for Incubator Hood TH 15



① Shelf TH 15 - No. 0052 126

Incubator hood TH 30

Excellent system thanks to reproducible ambient temperatures and useful accessories
Can be combined with all universal shakers or used as independent unit

Basic Equipment: Coated metal housing with 3 windows, the front door opens upwards
Safety cut-off of UV light and stop of the shaker, when the front door is opened, if combined with a Bühler Universal Shaker
Standard temperature controller 5°C above ambient to +50°C



Accessories
① Shelf TH 30
② 2-storey top frame SM / TH
Description see page 20

No. 0052 123
No. 0052 117

Temperature range: 5°C above ambient to +50°C
: to +60°C option
cooling possibility available as option
Temperature accuracy: + / - 1% of the set value
Inner dimensions: 660 x 540 x 430 mm (w x d x h)
height sufficient for 2000 ml Erlenmeyer flasks
Lighting: standard light 15W
UV light 15W

Technical Data
Temperature sensor: PT 100
Air circulation: 240 m³/h
Electrical supply: 230 V or 115 V, 50 / 60 Hz
please indicate in case of order

External Dimensions (w x d x h): 680 x 610 x 560 mm
Weight: 33 kg

TH 30 No. 6162 000

Options
Extended temperature range up to +60°C No. 0052 105
Installation only possible in our works
Cooling device TH 30 No. 0051 581
Installation only possible in our works
Description see page 20

Accessories for Incubator Hood TH 30



① Shelf TH 30 - No. 0052 123



② 2-storey top frame SM / TH - No. 0052 117