

WOLF LABORATORIES

Colenso House, Deans Lane, Pocklington, York YO42 2PX

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk

EXQUISIT FRIDGES

Model	F45	F60	F62G	KB60	C65	KS110	KS1401 (KS15)	KS1401RV (KS17RV)	R290
Gross Volume (L)	45	62	62	56	66	113	142 (138)	142 (138)	290
Gross Volume (cu ft)	1.6	2.2	2.2	2.0	2.4	4.0	5.0 (5.0)	5.0 (5.0)	10.0
Width (mm)	578	575	606	515	560	470	550 (550)	550 (550)	600
Height (mm)	415	470	522	510	600	826	860 (850)	860 (850)	1450
Depth (mm)	421	477	525	510	595	491	600 (570)	600 (570)	600
Defrost	manual	manual	manual	manual	auto	manual	manual	auto	auto
Energy Class	absorption	absorption	absorption	A	B	B	B (A)	B (A)	A

