


GX Wolf Microscopes

product datasheet

L Series

Illumination Systems for Stereo Microscopes



Specifications

Model	Lamp	Length mm	Fibre Diameter mm	Fibre Connector Diam mm
L1	24V 150W			16
L2		550	5.6	16
L3		550	5.6	16
L4		800	8	16
RFL4	220V 9W			
RFL8	220V 8W			60
RFL10	220V 10W			70

The L Series of illumination systems for stereo microscopes ensure that your samples can be viewed in bright white light showing every detail.

The stabilised lamphouse provides adjustable light intensity to the fibre optic light guides.

The fluorescence ring lights are ideal for those applications where more diffuse illumination is needed to reduce reflection effects.

Wolf Laboratories Limited

Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

June 04