


GX Wolf Microscopes

XDS-1 Optional Accessories

1. WF Eyepieces: 16X, 20x
2. **Long Distance** plan achromatic phase contrast objectives (cover glass: 1.2mm)
PLL25X/0.4PHP ----- must also order a 25X annular diaphragm plate
PLL40X/0.6PHP ----- must also order a 40X annular diaphragm plate

REMARK: XDS-1 already has PLL10X objectives, long distance condenser and centering telescope in the standard configuration



3. **Ultra-Long** working distance condenser, working distance 70mm



Ultra-Long working distance condenser does **NOT** work with phase contrast

4. **Special-Long** working distance phase contrast condenser, working distance 50mm
(Sliding, include a 2 hole adjustable plate)
OR 3 hole adjustable plate

Special-Long Distance plan achromatic phase contrast objectives (cover glass: 1.2mm)

PLL10X/0.25PHP2

PLL25X/0.4PHP2

PLL40X/0.6PHP2



Objectives (Optional)

3 Hole adjustable plate (Optional)

Special-Long
Distance
Condenser
(Optional)

2 hole
adjustable
plate
(Optional)

NOTE:

The difference between XDS-2 and XDS-1

	XDS-1	XDS-2
Long Distance plan achromatic phase contrast	Standard	Optional
Special-Long Distance plan achromatic phase contrast	Optional	Standard

Wolf Laboratories Limited

Colenso House

1 Deans Lane

Pocklington

York YO42 2PX

Tel 01759 301142

Fax 01759 301143

Email sales@wolflabs.co.uk

Web www.wolflabs.co.uk