

GX Wolf Microscopes

product datasheet

DCE-1

Video Eyepiece
Camera



DCE-1 on a Microscope



DCE-1 on a Macroscope

Features and Specifications

The DCE-1 is a versatile eyepiece video camera with a video connection for a PC or VCR. It provides a sharp, live image of the subject mounted under your microscope.

The eyepiece will fit most GX Microscopes simply by replacing one of the eyepieces. This not only saves the expense of a trinocular head but also no c-mount or bayonet mount camera adaptor is needed.

This popular camera is widely used in schools, colleges, universities, industry and research.

Specifications

CCD element:	640(H)x480(V)lines
CCD type:	CMOS
Image Type	Colour
Min illumination:	3lux
Interface:	Video output
Requires:	PC with Video Input
Software:	Not included requires video capture software supplied by video capture hardware vendor

Wolf Laboratories Limited
Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

June 04