


# GX Wolf Microscopes

product datasheet

## Eyepieces and Reticules



Huygens Eyepieces			
Multiple	Focal distance(mm)	Exit pupil distance (mm)	Field Range (mm)
5X	50	15.2	f 18
6X	41.7	6.6	f 16
10X	25	6	f 14
15X	16.6	6.11	f 11
Eyepieces with built-in 10mm scale reticules 23.2mm			
Multiple	Focal distance(mm)	Eyepoint Height (mm)	Field range (mm)
10X	25	16.1	f 17
Widefield Eyepieces 23.2mm			
Multiple	Focal distance(mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	19.8	f 20
10X	25	19.8	f 18
10X	25	19.2	f 18
10X	25	16.1	f 18
10X(point)	25	16.1	f 18
10X	25	15.8	f 18
10X	25	16.1	f 16
10X	25	18.5	f 18
12.5X	20	16.8	f 14
15X	16.6	15	f 11
16X	15.6	14.7	f 11
20X	12.5	15.3	f 11
20X	12.5	8.8	f 11
Eyepieces with built-in 10mm scale reticules 30mm			
Multiple	Focal Distance (mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	17.9	f 17
Wide Field Eyepieces 30mm			
Multiple	Focal Distance (mm)	Eyepoint Height (mm)	Field Range (mm)
10X	25	14.9	f 20
10X	25	17.9	f 20
15X	16.6	9.3	f 15
20X	12.5	11.1	f 12
25X	10	8.3	f 9

Also available 5mm scale and crosshair reticules.

Wolf Laboratories Limited  
Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

June 05