

# GX Wolf Microscopes

product datasheet

## GXHM Series Pocket Microscopes

A real breakthrough in affordable pocket microscopy.


Consisting of a durable ABS plastic body and glass optics, these the pocket sized microscopes produce remarkably clear images.

The models, magnification ranges and dimensions are:

GXHM30	30X	392x335x360mm
GXHMD60100	60X-100X	392x335x360mm
GXHMD60	60X	392x335x360mm

You can focus on the specimen by placing the microscope onto the specimen and then turn the focus control wheel. Built-in illumination is supplied in all models including a daylight filter (batteries are not included).

A specimen stage/stand with slide clips is available which is ideal for slide viewing.



GXHM30

GXHMD6010

GXHMD60  
with option  
stand

Wolf Laboratories Limited  
Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

Nov 05