



# GX Wolf Microscopes

product datasheet



**GXMZC-995CP**

## Hair & Skin Digital Imaging System

An extremely useful, modern tool for the examination and recording of images of skin and hair at a range of magnifications.

The system presents a live image direct to a PC via a USB link of the subject under examination. The built-in LED illumination provides a crisp bright image and the user can select from a range of magnifications by attaching the relevant lens (shown on right).

Magnifications:

X1  
X50  
X200  
X700

The imaging device is lightweight, simple to use and easily cleaned. It needs no separate power as it draws a current direct from the computer's USB port.

Wolf Laboratories Limited  
Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

Aug 07