

# GX Wolf Microscopes

## product datasheet

### Macroscope 8x30

Portable  
Measuring  
Microscope



The GXM Macroscope is an extremely useful field tool. It consists of three optical systems in one. Complete, the Macroscope is a x25 hand microscope with micrometer scale that measures both in inches and millimetres.

The scale and the microscope focus independently, and gives a field of view of 9.5mm, with the clear skirt allowing excellent illumination. The macroscope divides into two components - a x4 hand lens (field of view 70mm) and a telescope with variable focus which gives excellent clarity and is ideal as a spotting scope.

The complete unit weighs 290gms and measures 210mm in height.

Parameter	GXMMACROSCOPE8X30
Total Magnification	X25
Scale	Inches and mm focusable
Field of view	9.5mm
Hand Lens Mag	4X
Hand Lens Field of View	70mm
Telescope	Variable Focus
Weight	0.29kg
Height	210mm
Eyepiece Camera	optional

#### **Wolf Laboratories Limited**

Colenso House  
1 Deans Lane  
Pocklington  
York YO42 2PX

Tel 01759 301142  
Fax 01759 301143  
Email [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Web [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Due to our policy of continuous development specifications may change without notice

May 06