

GX Wolf Microscopes

product datasheet

Post Stands

Stand and illumination systems for SZ, PZ and MZ series microscopes



GX Microscopes Post Stands for Z series Microscopes	MA211101	MA211103	MA212111	MA221102	MA212112
Max Weight	4kg	4kg	4kg	4kg	4kg
Max Size of Base	320x305x26mm	320x305x26mm	320x305x26mm	320x305x26mm	320x305x26mm
Maximum Distance from vertical post to the front of base	273mm	273mm	273mm	273mm	273mm
Diameter of Post	32mm	32mm	32mm	32mm	32mm
Height of Post	305mm	305mm	305mm	305mm	305mm
Type and Diameter of Microscope holder ring/focus mount	For SZ series 76mm	For MZ series 39mm	For SZ series 76mm	For PZ series	For SZ series 76mm
Focus Range of Focus Mechanism	25mm	25mm	25mm	25mm	25mm
Illumination	None – requires accessory illumination	None – requires accessory illumination	Top and base illumination	None – requires accessory illumination	Top and base LED illumination
Power	None	None	85-265V 50/60Hz	85-265V 50/60Hz	85-265V 50/60Hz

Wolf Laboratories Limited

Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

June 04