

GX Wolf Microscopes

product datasheet

SX-C Prism Microscope



Extremely robust microscope for use in schools, colleges and in field studies. This simple to use instrument uses a tough prism to reflect the light through the slide into the achromatic objectives.

Supported on a metal stand the steel body houses the objective and eyepiece, securely locked into position.

The stand allows inclination of the viewing tube for comfort.

There are a range of models available giving a range of magnifications of up to x100.

Models

- SX-C-20 – x5eyepiece with x4 objective
- SX-C-40 – x10 eyepiece with x4 objective
- SX-C-50 – x5eyepiece with x10 objective
- SX-C-100 – x10eyepiece with x10 objective

Wolf Laboratories Limited

Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

July 04