


GX Wolf Microscopes

product datasheet



SX Series

Modern classroom
biological
microscopes



Specifications			
Model:		SX-A	SX-AL
Head	Monocular, 45deg inclination	STD	STD
Eyepiece	10X Wide-Field	STD	STD
Focusing	Smooth, stage-moving rack and pinion	STD	STD
Objectives	x4, x10, x40 achromatic, parfocal	STD	STD
Nosepiece	Triple position, rotating	STD	STD
Stage	Slide mounting clips with smooth stage surface	STD	STD
Illumination	Transmitted light 15W 110V or 220V		STD
	Concave Mirror	STD	

Features and Specifications

- Equipped with achromatic objectives and wide-field eyepieces, providing a clear, bright view.
- Monocular eyepiece tube with an inclination of 45deg for comfortable viewing and simple operation.
- Illuminated with 110/220V lamp (SX-AL) or concave mirror (SX-A)
- Objectives 4X, 10X and 40X mounted in rotating triple position nosepiece.
- Anti-Tamper features such as locked-in eyepiece.
- Modern, stage-moving focus mechanism
- Modern clean design.
- Great value for money.
- Tough metal construction with durable, easy to clean paint finish

Wolf Laboratories Limited

Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142

Fax 01759 301143

Email sales@wolflabs.co.uk

Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

June 04