

GX Wolf Microscopes

product datasheet

Slide Boxes

Quality, Durable Plastic
Slide Boxes

Features and Specifications

Wolf Laboratories range of high quality, incredibly good value slide boxes. Constructed from durable, shock resistant plastic, these containers have clip-lock mechanisms to ensure a long life for your precious specimen slides.

All except the 6 slide box include a slide listing sheet.

These containers are suitable for sending slides by post but it is a good idea to reduce rattling of the slides by applying tough tape over the tops of the slides before shipment. Also check that the foam is still in a good condition before sending by post.

All slide boxes are supplied in white. Other colours are available for large volume customers.



Wolf Laboratories Limited

Colenso House
1 Deans Lane
Pocklington
York YO42 2PX

Tel 01759 301142
Fax 01759 301143
Email sales@wolflabs.co.uk
Web www.wolflabs.co.uk

Due to our policy of continuous development specifications may change without notice

Jan 06